



## Artisan Technology Group is your source for quality new and certified-used/pre-owned equipment

- FAST SHIPPING AND DELIVERY
- TENS OF THOUSANDS OF IN-STOCK ITEMS
- EQUIPMENT DEMOS
- HUNDREDS OF MANUFACTURERS SUPPORTED
- LEASING/MONTHLY RENTALS
- ITAR CERTIFIED SECURE ASSET SOLUTIONS

### SERVICE CENTER REPAIRS

Experienced engineers and technicians on staff at our full-service, in-house repair center

### *InstraView*<sup>SM</sup> REMOTE INSPECTION

Remotely inspect equipment before purchasing with our interactive website at [www.instraview.com](http://www.instraview.com) ↗

### WE BUY USED EQUIPMENT

Sell your excess, underutilized, and idle used equipment. We also offer credit for buy-backs and trade-ins. [www.artisanng.com/WeBuyEquipment](http://www.artisanng.com/WeBuyEquipment) ↗

### LOOKING FOR MORE INFORMATION?

Visit us on the web at [www.artisanng.com](http://www.artisanng.com) ↗ for more information on price quotations, drivers, technical specifications, manuals, and documentation

**Contact us:** (888) 88-SOURCE | [sales@artisanng.com](mailto:sales@artisanng.com) | [www.artisanng.com](http://www.artisanng.com)

The choice of the critics

## DOG VENTILATOR

Cat. No. 5 0 2 5

New look Ventilator, featuring:-

- Compact valve block with negligible "dead spaces"
- Original variable stroke linkage, without sliding components
- No protruding elements: recessed controls and hook-up connections to the animal
- Attractive trolley which moves swiftly, thanks to its heavy-duty casters



### GENERAL

The **5025 Dog Ventilator** is a positive pressure pump (according to Starling's ventilation method), designed for use with dogs.

It basically consists of an asynchronous motor operated by a solid state speed control, coupled to a 700 (or 350) ml cylinder/piston assembly via an original variable stroke linkage (see picture) which features:-

- 1) The piston almost touches the cylinder end with each

stroke, regardless of the preset volume. This insures that all air taken into the pump is expelled with each stroke.

- 2) The stroke volume is adjustable by means of a thumb-wheel while the pump is either running or at standstill.

Volume is clearly indicated on a **stationary dial** and not on a mobile slotted link or by graduation marks on the cylinder, both systems leading to uncomfortable volume setting, the latter adjustable with the pump in motion only.

- 3) The mechanism does not incorporate any slotted link as all other ventilators do.

The reciprocating motion is generated, adjusted and transmitted to the piston by rods and articulated joints only. The lack of sliding friction leads to:-

- a) practically no wear
- b) no backlash and hence silent operation and exact stroke reproducibility.

A crank-link combination operates a slide valve, located perpendicularly to the cylinder axis, ensuring very short cylinder/valve connections and hence minimum dead space.

### Hook-up to Animal

A "Y" connection (not supplied) should be inserted between the connection to the animal (cannula etc.) and the two tubes coming from the pump. Dead space is minimized if the "Y" is as close to the animal as

Four ports (*Intake, To Animal, From Animal* and *Exhaust*) allow flexibility in air channelling. The input may be room air or any non-explosive gas mixture. The exhaust air may be partially or totally recycled or collected for analysis.

### Ventilator Controls

The speed control knob adjusts the geared motor to the desired speed, which is indicated on the 2-digit LED display labelled STROKES PER MINUTE.

The digital display assures an excellent reproducibility: you adjust the actual speed of the geared motor against a displayed number, in lieu of trying to overlap a pointer and a scale mark.

The stroke volume may be adjusted via the STROKE VOLUME thumb-wheel on the front panel, and read on the scale in the top panel.

**SPECIFICATIONS**

Rate: 10 to 50 strokes for minute  
Rate Read-out : digital display  
Stroke Volume : 0 to 700 or 0 to 350 de-  
pen-ding on cylinder/piston in-  
stalled  
Stroke Volume Scale : 0-700 ml  
Stroke Volume  
Reproducibility :  $\pm 2\%$   
Power Requirements : 115 or 230v, 50/60 Hz  
300 VA max.  
Operating Temperature : 10°C to 40°C  
Overall Dimensions : 45 x 38 x 83 (h) cm  
Net weight : Kg 30  
Shipping Weight : Kg 55 approx.

**ORDERING INFORMATION**

**5025** Dog Ventilator, complete with standard accessor-  
ies and 350 or 700 ml cylinder/piston assembly  
**5026** 700 ml cylinder/piston assembly, complete with  
gudgeon  
**5027** 350 ml cylinder/piston assembly, complete with  
gudgeon  
**5028** Set of lip-seal rings for 350 ml piston  
**5029** Set of lip-seal rings for 700 ml piston  
**5030\*** Valve Block, complete with spindle, gudgeon and 2  
open washers  
**5031\*** Perspex Hood  
**5032\*** Oiler (Oil for Bearings)  
**5033\*** Oiler (Oil for Cylinder)  
**7033\*** Grease Tube (Lithium Grease for slide valve)  
**5034\*** Set of 5 hex. wrenches (2, 2.5, 3, 4, 5)  
**5035\*** Instruction Manual  
**5037\*** M17 wrench (to tighten the trolley casters)

Set of two fuses (T 2.5 A) for either 115 or 220 V operation

\* Standard accessories, supplied with the instrument



## Artisan Technology Group is your source for quality new and certified-used/pre-owned equipment

- FAST SHIPPING AND DELIVERY
- TENS OF THOUSANDS OF IN-STOCK ITEMS
- EQUIPMENT DEMOS
- HUNDREDS OF MANUFACTURERS SUPPORTED
- LEASING/MONTHLY RENTALS
- ITAR CERTIFIED SECURE ASSET SOLUTIONS

### SERVICE CENTER REPAIRS

Experienced engineers and technicians on staff at our full-service, in-house repair center

### *InstraView*<sup>SM</sup> REMOTE INSPECTION

Remotely inspect equipment before purchasing with our interactive website at [www.instraview.com](http://www.instraview.com) ↗

### WE BUY USED EQUIPMENT

Sell your excess, underutilized, and idle used equipment. We also offer credit for buy-backs and trade-ins. [www.artisanng.com/WeBuyEquipment](http://www.artisanng.com/WeBuyEquipment) ↗

### LOOKING FOR MORE INFORMATION?

Visit us on the web at [www.artisanng.com](http://www.artisanng.com) ↗ for more information on price quotations, drivers, technical specifications, manuals, and documentation

**Contact us:** (888) 88-SOURCE | [sales@artisanng.com](mailto:sales@artisanng.com) | [www.artisanng.com](http://www.artisanng.com)