



Artisan Technology Group is your source for quality new and certified-used/pre-owned equipment

- FAST SHIPPING AND DELIVERY
- TENS OF THOUSANDS OF IN-STOCK ITEMS
- EQUIPMENT DEMOS
- HUNDREDS OF MANUFACTURERS SUPPORTED
- LEASING/MONTHLY RENTALS
- ITAR CERTIFIED SECURE ASSET SOLUTIONS

SERVICE CENTER REPAIRS

Experienced engineers and technicians on staff at our full-service, in-house repair center

*InstraView*SM REMOTE INSPECTION

Remotely inspect equipment before purchasing with our interactive website at www.instraview.com ↗

WE BUY USED EQUIPMENT

Sell your excess, underutilized, and idle used equipment. We also offer credit for buy-backs and trade-ins. www.artisanng.com/WeBuyEquipment ↗

LOOKING FOR MORE INFORMATION?

Visit us on the web at www.artisanng.com ↗ for more information on price quotations, drivers, technical specifications, manuals, and documentation

Contact us: (888) 88-SOURCE | sales@artisanng.com | www.artisanng.com

Lumina Hollow Cathode Lamp (2in) - Single Element

PerkinElmer's expertise in the design and manufacturing of light sources for atomic absorption spectrometers traces back to our original hollow cathode lamps developed nearly 40 years ago. The use of special light sources and careful wavelength selection permits the specific determination of one element in the presence of others. The quality of the light source is of primary importance in obtaining the best linearity, sensitivity, and precision in atomic absorption.

The use of special light sources and careful wavelength selection permits the specific determination of one element in the presence of others. The quality of the light source is of primary importance in obtaining the best linearity, sensitivity, and precision in atomic absorption.



For the Widest Selection of Elements

[View Larger Image](#)

The basic components of an HCL are a sealed glass cylinder filled with argon or neon, an anode, a cathode (a hollow cylinder of the element whose spectrum is to be produced), and either a glass or quartz end window. HCLs emit light by ionizing the fill gas using an electrical potential applied between the anode and the cathode. The positively charged ions collide with the negatively charged cathode and dislodge individual metal atoms in a process known as "sputtering." These "sputtered" metal atoms are excited through impact with fill gas ions, and light of the specific wavelengths for that element is emitted when the atom "decays" from the excited state back to the more stable state.

When you select a PerkinElmer lamp for your laboratory, you'll receive our pledge of satisfaction:

- **Long Life** – PerkinElmer's larger lamp volume results in longer lamp lifetime.
- **Automation** – Our coded lamps allow automatic setup of the correct wavelength, slit setting, and lamp current.
- **Complete Compatibility** – Lumina hollow cathode lamps can be used with PerkinElmer's entire range of atomic absorption spectrometers. Cables are available that allow Lumina lamps to be recognized by our instruments that recognize lamp coding.

Lumina single-element coded hollow cathode lamps are available in two series: with and without an Intensitron adapter cable. Series N30501XX Lumina lamps, shown below, do not include an Intensitron adapter cable and are designed to be used with PerkinElmer's AAnalyst series of instruments.



Artisan Technology Group is your source for quality new and certified-used/pre-owned equipment

- FAST SHIPPING AND DELIVERY
- TENS OF THOUSANDS OF IN-STOCK ITEMS
- EQUIPMENT DEMOS
- HUNDREDS OF MANUFACTURERS SUPPORTED
- LEASING/MONTHLY RENTALS
- ITAR CERTIFIED SECURE ASSET SOLUTIONS

SERVICE CENTER REPAIRS

Experienced engineers and technicians on staff at our full-service, in-house repair center

*InstraView*SM REMOTE INSPECTION

Remotely inspect equipment before purchasing with our interactive website at www.instraview.com ↗

WE BUY USED EQUIPMENT

Sell your excess, underutilized, and idle used equipment. We also offer credit for buy-backs and trade-ins. www.artisanng.com/WeBuyEquipment ↗

LOOKING FOR MORE INFORMATION?

Visit us on the web at www.artisanng.com ↗ for more information on price quotations, drivers, technical specifications, manuals, and documentation

Contact us: (888) 88-SOURCE | sales@artisanng.com | www.artisanng.com