Pacific Scientific ABS-2

Automatic Bottle Sampler

Limited Availability
Used and in Excellent Condition

Open Web Page

<https://www.artisantg.com/85816-1>

All trademarks, brandnames, and brands appearing herein are the property of their respective owners.

Your **definitive** source
for quality pre-owned
equipment.

Artisan Technology Group

(217) 352-9330 | sales@artisantg.com | artisantg.com

- Critical and expedited services
- In stock / Ready-to-ship

- We buy your excess, underutilized, and idle equipment
- Full-service, independent repair center


Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

DOC026.97.80280

ABS-2

02/2012, Edition 1

User Manual
Manuel d'utilisation

English.....	3
Français.....	20

Specifications

Specifications are subject to change without notice.

Dimensions (W x L x H)	32.4 x 45.7 x 68.6 cm (12.75 x 18 x 27 in.)
Weight	20 kg (44 lb)
Installation category	II
Pollution degree	2
Protection Class	I
Sample flow rate	1–200 mL/min (user adjustable, fluid-viscosity dependent)
Sample volume	5–100 mL (user adjustable, 5 mL increments)
Dead (Tare) volume	16 mL
Sample accuracy	±5% of selected volume
Repeatability	±1%
Fluid compatibility	Fluids compatible with stainless steel, glass, Teflon™ and Tygon™
Operating pressure range	Minimum: -20 PSI (-34 kPa) recommended Maximum: 50 PSI (344.7 kPa)
Temperature	Operating: 10 to 50 °C (50 to 122 °F)
	Non-operating: -40 to +65 °C (-40 to +149 °F)
Maximum altitude	Operating: 3000 m (9840 ft)
	Non-operating: 15240 m (50000 ft)
Relative humidity	10 to 80% RH (non-condensing)
Host I/O	RS232 (DB-25 female)
Certifications	CE

Power requirements	
Instrument	100–230 VAC, 50/60 Hz, 0.5 A
With optional pump	115 VAC, 60 Hz, 4.5 A
	100 VAC, 50 Hz, 4.5 A
	230 VAC, 50 Hz, 1.9 A

General information

In no event will the manufacturer be liable for direct, indirect, special, incidental or consequential damages resulting from any defect or omission in this manual. The manufacturer reserves the right to make changes in this manual and the products it describes at any time, without notice or obligation. Revised editions are found on the manufacturer's website.

Safety information

NOTICE

The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

Please read this entire manual before unpacking, setting up or operating this equipment. Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

Make sure that the protection provided by this equipment is not impaired. Do not use or install this equipment in any manner other than that specified in this manual.

Use of hazard information

⚠ DANGER
Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.
⚠ WARNING
Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.
⚠ CAUTION
Indicates a potentially hazardous situation that may result in minor or moderate injury.
NOTICE
Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

Precautionary labels

	This is the safety alert symbol. Obey all safety messages that follow this symbol to avoid potential injury. If on the instrument, refer to the instruction manual for operation or safety information.
	This symbol indicates that a risk of electrical shock and/or electrocution exists.
	This symbol indicates the presence of devices sensitive to Electrostatic Discharge (ESD) and indicated that care must be taken to prevent damage with the equipment.
	This symbol, when noted on the product, identifies the location of a fuse or current limiting device.

	This symbol identifies a risk of chemical harm and indicates that only individuals qualified and trained to work with chemicals should handle chemicals or perform maintenance on chemical delivery systems associated with the equipment.
	This symbol indicates the need for protective eye wear.
	Electrical equipment marked with this symbol may not be disposed of in European public disposal systems after 12 August of 2005. In conformity with European local and national regulations (EU Directive 2002/98/EC), European electrical equipment users must now return old or end-of-life equipment to the Producer for disposal at no charge to the user. <i>Note: For return for recycling, please contact the equipment producer or supplier for instructions on how to return end-of-life equipment, producer-supplied electrical accessories, and all auxiliary items for proper disposal.</i>

Certification

Canadian Radio Interference-Causing Equipment Regulation, IECS-003, Class A:

Supporting test records reside with the manufacturer.

This Class A digital apparatus meets all requirements of the Canadian Interference-Causing Equipment Regulations.

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FCC Part 15, Class "A" Limits

Supporting test records reside with the manufacturer. The device complies with Part 15 of the FCC Rules. Operation is subject to the following conditions:

1. The equipment may not cause harmful interference.
2. The equipment must accept any interference received, including interference that may cause undesired operation.

Changes or modifications to this equipment not expressly approved by the party responsible for compliance could void the user's authority to

operate the equipment. This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference, in which case the user will be required to correct the interference at their expense. The following techniques can be used to reduce interference problems:

1. Disconnect the equipment from its power source to verify that it is or is not the source of the interference.
2. If the equipment is connected to the same outlet as the device experiencing interference, connect the equipment to a different outlet.
3. Move the equipment away from the device receiving the interference.
4. Reposition the receiving antenna for the device receiving the interference.
5. Try combinations of the above.

Product overview

The instrument sends and measures a specified volume of liquid sample through a sensor. The instrument sends a start signal and a stop signal to a counter as the sample moves through the sensor.


This instrument is part of a batch fluid analysis system. For details on how to include the instrument in a specific application or system, refer to the manuals for the peripheral equipment. Contact a technical representative for compatibility or suitability questions.

Instrument connections

[Figure 1](#) shows the main connections on the back of the instrument.


[Figure 2](#) shows the layout and connections of components in the top of the instrument.

Figure 1 Rear panel connections

1 Vent outlet	7 Line power connection
2 Pressure/vacuum vent	8 AC fuse compartment
3 Pressure port	9 On/Off switch
4 Vacuum port	10 Sampler communication port (to counter)
5 Sample drain port	11 Sensor cable port (to counter)
6 Pressure/vacuum pump AC outlet	

Figure 2 Top compartment components

1 Micrometer valve inlet	5 Check valve
2 Micrometer valve outlet	6 Sight tube with ball valve
3 Micrometer valve bypass	7 Desiccant chamber
4 Check valve inlet	8 Line filter


Product components

⚠ WARNING

Personal injury hazard. Instruments or components are heavy. Use assistance to install or move.


Make sure that all components have been received (Figure 3 and Figure 4). If any items are missing or damaged, contact the manufacturer or a sales representative immediately.

Figure 3 Instrument components

1 ABS-2	7 Drain line tubing, 3/8-in.
2 Pressure/Vacuum port cap	8 Stir bars
3 O-rings	9 Communication cable
4 Line fuses (T, 6.3 A, 250 V) (500-050-0002)	10 Power cord (USA)
5 Ball valve	11 Power cord (European)
6 Ball valve warning label	

Figure 4 Sensor components

1 MC05 sensor

2 HRLD sensor

3 HRLD sensor cable

Installation

⚠ DANGER

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

Recommendations for safe installation

⚠ DANGER

Fire/explosion hazard. Do not use the instrument with samples or liquids that have an ignition point below 90 °C (194 °F).

⚠ DANGER

Fire and inhalation hazards. Potentially dangerous gases can be made when samples are moved or analyzed. If there is any doubt about a fire hazard or an inhalation hazard, keep the equipment under a fume hood during the analysis.

⚠ CAUTION

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current material safety data sheets (MSDS) for safety protocols.

⚠ CAUTION

Chemical exposure hazard. Dispose of chemicals and wastes in accordance with local, regional and national regulations.



Make sure that hazardous fluids or fluids that may produce hazardous vapors are handled safely. These fluids must be handled in accordance with local regulatory agency requirements on permissible exposure limits.


Recommendations for safe handling include:

- Connect the vent outlet to an approved fume exhaust.
- Install the ball valve to keep fumes out of the environment.
- Connect the sample drain port to an approved container.



Install the sensor

The illustrated steps show how to install the sensor.

Install the ball valve

Micrometer valve bypass

To prevent blockage of the micrometer valve when highly viscous or dirty fluids are sampled, bypass the micrometer valve. Refer to [Figure 2](#) on page 6.

1. Remove the tube and fittings between the micrometer valve outlet and the check valve inlet.
2. Remove the micrometer valve inlet tube and fittings from the micrometer valve.
3. Attach the micrometer valve inlet tube to the check valve inlet. Attach and tighten the fittings.

Connect the vent outlet

⚠ DANGER	
	Fire and inhalation hazards. Potentially dangerous gases can be made when samples are moved or analyzed. If there is any doubt about a fire hazard or an inhalation hazard, keep the equipment under a fume hood during the analysis.
⚠ CAUTION	
	Chemical exposure hazard. Dispose of chemicals and wastes in accordance with local, regional and national regulations.

Required items:


- 6.35 mm (¼ in.) O.D. vent service line
 - User defined fume exhaust
 - Nut and ferrule
1. Connect the vent service line to the fume exhaust with the nut and ferrule.
 2. Make sure that the connection is tight and fix all leaks.

Connect the sample drain

⚠ DANGER	
	Fire/explosion hazard. Do not use the instrument with samples or liquids that have an ignition point below 90 °C (194 °F).

The drain line is a 9.5 mm (3/8 in.) connection made from Tygon™ tubing. 122 cm (4 ft) of tubing is supplied with the instrument. Connect the sample drain as shown in [Figure 5](#).

Figure 5 Connect the sample drain


Connect the pressure inlet and vacuum exhaust

The pressure inlet is a 6.35 mm (¼-in.) connection. The connection keeps the ABS 2 line pressure at or below 60 PSI (413.7 kPa).

The pressure line is installed from either the house air (clean dry air) or from the pressure port of the pressure/vacuum pump. For house air and vacuum, install the 6.35 mm (¼-in.) cap supplied with the instrument on the pressure/vacuum exhaust port.

Connect the inlet and exhaust as shown in [Figure 6](#).

Figure 6 Connect pressure inlet and vacuum exhaust

Connect the power cord

⚠ DANGER

Electrical shock and fire hazards. Replacement power cords must:

- Have the correct plug style for the outlet connection
- Have a rating sufficient for the supply voltage and current. Refer to the requirements in the Specifications section.
- Meet or exceed local electrical code requirements

The ABS 2 must be located within 2.1 m (7 ft) of the grounded electrical power receptacle.

Do not use extension cords with the ABS 2.

Put the ABS 2 in a location to let the power switch be easily operated.


Do not connect the power cord to the power receptacle until all connections are made and checked and the system is ready for use.

If an optional pressure/vacuum pump is used, connect the power cord from the pump to the air pump power outlet.

1. Connect the power cord to the power receptacle on the back of the instrument.

2. Connect power to the instrument. Refer to [Figure 7](#).

Figure 7 Connect the power cords

Validate the instrument

Connect external equipment (according to the equipment documentation), then complete the steps in this procedure.

1. Make sure that sensor power and communication cables are correctly connected.
2. Make sure that the sample, drain, pressure, vent, vacuum and exhaust ports are connected and capped, and there are no leaks.
3. Make sure that the line pressure to the ABS 2 is not greater than 60 PSI (413.7 kPa).
4. Make sure that the zero level detector is set to 0 on the sight tube.
5. Before sample fluid is run, flush clean (particle free) water through the instrument and fix any leaks.

Operation

⚠ DANGER	
	Fire/explosion hazard. Do not use the instrument with samples or liquids that have an ignition point below 90 °C (194 °F).
⚠ WARNING	
 	Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current material safety data sheets (MSDS) for safety protocols.
⚠ WARNING	
	Personal injury hazard. Enclosed systems contain high pressure. Qualified personnel must remove pressure from the system before the sensor can be installed or removed.

The instrument can operate in Automatic, Manual, Pressure and Vacuum modes. In Manual mode, set the instrument to start and end the sampling procedure. In Automatic mode, a host particle counter sets the instrument when to start and when to end the sampling procedure. Automatic mode allows batch sampling. The Pressure mode pressurizes the sample liquid for transport through the sensor. The Vacuum mode evacuates the sample chamber and degasses the sample before analysis.


While in automatic or manual mode, if the sight tube is almost in an overflow condition, immediately stop the sample process and remove power from the instrument. Do not move the chamber release/lock handle while the instrument chamber is pressurized. Depressurize the pressure chamber before the handle is tightened or released.

Monitor the process at all times when the instrument operates in manual mode.

If an overflow condition occurs, immediately remove power from the instrument and do the cleanup procedures described in the maintenance section.

Figure 8 shows the operator controls. Table 1 shows control descriptions.

Figure 8 Operator controls

1 Flow control valve	5 Pump control switch
2 Stir bar speed control	6 Vacuum control valve
3 Vacuum/Pressure mode switch	7 Pressure control valve
4 Manual/Auto mode switch	8 Pressure/vacuum gauge

Table 1 Control descriptions

Control	Description
Flow control valve	Adjusts the rate at which the sample flows through the sensor.
Stir bar speed controller	Adjusts the rotation speed of the stir bar.
Vacuum/Pressure mode switch	Puts the instrument in Vacuum or Pressure operation mode.
Manual/Auto mode switch	Puts the instrument in Manual or Automatic operation mode.
Pump control switch	Turns the optional remote pressure/vacuum pump on or off.
Vacuum control valve	Adjusts the vacuum in the sample chamber. The pressure/vacuum gauge shows the amount of chamber vacuum.
Pressure control valve	Adjusts the line pressure in the sample chamber.
Pressure/vacuum gauge	The pressure/vacuum gauge shows the amount of chamber pressure.

Electromagnetic stir assembly

Particles in the sample liquid must stay suspended as a sampling operation occurs. A magnetic field makes the stir bar rotate in the sample. This action stirs the sample at a constant speed and keeps the particles suspended. The user adjusts a control on the front of the instrument to change the stir speed.

Make sure that the stir bars are:

- Clean (particle free)
- Sensitive to magnetic fields
- Covered in glass or another material that does not give off particles
- Non-magnetic

Adjust the flow rate

Note: Do not over-tighten the flow control valve as this can permanently damage the micro-valve assembly.

The initial settings for a sample operation are:

- Flow valve fully open
- Line pressure 20 psi (137.8 kPa)

Adjust the flow control valve and the line pressure so a specified volume fills the sight tube in a specified amount of time.

For best results, adjust the flow rate of the sample fluid to be the same as the calibrated specification for the sensor.

1. Set the volume level detector at the calibrated flow rate volume for the sensor.
2. Set the pressure/vacuum switch to Pressure.
3. Make sure that a clean (particle free) stir bar is in the sample bottle. Then put the sample bottle on the stir assembly and lock the stir assembly to the sample chamber.
4. Turn on the automatic bottle sampler.
5. Set the auto/manual switch to Automatic and configure the particle counter for volume mode operation.
6. On the counter, push START.
Sample fluid flows and becomes visible in the sight tube.
7. Start the timer when the sample is at the 0 in the sight tube.
8. Stop the timer when the counter stops the bottle sampler.
The remote pump (if used) turns off and the sample chamber fumes are removed.
9. Adjust the flow rate of the bottle sampler as necessary.
 - a. If the time necessary to get the expected sample volume is greater than 1 minute, turn the flow control valve clockwise or decrease the line pressure.
 - b. If the time necessary to get the expected sample volume is less than 1 minute, turn the flow control valve counterclockwise or increase the line pressure.

10. Do the procedure steps above again and again until the sample volume is completed in 1 minute.

Manual mode

In the Manual mode, the sampler operates as a fluid transport device and not as a volume sampler. The sight tube detectors are not used. The procedure below can be repeated as often as necessary.

To operate the sampler in Manual mode:

1. Set the vacuum/pressure mode switch to the pressure mode.
2. Set the manual/auto mode switch to automatic.
3. Turn on the automatic bottle sampler.
4. Configure the particle counter for manual mode operation.
5. Make sure that a clean (particle free) stir bar is in the sample bottle. Place the sample bottle on the stir assembly and lock the stir assembly to the sample chamber.
6. Push the automatic/manual switch.
The remote pump starts and there is fluid in the sight tube.
7. When the sample reaches the zero indicator on the sight tube, push **Start** on the counter.
8. When the sample reaches the correct indicator on the sight tube, or when the instrument has sent enough sample for data, push the Automatic/Manual switch.
The pump and counter stop.
9. After sample operations are complete, flush the sensor sample lines and the sight tube with clean (particle free) water or other clean liquid that can mix with the sample liquid.

Automatic mode

In the Automatic mode, the level detectors of the ABS 2 control the start and stop of the particle counting process.

1. Make sure that the volume level detector is set at the necessary sample volume indicator on the sight tube.

Note: If multiple sample runs are necessary, make sure that the total sample volume is sufficiently large. The total sample volume is equal to the sample volume and the instrument tare volume combined (16 mL).

2. Set the pressure/vacuum switch to pressure mode.
3. Set the automatic/manual switch to Automatic.
4. Turn on the automatic bottle sampler and the pump (if used).
5. Configure the particle counter for Volume mode operation.
6. Make sure that a clean (particle free) stir bar is in the sample bottle.
7. Put the sample bottle on the stir assembly and lock the stir assembly to the sample chamber.
8. Start the sampling.
The pump starts and fluid becomes visible in the sight tube. Counting begins when the sample level reaches the lower level detector. Counting ends when the sample level reaches the upper volume level detector. The sample drains from the sight tube and the instrument is ready for the next sample analysis.
9. After sample operations are complete, flush the sensor sample lines and the sight tube with clean (particle free) water or other clean liquid that can mix with the sample liquid.

Vacuum mode

⚠ DANGER	
	Personal injury hazard. Do not adjust the chamber release and the lock handle while the sample chamber is under pressure. Always depressurize the sample pressure chamber before the lock handle is tightened or released.

Put a sample in the vacuum chamber

To adjust the volume and zero volume detector, refer to [Figure 9](#). To put a sample in the vacuum chamber, refer to [Figure 10](#).

This diagram illustrates the vertical scale and its integration into the microscope's frame. On the left, a detailed view of the scale is shown, featuring a central vertical ruler with markings from 0 to 110. Two horizontal sliding mechanisms, labeled 1 and 2, are positioned on either side of the ruler. Label 1 points to the upper sliding mechanism, and label 2 points to the lower sliding mechanism. Label 3 points to the base of the scale assembly. On the right, a perspective view of the microscope shows the scale assembly (indicated by a black cone) mounted on the left side of the frame, adjacent to the objective and eyepiece lenses.

1 Volume detector	3 Zero volume detector
2 Sight tube	

The diagram illustrates the components of the HANNA HI9142 pH titrator. The main unit features a digital display at the top, a central titration chamber with a stirrer, and a sample container at the bottom. A control panel at the base includes a large analog meter and two smaller circular indicators. Five numbered callouts identify key parts: 1 points to the sample container, 2 to the titrant reservoir, 3 to the titrant delivery valve, 4 to the sample container lid, and 5 to the sample container. Two circular insets provide detailed views: the left inset shows the lid being locked (indicated by a padlock icon and a downward arrow) and unlocked (indicated by a padlock icon and an upward arrow); the right inset shows the titrant reservoir being filled with liquid from a graduated cylinder.

1 Stir assembly	4 Stir bar
2 Sample chamber	5 Sample bottle
3 Chamber release handle	

Degas a sample

In some situations, samples may develop air bubbles. To degas a sample before analysis, refer to steps 1-5.

Note: Some fluids contain volatile constituents that boil off under vacuum conditions. If a sample boils, stop the degassing operation.

1. Put a sample in the sample chamber (Figure 10 on page 15).
2. Set the pressure/vacuum switch to the vacuum mode.

3. Set the automatic/manual switch to automatic.
4. Turn on the automatic bottle sampler and the pump (if connected).
5. Look at the pressure/vacuum gauge and adjust the vacuum regulator until the pressure gauge shows a vacuum.
6. Give the sample at least 1 minute to degas.

Note: Degassing time is different for fluids with different viscosities.

Maintenance

⚠ DANGER	
	Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

General maintenance

⚠ DANGER	
	Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

⚠ DANGER	
	Fire and inhalation hazards. Potentially dangerous gases can be made when samples are moved or analyzed. If there is any doubt about a fire hazard or an inhalation hazard, keep the equipment under a fume hood during the analysis.

⚠ DANGER	
	Electrocution hazard. Do not disassemble the instrument for maintenance. If the internal components must be cleaned or repaired, contact the manufacturer.

⚠ WARNING

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current material safety data sheets (MSDS) for safety protocols.

Use a fume hood when the instrument is maintained or cleaned.

Note: Always use clean (particle free) water or other appropriate liquid to clean the instrument sampling system.

Flush the sampling system of the instrument at these times:

- When foreign material has entered or is suspected to have entered the sampling system.
- Before any instrument component is opened for maintenance or repair.
- Between the completion of one sample analysis and the start of another.

How to clean spills

⚠ CAUTION

Chemical exposure hazard. Dispose of chemicals and wastes in accordance with local, regional and national regulations.

1. Obey all facility safety protocols for spill control.
2. Discard the waste according to applicable regulations.

Flush the sampling system

Flush the system after sampling operations have been completed.

To switch from water based fluids to oil based fluids, flush the system with three 20 mL samples of Reagent Grade IPA (purity 99% or greater).

To switch from oil based to water based fluids, flush the system with three 20mL samples of a filtered solvent (QED, Mineral Spirits, etc.). Then flush the system with Reagent Grade IPA as above.

To flush the sampling system:

1. Make sure that the volume level detector is set at the correct sample volume indicator on the sight tube.
2. Set the pressure/vacuum switch to the pressure mode.
3. Set the automatic/manual switch to automatic.
4. Turn on the instrument and the pump (if connected).
5. Place a bottle of clean (particle free) water, water mixed with a mild detergent, or liquid as specified above on the stir assembly. Lock the assembly to the sample chamber.
6. On the counter, push START.
The pump will start (if connected to a remote pump) and the sight tube shows fluid. When the liquid drains from the sight tube the instrument is ready for the next cycle.

Clean the instrument

Sample spills and overflow can result in potential electrical and fire hazards. In the event of an overflow or spill, disconnect power and refer to [How to clean spills](#) on page 16.

Clean the exterior of the instrument with a moist cloth and a mild soap solution.

Clean the sight tube

Use an apparatus brush and a mild soap solution to clean the sight tube.

1. Turn off the instrument and disconnect the line power cord.
2. Push the release buttons on either side of the top assembly and open the top of the instrument.
3. Remove the sight tube ball valve.
4. Gently put the apparatus brush in to the sight tube assembly. Gently move the brush up and down five to ten times. Do not allow liquid to exit the sight tube or to spill on the rear panel.

5. Assemble all instrument components.

6. Flush the sample system.

Replace the line filter

A $\frac{9}{16}$ -inch or adjustable wrench is necessary to do this procedure.

Remove instrument power, disconnect the line power cord, secure all pressure lines and flush the system before doing the procedure below.

Replace the line filter when the:

- filter is wet
 - filter is contaminated with foreign material
 - desiccant filter is replaced
1. Push the release buttons on the top of the instrument and open the top.
 2. Use the $\frac{9}{16}$ -inch wrench to remove the fittings from each side of the line filter. Record the direction of the flow arrow on the filter.
 3. Remove the old filter and install the new filter and tighten the fittings by hand. Install the new filter so the arrow points in the same direction.
 4. Use the wrench to tighten the fittings another $\frac{1}{4}$ turn. Do not over tighten the fittings as this can cause permanent damage to the fittings.

Replace the desiccant

A $\frac{9}{16}$ -inch or adjustable wrench is necessary for this procedure.

Turn off the instrument, secure all pressure lines, and disconnect the power cord before this procedure.

The desiccant removes moisture from the pressurized line gas. Desiccant beads are normally blue. Install a new desiccant chamber (390020) when the beads turn white or brown. Install a new line filter at the same time.

1. Push the release buttons on the sides of the top assembly and open the top of the instrument.

2. Loosen the sample line fittings on both sides of the desiccant chamber.
3. Remove the old chamber and install the new one.
4. Tighten the fittings.

Clean the check valve

A $\frac{7}{16}$ -inch or adjustable wrench is necessary for this procedure.

The check valve keeps line pressure through the sensor consistent. Clean the valve when the necessary sample flow is difficult or impossible to maintain.

1. Turn off the instrument power and disconnect the power cord.
2. Push the release buttons on both sides of the top assembly and open the top of the instrument.
3. Loosen the line connectors and remove the valve. Record the direction of the flow arrow. The arrow must point in the same direction when the valve is installed after it is cleaned.
4. Soak the valve for a few minutes in a solvent compatible with the sample solution.
5. Flush the valve with clean water and install the valve. Install the valve with the correct orientation.
6. Attach and hand tighten the sample line connectors. Tighten the connectors with the wrench. Do not over tighten the connectors.

Line pressure drier maintenance

There is no easy access to the line pressure drier. If liquid is found under the instrument or if the line pressure drier fills with moisture, contact technical support.

Replace the line fuse

A small flat head screwdriver is necessary for this procedure.

⚠ DANGER


Fire hazard. Use the same type and current rating to replace fuses.

⚠ WARNING

Potential Electrocution Hazard. Always disconnect power to the instrument when making electrical connections.

Troubleshooting

To troubleshoot the performance of the ABS 2:

1. Examine the power and electrical parts of the unit. Make sure that the ABS 2 power switch is set to ON. Make sure that the fuses have not failed.
2. If the instrument still will not sample fluids, make sure that the pressure and vacuum hoses are connected to the correct ports on the rear of the sampler.
3. Make sure that the micrometer valve is fully open (set to around 20 on the micrometer scale).
4. Make sure that the sampler is configured to operate in the correct mode. Make sure that the unit is not in Manual mode if Automatic mode is necessary for the application.
5. If fluid runs through the sensor and the instrument does not give count data, there may be a problem with the sensor. Refer to the sensor operation manual.

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Tel. +41 22 594 6400

Fax +41 22 594 6499

Artisan Technology Group is an independent supplier of quality pre-owned equipment

Gold-standard solutions

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

Learn more!

Visit us at [artisan^{tg}.com](https://www.artisantg.com) for more info on price quotes, drivers, technical specifications, manuals, and documentation.

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We're here to make your life easier. How can we help you today?

(217) 352-9330 | sales@artisan^{tg}.com | [artisan^{tg}.com](https://www.artisantg.com)

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P.O. Box 389, Loveland, CO 80539-0389 U.S.A.
(800) 227-4224 (U.S.A. only)