

Genomic Solutions HiGro

## High-Capacity Microwell Plate Growth System



Limited Availability  
Used and in Excellent Condition

Open Web Page

<https://www.artisanng.com/76883-1>

All trademarks, brandnames, and brands appearing herein are the property of their respective owners.



Your **definitive** source  
for quality pre-owned  
equipment.

**Artisan Technology Group**

(217) 352-9330 | [sales@artisanng.com](mailto:sales@artisanng.com) | [artisanng.com](http://artisanng.com)

- Critical and expedited services
- In stock / Ready-to-ship

- We buy your excess, underutilized, and idle equipment
- Full-service, independent repair center

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

## HiGro®

A high-capacity microwell plate growth system.



The HiGro is a high-capacity, incubating shaker designed to accelerate sample growth in microwell plates. This uniquely engineered instrument combines a small shaking orbital, gas flow and temperature control system in a **user-friendly, compact layout**. Plate-holding cassettes allow for **easy loading** and unloading.

The convenience and control of the HiGro system make it the ideal solution for facilities with high-throughput growth needs as well as for laboratories using applications that are best optimized in a microwell format.

### Key Benefits

- Space-saving design
- Optimizes growth in 96- and 384-microwell plate formats
- Easy to load and unload
- Superior control over growth environment

## Design Specifications


<b>Dimensions</b>	26" W x 21" D x 22" H (0.66m x 0.53m x .56m) (small footprint allows placement on lab bench)
<b>Capacity</b>	4 growth chambers, each holding 1 cassette with 12 shallow-well or 6 deep-well plates. Total system capacity is 48 shallow-well plates or 24 deep-well plates.
<b>Temperature</b>	Independent control for each chamber. Range is ambient to 55°C. Resolution is 1°C.
<b>Shaking</b>	Rotational orbital is 8.0 mm (ideally suited for 96-well growth) or 4.0 mm (ideally suited for 384-well growth)*. Maximum speed is 650 rpm. Resolution is 1 rpm.
<b>Gas Flow</b>	Oxygen- or air-ready*. Gas flow rate and timing controls located on instrument front panel.

\*Depending on Model

## Data

**HiGro & Floor Shaker Comparison for Plasmid Growth**

**Figure 1**



**Figures 1 & 2**

**Conditions:** E.coli strain is DH10B carrying pUC 19; 1x TB media; inoculum: 5 µl glycerol stock culture; selection: Ampicillin; 37°C.


**HiGro:** 300 µl/well in 96-shallow-well, flat-bottom plate; 550 rpm; Air flow: no delay, 2 sec on, 30 sec off, 6 standard liters per minute; no plate lids.

**Figure 1:**

**Floor Shaker:** 1000 µl/well in 96-well Beckman deep-well block, 250 rpm in New Brunswick floor shaker (plates were taped down); aluminum foil placed loosely over each plate; wetted paper towel placed in shaker to maintain ambient humidity.


**Cross-Contamination Growth Test with No Plate Lids**

**Figure 2**



**Plasmid Yield Time Course from 384-Well Growth\***

**Figure 3**



>1µg plasmid after only 8 hours!  
(total yield from 90 µl growth volume)

DNA eluted into 50 µl EB

Each band in the ladder represents 25 ng DNA

**Figure 3:**

**Growth Conditions:** 90 µl/well Magnificent Media, 384-shallow-well, flat-bottom plate, rounded-square wells; inoculation (1 µl) from O/N culture; high copy plasmid, 3 kb total size; selection: Carbapenem.

**HiGro:** 37°C; 500 rpm; Air flow: no delay, 0.5 sec on, 0.5 min off, 1.5 standard liters per minute; no plate lids.

**Plasmid Prep Conditions:** 90 µl (total culture) pipetted into alkaline lysis reagents (100 µl). Neutralization reagents added (100 µl). Pelleted sample. Isolated DNA on commercially available column. Eluted into 50 µl EB.

**Gel Conditions:** 2 µl & 4 µl of each (50 µl) prep were run on 1% agarose gel. Ladder: two 10-band mass ladders combined (250 ng of each ladder).

Request the latest HiGro Application & Technical Notes.

\*Data obtained through customer collaboration

## Worldwide Headquarters

Digilab, Inc.  
84 October Hill Road  
Holliston, MA 01746  
USA

Phone: (508) 893-3130  
Toll Free: (800) 935-8007  
Fax: (508) 893-8011

E-Mail: [info@digilabglobal.com](mailto:info@digilabglobal.com)

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | [www.artisan-tg.com](http://www.artisan-tg.com)

**DIGILAB®**

# Artisan Technology Group is an independent supplier of quality pre-owned equipment

## Gold-standard solutions

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

## We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

## Learn more!

Visit us at [artisanng.com](https://www.artisanng.com) for more info on price quotes, drivers, technical specifications, manuals, and documentation.

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

**We're here to make your life easier. How can we help you today?**

(217) 352-9330 | [sales@artisanng.com](mailto:sales@artisanng.com) | [artisanng.com](https://www.artisanng.com)

