

Spex 1934C
Phosphorimeter



\$995.00

In Stock

Qty Available: 1

Used and in Excellent Condition

Open Web Page

<https://www.artisanng.com/74604-1>

All trademarks, brandnames, and brands appearing herein are the property of their respective owners.



Your **definitive** source
for quality pre-owned
equipment.

Artisan Technology Group

(217) 352-9330 | sales@artisanng.com | artisanng.com

- Critical and expedited services
- In stock / Ready-to-ship

- We buy your excess, underutilized, and idle equipment
- Full-service, independent repair center

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

Instruction Manual
1934D PHOSPHORIMETER



INDUSTRIES, INC. 3880 PARK AVE. EDISON, N.J. 08620

ELEX: 076341





1. 1922-1923



THE UNIVERSITY OF CHICAGO

Contents

1.0 Description	2
2.0 Installation	5
3.0 Operation	6
3.1 Start-up Procedure	6
3.2 Reference Detector	7
3.3 Controlling the Phosphorimeter from the DM3000.....	7
3.3a Acquisition Mode	8
3.3b Delay After Flash.....	8
3.3c Sample Window	8
3.3d Number of Flashes	9
3.3e Flash Rate	9
3.4 System Status	9
3.5 Emission Scans and Time Discrimination.....	9
3.5a Scanning a Emission Spectrum	9
3.5b An Example: Obtaining the Emission Spectrum of a Sample	12
3.6 Excitation Scans	12
3.7 Synchronous Scans	12
4.0 Determination of Error Codes	12
4.1 Error Reporting	12
4.2 Interpretation of Display.....	13
5.0 Troubleshooting Guide	14
6.0 Flash Lamp Replacement	14
6.1 Lamp Replacement Procedure	15
6.2 Lamp Alignment and Focusing Procedure	15
7.0 Specifications	18
Appendix A Error codes	19
Appendix B Modifying Configuration	21

1. The first part of the report deals with the general situation of the country and the progress of the work during the year.

2. The second part contains a detailed account of the work done in the various departments.

3. The third part gives a summary of the results of the work and a statement of the financial position.

4. The fourth part contains a list of the names of the persons who have been employed during the year.

5. The fifth part contains a list of the names of the persons who have been promoted during the year.

6. The sixth part contains a list of the names of the persons who have been dismissed during the year.

7. The seventh part contains a list of the names of the persons who have been transferred during the year.

8. The eighth part contains a list of the names of the persons who have been retired during the year.

9. The ninth part contains a list of the names of the persons who have been appointed during the year.

10. The tenth part contains a list of the names of the persons who have been promoted during the year.

1.0 DESCRIPTION

The SPEX Phosphorimeter is an accessory to the Fluorolog-2 Series of Spectrofluorometers that adds a programmable pulsed excitation source and selectable gating of signal from the photomultiplier tube (PMT). This provides time-discrimination capability to sort out competing luminescence emissions on the basis of their lifetimes.

Because the duration of each exciting pulse from the phosphorimeter is very short (about 3 microseconds), lamp interference during acquisition of decay curves is minimized. This allows the researcher to follow the decay of samples an order of magnitude faster than conventional systems that depend on mechanical choppers.

The phosphorimeter consists of three main modules: an illuminator, a phosphorimeter control module and a reference amplifier module. The reference channel circuitry is housed in the phosphorimeter control module.

The illuminator housing, or flash lamp, operates at rates up to 33 pulses per second. Each lamp pulse is triggered by a signal from the control module. When the start of the actual light output is sensed, a signal is sent to the phosphorimeter control module for timing purposes.

The phosphorimeter control module houses the signalgating circuitry that intercepts the sample signal from the pulse-counting emission photomultiplier tube (PMT), collects a selected, time-delimited portion of the signal, and later passes it to the DM3000. The reference amplifier module houses the fast amplifier for the reference channel.

A typical sequence of data-taking begins with a flash from the pulsed lamp, which is sensed by the control module as time zero (see Figure 1). The radiation enters the excitation spectrometer where it is dispersed. A selected wavelength band is passed on to excite the sample. After a delay, the sample emits luminescence that passes through the emission spectrometer to the PMT. The signal from the PMT is collected by the control module for a preset length of time. Any signal that arrives before or after this time is ignored. This sequence of excitation, delay, and sampling is repeated for each lamp flash and the signal is accumulated for a predetermined number of exciting pulses before the data is sent to the DM3000. The data can be manipulated by the DM3000's firmware routines or by user-designed programs.

1. The first part of the report deals with the general situation of the company and the results of the work done during the year. It is a summary of the main achievements and the progress made in the various departments.

2. The second part of the report deals with the financial results of the company. It shows the income and expenditure for the year and the position of the company's assets and liabilities.

3. The third part of the report deals with the personnel and the work done by the various departments. It shows the number of employees and the work done by each department.

4. The fourth part of the report deals with the future prospects of the company. It shows the plans for the next year and the prospects for the company's growth.

5. The fifth part of the report deals with the conclusions of the report. It shows the main findings of the report and the recommendations for the future.

6. The sixth part of the report deals with the appendix. It contains the detailed financial statements and the other documents referred to in the report.




Figure 1 Phosphorimeter event sequence

This sequence of events is controlled by the DM3000 through four phosphorimeter parameters - Delay After Flash, Sample Window, Flash Rate and Number of Flashes. These appear as additional Data Acquisition Parameters under the Define Experiment menu when the phosphorimeter hardware option is enabled. The signal collected during the sample time may be integrated over as many as 999 lamp flashes before the data is sent to the DM3000. The delay between excitation pulses of the sample and opening of Sample Window can be set from 0 to 10,000 milliseconds in 0.001 millisecond increments. The rate at which excitation pulses are delivered can be as high as 33 flashes per second.

The Delay After Flash parameter can be varied as a function of time to produce a decay curve (Figure 2). Spectra can be scanned to isolate different phosphorescing components based on the lifetime of the luminescence species present in the sample (Figure 3 and 4). These two techniques together can be used to produce three-dimensional plots. For example, you can plot successive scans with varying delay times, as shown in Figure 5. Each contour then isolates the species in a particular slice of time.




Figure 2 Decay curve produced by plotting intensity as a function of delay time

例 1 已知函数 $f(x) = x^2 + 2x + 1$ ，求 $f(1)$ 的值。

解 将 $x=1$ 代入函数表达式，得

$$f(1) = 1^2 + 2 \times 1 + 1 = 4$$

例 2 已知函数 $f(x) = x^2 + 2x + 1$ ，求 $f(x)$ 的表达式。

解 由已知条件，得 $f(x) = x^2 + 2x + 1$ 。将 $x=1$ 代入，得 $f(1) = 4$ 。将 $x=2$ 代入，得 $f(2) = 9$ 。将 $x=3$ 代入，得 $f(3) = 16$ 。将 $x=4$ 代入，得 $f(4) = 25$ 。将 $x=5$ 代入，得 $f(5) = 36$ 。将 $x=6$ 代入，得 $f(6) = 49$ 。将 $x=7$ 代入，得 $f(7) = 64$ 。将 $x=8$ 代入，得 $f(8) = 81$ 。将 $x=9$ 代入，得 $f(9) = 100$ 。将 $x=10$ 代入，得 $f(10) = 121$ 。将 $x=11$ 代入，得 $f(11) = 144$ 。将 $x=12$ 代入，得 $f(12) = 169$ 。将 $x=13$ 代入，得 $f(13) = 196$ 。将 $x=14$ 代入，得 $f(14) = 225$ 。将 $x=15$ 代入，得 $f(15) = 256$ 。将 $x=16$ 代入，得 $f(16) = 289$ 。将 $x=17$ 代入，得 $f(17) = 324$ 。将 $x=18$ 代入，得 $f(18) = 361$ 。将 $x=19$ 代入，得 $f(19) = 400$ 。将 $x=20$ 代入，得 $f(20) = 441$ 。将 $x=21$ 代入，得 $f(21) = 484$ 。将 $x=22$ 代入，得 $f(22) = 529$ 。将 $x=23$ 代入，得 $f(23) = 576$ 。将 $x=24$ 代入，得 $f(24) = 625$ 。将 $x=25$ 代入，得 $f(25) = 676$ 。将 $x=26$ 代入，得 $f(26) = 729$ 。将 $x=27$ 代入，得 $f(27) = 784$ 。将 $x=28$ 代入，得 $f(28) = 841$ 。将 $x=29$ 代入，得 $f(29) = 900$ 。将 $x=30$ 代入，得 $f(30) = 961$ 。将 $x=31$ 代入，得 $f(31) = 1024$ 。将 $x=32$ 代入，得 $f(32) = 1089$ 。将 $x=33$ 代入，得 $f(33) = 1156$ 。将 $x=34$ 代入，得 $f(34) = 1225$ 。将 $x=35$ 代入，得 $f(35) = 1296$ 。将 $x=36$ 代入，得 $f(36) = 1369$ 。将 $x=37$ 代入，得 $f(37) = 1444$ 。将 $x=38$ 代入，得 $f(38) = 1521$ 。将 $x=39$ 代入，得 $f(39) = 1600$ 。将 $x=40$ 代入，得 $f(40) = 1681$ 。将 $x=41$ 代入，得 $f(41) = 1764$ 。将 $x=42$ 代入，得 $f(42) = 1849$ 。将 $x=43$ 代入，得 $f(43) = 1936$ 。将 $x=44$ 代入，得 $f(44) = 2025$ 。将 $x=45$ 代入，得 $f(45) = 2116$ 。将 $x=46$ 代入，得 $f(46) = 2209$ 。将 $x=47$ 代入，得 $f(47) = 2304$ 。将 $x=48$ 代入，得 $f(48) = 2401$ 。将 $x=49$ 代入，得 $f(49) = 2500$ 。将 $x=50$ 代入，得 $f(50) = 2601$ 。将 $x=51$ 代入，得 $f(51) = 2704$ 。将 $x=52$ 代入，得 $f(52) = 2809$ 。将 $x=53$ 代入，得 $f(53) = 2916$ 。将 $x=54$ 代入，得 $f(54) = 3025$ 。将 $x=55$ 代入，得 $f(55) = 3136$ 。将 $x=56$ 代入，得 $f(56) = 3249$ 。将 $x=57$ 代入，得 $f(57) = 3364$ 。将 $x=58$ 代入，得 $f(58) = 3481$ 。将 $x=59$ 代入，得 $f(59) = 3600$ 。将 $x=60$ 代入，得 $f(60) = 3721$ 。将 $x=61$ 代入，得 $f(61) = 3844$ 。将 $x=62$ 代入，得 $f(62) = 3969$ 。将 $x=63$ 代入，得 $f(63) = 4096$ 。将 $x=64$ 代入，得 $f(64) = 4225$ 。将 $x=65$ 代入，得 $f(65) = 4356$ 。将 $x=66$ 代入，得 $f(66) = 4489$ 。将 $x=67$ 代入，得 $f(67) = 4624$ 。将 $x=68$ 代入，得 $f(68) = 4761$ 。将 $x=69$ 代入，得 $f(69) = 4900$ 。将 $x=70$ 代入，得 $f(70) = 5041$ 。将 $x=71$ 代入，得 $f(71) = 5184$ 。将 $x=72$ 代入，得 $f(72) = 5329$ 。将 $x=73$ 代入，得 $f(73) = 5476$ 。将 $x=74$ 代入，得 $f(74) = 5625$ 。将 $x=75$ 代入，得 $f(75) = 5776$ 。将 $x=76$ 代入，得 $f(76) = 5929$ 。将 $x=77$ 代入，得 $f(77) = 6084$ 。将 $x=78$ 代入，得 $f(78) = 6241$ 。将 $x=79$ 代入，得 $f(79) = 6400$ 。将 $x=80$ 代入，得 $f(80) = 6561$ 。将 $x=81$ 代入，得 $f(81) = 6724$ 。将 $x=82$ 代入，得 $f(82) = 6889$ 。将 $x=83$ 代入，得 $f(83) = 7056$ 。将 $x=84$ 代入，得 $f(84) = 7225$ 。将 $x=85$ 代入，得 $f(85) = 7396$ 。将 $x=86$ 代入，得 $f(86) = 7569$ 。将 $x=87$ 代入，得 $f(87) = 7744$ 。将 $x=88$ 代入，得 $f(88) = 7921$ 。将 $x=89$ 代入，得 $f(89) = 8100$ 。将 $x=90$ 代入，得 $f(90) = 8281$ 。将 $x=91$ 代入，得 $f(91) = 8464$ 。将 $x=92$ 代入，得 $f(92) = 8649$ 。将 $x=93$ 代入，得 $f(93) = 8836$ 。将 $x=94$ 代入，得 $f(94) = 9025$ 。将 $x=95$ 代入，得 $f(95) = 9216$ 。将 $x=96$ 代入，得 $f(96) = 9409$ 。将 $x=97$ 代入，得 $f(97) = 9604$ 。将 $x=98$ 代入，得 $f(98) = 9801$ 。将 $x=99$ 代入，得 $f(99) = 10000$ 。将 $x=100$ 代入，得 $f(100) = 10201$ 。将 $x=101$ 代入，得 $f(101) = 10404$ 。将 $x=102$ 代入，得 $f(102) = 10609$ 。将 $x=103$ 代入，得 $f(103) = 10816$ 。将 $x=104$ 代入，得 $f(104) = 11025$ 。将 $x=105$ 代入，得 $f(105) = 11236$ 。将 $x=106$ 代入，得 $f(106) = 11449$ 。将 $x=107$ 代入，得 $f(107) = 11664$ 。将 $x=108$ 代入，得 $f(108) = 11881$ 。将 $x=109$ 代入，得 $f(109) = 12100$ 。将 $x=110$ 代入，得 $f(110) = 12321$ 。将 $x=111$ 代入，得 $f(111) = 12544$ 。将 $x=112$ 代入，得 $f(112) = 12769$ 。将 $x=113$ 代入，得 $f(113) = 12996$ 。将 $x=114$ 代入，得 $f(114) = 13225$ 。将 $x=115$ 代入，得 $f(115) = 13456$ 。将 $x=116$ 代入，得 $f(116) = 13689$ 。将 $x=117$ 代入，得 $f(117) = 13924$ 。将 $x=118$ 代入，得 $f(118) = 14161$ 。将 $x=119$ 代入，得 $f(119) = 14400$ 。将 $x=120$ 代入，得 $f(120) = 14641$ 。将 $x=121$ 代入，得 $f(121) = 14884$ 。将 $x=122$ 代入，得 $f(122) = 15129$ 。将 $x=123$ 代入，得 $f(123) = 15376$ 。将 $x=124$ 代入，得 $f(124) = 15625$ 。将 $x=125$ 代入，得 $f(125) = 15876$ 。将 $x=126$ 代入，得 $f(126) = 16129$ 。将 $x=127$ 代入，得 $f(127) = 16384$ 。将 $x=128$ 代入，得 $f(128) = 16641$ 。将 $x=129$ 代入，得 $f(129) = 16900$ 。将 $x=130$ 代入，得 $f(130) = 17161$ 。将 $x=131$ 代入，得 $f(131) = 17424$ 。将 $x=132$ 代入，得 $f(132) = 17689$ 。将 $x=133$ 代入，得 $f(133) = 17956$ 。将 $x=134$ 代入，得 $f(134) = 18225$ 。将 $x=135$ 代入，得 $f(135) = 18496$ 。将 $x=136$ 代入，得 $f(136) = 18769$ 。将 $x=137$ 代入，得 $f(137) = 19044$ 。将 $x=138$ 代入，得 $f(138) = 19321$ 。将 $x=139$ 代入，得 $f(139) = 19600$ 。将 $x=140$ 代入，得 $f(140) = 19881$ 。将 $x=141$ 代入，得 $f(141) = 20164$ 。将 $x=142$ 代入，得 $f(142) = 20449$ 。将 $x=143$ 代入，得 $f(143) = 20736$ 。将 $x=144$ 代入，得 $f(144) = 21025$ 。将 $x=145$ 代入，得 $f(145) = 21316$ 。将 $x=146$ 代入，得 $f(146) = 21609$ 。将 $x=147$ 代入，得 $f(147) = 21904$ 。将 $x=148$ 代入，得 $f(148) = 22201$ 。将 $x=149$ 代入，得 $f(149) = 22500$ 。将 $x=150$ 代入，得 $f(150) = 22801$ 。将 $x=151$ 代入，得 $f(151) = 23104$ 。将 $x=152$ 代入，得 $f(152) = 23409$ 。将 $x=153$ 代入，得 $f(153) = 23716$ 。将 $x=154$ 代入，得 $f(154) = 24025$ 。将 $x=155$ 代入，得 $f(155) = 24336$ 。将 $x=156$ 代入，得 $f(156) = 24649$ 。将 $x=157$ 代入，得 $f(157) = 24964$ 。将 $x=158$ 代入，得 $f(158) = 25281$ 。将 $x=159$ 代入，得 $f(159) = 25600$ 。将 $x=160$ 代入，得 $f(160) = 25921$ 。将 $x=161$ 代入，得 $f(161) = 26244$ 。将 $x=162$ 代入，得 $f(162) = 26569$ 。将 $x=163$ 代入，得 $f(163) = 26896$ 。将 $x=164$ 代入，得 $f(164) = 27225$ 。将 $x=165$ 代入，得 $f(165) = 27556$ 。将 $x=166$ 代入，得 $f(166) = 27889$ 。将 $x=167$ 代入，得 $f(167) = 28224$ 。将 $x=168$ 代入，得 $f(168) = 28561$ 。将 $x=169$ 代入，得 $f(169) = 28900$ 。将 $x=170$ 代入，得 $f(170) = 29241$ 。将 $x=171$ 代入，得 $f(171) = 29584$ 。将 $x=172$ 代入，得 $f(172) = 29929$ 。将 $x=173$ 代入，得 $f(173) = 30276$ 。将 $x=174$ 代入，得 $f(174) = 30625$ 。将 $x=175$ 代入，得 $f(175) = 30976$ 。将 $x=176$ 代入，得 $f(176) = 31329$ 。将 $x=177$ 代入，得 $f(177) = 31684$ 。将 $x=178$ 代入，得 $f(178) = 32041$ 。将 $x=179$ 代入，得 $f(179) = 32400$ 。将 $x=180$ 代入，得 $f(180) = 32761$ 。将 $x=181$ 代入，得 $f(181) = 33124$ 。将 $x=182$ 代入，得 $f(182) = 33489$ 。将 $x=183$ 代入，得 $f(183) = 33856$ 。将 $x=184$ 代入，得 $f(184) = 34225$ 。将 $x=185$ 代入，得 $f(185) = 34596$ 。将 $x=186$ 代入，得 $f(186) = 34969$ 。将 $x=187$ 代入，得 $f(187) = 35344$ 。将 $x=188$ 代入，得 $f(188) = 35721$ 。将 $x=189$ 代入，得 $f(189) = 36100$ 。将 $x=190$ 代入，得 $f(190) = 36481$ 。将 $x=191$ 代入，得 $f(191) = 36864$ 。将 $x=192$ 代入，得 $f(192) = 37249$ 。将 $x=193$ 代入，得 $f(193) = 37636$ 。将 $x=194$ 代入，得 $f(194) = 38025$ 。将 $x=195$ 代入，得 $f(195) = 38416$ 。将 $x=196$ 代入，得 $f(196) = 38809$ 。将 $x=197$ 代入，得 $f(197) = 39204$ 。将 $x=198$ 代入，得 $f(198) = 39601$ 。将 $x=199$ 代入，得 $f(199) = 40000$ 。将 $x=200$ 代入，得 $f(200) = 40401$ 。将 $x=201$ 代入，得 $f(201) = 40804$ 。将 $x=202$ 代入，得 $f(202) = 41209$ 。将 $x=203$ 代入，得 $f(203) = 41616$ 。将 $x=204$ 代入，得 $f(204) = 42025$ 。将 $x=205$ 代入，得 $f(205) = 42436$ 。将 $x=206$ 代入，得 $f(206) = 42849$ 。将 $x=207$ 代入，得 $f(207) = 43264$ 。将 $x=208$ 代入，得 $f(208) = 43681$ 。将 $x=209$ 代入，得 $f(209) = 44100$ 。将 $x=210$ 代入，得 $f(210) = 44521$ 。将 $x=211$ 代入，得 $f(211) = 44944$ 。将 $x=212$ 代入，得 $f(212) = 45369$ 。将 $x=213$ 代入，得 $f(213) = 45796$ 。将 $x=214$ 代入，得 $f(214) = 46225$ 。将 $x=215$ 代入，得 $f(215) = 46656$ 。将 $x=216$ 代入，得 $f(216) = 47089$ 。将 $x=217$ 代入，得 $f(217) = 47524$ 。将 $x=218$ 代入，得 $f(218) = 47961$ 。将 $x=219$ 代入，得 $f(219) = 48400$ 。将 $x=220$ 代入，得 $f(220) = 48841$ 。将 $x=221$ 代入，得 $f(221) = 49284$ 。将 $x=222$ 代入，得 $f(222) = 49729$ 。将 $x=223$ 代入，得 $f(223) = 50176$ 。将 $x=224$ 代入，得 $f(224) = 50625$ 。将 $x=225$ 代入，得 $f(225) = 51076$ 。将 $x=226$ 代入，得 $f(226) = 51529$ 。将 $x=227$ 代入，得 $f(227) = 51984$ 。将 $x=228$ 代入，得 $f(228) = 52441$ 。将 $x=229$ 代入，得 $f(229) = 52900$ 。将 $x=230$ 代入，得 $f(230) = 53361$ 。将 $x=231$ 代入，得 $f(231) = 53824$ 。将 $x=232$ 代入，得 $f(232) = 54289$ 。将 $x=233$ 代入，得 $f(233) = 54756$ 。将 $x=234$ 代入，得 $f(234) = 55225$ 。将 $x=235$ 代入，得 $f(235) = 55696$ 。将 $x=236$ 代入，得 $f(236) = 56169$ 。将 $x=237$ 代入，得 $f(237) = 56644$ 。将 $x=238$ 代入，得 $f(238) = 57121$ 。将 $x=239$ 代入，得 $f(239) = 57600$ 。将 $x=240$ 代入，得 $f(240) = 58081$ 。将 $x=241$ 代入，得 $f(241) = 58564$ 。将 $x=242$ 代入，得 $f(242) = 59049$ 。将 $x=243$ 代入，得 $f(243) = 59536$ 。将 $x=244$ 代入，得 $f(244) = 60025$ 。将 $x=245$ 代入，得 $f(245) = 60516$ 。将 $x=246$ 代入，得 $f(246) = 61009$ 。将 $x=247$ 代入，得 $f(247) = 61504$ 。将 $x=248$ 代入，得 $f(248) = 62001$ 。将 $x=249$ 代入，得 $f(249) = 62500$ 。将 $x=250$ 代入，得 $f(250) = 63001$ 。将 $x=251$ 代入，得 $f(251) = 63504$ 。将 $x=252$ 代入，得 $f(252) = 64009$ 。将 $x=253$ 代入，得 $f(253) = 64516$ 。将 $x=254$ 代入，得 $f(254) = 65025$ 。将 $x=255$ 代入，得 $f(255) = 65536$ 。将 $x=256$ 代入，得 $f(256) = 66049$ 。将 $x=257$ 代入，得 $f(257) = 66564$ 。将 $x=258$ 代入，得 $f(258) = 67081$ 。将 $x=259$ 代入，得 $f(259) = 67600$ 。将 $x=260$ 代入，得 $f(260) = 68121$ 。将 $x=261$ 代入，得 $f(261) = 68644$ 。将 $x=262$ 代入，得 $f(262) = 69169$ 。将 $x=263$ 代入，得 $f(263) = 69696$ 。将 $x=264$ 代入，得 $f(264) = 70225$ 。将 $x=265$ 代入，得 $f(265) = 70756$ 。将 $x=266$ 代入，得 $f(266) = 71289$ 。将 $x=267$ 代入，得 $f(267) = 71824$ 。将 $x=268$ 代入，得 $f(268) = 72361$ 。将 $x=269$ 代入，得 $f(269) = 72900$ 。将 $x=270$ 代入，得 $f(270) = 73441$ 。将 $x=271$ 代入，得 $f(271) = 73984$ 。将 $x=272$ 代入，得 $f(272) = 74529$ 。将 $x=273$ 代入，得 $f(273) = 75076$ 。将 $x=274$ 代入，得 $f(274) = 75625$ 。将 $x=275$ 代入，得 $f(275) = 76176$ 。将 $x=276$ 代入，得 $f(276) = 76729$ 。将 $x=277$ 代入，得 $f(277) = 77284$ 。将 $x=278$ 代入，得 $f(278) = 77841$ 。将 $x=279$ 代入，得 $f(279) = 78400$ 。将 $x=280$ 代入，得 $f(280) = 78961$ 。将 $x=281$ 代入，得 $f(281) = 79524$ 。将 $x=282$ 代入，得 $f(282) = 80089$ 。将 $x=283$ 代入，得 $f(283) = 80656$ 。将 $x=284$ 代入，得 $f(284) = 81225$ 。将 $x=285$ 代入，得 $f(285) = 81796$ 。将 $x=286$ 代入，得 $f(286) = 82369$ 。将 $x=287$ 代入，得 $f(287) = 82944$ 。将 $x=288$ 代入，得 $f(288) = 83521$ 。将 $x=289$ 代入，得 $f(289) = 84100$ 。将 $x=290$ 代入，得 $f(290) = 84681$ 。将 $x=291$ 代入，得 $f(291) = 85264$ 。将 $x=292$ 代入，得 $f(292) = 85849$ 。将 $x=293$ 代入，得 $f(293) = 86436$ 。将 $x=294$ 代入，得 $f(294) = 87025$ 。将 $x=295$ 代入，得 $f(295) = 87616$ 。将 $x=296$ 代入，得 $f(296) = 88209$ 。将 $x=297$ 代入，得 $f(297) = 88804$ 。将 $x=298$ 代入，得 $f(298) = 89401$ 。将 $x=299$ 代入，得 $f(299) = 90000$ 。将 $x=300$ 代入，得 $f(300) = 90601$ 。将 $x=301$ 代入，得 $f(301) = 91204$ 。将 $x=302$ 代入，得 $f(302) = 91809$ 。将 $x=303$ 代入，得 $f(303) = 92416$ 。将 $x=304$ 代入，得 $f(304) = 93025$ 。将 $x=305$ 代入，得 $f(305) = 93636$ 。将 $x=306$ 代入，得 $f(306) = 94249$ 。将 $x=307$ 代入，得 $f(307) = 94864$ 。将 $x=308$ 代入，得 $f(308) = 95481$ 。将 $x=309$ 代入，得 $f(309) = 96100$ 。将 $x=310$ 代入，得 $f(310) = 96721$ 。将 $x=311$ 代入，得 $f(311) = 97344$ 。将 $x=312$ 代入，得 $f(312) = 97969$ 。将 $x=313$ 代入，得 $f(313) = 98596$ 。将 $x=314$ 代入，得 $f(314) = 99225$ 。将 $x=315$ 代入，得 $f(315) = 99856$ 。将 $x=316$ 代入，得 $f(316) = 100489$ 。将 $x=317$ 代入，得 $f(317) = 101124$ 。将 $x=318$ 代入，得 $f(318) = 101761$ 。将 $x=319$ 代入，得 $f(319) = 102400$ 。将 $x=320$ 代入，得 $f(320) = 103041$ 。将 $x=321$ 代入，得 $f(321) = 103684$ 。将 $x=322$ 代入，得 $f(322) = 104329$ 。将 $x=323$ 代入，得 $f(323) = 104976$ 。将 $x=324$ 代入，得 $f(324) = 105625$ 。将 $x=325$ 代入，得 $f(325) = 106276$ 。将 $x=326$ 代入，得 $f(326) = 106929$ 。将 $x=327$ 代入，得 $f(327) = 107584$ 。将 $x=328$ 代入，得 $f(328) = 108241$ 。将 $x=329$ 代入，得 $f(329) = 108900$ 。将 $x=330$ 代入，得 $f(330) = 109561$ 。将 $x=331$ 代入，得 $f(331) = 110224$ 。将 $x=332$ 代入，得 $f(332) = 110889$ 。将 $x=333$ 代入，得 $f(333) = 111556$ 。将 $x=334$ 代入，得 $f(334) = 112225$ 。将 $x=335$ 代入，得 $f(335) = 112896$ 。将 $x=336$ 代入，得 $f(336) = 113569$ 。将 $x=337$ 代入，得 $f(337) = 114244$ 。将 $x=338$ 代入，得 $f(338) = 114921$ 。将 $x=339$ 代入，得 $f(339) = 115600$ 。将 $x=340$ 代入，得 $f(340) = 116281$ 。将 $x=341$ 代入，得 $f(341) = 116964$ 。将 $x=342$ 代入，得 $f(342) = 117649$ 。将 $x=343$ 代入，得 $f(343) = 118336$ 。将 $x=344$ 代入，得 $f(344) = 119025$ 。将 $x=345$ 代入，得 $f(345) = 119716$ 。将 $x=346$ 代入，得 $f(346) = 120409$ 。将 $x=347$ 代入，得 $f(347) = 121104$ 。将 $x=348$ 代入，得 $f(348) = 121801$ 。将 $x=349$ 代入，得 $f(349) = 122500$ 。将 $x=350$ 代入，得 $f(350) = 123201$ 。将 $x=351$ 代入，得 $f(351) = 123904$ 。将 $x=352$ 代入，得 $f(352) = 124609$ 。将 $x=353$ 代入，得 $f(353) = 125316$ 。将 $x=354$ 代入，得 $f(354) = 126025$ 。将 $x=355$ 代入，得 $f(355) = 126736$ 。将 $x=356$ 代入，得 $f(356) = 127449$ 。将 $x=357$ 代入，得 $f(357) = 128164$ 。将 $x=358$ 代入，得 $f(358) = 128881$ 。将 $x=359$ 代入，得 $f(359) = 129600$ 。将 $x=360$ 代入，得 $f(360) = 130321$ 。将 $x=361$ 代入，得 $f(361) = 131044$ 。将 $x=362$ 代入，得 $f(362) = 131769$ 。将 $x=363$ 代入，得 $f(363) = 132496$ 。将 $x=364$ 代入，得 $f(364) = 133225$ 。将 $x=36$




Figure 3 Emission spectrum from a mixture of EuCl_3 and TbCl_3 with a 10 μs delay time




Figure 4 Emission spectrum from a mixture of EuCl_3 and TbCl_3 with a 40 μs delay time




Figure 5 3-D plot of phosphorescence spectra for a mixture of terbium and uranyl with varying delay times



[Faint, illegible text, likely bleed-through from the reverse side of the page.]

7. Connect the rectangular end of the 8-ft. round power cable with the 9-pin D connector to the REFERENCE PREAMP POWER connector phosphorimeter control module. Connect the other end of the PREAMP POWER connector on the reference amplifier module.

Connecting The Sample Channel

8. Disconnect the BNC connector from the DM302. Install the BNC "T" connector. Reattach BNC connector to one side of the T. Reattach BNC connector end of 8' cable to other side of "T" and 15-pin end to sample power connector on the control module.

Connecting The Phosphorimeter To The DM3000

9. Connect one end of the RS232 cable (25 pin 'D') to the RS232 communication connector on the Phosphorimeter Control Module. Connect the other end of the RS232 cable to COM 1 (or COM 2) on the D3000. (Note 1) it may be necessary to use a 9-pin to 25-pin adapter cable 2) refer to Appendix B if it is necessary to use COM 2 instead of COM 1).

3.0 OPERATION

3.1 Start-up Procedure

1. Before turning on the DM3000, turn on the illuminator, the control module and all peripheral devices. Then on the DM3000.
2. When the phosphorimeter control module is first turned on, it will display its firmware revision level in the front panel LED's for a few seconds. All LED'S will then go out and the phosphorimeter will wait for the DM3000 to initiate self-testing. Hardware errors will not be reported (Appendix A) until the DM3000 attempts to initialize the RS232 link.
3. To use the DM3000 with a cw source without disconnecting the phosphorimeter control module from the DM302, the phosphorimeter control module must be turned on.
4. To activate the phosphorimeter, go to the Set Defaults, Optional Hardware menu. Use the space bar to show Yes on the Phosphorimeter line, then press the Enter key to actually make the selection. Move the scroll bar to Go and press then Enter key to initiate the installation. (Note: to always have the DM3000 start with Phosphorimeter enabled, move the scroll bar to Save Defaults (on the Set Defaults menu) the press Enter.

1. The first part of the document discusses the importance of maintaining accurate records of all transactions and activities. It emphasizes the need for transparency and accountability in financial reporting.

2. The second part of the document outlines the various methods and techniques used to collect and analyze data. It includes a detailed description of the experimental procedures and the statistical analysis performed.

3. The third part of the document presents the results of the study. It includes a series of tables and graphs that illustrate the findings. The results show a clear correlation between the variables studied, supporting the hypothesis.

4. The fourth part of the document discusses the implications of the findings. It highlights the potential applications of the research in the field of finance and economics. The authors conclude that the study provides valuable insights into the relationship between the variables.

5. The fifth part of the document provides a summary of the key findings and conclusions. It reiterates the importance of the research and the need for further studies in this area. The authors express their gratitude to the funding agency and the participants.

6. The sixth part of the document includes a list of references and a bibliography. It cites the works of other researchers in the field, providing a context for the current study.

7. The seventh part of the document contains a list of appendices and supplementary materials. These include additional data, charts, and detailed descriptions of the experimental setup. The authors provide contact information for further inquiries.

The SPEX phosphorimeter also has reference channel to monitor the output of the pulsed lamp. When used in the ratio mode, this reference signal will cancel out variations in lamp intensity to produce corrected excitation spectra or to compensate for time-dependent variations.

2.0 INSTALLATION

The phosphorimeter has three main modules: the illuminator, the phosphorimeter control module and the reference amplifier module.

CAUTION: Before installing these units, turn off the power to the DM3000 and all peripherals.

The pulsed-lamp illuminator contains the xenon flash lamp and the lamp power supply. It is installed in place of the normal continuous-wave illuminator at the entrance slit of the excitation spectrometer.

Connecting The Illuminator

1. Slide the xenon lamp illuminator housing aside.
2. Carefully slide the illuminator opening over the entrance slit of the excitation spectrometer. Adjust the feet on the bottom of the illuminator so that no strain is exerted on the slit. Replace the cover on the flash lamp.
3. Connect the round cable with the 9-pin D connector between the LAMP TRIGGER connector on the phosphorimeter control module and LAMP TRIGGER connector on the illuminator module.
4. Connect the illuminator power cable to the proper line voltage.

Connecting The Reference Channel

Note: For the reference channel to work properly, a cuvette filled with rhodamine (reference quantum counter) must be inserted into the reference holder in the sample module (see the Fluorolog manual for information on how to prepare the rhodamine solution)

5. Disconnect the BNC (to the reference detector) from the INPUT connector of the DM303. Connect this cable to the IN connector on the reference amplifier module.
6. Connect one end of an 8-ft BNC cable to the OUT connector on the reference amplifier module. Connect the other end of the REFERENCE INPUT connector on the phosphorimeter control module.

1. The first part of the document discusses the importance of maintaining accurate records of all transactions and activities. It emphasizes the need for transparency and accountability in financial reporting.


2. The second part of the document outlines the various methods and techniques used to collect and analyze data. It includes a detailed description of the experimental procedures and the statistical analysis performed.

3. The third part of the document presents the results of the study. It includes a series of tables and graphs that illustrate the findings of the research. The data shows a clear trend in the relationship between the variables studied.

4. The fourth part of the document discusses the implications of the findings. It highlights the potential applications of the research in various fields and the need for further investigation in this area.

5. The fifth part of the document provides a conclusion and summarizes the key points of the study. It reiterates the importance of the research and the need for continued efforts in this field.

6. The sixth part of the document includes a list of references and a bibliography. It cites the various sources used in the research and provides a comprehensive overview of the literature in this area.



Cable Connections for 1934D Phosphorimeter

TM9

3.2 REFERENCE DETECTOR

1. The reference detector, located in the sample module, uses the Fluorolog 2 reference quantum counter. The DC Data Acquisition Module (DM303) is replaced by an analog detector module which contains a reference amplifier module. The output of this module is scaled in the phosphorimeter control module, then reported to the DM3000.

The reference detector monitors the excitation radiation. It is integrated for fixed window independent of parameters governing acquisition of the emission (sample) signal. This window is opened before the lamp flash is triggered and remains open for a time sufficient to include the entire flash. Timing of the window raises less than one micro second from flash to flash.

2. If the high voltage applied to the reference detector is too large, both the photomultiplier and reference amplifier can be saturated, leading to an apparent 0 signal. To avoid this problem, adjust the reference high voltage (HV2) while monitoring Reference signal level in the status window (section 3.4)

Set the excitation monochromator to the wavelength of maximum lamp intensity expected for excitation and synchronous scans (usually about 390nm), or to the excitation position for emission scans. (To set the excitation monochromator, enter the desired position in monochromator 1: Move to X, then press the Enter key to start monochromator movement.)

3. Once the monochromator is positioned, set HV2 to 150 to 200 volts, then slowly increase it until the Reference Input Value is from 12,000 to 14,000. The maximum possible signal on the reference channel is 16,128, so a value of 12-14,000 should maximize dynamic range of the result. (After determining an optimum high voltage, remember to make sure that this value is also entered into the Define Experiment Windows).

3.3 Controlling the Phosphorimeter from the DM3000

Operation of a DM3000F with a phosphorimeter is similar to operation with a continuous lamp. Excitation, emission, and synchronous scans are all available. The only differences are additional parameters required to control the phosphorimeter.

1. The first part of the report deals with the general situation of the country and the position of the various groups. It is a very good introduction to the subject and gives a clear picture of the situation as it is at present.

2. The second part of the report deals with the various groups and their position. It is a very good introduction to the subject and gives a clear picture of the situation as it is at present.

3. The third part of the report deals with the various groups and their position. It is a very good introduction to the subject and gives a clear picture of the situation as it is at present.

4. The fourth part of the report deals with the various groups and their position. It is a very good introduction to the subject and gives a clear picture of the situation as it is at present.

5. The fifth part of the report deals with the various groups and their position. It is a very good introduction to the subject and gives a clear picture of the situation as it is at present.

6. The sixth part of the report deals with the various groups and their position. It is a very good introduction to the subject and gives a clear picture of the situation as it is at present.

3.3a Acquisition Modes

Available acquisition modes can be viewed by pressing the Help (F1) key while on the Acquisition Mode line in a Data Acquisition Parameters window. Normally, the choices will be:

s	- sample signal only
r	- reference signal only
s & r	- both sample and reference signals; two result files are generated
s/r	- ratio of sample to reference signal

Modes can be selected from the Help window by pressing Enter when the scroll bar is on an appropriate mode, from the Acquisition Mode line by pressing spacebar until the desired mode is displayed then pressing Enter, or by directly typing in the desired mode.

NOTE: The phosphorimeter will not operate in a chopped mode (F222 systems).

3.3b Delay After Flash

The Delay After Flash parameter sets the time, in milliseconds, between the start of the lamp flash and the beginning of data acquisition (opening of the SAMPLE WINDOW). This delay can be between 0 and 10,000 milliseconds, in increments of 0.001 millisecond. The sample delay from one flash to the next will vary by less than 0.001 millisecond.

The Delay After Flash should be long enough so that fluorescence emission and lamp decay have occurred. The result spectrum will represent the phosphorescence spectrum.

3.3c Sample Window

The SAMPLE WINDOW parameter sets the duration of signal acquisition in milliseconds. The sample window can be set between 0.01 and 10,000 milliseconds, in increments of .001 millisecond.

The SAMPLE WINDOW opens when the Delay After Flash ends (see Figure 1). When the SAMPLE WINDOW is open, the signal that enters the control module is counted and integrated for subsequent output to the DM3000. After the preset time, the SAMPLE WINDOW closes and the signal is ignored.

If you want to record fluorescence emissions, set Delay After Flash to zero.

1. The first part of the report deals with the general situation of the country and the position of the various groups.

2. The second part of the report deals with the economic situation and the measures taken to improve it.

3. The third part of the report deals with the social situation and the measures taken to improve it.

4. The fourth part of the report deals with the political situation and the measures taken to improve it.

5. The fifth part of the report deals with the cultural situation and the measures taken to improve it.

6. The sixth part of the report deals with the environmental situation and the measures taken to improve it.

7. The seventh part of the report deals with the international situation and the measures taken to improve it.

8. The eighth part of the report deals with the future prospects of the country and the measures taken to improve it.

9. The ninth part of the report deals with the conclusion of the report and the measures taken to improve it.

3.3d Number of Flashes

The Number of Flashes parameter sets the number of flashes (lamp pulses) that will contribute to each data point. Number of Flashes can be set from 1 to 999. The signal collected during the sampling time is integrated over the number of lamp pulse before the data is passed to the DM3000.

3.3e Flash Rate

Flash Rate sets the repetition rate of the pulse lamp. Flash rate must be between 0.05 and 33 flashes/second, and in addition, must be slow enough that the Sample Window will be completed before another flash begins. Alternatively, Flash Rate can be set to 0 and the DM3000 will determine the maximum allowable rate. Accuracy in timing of Flash Rate is +/- 1 millisecond.

NOTE: Maximum flash rate can be calculated as:

$$\text{Rate} < \frac{1000}{\text{Delay After Flash} + \text{Sample Window} + 10}$$

where Delay After Flash and Sample Window are in milliseconds

3.4 System Status

The effect of possible phosphorimeter parameters can be viewed in the System Status Window (F9) without having to run a scan. All results displayed are for single flashes of the lamp. Delay After Flash and Sample Window have the same range as described in Section 3.3. Time Between Flashes is controlled by the system clock in the DM3000 rather than the more accurate clock in the phosphorimeter control module, and is accurate to 0.1 second. Allowable values for Time Between Flashes vary from 0.1 to 20 seconds in increments of 0.1 seconds. Time Between Flashes must also be greater than Delay After Flash plus Flash Window (in seconds).

System Status is particularly useful for determining the best setting for reference detector high voltage (HV2).

3.5 Emission Scans and Time Discrimination

3.5a Scanning an Emission Spectrum

To scan an emission spectrum with the phosphorimeter, use the following procedure.

1. The first part of the report deals with the general situation of the country and the position of the various groups of the population. It is a very interesting and informative study of the social and economic conditions of the country.

2. The second part of the report deals with the political situation of the country and the position of the various groups of the population. It is a very interesting and informative study of the political conditions of the country.

3. The third part of the report deals with the economic situation of the country and the position of the various groups of the population. It is a very interesting and informative study of the economic conditions of the country.

4. The fourth part of the report deals with the cultural situation of the country and the position of the various groups of the population. It is a very interesting and informative study of the cultural conditions of the country.

5. The fifth part of the report deals with the social situation of the country and the position of the various groups of the population. It is a very interesting and informative study of the social conditions of the country.

6. The sixth part of the report deals with the legal situation of the country and the position of the various groups of the population. It is a very interesting and informative study of the legal conditions of the country.

7. The seventh part of the report deals with the administrative situation of the country and the position of the various groups of the population. It is a very interesting and informative study of the administrative conditions of the country.

8. The eighth part of the report deals with the military situation of the country and the position of the various groups of the population. It is a very interesting and informative study of the military conditions of the country.

9. The ninth part of the report deals with the foreign relations of the country and the position of the various groups of the population. It is a very interesting and informative study of the foreign relations of the country.

10. The tenth part of the report deals with the conclusion of the report and the position of the various groups of the population. It is a very interesting and informative study of the conclusion of the report.

11. The eleventh part of the report deals with the appendix of the report and the position of the various groups of the population. It is a very interesting and informative study of the appendix of the report.

1. Go to the Define Experiment, Emission Scan, Data Acquisition Parameters Window.
2. Enter the number of times the scan is to be repeated and averaged (Number of Scans), the Start, End, and Increment wavelengths.
3. Enter Delay After Flash. The type of results you can expect for varying delay are shown in Figure 3 and 4. In general, when viewing phosphorescence, a delay of at least 0.05 milliseconds (50 microseconds) will discriminate against lamp pulse decay. As shown in Figure 6, the main component of the lamp pulse has a lifetime in the 3 microsecond range, but there is a long tail on the curve which does not fall below 1% of the peak height until after 20 microseconds. Setting a delay of 50 microseconds will discriminate against lamp pulse decay, fluorescence and stray light.




Figure 6 Pulse-lamp output

4. Set the SAMPLE WINDOW to cover about 5 to 10 lifetimes of the luminescence you wish to record. (An alternative technique is to make the window a fraction (e.g. 1/2 or 1/3) of the lifetime and vary the delay to take individual slices. If the lifetime is unknown, start with a small SAMPLE WINDOW, then increase it until good spectra are obtained.

If the SAMPLE WINDOW is too long, the background signal will continue to be recorded after the luminescence has decayed. This will add unnecessary noise to the signal. If the sample window is too short, second and third lifetime components may be missed. To obtain the best spectra when two or more long-lived species are present, it may be necessary to vary the Delay After Flash and the SAMPLE WINDOW.

The first part of the paper is devoted to a review of the literature on the topic. It is found that there is a general consensus that the model is valid for the range of parameters considered. However, there is some disagreement as to the exact form of the model.

In the second part of the paper, the model is applied to the case of a specific system. It is found that the model predicts the behavior of the system very well. This is a strong indication that the model is valid for this system.



The third part of the paper discusses the implications of the results. It is concluded that the model is a good approximation of the system's behavior. This is important because it allows us to predict the system's behavior without having to solve the full equations of motion.

Finally, the paper discusses some of the limitations of the model. It is noted that the model is only valid for a certain range of parameters. Outside of this range, the model may not accurately describe the system's behavior.

5. Enter the Flash Rate. In general the value that will give the maximum repetition rate to accumulate data quickly can be calculated using the formula given in section 3.3e. Alternatively, Flash Rate can be set to 0 and the DM3000 will calculate the maximum rate.

NOTE: Flash Rate must be 0 or between 0.05 and 33 Hz.

6. Enter the Number of Flashes you wish to record. The more flashes accumulated, the more events that will contribute to your data and therefore increase the signal-to-noise ratio. The primary effect is similar to increasing the integration time under normal operation.

7. Enter Excitation Position and select an Acquisition Mode.

NOTES: a) All fields have default values, so you only need only enter those values which will be changed.

b) Auto zero does not work in release 2.40.

8. Press the Esc key to return to Emission Scan menu. Select the High Voltage menu to check that the high voltages specified match those previously determined (section 3.2.2).
9. Press the Esc key to return to Emission Scan menu. Select Specify Axis Titles. The titles shown will appear on any plot or graph of the data.
10. Press the Esc key to return to Emission Scan menu. Select Run Experiment. Enter all File Name for the result of the experiment. Optionally, enter a File Title to help describe the file (for viewing in directory displays).

NOTE: Step 8 and 9 are optional.

11. Select GO and press the Enter key to begin acquiring data. Data are automatically displayed as the scan progresses.

1. The first part of the report is a general introduction to the subject of the study. It discusses the importance of the study and the objectives of the research.

2. The second part of the report is a detailed description of the methodology used in the study. It includes information about the sample, the data collection methods, and the statistical analysis.

3. The third part of the report is a presentation of the results of the study. It includes a summary of the findings and a discussion of their implications. The results are presented in a clear and concise manner, using tables and figures where appropriate.

4. The fourth part of the report is a conclusion and a list of references. The conclusion summarizes the main findings of the study and provides recommendations for future research. The references list the sources of information used in the study.

5. The fifth part of the report is an appendix containing additional information that is not included in the main body of the report. This may include raw data, detailed calculations, or other supporting materials.

6. The sixth part of the report is a bibliography listing the sources of information used in the study.

7. The seventh part of the report is a list of figures and tables. Each figure and table is accompanied by a brief description of its content and a reference to the text where it is discussed.

8. The eighth part of the report is a list of abbreviations and symbols used in the study. This helps to ensure that the reader can understand the terminology used in the report.

9. The ninth part of the report is a list of acknowledgments. This is a place to thank the people and organizations that have helped to make the study possible. It may include the sponsor of the study, the advisor, and other individuals who have provided support or assistance.

10. The tenth part of the report is a list of appendices. This is a place to list the additional information that is included in the report, such as raw data, detailed calculations, or other supporting materials.

11. The eleventh part of the report is a list of references. This is a list of the sources of information used in the study, including books, articles, and other documents. It is important to provide a complete and accurate list of references to allow the reader to locate the sources of information used in the study.

3.5b An Example: Obtaining the Emission Spectrum of a Sample

To familiarize yourself with the operation of the phosphorimeter, run an emission scan on a sample of 100 ppm uranyl nitrate in 10% phosphoric acid (H_3PO_4) by volume. Set the excitation wavelength to 295 nm by entering 295 for the excitation position parameter in Data Acquisition Parameters window. Set start to 450 nm, End to 650 nm, and Increment to 1 nm. Suggested slit settings are 1.5 nm for the emission monochromator and 5mm for the excitation monochromator.

Set:

Delay After Flash to 0.05 msec.,
Sample Window to 0.1 msec.,
Flash Rate to 33, and
Number of Flashes to 100.

3.6 Excitation Scans

Scanning an excitation spectrum is analogous to scanning an emission spectrum (Section 3.5), but now the emission spectrometer is kept stationary on an emitting peak while the excitation spectrometer is scanned. The reference detector enables the excitation scan to be corrected automatically for the wavelength-dependent variations of the pulse lamp and the gratings in the excitation spectrometer.

3.7 Synchronous Scans

Scanning a synchronous spectrum is analogous to scanning excitation or emission spectra (Sections 3.5 and 3.6) but now both spectrometers are scanned with a constant offset between them. Start and End parameters refer to the excitation spectrometer. The reference detector enables the synchronous scan to be automatically corrected for wavelength dependent variations of the pulse lamp and grating(s) in the excitation spectrometer.

4.0 DETERMINATION OF ERROR CODES

4.1 Error Reporting

Phosphorimeter errors are reported by both the DM3000 and the phosphorimeter control module. Errors reported by the DM3000 are a limited set and are explained by the text in the error window which appears. When an error window appears, check the LED's on the front of the phosphorimeter control module to determine the phosphorimeter's error state.

NOTE: If an error occurs during a scan, the DM3000 will attempt to continue the scan by restarting the phosphorimeter.

1. The first part of the report discusses the general situation of the company and the results of the work done during the year. It also mentions the financial position and the state of the company's affairs.

2. The second part of the report deals with the specific results of the work done during the year. It mentions the progress made in the various departments and the results of the different projects.

3. The third part of the report discusses the financial position of the company and the state of the company's affairs. It mentions the income and expenses for the year and the balance sheet at the end of the year.

4. The fourth part of the report discusses the results of the work done during the year. It mentions the progress made in the various departments and the results of the different projects.

5. The fifth part of the report discusses the financial position of the company and the state of the company's affairs. It mentions the income and expenses for the year and the balance sheet at the end of the year.

6. The sixth part of the report discusses the results of the work done during the year. It mentions the progress made in the various departments and the results of the different projects.

7. The seventh part of the report discusses the financial position of the company and the state of the company's affairs. It mentions the income and expenses for the year and the balance sheet at the end of the year.

4.2 Interpretation of Display

The eight LED displays on the front panel of the phosphorimeter control module will light in a particular combination when hardware, software, or operator errors are detected. The source of the error is identified in hexadecimal code. The eight LEDs are divided into two sections, each section revealing one digit of the code. Individual codes are decoded as follows:

LED Display	o	o	o	o		o	o	o	o
light is on	8	4	2	1		8	4	2	1
light is off	0	0	0	0		0	0	0	0

To determine the first digit of the error code, add the values of the lit LEDs on the left side. Find the corresponding hexadecimal equivalent from Table 1. Follow the same procedure with the four rightmost LEDs to determine the second digit.

Example:

(*=lit, o=not lit)	*	*	o	*		*	o	*	o
Decimal values	8	4	0	1		8	0	2	0
	-----					-----			
Total	13					10			
Hexadecimal equivalents (from Table 1)	D					A			

Decimal	Hex	Decimal	Hex
0	0	8	8
1	1	9	9
2	2	10	A
3	3	11	B
4	4	12	C
5	5	13	D
6	6	14	E
7	7	15	F

Table 1 - Hexadecimal equivalents of decimal numbers

The errors that the codes represent are listed in Appendix A. Most of the errors are hardware errors (for example, a board must be replaced) that cannot be corrected by the operator.

1. The first part of the report deals with the general situation of the country. It is a very interesting and informative study of the country's development. The author has done a great deal of research and has gathered a wealth of material. The report is well written and is easy to read. It is a valuable contribution to the study of the country's development.

2. The second part of the report deals with the economic situation of the country. It is a very interesting and informative study of the country's economic development. The author has done a great deal of research and has gathered a wealth of material. The report is well written and is easy to read. It is a valuable contribution to the study of the country's economic development.

3. The third part of the report deals with the social situation of the country. It is a very interesting and informative study of the country's social development. The author has done a great deal of research and has gathered a wealth of material. The report is well written and is easy to read. It is a valuable contribution to the study of the country's social development.

4. The fourth part of the report deals with the political situation of the country. It is a very interesting and informative study of the country's political development. The author has done a great deal of research and has gathered a wealth of material. The report is well written and is easy to read. It is a valuable contribution to the study of the country's political development.

After determining the error code, call the SPEX service department for assistance if the problem can be corrected.

The following error codes indicate problems which can sometimes be corrected by the operator.

Error Code	Problem	Possible Remedy
40	No sample board	Check that a sample board is present and properly installed.
60	No reference board	Check that a reference board is present properly installed.
92	No flash detected	Check that the illuminator is plugged in and turned on. Check that the illuminator cover is on tight (there is a safety interlock)

5.0 TROUBLESHOOTING GUIDE

Problem	Possible Remedy
There is no signal or the signal level is low	<p>Check that the High Voltage is on and at the proper level (when the high voltage is on the background of the HV (F10) key is on.)</p> <p>Check that the cables are connected properly.</p> <p>Check that the illuminator is aligned.</p>
The reference signal saturates the reference detector	Check that the cuvette with the rhodamine B solution is present

6.0 FLASH LAMP REPLACEMENT

The flash lamp in the phosphorimeter illuminator typically has a half-intensity life of more than 10 million flashes. When it must be replaced, follow the procedure outlined below.

The following information is being furnished to you for your information and is not to be used for any other purpose.

It is requested that you keep this information confidential and not disclose it to any other person.

The information is being furnished to you for your information and is not to be used for any other purpose.

It is requested that you keep this information confidential and not disclose it to any other person.

The information is being furnished to you for your information and is not to be used for any other purpose.

IMPORTANT: THE ILLUMINATOR CONTAINS SEVERAL HIGH VOLTAGE SOURCES.

TO AVOID POSSIBLE FATAL SHOCKS, THE FOLLOWING SAFETY PRECAUTIONS SHOULD BE OBSERVED.

BEFORE REMOVING THE LAMP COVER, UNPLUG THE POWER CORD AND WAIT FOR AT LEAST ONE MINUTE WHILE THE INTERNAL CAPACITORS DISCHARGE. FAILURE TO DO THIS COULD RESULT IN SEVERE INJURY CAUSED BY ELECTRIC SHOCKS.

NEVER OPERATE THE LAMP WITH THE COVER REMOVED.

6.1 Lamp Replacement Procedure

Only handle the lamp when wearing clean gloves or using a paper towel. The skin oils from your hands can degrade the UV performance of the lamp and shorten the life of the lamp.

1. Unplug the power cable from the illuminator and wait at least one minute.
2. Remove the screws holding the lamp cover in place, then remove the lamp cover.
3. Remove the lamp from the socket by pulling upward with a steady motion.
4. Align the pins of the new lamp with the lamp socket then push the lamp firmly into place. If the lamp tilts, realign the pins.
5. Replace the cover on the lamp. Never operate it without the cover.
6. For optimum results, you should check the alignment and focus of the the lamp.

6.2 Lamp Alignment and Focusing Procedure

To align and focus the lamp, you will need a flathead screwdriver and a focusing target. You can make the target by tracing the outlines in Figure 7. After cutting them out, roll the rectangular one (A) from the short end to make a tube with a 3.5-inch diameter. Tape the circular cut-out (B) to one end of this tube. The focusing target is now ready to use.

The objective of the alignment procedure is to get the light image as close as possible to the center of the target (marked with an X). The objective of the focusing procedure is to obtain the sharpest light image possible.

During parts of this procedure, the lamp must flash rapidly.

Select the system Status Window (F9) and set the Time Between Flashes to 0.1

1. Tape the target to the exit port of the illuminator.
2. Turn on the illuminator. Observe where the light image is on the target. If the image is on the X, no alignment is necessary and you can go directly to step 8.
3. Unplug the illuminator, then remove the cover.
4. Loosen the two screws that hold the lamp socket in place so that the socket can move easily in circular motion in the horizontal plane.
5. If the light image was high on the target, move the lamp socket forward toward the exit slit. If the lamp image was low on the target, move the lamp away from the exit slit. Replace the illuminator cover and plug in the power cord.
6. Turn on the illuminator and observe where the image is located on the target.
7. If the image is centered, proceed to step 8. If the light image is not yet centered, unplug the illuminator and wait at least one minute. Remove the cover. Repeat steps 5 and 6 until the image is centered. Note the sharpness of the image.
8. After the image is centered, turn off the illuminator and wait for one minute. Tighten the two screws that hold the lamp socket in place. Plug in the illuminator.
9. If the image was not sharp, the lamp must be focused. Insert a straight blade screwdriver into the opening in the back panel of the illuminator. Turn the screw to move the lamp up and down until the sharpest image is obtained.

1. The first part of the report deals with the general situation of the country and the position of the various groups of the population. It is a very interesting and informative study of the social and economic conditions of the country.

2. The second part of the report deals with the political situation of the country. It is a very interesting and informative study of the political conditions of the country.

3. The third part of the report deals with the cultural situation of the country. It is a very interesting and informative study of the cultural conditions of the country.

4. The fourth part of the report deals with the economic situation of the country. It is a very interesting and informative study of the economic conditions of the country.

5. The fifth part of the report deals with the social situation of the country. It is a very interesting and informative study of the social conditions of the country.

6. The sixth part of the report deals with the legal situation of the country. It is a very interesting and informative study of the legal conditions of the country.


7. The seventh part of the report deals with the educational situation of the country. It is a very interesting and informative study of the educational conditions of the country.

8. The eighth part of the report deals with the health situation of the country. It is a very interesting and informative study of the health conditions of the country.

9. The ninth part of the report deals with the environmental situation of the country. It is a very interesting and informative study of the environmental conditions of the country.

10. The tenth part of the report deals with the international situation of the country. It is a very interesting and informative study of the international conditions of the country.

A



B




Figure 11
Focusing Target

7.0 SPECIFICATIONS

Lamp - UV xenon flash tube

Flash Rate - 0.05 to 33 flashes/second

Flash Duration - 3 microseconds, typical, measured at 1/2 peak amplitude

Delay After Flash - variable from 0 to 10,000 milliseconds in 1 microsecond increments, 4 significant digits

Sample Window - variable from 0.01 to 10,000 milliseconds in 1 microsecond increments, 4 significant digits.

Flashes/Data Point - variable from 1 to 999 flashes

Power Requirements - 115/230V, 50/60 Hz, 50W

Precision - 32 bits pulse-counting

1947, 1948

1949, 1950, 1951

1952, 1953, 1954

1955, 1956, 1957, 1958, 1959

1960, 1961

1962, 1963, 1964, 1965, 1966

1967, 1968, 1969, 1970, 1971

1972, 1973, 1974, 1975, 1976

1977, 1978, 1979, 1980, 1981

1982, 1983, 1984, 1985, 1986

1987, 1988, 1989, 1990, 1991

1992, 1993, 1994, 1995, 1996

Appendix A - Error Codes

Communications Board Errors

21	Reset is bad	2A	Str relay or sense bad
22	Flops are bad	2B	Select stuck off or busy high short at interface
23	Irq flop on	2C	Busy low short
24	Ready flop on	2D	Str relay off
25	Select flop on	2E	Str relay on
26	Busy flop on	2F	Str sense bad
27	Strobe relay closed or sense shorted low	30	Irq flop off
28	Str relay off or busy low short at interface	31	Ready flop off
29	Set relay off or busy low short at interface	32	Select flop off
		33	Busy flop off

Sample Board Errors

40	No sample board	4D	Test gate bad (no stop)
41	Control regs bad	4E	Irq stuck low
42	Too many irq's	4F	AD chip busy at PU
43	B side won't clear	50	AD chip too busy
44	A side won't clear	51	AD ready too soon
45	Timer gate on	52	AD never ready
46	Test count error	53	AD data non zero high byte
47	Timer out of range	54	AD data not within 1 count
48	A counter flaky	55	AD data not within 1 count
49	B counter flaky	56	HAD chip too busy
4A	Counters won't stop	57	HAD ready too soon
4B	Clear is bad	58	HAD never ready
4C	Timer gate is on	59	AD data not all high

Reference Board Errors

60	No reference board	6D	Test gate bad (can't stop)
61	Control regs bad	6E	Irq stuck low
62	Too many irq's	6F	RAD chip busy at PU
63	B side won't clear	70	AD chip too busy
64	A side won't clear	71	AD ready too soon
65	Timer gate on	72	AD never ready
66	Test count on	73	AD data not zero high byte
67	Test count error	74	AD data not zero low byte
68	A counter flaky	75	AD data not within 1 count
69	B counter flaky	76	HAD chip too busy
6A	Counters won't stop	77	HAD chip ready too soon
6B	Clear bad	78	HAD never ready
6C	Timer gate is on	79	AD data not high

Table 1: Summary of Data

Year	Q1	Q2	Q3	Q4	Annual Total
2010	100	120	150	180	550
2011	110	130	160	190	590
2012	120	140	170	200	630
2013	130	150	180	210	670
2014	140	160	190	220	710
2015	150	170	200	230	750
2016	160	180	210	240	790
2017	170	190	220	250	830
2018	180	200	230	260	870
2019	190	210	240	270	910
2020	200	220	250	280	950

Table 2: Detailed Data

Year	Q1	Q2	Q3	Q4	Annual Total
2010	100	120	150	180	550
2011	110	130	160	190	590
2012	120	140	170	200	630
2013	130	150	180	210	670
2014	140	160	190	220	710
2015	150	170	200	230	750
2016	160	180	210	240	790
2017	170	190	220	250	830
2018	180	200	230	260	870
2019	190	210	240	270	910
2020	200	220	250	280	950

Table 3: Summary of Data

Year	Q1	Q2	Q3	Q4	Annual Total
2010	100	120	150	180	550
2011	110	130	160	190	590
2012	120	140	170	200	630
2013	130	150	180	210	670
2014	140	160	190	220	710
2015	150	170	200	230	750
2016	160	180	210	240	790
2017	170	190	220	250	830
2018	180	200	230	260	870
2019	190	210	240	270	910
2020	200	220	250	280	950

Timer Control Board Errors

80	Error output not clear; circuit board is missing	95	Bin strt ctr3
81	Side 1 of U11 won't clear	96	Bin strt ctr2
82	U6 output won't set	97	Bin strt ctrl
83	Ctr -5 won't clear (get back to zero)	98	E bin ct ctr5
84	Ctr -4 won't clear (get back to zero)	99	E bin ct ctr4
85	Ctr -3 won't clear	9A	E bin ct ctr3
86	Ctr -2 won't clear	9B	E bin ct ctr2
87	Ctr -1 won't clear	9C	E bin ct ctrl
88	Side 1 of U11 output can't be set or clocked	9D	E binary rate
89	Side 2 of U9 won't clear	9E	E 1 MHz rate
8A	Side 2 of U11 won't set	9F	E 100 KHz rate
8B	U6 output -2 not functional	A0	E 10 KHz rate
8C	Samp not clear	A1	E 1 KHz rate
8D	Delay not set	A2	E 100 Hz rate
8E	Samp not set	A3	E out inactive
8F	Samp not clock	A4	Too many flashes
90	Flash not set	A5	9513 busy too long
91	Flash not clear	F7	IRQ full
92	No flash check illuminator power	F8	FIRQ full
93	Bin strt ctr5	F9	RES error
94	Bin strt crt4	FA	Nmi error
		FB	Sw3 error
		FC	Sw2 error
		FD	Sw1 error
		FE	IRQ error
		FF	FIRQ error

1	1994-1995	1994-1995
2	1994-1995	1994-1995
3	1994-1995	1994-1995
4	1994-1995	1994-1995
5	1994-1995	1994-1995
6	1994-1995	1994-1995
7	1994-1995	1994-1995
8	1994-1995	1994-1995
9	1994-1995	1994-1995
10	1994-1995	1994-1995
11	1994-1995	1994-1995
12	1994-1995	1994-1995
13	1994-1995	1994-1995
14	1994-1995	1994-1995
15	1994-1995	1994-1995
16	1994-1995	1994-1995
17	1994-1995	1994-1995
18	1994-1995	1994-1995
19	1994-1995	1994-1995
20	1994-1995	1994-1995
21	1994-1995	1994-1995
22	1994-1995	1994-1995
23	1994-1995	1994-1995
24	1994-1995	1994-1995
25	1994-1995	1994-1995
26	1994-1995	1994-1995
27	1994-1995	1994-1995
28	1994-1995	1994-1995
29	1994-1995	1994-1995
30	1994-1995	1994-1995
31	1994-1995	1994-1995
32	1994-1995	1994-1995
33	1994-1995	1994-1995
34	1994-1995	1994-1995
35	1994-1995	1994-1995
36	1994-1995	1994-1995
37	1994-1995	1994-1995
38	1994-1995	1994-1995
39	1994-1995	1994-1995
40	1994-1995	1994-1995
41	1994-1995	1994-1995
42	1994-1995	1994-1995
43	1994-1995	1994-1995
44	1994-1995	1994-1995
45	1994-1995	1994-1995
46	1994-1995	1994-1995
47	1994-1995	1994-1995
48	1994-1995	1994-1995
49	1994-1995	1994-1995
50	1994-1995	1994-1995
51	1994-1995	1994-1995
52	1994-1995	1994-1995
53	1994-1995	1994-1995
54	1994-1995	1994-1995
55	1994-1995	1994-1995
56	1994-1995	1994-1995
57	1994-1995	1994-1995
58	1994-1995	1994-1995
59	1994-1995	1994-1995
60	1994-1995	1994-1995
61	1994-1995	1994-1995
62	1994-1995	1994-1995
63	1994-1995	1994-1995
64	1994-1995	1994-1995
65	1994-1995	1994-1995
66	1994-1995	1994-1995
67	1994-1995	1994-1995
68	1994-1995	1994-1995
69	1994-1995	1994-1995
70	1994-1995	1994-1995
71	1994-1995	1994-1995
72	1994-1995	1994-1995
73	1994-1995	1994-1995
74	1994-1995	1994-1995
75	1994-1995	1994-1995
76	1994-1995	1994-1995
77	1994-1995	1994-1995
78	1994-1995	1994-1995
79	1994-1995	1994-1995
80	1994-1995	1994-1995
81	1994-1995	1994-1995
82	1994-1995	1994-1995
83	1994-1995	1994-1995
84	1994-1995	1994-1995
85	1994-1995	1994-1995
86	1994-1995	1994-1995
87	1994-1995	1994-1995
88	1994-1995	1994-1995
89	1994-1995	1994-1995
90	1994-1995	1994-1995
91	1994-1995	1994-1995
92	1994-1995	1994-1995
93	1994-1995	1994-1995
94	1994-1995	1994-1995
95	1994-1995	1994-1995
96	1994-1995	1994-1995
97	1994-1995	1994-1995
98	1994-1995	1994-1995
99	1994-1995	1994-1995
100	1994-1995	1994-1995

Appendix B - Modification of Phosphorimeter Parameters

Some parameters affecting phosphorimeter operation can be changed by editing the text file ARRAY.INI, which is located in directory \DM3000F\RUN. This file is shared by the phosphorimeters and array detectors (available only with DM3000R software). In particular the array attached flag must be 0, and phosphorimeter attached flag must be 1.

The only parameters whose modification will affect operation are com port and baud rate.

Valid values for com port are (all upper case letters):

COM 1
COM 2

Valid values for baud rate are 9600 and 19200.

```
;all ver >= 2.4    4/12/88      rdw
'NOTE all lines must be 80 characters or less !!!!!!!
;NOTE semi-colon is the comment delimiter in the spirit of
  assembler,
;ie skip to end of line

-5 -5 ; extra variables extra_x, extra_y (normal values are 0,0)
      ; extra_x (-5, 90), extra_y (-5, 26)

; added 7/29/87 array detector information
; unit we have is 4.032mhz, 1024 pixels, db25 rs232  only pins
  2,3,7

0      ;array attached flag (0-> not there, 1-> is there, -1->
      fake it)

;added 4/12/88 phosphorimeter information
;db25 rs232 only pins 2,3,7 required, others optional.

1      ;phos attached flag (0-> not there, 1-> is there, -1->
      fake it)
      ;NOTE either array or phos may be attached, NOT both.

COM1    ; com port for either detector, use COM1  or COM2
19200   ; phos <= 19.2K standard baud 9600, could be higher

ST-100 ;model number of detector EXPECTED
4.032  ;or 8.064 Mhz speed of detector
1024   ;or 512 number of pixels (should be less than 4097)
0.025  ;pixel spacing in mm
```

Typical contents of ARRAY.INI

THE SECRETARY OF THE BOARD OF DIRECTORS
OF THE AMERICAN RED CROSS
WASHINGTON, D. C.

DEAR MR. SECRETARY:

I have the honor to acknowledge the receipt of your letter of the 15th inst.

and in reply to inform you that the same has been forwarded to the proper authorities for their consideration.

I am, Sir, very respectfully,
Yours truly,
J. H. HARRIS

Enclosed for you are two copies of a report of the Committee on the Administration of the American Red Cross, which was presented to the Board of Directors at its meeting on the 10th inst.

I am, Sir, very respectfully,
Yours truly,
J. H. HARRIS

Very truly,
J. H. HARRIS

Very truly,
J. H. HARRIS

Very truly,
J. H. HARRIS

Very truly,
J. H. HARRIS

Very truly,
J. H. HARRIS

Very truly,
J. H. HARRIS

Very truly,
J. H. HARRIS

Very truly,
J. H. HARRIS

Artisan Technology Group is an independent supplier of quality pre-owned equipment

Gold-standard solutions

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

Learn more!

Visit us at [artisanng.com](https://www.artisanng.com) for more info on price quotes, drivers, technical specifications, manuals, and documentation.

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

We're here to make your life easier. How can we help you today?

(217) 352-9330 | sales@artisanng.com | [artisanng.com](https://www.artisanng.com)

