

Multi-Channel Potentiostat / Galvanostat



Limited Availability
Used and in Excellent Condition

Open Web Page

<https://www.artisanTG.com/80699-1>

All trademarks, brandnames, and brands appearing herein are the property of their respective owners.



Your **definitive** source
for quality pre-owned
equipment.


Artisan Technology Group

(217) 352-9330 | sales@artisanTG.com | artisanTG.com

- Critical and expedited services
- In stock / Ready-to-ship

- We buy your excess, underutilized, and idle equipment
- Full-service, independent repair center

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.



Metrohm Autolab

Instruments for Electrochemistry

History of Metrohm Autolab


- Founded in 1986 as Eco Chemie in Utrecht
- Develops state of the art instruments for electrochemistry
- Joined the Metrohm group in 1999
- Name changed to Metrohm Autolab in 2009




Metrohm Group



International presence





Instruments for Electrochemistry

Four instrument ranges





Single channel, modular

- N series of Autolab instruments
- Potentiostat/Galvanostat with USB
- Module bays for instrument extensions



The N series of PGSTATs

- Three models available



PGSTAT128N
12 V / 800 mA




PGSTAT302N
30 V / 2 A



PGSTAT100N
100 V / 250 mA

N Series Summary

	PGSTAT128N	PGSTAT302N	PGSTAT100N
Compliance voltage	12 V	30 V	100 V
Maximum current	800 mA	2 A	250 mA
Bandwidth	500 kHz	1 MHz	400 kHz




Multi Channel, Modular

Twelve Channel Instruments

Multi channel, modular

- 12 channel cabinet with USB
- Fitted with M101 PGSTAT modules
- Fitted with optional modules




M101 module




- Single channel PGSTAT module
 - Output power: 10 V / 100 mA
 - Potential range: 10 V
 - Current range: 10 nA to 10 mA in 7 decades
- Up to 12 modules can be installed

Multi channel, modular

- Combination with optional modules



Flexible System



PGSTAT only
12 M101 modules per
frame



PGSTAT/Module combination
1-11 M101 modules per frame
1 optional module / M101

Compatible Modules

- Optional modules can be installed alongside each M101 PGSTAT module
 - pX1000
 - MUX
 - FRA32M
 - BA
 - EQCM

Unique Features

- Built-in hub allows control from 3 different computers
- Synchronization between channels possible



Compact design

- Potentiostat/Galvanostat with USB
- Compact design at entry level price



PGSTAT101
10 V / 100 mA



PGSTAT204
20 V / 400 mA

New product!

PGSTAT204



PGSTAT101 vs PGSTAT204

	PGSTAT101	PGSTAT204
Compliance voltage	10 V	20 V
Potential range	10 V	10 V
Maximum current	100 mA	400 mA
Modularity	No	Yes
Impedance option	No	Yes



PGSTAT302F

- Switchable floating mode instrument
- Suitable for working on grounded cells or setups requiring floating instrument
- Can be fitted with impedance spectroscopy module

PGSTAT302F




- Entry level modular PGSTAT
 - Output power: 30 V / 2 A
 - Potential range: 10 V
 - Current range: 10 nA to 1 A in 9 decades
 - Switchable from grounded to floating

MultiBA Systems

- Multi working electrode system

Main working electrode


Additional working electrodes



MultiBA Systems



- Multi working electrode system
- Shares the same Reference and same Counter electrode with up to 6 working electrodes
- Independent control on all electrodes



Module Extensions

Customizing the Modular Autolab PGSTAT

Modules

- Extension cards for the modular PGSTATs
- Can be added at any time
- Reverse compatible
- Can be installed onsite

Modules

Module	Description	Main particularities
FRA32M	Electrochemical Impedance Spectroscopy	10 μ Hz-32 MHz*, single sine and multi sine, fit and simulation
SCAN250	True Analog Scan and Fast Scan	250 kV/s, analog scan generator
ADC10M	Fast Data Acquisition	10 MHz data acquisition card for very fast measurements
BA, MultiBA	Bi-Potentiostat	Bipot and Array, up to 6 WE
ECD	Low Current Amplifier	300 aA current resolution, RC filter included

* 1 MHz when used with a PGstat

Modules

Module	Description	Main particularities
pX1000	pH, pX, Voltage and T	Additional diffamp for voltage measurements
EQCM	Electrochemical Quartz Crystal Microbalance	6 MHz, < $\mu\text{g}/\text{cm}^2$ possible
FI20	Yes	Real time Filter & Integrator
MUX	Multiplexer	Up to 64 cells or 255 WEs
ECN	Electrochemical Noise	800 nV voltage resolution
Booster10A, 20A	Current Boosters	Apply and measure up to 20A

Impedance Modules

- Frequency Response Analysis module (FRA) for impedance spectroscopy
 - FRA32M (2011)
 - FRA2 (→ 2011)

Impedance Modules



- Frequency Response Analysis module (FRA) for impedance spectroscopy
 - FRA32M (2011)
 - FRA2 (→ 2011)

FRA32M Module

- Module for electrochemical impedance spectroscopy
 - 32 MHz to 10 μ Hz
 - 1 MHz to 10 μ Hz in combination with the PGSTAT
 - DC range: up to 10 V
 - Single sine and multi-sine



FRA32M/FRA2 Module

- Electrochemical impedance spectroscopy



FRA32M/FRA2 Module


- Electrochemical impedance spectroscopy



Advanced plots for raw measured data

FRA32M / FRA2 Module

- Software integrates powerful data analysis



The screenshot displays the 'Equivalent Circuit Editor' software interface. The main window shows a circuit diagram with three components: a resistor R_s (value $R = 100 \Omega$), a resistor R_p (value $R = 999 \Omega$), and a capacitor C_d (value $C = 981 \text{ nF}$). The circuit is connected to an input and output. A 'Library' window on the left lists pre-defined circuits, with $R(RC)$ selected. A 'Properties' window on the right shows the 'Capacitance (C)' settings, with 'Name' set to C_d , 'Start' at $1E-06$, and 'Fitted' at $9.806E-07$ (0.04%). A 'Circuit Report' window on the bottom right displays the following data:

Element	Parameter	Value	Estimated Error (%)
R_s	R	100.08	0.040
R_p	R	999.09	0.022
C_d	C	9.806E-07	0.039
	χ^2	5.7258E-05	

FRA32M/FRA2

- Module compatibility

	FRA32M	FRA2
PGSTAT128N	Yes	Yes
PGSTAT302N	Yes	Yes
PGSTAT100N	Yes	Yes
PGSTAT101	No	No
M101	Yes	No
PGSTAT204	Yes	No
PGSTAT302F	Yes	Yes
MultiBA	Yes	Yes


SCAN250 Module

- True analog scan generator module
 - Scan rate: 10 mV/s to 250 kV/s
 - Scan range: ± 5 V vs DC offset
 - Linear scan (no staircase)

SCAN250 Module



- Critical difference between staircase and linear scan
- Significant reduction in currents with staircase methods
- Can lead to the wrong conclusions

SCAN250 Module




staircase

linear




SCAN250 Module

staircase



linear



Step: $\frac{\partial V}{\partial t} \rightarrow \infty$

Current spike:

$$i_c = C_{dl} \frac{\partial V}{\partial t}$$

Plateau: $\frac{\partial V}{\partial t} = 0$

Current decay:


$$i = \frac{E_{step}}{R_u} e^{-\frac{t}{R_u C_{dl}}}$$

Ramp: $\frac{\partial V}{\partial t} = const \neq 0$

Const charging current: $i_c = C_{dl} \frac{\partial V}{\partial t}$

SCAN250 Module

- Linear scan vs Staircase



SCAN250 Module

- Module compatibility

	SCAN250
PGSTAT128N	Yes
PGSTAT302N	Yes
PGSTAT100N	Yes
PGSTAT101	No
M101	No
PGSTAT204	No
PGSTAT302F	No
MultiBA	No

ADC10M Module


- Ultra fast sampling synchronous A/D converter
- Sampling rate: 10 M samples/s
- Lowest interval time: 100 ns

ADC10M Module

- Sampling rate: 10 M samples/s
- Lowest interval time: 100 ns
- Maximum scan rate in combination with SCAN250: 250 kV/s
- Maximum number of points: 1 M

ADC10M Module

- 100 ns interval time



ADC10M Module

- Module compatibility

	ADC10M
PGSTAT128N	Yes
PGSTAT302N	Yes
PGSTAT100N	Yes
PGSTAT101	No
M101	No
PGSTAT204	No
PGSTAT302F	No
MultiBA	No

BA Module


- Dual mode Bipotentiostat module
 - Bipot (vs RE) or Scanning Bipot (vs WE)
 - Potential range: 10 V
 - Maximum current: 50 mA
 - Current ranges: 10 nA to 10 mA in 7 decades

BA module

- Measurements on 2 working electrodes

RRDE

Rotating Ring Disc
Electrode
measurements



BA Module

- Module compatibility


	BA
PGSTAT128N	Yes
PGSTAT302N	Yes
PGSTAT100N	Yes
PGSTAT101	No
M101	Yes
PGSTAT204	Yes
PGSTAT302F	No
Multi BA	Yes (max. 5)

ECD Module

- Low current amplifier module
- Current ranges: 100 μ A to 100 pA in 7 decades
- Integrated 3rd order Sallen-Key filter (0.1 s, 1 s & 5 s)
- Resolution: 0.003 % of current range
- Maximum resolution: 300 aA

ECD Module

- Built-in offset compensation for low current measurements




ECD Module

- Module compatibility

	ECD
PGSTAT128N	Yes
PGSTAT302N	Yes
PGSTAT100N	Yes
PGSTAT101	No
M101	No
PGSTAT204	No
PGSTAT302F	No
MultiBA	No

Booster10A/20A



External high current amplifiers

Booster10A/20A

- High current amplifier
 - 20 A maximum current (Booster20A)
 - 10 A maximum current (Booster10A)

Booster10A/20A

	Booster10A	Booster20A
Compliance voltage	20 V	20 V
Maximum voltage	10 V	10 V
Maximum current	10 A	20 A
Bandwidth	10 kHz	20 kHz
Compatibility	PGSTAT128N/302N/100N	PGSTAT302N

pX1000 Module

- pH and T measurements
- Input range: 10 V
- Potential resolution: 30 μV
- Temperature resolution: 0.1°C

pX1000 Module

- pH sensor connection
- Automatic pH correction with temperature
- Built-in calibration application



pX1000 Module

- Module compatibility


	pX1000
PGSTAT128N	Yes
PGSTAT302N	Yes
PGSTAT100N	No
PGSTAT101	No
M101	Yes
PGSTAT204	Yes
PGSTAT302F	No
MultiBA	No

ECN module

- Electrochemical noise and galvanic coupling
 - Input range: 2.5 V
 - Maximum resolution: 760 nV
 - DC offset compensation circuit

ECN module

- Current and potential noise measurements



ECN Module

- Module compatibility

	ECN
PGSTAT128N	Yes
PGSTAT302N	Yes
PGSTAT100N	No
PGSTAT101	No
M101	No
PGSTAT204	No
PGSTAT302F	No
MultiBA	No

FI20 module

- Analog filter and integrator
 - Integrated 3rd order Sallen-Key filter: 0.1 s, 1 s & 5 s)
 - Available integration time constants: 0.01 s, 0.1 s, 1 s & 10 s

FI20 Module

- Module compatibility

	FI20 - Filter	FI20 - Integrator
PGSTAT128N	Yes	Yes
PGSTAT302N	Yes	Yes
PGSTAT100N	Yes	Yes
PGSTAT101	No	Built-in
M101	No	Built-in
PGSTAT204	No	Built-in
PGSTAT302F	No	No
MultiBA	No	No

MUX module

- Multiplexer device for single or multiple electrodes
 - **MUX-MULTI4:** monitor up to **64** electrochemical cells in sequence
 - **MUX-SCANNER16:** monitor up to **255** working electrodes in the same cell in sequence

MUX Module



Sequential measurements on independent electrochemical cells
or on several working electrode

MUX Module

- Module compatibility


	MUX
PGSTAT128N	Yes
PGSTAT302N	Yes
PGSTAT100N	No
PGSTAT101	No
M101	Yes
PGSTAT204	Yes
PGSTAT302F	No
MultiBA	No

EQCM module

- Quartz crystal microbalance module
 - 6 MHz crystal oscillator
 - Dynamic range: 80 kHz
 - Integrated temperature sensor
 - Detection limit: $< 10 \text{ ng/cm}^2$

EQCM module

- Electrogravimetric measurements



EQCM Module

- Module compatibility

	EQCM
PGSTAT128N	Yes
PGSTAT302N	Yes
PGSTAT100N	No
PGSTAT101	No
M101	Yes
PGSTAT204	Yes
PGSTAT302F	No
MultiBA	No



Software

Advanced NOVA Electrochemical Software

Software Overview

- Two software packages available:
 - NOVA – Advanced electrochemical software
 - SDK – Software Development Kit

Autolab NOVA


- One application for any electrochemical experiment
- New software for the Autolab Potentiostat/Galvanostat with USB
- From simple to advanced experiments
- Real time data handling for dynamic procedures

Autolab SDK

- Software Development Kit
- Provides controls for the Autolab instruments through 3rd party application like LabVIEW, Matlab, ...

Autolab SDK

Library of LabVIEW
VIs provided with
the SDK





Rotating Disc Electrode (RDE)

- Supplied with external motor controller



Rotating Disc Electrode

- Mercury contact for friction-less measurements
- 0 – 10000 RPM (manual or software)
- 4000 RPM/s acceleration
- Compatible with 3 and 5 mm electrode tips

Autolab Optical Bench



Autolab Optical Bench

- For studying the optical behaviour of photovoltaic cells
- Software integrated control of light intensity
- Up to 10 interchangeable LEDs


Autolab Optical Bench



Autolab Optical Bench

Article code	Color	Wavelength (nm)
AUT.LED.LDC655	Deep-Red	655
AUT.LED.LDC627	Red	627
AUT.LED.LDC617	Red-Orange	617
AUT.LED.LDC590	Amber	590
AUT.LED.LDC530	Green	530
AUT.LED.LDC505	Cyan	505
AUT.LED.LDC470	Blue	470
AUT.LED.LDC	Cool white	-
AUT.LED.LDC	Warm white	-
AUT.LED.LDC	Neutral white	-

Voltammetric analysis stand



Hg drop electrode stand (663 VA Stand) and
IME663 interface

Voltammetric analysis stand

- Remote control through the IME663 interface
- Stirrer, purging and new drop control
- DME, SDME, HDME control
- Voltammetric analysis methods

Multi Mode Electrode

- Hg drop
 - DME
 - SMDE
 - HDME



Electrochemical Cells

- General purpose cells and stands



Corrosion Cells

- Dedicated cells for corrosion measurements
 - 1 L cell
 - 400 mL cell
 - Flat cell

Corrosion Cell (1 L)

- ASTM compliant
- 1 cm² sample holder
- 6 conical fittings
- Luggin capillary



Corrosion Cell (400 mL)

- 5 conical fitting
- Sample area 0.785 cm²
- Luggin capillary
- Thermostatic jacket



Flat Cell

- Corrosion cell for large surface areas
- Surface area $> 15 \text{ cm}^2$



Working Electrodes

- 3 mm and 5 mm electrode tips
- Au, Pt, GC, Ag
- Empty tip for custom material



Microelectrodes

- Different μm range diameters available

Material	Diameter (μm)
Platinum	10, 20, 25, 50, 100, 200, 300, 500
Gold	10, 25, 40, 100, 200, 300, 500
Palladium	25, 60, 100, 200, 300, 500
Silver	25, 30, 60, 100, 200, 300, 500
Iridium	75



Reference Electrodes

- Different reference electrodes available
- Single and double junction



Counter Electrodes

- Range of materials available
- Pt, GC, Stainless steel

Faraday cage

- Provides shielding for the electrochemical cell
- Compatible with Metrohm cell stands



www.metrohm-autolab.com

Logo: Metrohm Autolab s.v. | Navigation: Home | Autolab Products | Applications | Company | Support | Login



News

- 03.10.2011, 13:00
[Metrohm Analytics Romania to distribute Autolab products in Romania](#)
- 02.10.2011, 06:00
[New website content](#)
- 01.10.2011, 12:00
[SPR Product line discontinued](#)
- 23.09.2011, 17:32
[Metrohm Autolab Testimonials](#)
- [More news...](#)

Shows and conferences

- 09.10.2011 - 14.10.2011
[ECS Fall Meeting](#)
Metrohm USA
- 04.11.2011
Annual Symposium of the KNCV Working Group on Electrochemistry
Metrohm Belgium
[Program](#)
- 28.11.2011 - 02.12.2011
[MRS Fall Meeting](#)
Metrohm USA
- [More...](#)

New software for Autolab Potentiostats is available: check our latest release of [NOVA](#).



facebook [Join our group](#)

Metrohm worldwide
Please select your Country:

Contact form
Fill in the contact form and receive a brochure on your mail.

Logo: Metrohm Autolab s.v. | Navigation: Home | Autolab Products | Applications | Company | Support

Application notes

The Autolab application notes are available in PDF format for download.

General Electrochemistry

- EC01 - Study of the Mass Transport Characteristics of $K_3[Fe(CN)_6]/K_4[Fe(CN)_6]$ Oxidation and Reduction Reaction Using AUTOLAB RDE
- EC02 - Reference Electrodes and Their Usage
- EC03 - Ohmic Drop: Part 1 - Basic Principles
- EC04 - Ohmic Drop: Part 2 - Measurement
- EC05 - EQCM Study of Underpotentially Deposited (UPD) Lead Adlayer on Gold
- EC06 - A Simple Experimental Protocol for Pt Deposition on Au Electrodes
- EC07 - Comparison between Staircase Cyclic Voltammetry and Cyclic Voltammetry Linear Scan

Electrochemical Impedance Spectroscopy


- EIS01 - Electrochemical Impedance Spectroscopy (EIS): Part 1 - Basic Principles
- EIS02 - Electrochemical Impedance Spectroscopy (EIS): Part 2 - Experimental

Service and Support

- All instruments covered by 3 years factory warranty

A large, teal, curved graphic element that frames the text '3 year instrument warranty'. The text is in a large, bold, sans-serif font, with '3 year' on the top line and 'instrument warranty' on the bottom line.

3 year
instrument warranty



Metrohm Autolab

Instruments for Electrochemistry

Artisan Technology Group is an independent supplier of quality pre-owned equipment

Gold-standard solutions

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

Learn more!

Visit us at [artisanng.com](https://www.artisanng.com) for more info on price quotes, drivers, technical specifications, manuals, and documentation.

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

We're here to make your life easier. How can we help you today?

(217) 352-9330 | sales@artisanng.com | [artisanng.com](https://www.artisanng.com)

