Mettler Toledo Aquastar 51108753

Generator Electrode without Diaphragm

Limited Availability
Parts Grade / As-Is

Open Web Page

<https://www.artisanng.com/63576-9>

All trademarks, brandnames, and brands appearing herein are the property of their respective owners.

Your **definitive** source
for quality pre-owned
equipment.

Artisan Technology Group

(217) 352-9330 | sales@artisanng.com | artisanng.com

- Critical and expedited services
- In stock / Ready-to-ship


- We buy your excess, underutilized, and idle equipment
- Full-service, independent repair center

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

Operating Instructions

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Mode d'emploi

Drying Oven
Trocknungs-ofen
Four de séchage

D0308

METTLER TOLEDO

Français

Deutsch

English

Contents

1. Introduction	3
2. Safety measures	3
3. Overview	5
4. Installation	6
4.1 Assembly of the oven	6
4.2 Installing the sample transport assembly	7
4.3 Connecting the tubing	8
4.3.1 Gas drying unit	8
4.3.2 Gas outlet tubing "cold zone": connection to the stop cock	9
4.3.3 Transfer tube: connection to the titrator	9
5. Settings and Programming	10
5.1 The menu structure	10
5.2 The display	11
5.3 The basic operating concept	11
5.4 Operating the DO308 in manual mode	12
5.4.1 Switching to manual mode	12
5.4.2 Changing the target temperature	12
5.4.3 Edit, create or delete a program in manual mode	12
5.5 Operating the DO308 in program mode	14
5.5.1 Switching to program mode	14
5.5.2 Load a program	14
5.5.3 Edit or create a program	14
5.5.4 View a program	16
5.5.5 Delete a program	16
6. Putting into operation	17
6.1 Purge gas flow	17
6.2 Initial drying out of the sample tube and glass boat	17
6.3 Heating a sample to dryness	18
6.3.1 Introducing sample with the transport assembly	18
6.3.2 Introducing sample via the tapered joint	19
6.4 Notes on the moisture determination	20
7. Options	21
7.1 Air pump box	21
8. Error messages	22

Contents

9. Maintenance		23
9.1	Changing the fuses	23
9.2	Cleaning the apparatus	23
9.3	Cleaning the sample tube	23
9.4	Changing the silica gel and molecular sieves	23
9.5	Replacing the sample tube	24
10. Accessories		25
10.1	Standard equipment	25
10.2	Optional equipment	26
11. Technical Data		27
12. Heating temperatures for different substances		28
13 Disposal		29

1. Introduction

The METTLER TOLEDO DO308 Drying Oven is used for moisture determination by the Karl Fischer method for samples where direct moisture determination is not possible. Examples are insoluble samples or those whose determination yield false results due to side reactions with the KF reagent. The oven can also be an alternative to external extraction. The drying oven can be used together with both the volumetric and coulometric METTLER TOLEDO Karl Fischer titrators.

The sample in a glass boat, or in an aluminium foil insert, is slid into a glass tube, where it is heated. This heating drives off any moisture from the sample in the form of vapor, which is then fed from the heated glass tube to the KF titration cell by a stream of dry purge gas, such as dry air or nitrogen.

In these operating instructions we explain the functions of the oven; we assume familiarity with the titrator used.

2. Safety measures

The DO308 has been tested for the experiments and intended purposes documented in these operating instructions. However, this does not absolve you from the responsibility of performing your own tests of the product supplied by us regarding its suitability for the methods and purposes you intend to use it for. You should therefore observe the following safety measures.

Measures for your protection

Risk of fire

- Never place or install objects on or above the apparatus! These could catch fire.
- Avoid heat buildup above the apparatus by leaving a free space of approx. 1 meter! Readily flammable materials could ignite.
- For safety reasons, there must be at least 30 cm between the wall or other equipment and the rear and sides of the device! No containers, chemicals or other devices should be placed behind the device!

Risk of burns

- Never touch the guard of the oven as well as the flange of the sample inlet which is attached to the flange ring of the heating tube during or after operation! You could burn yourself.

Risk of electric shock

- Install the power cable so that it doesn't come into contact with the hot parts of the oven! Otherwise the cable could melt and cause a lethal technical fault.
- Ensure that you plug the power cable supplied into a receptacle outlet that is grounded! In the absence of grounding, a technical fault could be lethal.
- Switch the apparatus off and disconnect the power cable before you change blown fuses! An electric shock could be lethal.

Risk of explosion

Risk of burns

Risk of fire

- Never work in an environment subject to explosion hazards! The housing of the apparatus is not gas tight (explosion hazard due to spark formation, corrosion caused by the ingress of gases).
- Never use a reactive gas as a purge gas! It could lead to an explosion.
- In the determination of water contents of organic substances and when using ambient air as purge gas, ensure that no explosive gas mixtures can be formed! Risk of explosion!
- Check the gas paths and drying tubes for free gas flow ! This prevents buildup of an overpressure in the titration cell, which could lead to an explosion.
- Open the stop cock before you stop the gas flow, switch off the apparatus or lower the temperature! This prevents siphoning back of the titration solution into the hot drying oven. Failure to do this could result in an explosion hazard.
- Don't use flammable liquids with a flashpoint below 325°C in conjunction with the oven DO308! Otherwise they could catch fire.
- Use flammable liquids only in quantities, too small to cause spread of fire!
- When using chemicals and solvents, comply with the instructions of the producer and the general lab safety rules!
- When working with chemicals or solvents of unknown composition, it is the operator's responsibility to obtain information from appropriate sources about possible hazards associated with the use of such chemicals or solvents!
- The drying oven DO308 has been designed and built as a laboratory device. It is only intended to be used to drive out water from a small amount of sample at temperatures from room temperature up to 300°C. Using the device for any purpose other than that stated above, or using it in any application that does not correspond to the technical data, is considered to be improper usage. Any damage resulting from such use is the sole responsibility of the operator!


Caution

Measures for operational safety


- Use only fuses of the type specified in the Operating Instructions!
- Have the apparatus opened and serviced only by METTLER TOLEDO service!
- Exclude the following environmental influences:
 - powerful vibrations,
 - direct sunlight,
 - atmospheric humidity greater than 75%,
 - temperatures below 5 °C and above 40 °C,
 - powerful electric or magnetic fields!

3. Overview

D0308 with connection to a titrator and to a pump

a	Transfer tube to titration cell	i	Pump (optional)
b	Touch guard	j	Multifunctional adjusting knob
c	Jacket tube	k	On/off switch
d	Flanged ring of the sample inlet	l	Display
e	Sample tube fixing	m	Stop button
f	Sample tube	n	Start button
g	Screw cap	o	Stop cock
h	Drying unit		

Rear view


a	Jacket tube	e	Power supply connection
b	Touch guard	f	Fuse holder
c	Heating tube	g	Warning sign
d	Cable of the oven control		

4. Installation**4.1 Assembly of the oven**

The DO308 is delivered disassembled in two parts, the oven unit and the sample tube unit.


To mount the oven, proceed as follows:

- 1) Remove the three screws from the oven unit (position 1, 2 and 3).
- 2) Carefully slide the sample tube unit into the oven unit, as shown in the picture below.
- 3) Tighten the sample tube unit to the oven unit with the three screws (A, B and C), removed before.


4.2 Installing the sample transport assembly

This unit is used to insert the sample into the heated zone of the sample tube.

- Place an aluminum insert into the glass boat. (Only necessary if you have a sample that sticks to the glass boat.)
- Insert the hook of the guide rod into the eye of the glass boat.

- Remove the Screw cap from the Sample tube.
- Fully slide back the slider towards the opening of the tube (as illustrated).
- Carefully introduce the transport assembly into the sample tube until the position of the guide rod is fixed by the attraction of the slider.

The boat is now positioned underneath the tapered joint. You can introduce samples through the tapered joint.


- Insert a silicon stopper (septum) or NS stopper into the tapered joint.
- Close the tube with the screw cap as shown.

If you move the slider in towards the oven, the glass boat will be slid into the heating zone of the sample tube.


4.3 Connecting the tubing


4.3.1 Gas drying unit

As the purge gas used must be 100% dry, we recommend to dry it before it flows into the oven using the gas drying unit (standard equipment). The drying agents silica gel and molecular sieves shown in the diagram are not part of the standard equipment; you must connect the connection tubing yourself.


- First fill the bottles: Fill the first bottle with silica gel whose capacity for water uptake is higher; fill the second bottle with 3 Å molecular sieve to trap the remaining water.
- Connect the flow meter to the first bottle with the 27 cm tube: Press the ends into the quick fit connectors. In the case of the bottle; make sure to use the connection which leads the gas to the base of the bottle (IN).
- Connect the first bottle and the second bottle with the 21 cm tube: In the case of the second bottle insert the tubing in the connection which leads the gas to the base of the bottle (IN).
(to loosen the tubing, press on the clamping rings with your fingertips and take out the tubing.)
- Connect the hose connector with barbed fitting to the gas source using tubing.
- Set the gas flow rate on the flow meter (30 - 240 mL/min).
- To detach the hose connector, press the quick acting coupler in the direction of the arrow.


- Connect the second bottle (OUT) of the drying unit to the gas inlet of the sample tube using a 37 cm tube.
(Push the coupling of the tube over the gas inlet of the sample tube in direction of the arrow.)

4.3.2 Gas outlet tubing "cold zone": connection to the stop cock

- Connect the stop cock with the gas outlet of the "cold zone" (see Chapter "6.1 Purge gas flow") using a 37 cm tube.
(To remove the tube from the stop cock, press on the clamping rings with your fingertips and pull out the tubing.)


4.3.3 Transfer tube: connection to the titrator

- Connect the transfer tube with the gas outlet of the sample tube.
- Connect the other end of the transfer tube with the gas inlet of the titrator.

5. Settings and Programming

5.1 The menu structure

The menu structure of the DO308 consists of three submenus:





- Mode
- Programs
- Options

In the menu "Mode" the operation mode of the oven, either "Manual" or "Program" can be chosen.

In the manual mode it is only possible to set a desired-value for the temperature and to switch the heating on or off. In the program mode you have the opportunity to create and run programs with up to three different temperature values and heating durations. These programs can be created, loaded, edited and viewed in the menu "program".

In the "Options" menu the instrument language, the contrast of the display and the temperature unit can be adjusted. (The menu "Vakuum Control" has no function.)


5.2 The display

a	Operating condition (On or Off)
b	Information about the operation mode - in the manual mode "Man" is printed - In the program mode the number of the running programm and below the current stage (I, II oder III) is printed
c	Display of the elapsed time - in the manual mode the heating duration is shown - in the program mode the remaining time for the program and, below, for the current stage is shown
d	the choosen temperature unit
e	the current temperature
f	graphical display of the heating intensity
g	Setpoint adjustment
h	Menu selection
i	name of the current menu
j	menu structure, indicates the present level as a black square
k	Scrollbar of the present menu

1	Information screen
2	Working menu

5.3 The basic operating concept

The DO308 is operated using the multifunction button (m) and the start and stop buttons (k and l) only.

By pressing the multifunction button (m button) you confirm a selection or an adjusted value.


By turning the m button clockwise or counter clockwise you can move the black scroll bar up or down within a menu to select a menu item.

All Settings are done in the working menu (2).

After a short time without any action taken by the user the display switches back to the information screen (1). In the program mode a program is started with the green "Start" button and stopped with the red "Stop" button. In the manual mode the heating of the oven is switched on with the "Start" button and off with the "Stop" button.


5.4 Operating the D0308 in manual mode

5.4.1 Switching to manual mode



- Press or turn the multifunction button to open the working menu.
- Turn the m button counter clockwise to choose the menu item "Menu" and confirm your selection by pressing the m button.
- Choose the menu item "Mode" and confirm your selection by pressing the m button.
- Choose the option "Manual" and confirm your selection.
- By choosing "Back" you get back to the working menu.

5.4.2 Changing the target temperature

- Press or turn the multifunction button to open the working menu.
- Press the m button once again to choose the menu item "Setpoint".
- By turning the m button, select the desired temperature value.
- For room temperature select "RT".
- Press the m button to confirm your selection.

5.4.3 Edit, create or delete a program in manual mode

Even in the manual mode and during the oven is heating it is possible to view, create, edit or delete programs.



- Press or turn the multifunction button to open the working menu.
- Turn the m button counter clockwise to choose the menu item "Menu" and confirm your selection by pressing the m button.
- Choose the menu item "Programs" and confirm your selection.
- Choose "View" to view a program and confirm your selection.
- Choose "Select:" and confirm your selection.
- By turning the m button you can choose the appropriate program.
- After the confirmation of your selection, the program is shown on the display.
- With "Back" you get back to the "Programs" menu.

- Choose "Edit" to change a program. To create a new program choose "New" instead of the name of an existing program and save the new created program later on one of the ten possible storage locations (P01-P10). By choosing "Back" you get back to the working menu.
- Choose "Delete" to delete an existing program. Select the name of the program you want to delete by turning the m button. After you confirm your selection you will be asked if you want to delete the selected program. Confirm this question with "Yes" or choose "No" if you are not sure you want to delete the program.

For a more detailed description of the program handling, see "5.5 Operating the DO308 in program mode".


Notes:

1. As long as the heating of the oven is in use (Operating Condition: "On"), it is not possible to choose the menu items "Mode" and "Load".
2. If you load a program in the manual mode, the DO308 automatically changes to the program mode.


5.5 Operating the D0308 in program mode

5.5.1 Switching to program mode


- Press or turn the multifunction button to open the working menu.
- Turn the m button counter clockwise to choose the menu item "Menu" and confirm your selection by pressing the m button.
- Choose the menu item "Mode" and confirm your selection by pressing the m button.
- Choose the option "Program" and confirm your selection.
- By choosing "Back" you get back to the working menu.

5.5.2 Load a program

- Press or turn the multifunction button to open the working menu.
- Turn the m button counter clockwise to choose the menu item "Menu" and confirm your selection by pressing the m button.
- Choose the menu item "Programs" and confirm your selection by pressing the m button.
- Choose "Load" and confirm your selection.
- Choose an existing program (P01-P10) by turning the m button and confirm your selection.
- After the display switches back to the working menu you can start the program by pressing the green "Start" button.

5.5.3 Edit or create a program

- Press or turn the multifunction button to open the working menu.
- Turn the m button counter clockwise to choose the menu item "Menu" and confirm your selection by pressing the m button.
- Choose the menu item "Programs" and confirm your selection by pressing the m button.
- Choose "Edit" and confirm your selection.
- Choose "Select:" and confirm your selection.

Edit		
Select:		New
Step 1	100°C	00:10
Step 2	150°C	00:20

- By turning the m button, choose an existing program or choose "New", if you want to create a new program. Confirm your selection.
- The scroll bar jumps to "Step 1".
- Enter one after another for steps 1 to 3 the desired heating temperature and heating duration. By pressing the m button you can navigate through the input fields for temperature, heating duration hours and heating duration minutes. By turning the m button you choose the appropriate values and by pressing the m button you confirm them.
- Choose the menu item "Save" and confirm your selection.
- Choose the name of the program under which you want to save the edited program by turning the m button and confirm your selection by pressing the m button. (If you created a new program the name of an empty program storage place is offered. If already all ten program storage places are occupied, an existing program has to be overwritten.)
- If you want to overwrite an existing program, you have to answer the subsequent question with "Yes".
- By choosing "Back" twice, you get back to the working menu.

Notes

1. When you enter the heating duration, keep in mind, that the selected time is from the start of heating and not from once the set temperature is reached.
2. Enter a heating duration of 00:00 for unneeded steps in programs. Those steps are then ignored.

5.5.4 View a program


- Press or turn the multifunction button to open the working menu.
- Turn the m button counter clockwise to choose the menu item "Menu" and confirm your selection by pressing the m button.
- Choose the menu item "Programs" and confirm your selection by pressing the m button.
- Choose the menu item "View" and confirm your selection.
- Choose the menu item "Select" and confirm your selection.
- Choose the name of the program you want to view by turning the m button and confirm your selection.
- The scroll bar jumps to "Step 1".
- You can now view the program but can not change it.
- With three times "Back" you get back to the working menu.

5.5.5 Delete a program


- Press or turn the multifunction button to open the working menu.
- Turn the m button counter clockwise to choose the menu item "Menu" and confirm your selection by pressing the m button.
- Choose the menu item "Programs" and confirm your selection by pressing the m button.
- Choose the menu item "Delete" and confirm your selection.
- Choose the name of the program you want to delete by turning the m button and confirm your selection.
- Answer the subsequent question with "Yes".
- The program is deleted.
- With two times "Back" you get back to the working menu.

6. Putting into operation

6.1 Purge gas flow

Gas flow with stop cock in **open** position. This allows purging of the "cold zone" of the sample tube (back purge).

Gas flow with stop cock in **closed** position. This allows purging of the "hot zone" of the sample tube (cell purge) for moisture transfer.

6.2 Initial drying out of the sample tube and glass boat

- Connect the apparatus to the power supply with the power cable.
- Switch the instrument on: The operating condition should be "Off" and the stop cock "open".
- Set the gas flow rate on the flow meter.
- Purge the "cold zone" of the sample tube for approx. 5 minutes. Meanwhile prepare the titrator so that the drift can be measured.
- Baking out the oven and the glass boat:
 - Close the stop cock to allow the purging of the "hot zone".
 - Slide the boat to the end of the sample tube and set a temperature of 300 °C. The moisture in the sample tube is transferred to the titration cell.
 - Start the heating of the oven (Press Start-button, operating condition "On")
 - Wait until the drift has dropped to below 50 µg/min (time required: approx. 15 minutes).

The system is now free from moisture and ready for measurement.


Notes

1. If the DO308 has been switched off for a relatively long time, you should always predry the sample tube and glass boat before the first measurement.
2. If the DO308 is in daily use, we advise you not to shut off the purge gas, but simply reduce the flow rate to approx. 60 mL/min.
3. Open the stop cock before you stop the gas flow, switch off the apparatus or lower the temperature! This prevents siphoning back of the titration solution into the hot drying oven.


6.3 Heating a sample to dryness

6.3.1 Introducing sample with the transport assembly


- Set the oven temperature appropriate for the sample. Open the stop cock to allow purging of the "cold zone" of the sample tube during heating.
- After the set temperature has been reached, close the stop cock to allow gas purge via the "hot zone" into the titration cell of the connected Karl Fischer titrator. Wait until a stable drift has been reached.
- Place the aluminium insert into the glass boat and weigh the sample therein.
- Open the stop cock to purge the "cold zone" of the sample tube and remove the screw cap.
- Insert the sample with the transport assembly into the sample tube, close the tube with the screw cap and purge the cold zone of the sample tube for approximately 2...3 minutes.
- Close the stop cock to purge via the hot zone into the titration cell.
- After a stable drift has been reached, start the titration at the titrator and slide the sample boat to the end of the sample tube. The moisture is then transferred into the titration cell.
- After completion of the titration open the stop cock. Pull back the sample boat and allow it to cool before unscrewing the cap.
- Remove the transport assembly and close the sample tube. Leave the stop cock open.

Stop cock closed

Stop cock open


Note

This way of introducing the sample is recommended mainly for samples that are difficult to handle and can not be introduced easily via the tapered joint. (See 6.3.2 Introducing sample via the tapered joint)

6.3.2 Introducing sample via the tapered joint

If you add your sample via the tapered joint, proceed as follows:

- Open the stop cock to purge the "cold zone" of the sample tube.
- Weigh the sample and tare the balance with the the filled syringe or weighing boat. Either inject the sample through the silicon stopper (septum) or take out the stopper and add the sample to a sample boat positioned below the inlet.

- Back weigh the empty syringe or weighing boat to get actual sample weight.
- Close the stop cock to purge via the "hot zone" into the titration cell.
- After a stable drift has been reached, start the titration at the titrator and slide the sample boat to the end of the sample tube. The moisture is then transferred into the titration cell.
- After completion of the titration open the stop cock. Pull back the sample boat and allow it to cool before unscrewing the cap.
- Remove the transport assembly and close the sample tube. Leave the stop cock open.

Notes

1. This way of introducing the sample is recommended mainly for samples that contain a certain quantity of surface humidity. This amount of water from the surface would get lost for the determination during the purging of the cold zone, which would be necessary, if the sample is introduced into the sample tube together with the transport assembly. (See 6.3.1 Introducing sample with the transport assembly)
2. You can heat several samples in the same aluminum insert to dryness if the sample residues do not have an adverse effect on the subsequent determinations.
3. The oven temperature should be as high as possible to liberate moisture quickly from the sample. However, decomposition, vaporization or sublimation of the sample should be avoided. Heating temperatures for several samples are listed in Section 12.

6.4 Notes on the moisture determination

1. Enter the time required to liberate the moisture from the sample on the titrator as a stirring time to avoid premature termination of the titration (applies to volumetric titration). The liberation rate depends on the sample.
2. You should avoid high flow rates (greater than 200 mL/min). If you do not, the moisture introduced as vapor would not be completely "trapped" by the titration solution, but would be carried away from the titration cell by the escaping gas.
3. The tubing transferring the water to the titration cell must be immersed in the titration solution by at least 4 mm! Introduction of the "hot" vapor into the titration cell causes evaporation of the volatile fraction of the titration solution from the cell (up to 30 mL within 8 hours at an oven temperature of 200 °C and a gas flow rate of 200 mL/min).
 - With coulometric titrations, you should thus always top up the anolyte solution, with volumetric titrations with the particular solvent used.
4. For the DL32/DL39 Coulometer we advise replacing the straight drying tube of the generator electrode by the bent drying tube. This prevents any solvent condensed on the drying tube from dripping into the catholyte.


7. Options

7.1 Air pump box

If you wish to use air for purging instead of an inert gas, you can attach an air pump box (see Section 10.2: Optional equipment) to the standard gas drying unit .

- Switch off the DO308.
- Connect the pump and the flow meter using the supplied silicon tubing.
- Connect the pump to the oven using the 5-pin cable supplied.
- Switch on the DO308 and set the flow rate on the flow meter (200 ± 20 mL/min).

Switching on the oven also starts the air pump. The air is now pumped through the drying unit into the oven.

Note

Air contains oxygen which at elevated temperatures and above all with organic compounds, e.g. polymers and carbohydrates can lead to decomposition!

8. Error messages

Coded error messages can appear on the display of the oven; these and other malfunctions are listed in the following table:

Malfunctions	Cause / Possible Cause	Corrective action
Error message E1	Internal temperature sensor: defect	Contact METTLER TOLEDO Service.
Error message E2	Internal temperature sensor: defect	Contact METTLER TOLEDO Service.
Error message E3	Internal temperature sensor: overheating	Contact METTLER TOLEDO Service.
Error message E4...En		Contact METTLER TOLEDO Service.
Main switch does not light up	Device is turned off.	Switch on the device.
	Power supply voltage is missing.	Check power supply voltage.
	Fuse blown.	Check fuses.
No temperature display	LED display or board faulty.	Check line voltage and fuses, replace if need be, otherwise call METTLER TOLEDO service.
Oven does not heat up	Cable not connected.	Check connections.
	Blown fuse.	Replace fuse (see "9.1 Changing the fuses").
	Temperature set too low.	Set correct temperature.
	Heating faulty.	Contact METTLER TOLEDO Service.
No gas flow in titration cell at 200 ml/min	Gas inlet / gas outlet connections not tight.	Check connections.
	Screw cap not tight.	tighten screw cap, check sealing.
Drift too high	Titration cell not tight.	See Operating Instructions of the titrator.
	Gas inlet / gas outlet connections not tight.	Check connections.
	Drying agent for purge gas exhausted.	Replace drying agent in gas drying unit.
Moisture content found less than expected	No purge gas.	Check purge gas supply.
	Discharge tip of the gas outlet tubing blocked.	Clean discharge tip with solvent.
Subsequent consumption of titrant too high	Titration terminated too early.	Increase stirring time.
	Temperature too low.	Set higher temperature.

9. Maintenance


9.1 Changing the fuses

Risk of
electric shock

Switch off the DO308 and disconnect the power cable!

- Loosen the fuse holder from above with a flat screwdriver and take it out.
- Replace blown fuses with new ones with the same rating and push the fuse holder back in.

Caution

Use only fuses
of the following type:
T 3.15A L 250V

9.2 Cleaning the apparatus

The housing of the oven is not splash proof. In case of contamination, it can be cleaned with a cloth moistened with alcohol.

9.3 Cleaning the sample tube


You can clean the sample tube with chromic acid or other glass cleaning agents.

- To dismantle it, see description in Section 9.5.
- Dry the tube before installation and bake out following installation (see Section 6.2).

9.4 Changing the silica gel and molecular sieves

- You should change the silica gel of the gas drying unit as soon as it starts to change its color.
- You should change the molecular sieves when the drift value of the titration solution is two to three times higher than that with fresh molecular sieves.

9.5 Replacing the sample tube

- Switch off the DO308, disconnect the power cable and let the oven cool down.
- Unscrew the screw cap (a) and remove the sample transport assembly (guide rod and glass boat).
- Disconnect all tubing from the sample tube.
- Unscrew the two screws on the holder (b) with the provided screwdriver and remove the spacing plate (i).
- Remove the holder from the guide rods (e) and take it away from the sample tube.
- Remove the slider (h) from the square bar (g).
- Unscrew the two screws from the sample tube fixing (c) in the middle of the flanged ring (f).
- Carefully pull the sample tube, together with the sample tube fixing (c) and the gasket (d), out of the oven.
- Push the gasket on the new sample tube and place the sample tube fixing between the upper gas outlet and the thickening of the sample tube, as shown in the picture above. (Take care, that the counter sinkings for the screws have to point frontwards.)
- Carefully slide the sample tube into the oven, with the guide rods sliding to the lower holes of the sample tube fixing.
- Tighten the sample tube fixing.

Do not tighten the screws too tightly!

- Reinstall the slider to the square bar.
- Slide the holder from above on the sample tube and attach it to the guide rods.
- Attach the spacing plate to the square bar and fix it on the holder with the two screws.


Do not tighten the screws too tightly!

- Insert the sample transport assembly and reconnect the tubing.


10. Accessories

10.1 Standard equipment


You can order each part listed with an order number from METTLER TOLEDO


Drying oven with sample tube and screw cap
D0308


Guide rod
51371101


Gas drying unit
(without tubing set)
51108660


Sample boat (glass) (2 Pcs.)
51108657


Operating instructions in three languages (English, German, French)
51710380


Aluminum insert
51108649
(Set of 50)


Power cable
according to order


Silicon stopper (septum), NS 14,5
23276
(Set of 10)


TORX screwdriver 10
51191659


Three-hole adapter
(with plug stoppers for volumetric titrations)
23982


TORX screwdriver 20
51191660


Stop cock
51191130


Tubing set
51108676

Sample plug (TEFLON)
51108668


Transfer tube (insulated)
51108837


Gas inlet plug for DL32/DL39
51108761


Accessories


Gas inlet DL38
105091


Spare fuses T3.15L250V
54286
(Set of 3)


Gas inlet DL39
51108669

Guideway oil (10 ml)


Glass wool (2g)


Drying tube bent for
DL32/DL39
51108639


10.2 Optional equipment


Air pump with cable
51108673


Male stud fitting with tail piece
adapter
Connection to gas source
51107375


Sample tube
51371206


Silica gel (1 kg)
105080


Screw cap for sample tube
51191129

Molecular sieves 3 Å (250 g)
71478

Gasket for sample tube
51371207

Silicon grease (1 Tube)
71300

11. Technical Data

Power supply connection

- Voltage 100 - 230 VAC
- Admissible voltage fluctuations $\pm 10\%$
- Frequency 50/60 Hz
- Power consumption 450 W
- Fuses T 3.15A L 250V

Temperature

- Control range 40 °C ... 300 °C
- Measurement deviation ± 5 °C in the middle of the sample tube at 300 °C

Interface

- 5-pin mini Din socket for air pump box 24 VDC regulated, max. 10 W

Dimensions

- Footprint housing width 256 mm, depth 180 mm
- Dimensions width 300 mm, height 300 mm, depth 550 mm
- Sample tube diameter 30 mm, length 390 mm
- Effective boat volume 15 cm³
- Weight 9.1 kg

Ambient conditions

Ambient temperature: 5...40 °C
 max. relative humidity (not condensing):
 80% for temperatures up to 30 °C, linear decreasing up to
 50% relative humidity at 40 °C
 Use indoors
 Height up to 2000 m
 Installation-/ Overvoltage category II
 Pollution degree 2

Materials


- Housing Sheet steel
- Sample tube, heating tube, jacket tube Borosilicate glass
- Flanged rings Aluminum
- Gas carrying tubing Polyamide
- Sample boat Borosilicate glass

12. Heating temperatures for different substances

Substance	Heating temperature [°C]	Purge gas
α-Polyoxymethylene (POM)	180	N ₂
ABS resin	160	N ₂
apura® Water Standard Oven 1%	170	air or N ₂
Calcium carbonate	100 / 300	air or N ₂
Coal, carbon black	180 - 200	N ₂
Copper sulfate	250	air or N ₂
Epoxide resin	230	N ₂
Ferrite	200	N ₂
Fluor	150	N ₂
HYDRANAL® Water Standard KF-Oven 5.55%	220	air or N ₂
Iron oxide	300	air or N ₂
Manganese carbonate	300	air or N ₂
Phenol resin	110	N ₂
Polyamide	140 - 200	N ₂
Polyacetal	200	N ₂
Polycarbonate (PC)	150 - 200	N ₂
Polyethylene (PE)	140 - 200	N ₂
Polyethylene terephthalate (PET)	180 - 200	N ₂
Polypropylene (PP)	160 - 180	N ₂
Polystyrene (foamed)	130	N ₂
Polystyrene (PS)	150	N ₂
Polyvinyl alcohol (PVA)	150	N ₂
Polyvinyl chlorid (PVC)	130 - 160	N ₂
Rice flour	150	N ₂
Rubber	150	N ₂
Sodium nitrite	150	N ₂
Sodium sulfite	150	N ₂
Sodium tartrate * 2 H ₂ O	200	N ₂
Starch (grain, rice, soja)	135 - 150	N ₂

13 Disposal

- If you need to dispose of your DO308, please contact METTLER TOLEDO Service.

- At the end of its life cycle the instrument and its parts should be disposed of in an environment friendly manner. Some of the components and accessories of the instrument contain materials that must be treated as toxic waste on disposal. The symbol at the back of the instrument reminds you that the instrument and its parts should not be treated as normal waste:
- Please contact your METTLER TOLEDO service support for details and advice on how to dispose of your DO308 and its accessories. Please ask your service engineer for more details.

EC - DECLARATION OF CONFORMITY

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KD-Nr.: 51371200

Doku-Nr.: 20050012

The undersigned, representing the following manufacturer

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Sonnenbergstrasse 74

CH-8603 Schwerzenbach, Switzerland

herewith declares that the product

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certified model:

--

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is in conformity with the provisions of the following EC directives (incl. all applicable amendments)

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Last two digits of the year in which the CE marking was affixed: **05**

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.....
Rolf Truttmann

Manager Business Area AnaChem

References of standards for this declaration of conformity, or parts thereof:

Harmonized standards of Europe and Switzerland:

Safety standards:

IEC/EN61010-1:2001

EMC standards:

EN61326:1997+A1:98+A2:01+A3:03 (class B)

EN61326:1997+A1:98+A2:01+A3:03 (Minimal requirements)

Metrological standards:

--

IP standards:

--

Standards for Canada, USA and Australia

CAN/CSA C22.2 No. 61010-1-04

UL Std No. 61010A-1

FCC, Part 15, class A

AS/NZS4251.1, AS/NZS4252.1

ISO 9001 certificate for METTLER TOLEDO

The Mettler-Toledo AG company, Greifensee, was examined and evaluated in 1991 by the Swiss Association for Quality and Management Systems (SQS), and was awarded the ISO 9001 certificate. This certifies that Mettler-Toledo AG, Greifensee, has a quality management system that conforms with the international standards of the ISO 9000 series.

Repeat audits are carried out by the SQS at intervals to check that the quality management system is operated in the proper manner and is continuously updated in relation to changes brought about.

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Quality certificate. Development, production and testing according to ISO9001 . Environmental management system according to ISO14001 .

Worldwide service. Our extensive service network is among the best in the world and ensures maximum availability and service life of your product.

European conformity. The CE conformity mark provides you with the assurance that our products comply with the most recent EU directives.

On the Internet. You will quickly find lots of essential information about our products, our services, and our company at <http://www.mt.com>

Printed on 100% chlorine-free paper, for the sake of our environment.
Gedruckt auf 100% chlorfrei hergestelltem Papier. Unserer Umwelt zuliebe.
Imprimé sur papier 100% exempt de chlore, par soucie d'écologie.

P51710380

Subject to technical changes
Technische Änderungen vorbehalten
Sous réserve de modifications techniques

Mettler-Toledo AG 2005, 2008 ME-51710380 Printed in Switzerland 0806/2.31

Mettler-Toledo AG, Analytical, Sonnenbergstrasse 74, CH-8603 Schwerzenbach
Tel. +41 44 806 77 11, Fax +41 44 806 73 50, Internet: <http://www.mt.com>

Artisan Technology Group is an independent supplier of quality pre-owned equipment

Gold-standard solutions

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

Learn more!

Visit us at [artisanng.com](https://www.artisanng.com) for more info on price quotes, drivers, technical specifications, manuals, and documentation.

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We're here to make your life easier. How can we help you today?

(217) 352-9330 | sales@artisanng.com | [artisanng.com](https://www.artisanng.com)