

NEOS FOAOM

Two Port Coupled Acousto-Optical Modulator



Limited Availability
Used and in Excellent Condition

Open Web Page

<https://www.artisanng.com/46135-11>

All trademarks, brandnames, and brands appearing herein are the property of their respective owners.



Your **definitive** source
for quality pre-owned
equipment.

Artisan Technology Group

(217) 352-9330 | sales@artisanng.com | artisanng.com

- Critical and expedited services
- In stock / Ready-to-ship

- We buy your excess, underutilized, and idle equipment
- Full-service, independent repair center

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

OPERATING MANUAL

DOCUMENT NUMBER: 51A02562A

JUNE 2002

27 MHz INFRARED AOM

MODEL NUMBER:

26027-2-1.55-LTD

DOC. APPROVAL: W SEALE DATE: 1/14/02

NEOS TECHNOLOGIES, INC.
4005 OPPORTUNITY DRIVE
MELBOURNE, FL 32934
(321) 242-7818
FAX (321) 242-1019

TABLE OF CONTENTS

<u>SECTION</u>	<u>DESCRIPTION</u>	<u>PAGE</u>
I.	INSPECTION PROCEDURE	3
II.	DESCRIPTION	4
III.	SPECIFICATIONS	5
IV.	OUTLINE DRAWING	6
V.	CALCULATIONS	7
VI.	OPERATING INSTRUCTIONS	8
VII.	CLEANING INSTRUCTIONS	9
	ACCEPTANCE TEST RESULTS	Attached at end

SECTION I

INSPECTION PROCEDURE

Examine the shipping carton for damage. If the shipping carton or packing material is damaged it should be kept for the carrier's inspection. Notify the carrier and NEOS Technologies. Check the contents of the shipment for completeness, mechanical damage, and then test the equipment electronically. Operating procedures are contained in Section VII. If the contents are incomplete, or the equipment does not pass the electrical testing please notify NEOS Technologies.

If there is any problem with the use of this equipment, or if the equipment fails to function as expected contact NEOS Technologies, do not try to trouble shoot or repair this equipment. Consult with a NEOS service engineer. If the equipment needs repair or replacement, contact NEOS Technologies, Inc for a Return Authorization Number.


SECTION II

DESCRIPTION

26027-2-1.55-LTD

The modulator AO device is an AMTIR ($\text{Ge}_{33}\text{Se}_{55}\text{As}_{12}$) substrate with a lithium niobate transducer and is designed to be used for wavelengths between 1 and 2 μm . The RF input should not exceed 1.0 watt CW. The input laser beam diameter can be varied to produce the desired rise time, which can be as much as 520 nanoseconds with a 2 mm beam diameter. The modulator assembly should be mounted on a fixture to provide sufficient adjustment to peak the modulator efficiency.

The modulator can be driven by any good driver with a nominal 50 ohm output, however, it is recommended that a NEOS driver be used to drive this modulator for the system to achieve optimum performance. Be extremely careful not to focus the laser beam on the gold bond wires. NEOS will not warranty any such damage.



53D2200


SECTION III
SPECIFICATIONS

26027-2-1.55-LTD

<u>PARAMETER</u>	<u>SPECIFICATION</u>
Interactive Material	AMTIR
Acoustic Mode	Longitudinal
Operating Wavelength	1060-2100nm
Window Configuration	AR Coated
Operating Frequency	27.12MHz
Diffraction Efficiency	85% @ 1060nm 85% @ 1550nm 55% @ 2100nm
Light Polarization	Random
Acoustic Aperture Size	2mm
Rise Time	260ns/mm Beam Diameter
Deflection Angle	11.2mrad @ 1060nm 15.9mrad @ 1550nm 22.3mrad @ 2100nm
RF Power Level	<1Watt
Impedance	50 Ohms
VSWR	1.2:1 @ 27.12MHz
Package:	53D2200
Acceptance Test Procedure:	42A12285
Acceptance Test Results form:	52A02488
Recommended Drivers:	
Analog Driver System: 21027-1AS	Digital Driver System: 21027-1DS
Analog Driver Module: 21027-1AM	Digital Driver Module: 21027-1DM

SECTION IV

OUTLINE DRAWING



53D2200

SECTION V

CALCULATIONS

The equations to determine the AOM rise time " t_r " are as follows:

First determine the waist size by the equation,
$$d_0 = \frac{4f\lambda}{\pi d_1}$$

Where: f = lens focal length in mm

λ = the optical wavelength in 10^{-6}m

d_1 = the input optical beam diameter in mm

d_0 = the waist diameter inside the modulator in 10^{-6}m

Knowing the waist size inside the modulator, then the modulator rise time can be calculated from the relationship:

$$t_r = \frac{1.3d_0}{2V}$$

Where: V = the acoustic velocity of the modulator material which is 2510 m/s

The focal length of the lens is the $F\#$ of the lens times the input spot diameter:

$$F\# d_1 = f_{\text{lens}}$$

The deflection angle " \varnothing_d " is defined as the acoustic drive frequency in megahertz times the wavelength, divided by the acoustic velocity of the material:

$$\varnothing_d = 2\theta_{\text{Bragg}} = \frac{27.12 \times 10^6 \lambda}{2510 \text{ m/sec}}$$

Where: θ_{Bragg} = Bragg angle of the modulator.

SECTION IV. OPERATING INSTRUCTIONS

The AOM system is not polarization sensitive, so any optical polarization can be used. Be extremely careful not to focus the laser beam on the gold bond wires on the acoustic transducer, which may vaporize the bond wires. NEOS will not warranty any such damage. The modulator mount assembly should have sufficient adjustments to peak the modulator efficiency. For optimum results, the Bragg angle must be precisely adjusted.

The initial Bragg alignment can be done as follows:

If using a NEOS driver, set the mode switch is in the CW position. Make sure that the RF power does not exceed 1 watt. NEOS will not warranty any failure resulting from the application of too much RF power. Using a 50 Ohm coaxial cable, connect the "RF out" of the driver to the modulator. Turn on the RF power.

Direct a laser beam onto the optical crystal, with the optical beam close to the transducer. View the output beams on a surface one meter away using an IR viewer. Adjust the Bragg angle, by rotating the AOM, to allow the diffracted first order beam away from the transducer (connector end) to be the most intense.

Install the input lens, one f away, and adjust the height of the modulator to achieve diffraction. Make changes in the Bragg adjustment to obtain optimum efficiency. Adjust the RF driver for power level to obtain maximum diffraction efficiency. Install the output lens, one f away, to collimate the output beam.

The modulator has been designed and verified to satisfy the specifications.

To operate the modulator use the first order diffracted beam with the NEOS driver mode switch set to normal. See the driver manual for information on signal needed for modulation.

SECTION VII

OPTICAL CLEANING

Periodic cleaning of the AO device is a normal part of maintaining an optical system. When the device is installed in an optical system, make sure that there is access to allow removal of the protective cover and room to clean the device. If removal from the system is necessary, then follow the alignment procedure in this manual to reinstall, realign and, adjust the AO device.

To clean the AO device, remove the screws that hold the cover to the mount. Caution must be used when placing a screw driver near the window opening in the cover, as it is best to protect the opening with tape or cover the opening with your finger (without touching the crystal) to protect it. NEOS will not warrant any damage or scratches caused by inserting the screwdriver into the window opening.

- Remove the protective cover.
- Blow off any visible dust with canned air. Do not use an air gun unless it is filtered and water and oil free!
- Fold (4 times) a new lens tissue into a triangle to make a cleaning tool.
- Dip the tip of the lens tissue into fresh acetone or spray fresh acetone from a squeeze bottle onto it. Then shake excess fluid out of the lens tissue. Do not handle the wet area of the tissue, as your finger oil will be absorbed and contaminate the optical surface of the crystal.
- Wipe (only once) across the crystal in an even motion, starting near the transducer and drawing the tissue across the optical surface toward the other end. Do not damage the bond wires! Do not reuse the tissue as the mounting silver epoxy may be spread onto the window of the crystal.
- Repeat with a new tissue each time and for each surface that needs cleaning.
- Replace the protective cover and screws.
- Realign the device in your system and adjust the Bragg angle for maximum diffraction efficiency.

Notes:

- The lens tissue must be lint free and the best grade available.
- Only use each tissue once, for only one surface. Do not reuse the tissue, as it will redistribute the removed dust or mounting silver epoxy.
- The acetone must be electronic grade. The acetone must be fresh from a new bottle, as the acetone will absorb water from the air and cause streaks. Discard any acetone, which has been exposed to the air for more than 4 hours. If the bottle is half- empty, do not use.

Artisan Technology Group is an independent supplier of quality pre-owned equipment

Gold-standard solutions

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

Learn more!

Visit us at [artisanTG.com](https://www.artisanTG.com) for more info on price quotes, drivers, technical specifications, manuals, and documentation.

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

We're here to make your life easier. How can we help you today?

(217) 352-9330 | sales@artisanTG.com | [artisanTG.com](https://www.artisanTG.com)

