

Qioptiq Optem Zoom 70XL
7:1 Optical System



\$250.00

In Stock

Qty Available: 1

Used and in Excellent Condition

Open Web Page

<https://www.artisanng.com/52076-1>

All trademarks, brandnames, and brands appearing herein are the property of their respective owners.



Your **definitive** source
for quality pre-owned
equipment.

Artisan Technology Group

(217) 352-9330 | sales@artisanng.com | artisanng.com

- Critical and expedited services
- In stock / Ready-to-ship

- We buy your excess, underutilized, and idle equipment
- Full-service, independent repair center

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.



INSPECTION | MACHINE VISION | RESEARCH



OPTEM ZOOM 70XL 7:1 OPTICAL SYSTEM

EXCEPTIONAL RELIABILITY
250,000 Zoom Cycles Guaranteed

QIOPTIQ
Optics with Intelligence

www.qioptiqimaging.com

Qioptiq Imaging Solutions

Table of Contents

New Features	2
How to Specify	3
Upper Zoom Modules	3
Lower Function Modules	3
Illumination Options	4
Motorization Options	4
Magnification Options	4
Coaxial Performance	5
System Diagram	6
Oblique Performance	8



ZOOM 70XL



ZOOM 70XL

OPTOMECHANICAL RELIABILITY

The Zoom 70XL is specifically designed to stand up to the grueling conditions of today's automated imaging.

"XL" Stands For Xtended Life

Backed by superior engineering and only the highest-quality components and materials, the Optem™ Zoom 70XL delivers a guaranteed 250,000 cycles without mechanical failure! This robust reliability makes the Zoom 70XL the ideal solution for demanding and mission critical applications.

**250,000
CYCLES
GUARANTEED!**



Zoom 70XL Specs

Zoom Range	7:1
Magnification Range	0.75 – 5.25X
Resolution	72 – 240 lp/mm
Numerical Aperture	0.0240 – 0.080
Depth-of-Field	0.98 – 0.09 mm
Field-of-View (Low-Mag.)	6.40 x 8.53 mm
Field-of-View (High-Mag.)	0.91 x 1.22 mm
Working Distance	89 mm

Nominal 1X/1X Configuration with 1/2" Camera

NEW FOR ZOOM 70XL

NEW LED Coaxial Illuminators

Programmable 1-Watt white-light Coaxial LED Illuminators feature compact design and lightweight cable for streamlined OEM integration. **PAGE 4**

NEW Dual-Magnification Accessory


Now you can simultaneously integrate two upper lens modules with independent cameras over the same subject. Mix and match Zoom 70 and FMOS upper modules for added versatility. **PAGE 4**

NEW Long-Working Right Angle Module

Integrate a 90° turn in the optical path above the Lower Function Module. Unlike the below system right-angle adapter, this unique accessory affords configuration flexibility to integrate all lower functions without sacrificing any working distance. **PAGE 4**

NEW Optem M-Plan APO Objectives

Get the most from your high magnification imaging. Specify the newly expanded line of Optem Long-Working Distance Objectives, including the new M-Plan APO family for exceptional color accuracy and field-flatness. **PAGE 4**



HOW TO SPECIFY YOUR ZOOM 70XL LENS SYSTEM

Follow the seven easy steps below to configure the Zoom 70XL Lens System that meets your optical, functional and physical configuration needs.

- Referring to the **Optical Performance Matrix for Substage or Oblique Illumination** on the back cover, locate your desired performance parameters (*typically field-of-view at high- and low-end zoom as relates to your intended camera format 1/4" – 1"*). If you desire coaxial illumination, work from the **Optical Performance Matrix for Coaxial illumination**. **PAGE 5**

1.0x TV Tube		Low Mag		High Mag	
NO LENS REQ'D		0.01		0.001	
1.0x Confa: WD: 89 mm	NA	0.019		0.10	
	LP/mm	57		300	
	Mag	0.52		6.5	
	DOF	1.57		0.057	
	FOV	5.31 x 7.08		0.42 x 0.57	
Cam Format	1/4"	8.92 x 9.23		0.55 x 0.74	
	1/2"	9.23 x 12.31		0.74 x 0.98	
	2/3"	11.00 x 14.67		1.02 x 1.35	
	1"	14.80 x 19.40		1.48 x 1.97	

- Cross reference your desired performance block and note the **TV Tube magnification (TOP axis)** and **Auxiliary Lens magnification (LEFT axis)** to arrive at the combination of magnification components needed to yield your desired field-of-view [FOV]... consider the varying resolution, depth of field and/or working distance factors that are dictated by the magnification components.

- Once you have determined the **TV Tube** and **Auxiliary Lens magnifications** required, navigate down through the System Diagram (**PAGE 6-7**) starting from the camera. Select the physical shape of TV Tube and/or camera mount type desired.

- Select the **Upper Zoom Module** with the functional features you require... Manual, Iris, detent and DC or Stepper motorized Zoom options are available.

- Combine the Upper Zoom Module with one of the **Lower Function Modules** that incorporates manual or motorized fine focus and/or coaxial illumination options. If you have chosen both motorized zoom and focus, the compact frame constraints of Zoom 70XL require that you order a Dual-Motorized Non-Modular Zoom 70XL (see System Diagram) **PAGE 6-7**

- Now Specify the correct **Auxiliary Lens magnification** to complete the optical components of your Zoom 70XL Lens System as dictated in step 1 above. By default, the Zoom 70XL yields a 1X configuration without an Auxiliary Lens.

- Finally outfit your Zoom 70XL Lens System with the appropriate **Coaxial or Ringlight Illumination accessories**. If you chose motorized zoom or focus functions, select the appropriate **DC or Stepper controller and power supply**. You may also wish to specify **mounting hardware** to ease integration.

Components of Functional Zoom 70XL

• TV Tube (Top Magnification)

• Upper Zoom Module

• Lower Function Module

• Auxiliary Lens (Bottom Magnification)

UPPER ZOOM MODULES

The specific Upper Zoom Module selected will determine the type of zoom operation for your Zoom 70XL System.



Manual Module – Provides basic hand-driven 7:1 zoom function.



Iris Diaphragm Module –

Provides manual 7:1 zoom while also allowing better illumination control. With Iris, you can increase depth-of-field and/or maintain more consistent image brightness across your entire field.

Detent Module – Obtain specific and repeatable magnification stops throughout the 7:1 zoom range without the complexity and cost of motorization. Detents are ideal for metrology applications where each position can be calibrated.



Motorized Zoom Module – Provides programmable automated zoom in either DC or Stepper Motorized Versions. All Stepper Motorized Models come complete with programmable homing sensors.



LOWER FUNCTION MODULES

When configuring your Zoom 70XL Lens System, you may choose between the following Manual or motorized Lower Function Module options:

- Standard Module
- Internal Focus Module
- Coaxial Illumination Module
- Internal Focus with Coaxial Module

Variable Working Distance Auxiliary Lens

- Intended for applications requiring working distances between 127mm (5") and 432mm (17"). Covers magnification factors 0.288X at 127mm, and 0.109X at 432mm. See "Variable Working Distance Auxiliary Lens" Chart available as a download at www.qioptiqimaging.com.



Objective Modules - Incorporate Infinity-Corrected Objectives utilizing one of four Objective Lower Function Modules. All are equipped with M26 x 36T threads to accept Optem and Mitutoyo Objectives. An RMS thread adapter is available for use with Nikon and Olympus Objectives.

DUAL MOTORIZED NON-MODULAR MODELS

With it's extremely space efficient design, the Zoom 70XL must be ordered in a non-modular configuration when both motorized zoom and focus functions are desired. Select from four single-bodied Dual-Motorized Zoom 70XLs in either DC or Stepper motorization with various illumination options.





ILLUMINATION OPTIONS

The Zoom 70XL Lens System offers a variety of illumination options to meet a variety of imaging requirements.

LED Coaxial Illuminators

- Compatible with all Zoom 70XL Coaxial Lower Function Modules, new Programmable 1-Watt LED Coaxial illuminators offer reduced power requirement and heat generation with substantial service life gains.



Available in Straight and Right-angle configurations, Optem Coaxial LEDs feature compact designs and simplified cable management. Driven by a compact single-channel programmable controller, Optem Coaxial LEDs emit brilliant cool light in the visible spectrum and deliver virtually identical optical performance to our traditional Fiber Optic Coaxial illuminator.

Fiber Optic Coaxial Illuminators - Available with 40 and 60-inch flexible fiber bundles powered by 110V or 220V Optem VSI Fiber Optic Illuminators.

Polarized Light - When imaging highly reflective subjects, Polarizer Modules with built-in Analyzers are available to introduce polarization to both LED and Fiber Optic Coaxial Illumination paths.

Fiber Optic Ring Lights - Oblique ringlight Illumination is ideal to better define features of dimension-rich subjects. Driven with the same VSI Fiber Optic Illuminator as conventional coaxial systems, Optem Ringlights feature optional horizontal, vertical and Objective configurations.

CONFIGURATION ACCESSORIES

Extend the versatility of your Zoom 70XL Lens System with two new innovations in form and function.

Long-Working, Right-Angle Module -

This accessory allows the introduction of a 90° of your optical axis at the mid-body point of your Zoom 70XL Lens System. Users can now integrate any Zoom 70XL Lower Function Module below the right-angle turn maintaining full working distance and affording greater configuration flexibility.



NEW Dual-Magnification Module -

Simultaneously integrate two Zoom 70 or Optem FMOS upper modules over one subject. Mix and match fixed and zoom magnifications, camera types/formats. Choose from the full range of Zoom 70XL Lower Function Modules or integrate compact Coaxial illumination through the coaxial block.



MOTORIZING ZOOM 70XL



The Zoom 70XL Optical System can be specified with motorized zoom and/or focus functions.

Stepper Motor Models - Stepper Motors feature Hall-Effect Sensors to ensure pinpoint mechanical accuracy and repeatability. Control and power are easily integrated through RS232 serial port interface with a Rocker Switch controller and Windows GUI. An OEM Board is available for more streamlined OEM integration.. A Stepper Motor VI library is available for LabVIEW.



DC Motor Models



DC Motors allow continuous movement throughout either the zoom or focus range. Two RS232 Serial ports in the back of the joystick DC Motor Controller accept both zoom and focus motorized functions.

OPTEM LONG-WORKING DISTANCE OBJECTIVES

NEW Achieve significantly higher magnifications and increased resolution. Combine your Zoom 70XL with the expanded line of Optem infinity-corrected objectives.



Optem M-Plan APO - Eliminate Chromatic aberration across exceptionally flat fields for the ultimate in high-magnification accuracy. Select from 2X, 5X, 10X, 20X and 50X Long-Working Distance Objectives. These objectives are exact replacements for Mitutoyo 378 series objectives and are ideal for metrology applications.

Optem High-Resolution -

Specifically designed to capture maximum resolution at the high-end magnifications of Optem Zoom Lenses, the 5X, 10X and 20X Optem HR Objectives are ideal for applications where distinguishing every finite detail is critical.



Coaxial Illumination Zoom 70XL Optical Performance Matrix

NOTE: 0.18x, 0.25x and 0.5x Aux Lenses NOT recommended w/ coax illumination		0.375x TV Tube		0.5x TV Tube		0.67x TV Tube		1.0x TV Tube		1.5x TV Tube		2.0x TV Tube	
		Low Mag	High Mag	Low Mag	High Mag	Low Mag	High Mag	Low Mag	High Mag	Low Mag	High Mag	Low Mag	High Mag
0.75X Aux Lens 29-20-39-000 WD: 114 mm	NA	0.018	0.060	0.018	0.060	0.018	0.060	0.018	0.060	0.018	0.060	0.018	0.060
	LP/mm	54	180	54	180	54	180	54	180	54	180	54	180
	Mag	0.21	1.47	0.28	1.96	0.38	2.66	0.56	3.92	0.85	5.95	1.10	7.70
	DOF	1.75	0.16	1.75	0.16	1.75	0.16	1.75	0.16	1.75	0.16	1.75	0.16
	FOV	7.70 x 10.26	1.88 x 2.50	7.01 x 9.35	1.41 x 1.88	6.10 x 8.14	1.04 x 1.38	4.91 x 6.54	0.70 x 0.94	3.27 x 4.36	0.46 x 0.62	2.45 x 3.27	0.36 x 0.48
	Cam Format	1/4"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"
	1/3"	7.76 x 10.34	2.45 x 3.27	8.00 x 10.67	1.84 x 2.45	6.82 x 9.10	1.35 x 1.80	5.82 x 7.76	0.92 x 1.22	4.27 x 5.69	0.61 x 0.81	3.20 x 4.27	0.47 x 0.62
	1/2"	-	-	-	-	7.76 x 10.34	2.45 x 3.27	7.96 x 10.61	1.80 x 2.41	5.57 x 6.90	0.81 x 1.08	3.27 x 5.69	0.62 x 0.83
	2/3"	-	-	-	-	-	-	7.96 x 10.61	2.48 x 3.31	7.82 x 10.43	1.68 x 2.24	6.52 x 6.66	1.11 x 1.48
	1"	-	-	-	-	-	-	-	-	7.76 x 10.34	2.45 x 3.27	7.59 x #####	1.61 x 2.15
1.0X Config. NO LENS REQ'D WD: 89 mm	NA	0.024	0.080	0.024	0.080	0.024	0.080	0.024	0.080	0.024	0.080	0.024	0.080
	LP/mm	72	240	72	240	72	240	72	240	72	240	72	240
	Mag	0.29	2.03	0.38	2.66	0.50	3.50	0.75	5.25	1.10	7.70	1.50	10.50
	DOF	0.98	0.089	0.98	0.089	0.98	0.089	0.98	0.089	0.98	0.089	0.98	0.089
	FOV	8.18 x 10.90	1.36 x 1.81	7.36 x 9.81	1.04 x 1.38	5.49 x 7.32	0.79 x 1.05	3.68 x 4.91	0.53 x 0.70	2.45 x 3.27	0.36 x 0.48	1.84 x 2.45	0.26 x 0.35
	Cam Format	1/4"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"
	1/4"	6.74 x 8.98	1.77 x 2.36	8.73 x 11.64	1.35 x 1.80	7.16 x 9.55	1.03 x 1.37	4.80 x 6.40	0.69 x 0.91	3.20 x 4.27	0.47 x 0.62	2.40 x 3.20	0.34 x 0.46
	1/2"	-	-	-	-	6.74 x 8.98	1.80 x 2.41	6.88 x 11.58	1.37 x 1.83	6.40 x 8.53	0.91 x 1.22	4.27 x 5.69	0.62 x 0.83
	2/3"	-	-	-	-	-	-	6.57 x 8.76	1.89 x 2.51	8.80 x 11.73	1.26 x 1.68	5.87 x 7.82	0.86 x 1.14
	1"	-	-	-	-	-	-	-	-	6.74 x 8.98	1.83 x 2.44	8.53 x #####	1.25 x 1.66
1.5X Aux Lens 29-20-40-000 WD: 52 mm	NA	0.036	0.12	0.036	0.12	0.036	0.12	0.036	0.12	0.036	0.12	0.036	0.12
	LP/mm	108	360	108	360	108	360	108	360	108	360	108	360
	Mag	0.43	3.00	0.56	3.90	0.75	5.25	1.10	7.70	1.70	11.90	2.30	16.10
	DOF	0.44	0.039	0.44	0.039	0.44	0.039	0.44	0.039	0.44	0.039	0.44	0.039
	FOV	5.03 x 6.71	0.92 x 1.23	4.91 x 6.54	0.71 x 0.94	3.66 x 4.88	0.53 x 0.70	2.45 x 3.27	0.36 x 0.48	1.64 x 2.18	0.23 x 0.31	1.23 x 1.64	0.17 x 0.23
	Cam Format	1/4"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"
	1/4"	4.06 x 5.42	1.20 x 1.60	4.92 x 6.56	0.92 x 1.23	4.78 x 6.37	0.69 x 0.91	3.20 x 4.27	0.47 x 0.62	2.13 x 2.84	0.30 x 0.40	1.60 x 2.13	0.22 x 0.30
	1/2"	-	-	-	-	4.06 x 5.42	1.23 x 1.64	4.90 x 6.53	0.91 x 1.22	2.84 x 3.79	0.40 x 0.54	2.13 x 2.84	0.30 x 0.40
	2/3"	-	-	-	-	-	-	3.98 x 5.31	1.26 x 1.68	5.33 x 7.11	0.86 x 1.14	2.93 x 3.91	0.41 x 0.55
	1"	-	-	-	-	-	-	-	-	4.06 x 5.42	1.25 x 1.66	5.69 x 7.59	0.81 x 1.08
2.0X Aux Lens 29-20-41-000 WD: 32 mm	NA	0.048	0.16	0.048	0.16	0.048	0.16	0.048	0.16	0.048	0.16	0.048	0.16
	LP/mm	144	480	144	480	144	480	144	480	144	480	144	480
	Mag	0.57	3.99	0.75	5.25	1.00	7.00	1.50	10.50	2.30	16.10	3.00	21.00
	DOF	0.25	0.022	0.25	0.022	0.25	0.022	0.25	0.022	0.25	0.022	0.25	0.022
	FOV	2.34 x 3.12	0.69 x 0.92	2.45 x 3.27	0.53 x 0.70	2.50 x 3.33	0.39 x 0.53	1.84 x 2.45	0.26 x 0.35	1.23 x 1.64	0.17 x 0.23	0.92 x 1.23	0.13 x 0.18
	Cam Format	1/4"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"
	1/4"	2.06 x 2.75	0.90 x 1.20	2.29 x 3.05	0.69 x 0.91	2.39 x 3.18	0.51 x 0.69	2.40 x 3.20	0.34 x 0.46	1.60 x 2.13	0.22 x 0.30	1.20 x 1.60	0.17 x 0.23
	1/2"	-	-	-	-	2.27 x 3.03	0.69 x 0.91	2.46 x 3.28	0.46 x 0.61	2.13 x 2.84	0.30 x 0.40	1.60 x 2.13	0.23 x 0.30
	2/3"	-	-	-	-	-	-	2.32 x 3.09	0.63 x 0.84	2.67 x 3.56	0.41 x 0.55	2.20 x 2.93	0.31 x 0.42
	1"	-	-	-	-	-	-	-	-	2.37 x 3.16	0.60 x 0.80	2.46 x 3.28	0.46 x 0.61
OPTM Infinity-Corrected Objectives													
OPTM 2X M Plan Apo 28-21-02-000 WD: 34 mm	NA	0.020	0.060	0.020	0.060	0.020	0.060	0.020	0.060	0.020	0.060	0.020	0.060
	LP/mm	60	180	60	180	60	180	60	180	60	180	60	180
	Mag	0.27	1.89	0.36	2.52	0.48	3.38	0.72	5.04	1.08	7.56	1.44	10.08
	DOF	1.01	0.19	1.01	0.19	1.01	0.19	1.01	0.19	1.01	0.19	1.01	0.19
	FOV	3.88 x 5.17	1.46 x 1.95	4.26 x 5.68	1.10 x 1.46	5.20 x 6.94	0.82 x 1.09	3.83 x 5.11	0.55 x 0.73	2.56 x 3.41	0.37 x 0.49	1.92 x 2.55	0.27 x 0.37
	Cam Format	1/4"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"
	1/4"	-	-	-	-	3.85 x 5.13	1.43 x 1.90	4.39 x 5.85	1.07 x 1.42	5.00 x 6.67	0.71 x 0.95	2.22 x 2.96	0.36 x 0.48
	1/2"	-	-	-	-	-	-	3.83 x 5.10	1.42 x 1.90	4.44 x 5.93	0.95 x 1.27	3.33 x 4.44	0.48 x 0.63
	2/3"	-	-	-	-	-	-	-	-	3.99 x 5.31	1.31 x 1.75	4.70 x 6.27	0.87 x 1.16
	1"	-	-	-	-	-	-	-	-	-	-	4.04 x 5.39	1.27 x 1.69
OPTM 5X M Plan Apo 28-21-05-000 WD: 32 mm	NA	0.06	0.14	0.06	0.14	0.06	0.14	0.06	0.14	0.06	0.14	0.06	0.14
	LP/mm	180	420	180	420	180	420	180	420	180	420	180	420
	Mag	0.68	4.88	0.89	6.41	1.19	8.60	1.78	12.83	2.66	19.24	3.55	25.66
	DOF	0.16	0.029	0.16	0.03	0.16	0.029	0.16	0.029	0.16	0.029	0.16	0.029
	FOV	2.24 x 2.99	0.57 x 0.75	2.56 x 3.41	0.43 x 0.57	2.29 x 3.05	0.32 x 0.43	1.53 x 2.04	0.22 x 0.29	1.02 x 1.36	0.14 x 0.19	0.77 x 1.02	0.11 x 0.14
	Cam Format	1/4"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"
	1/4"	1.88 x 2.51	0.74 x 0.98	2.22 x 2.96	0.56 x 0.75	2.71 x 3.62	0.42 x 0.56	2.00 x 2.67	0.28 x 0.37	1.33 x 1.78	0.19 x 0.25	0.88 x 1.19	0.14 x 0.19
	1/2"	-	-	-	-	1.84 x 2.45	0.75 x 1.00	2.34 x 3.12	0.56 x 0.74	1.78 x 2.37	0.25 x 0.33	1.33 x 1.78	0.19 x 0.25
	2/3"	-	-	-	-	-	-	-	-	2.44 x 3.26	0.34 x 0.46	1.83 x 2.44	0.26 x 0.34
	1"	-	-	-	-	-	-	-	-	-	-	2.37 x 3.16	0.50 x 0.67
OPTM 10X M Plan Apo 28-21-10-000 WD: 33 mm	NA	0.12	0.30	0.12	0.30	0.12	0.30	0.12	0.30	0.12	0.30	0.12	0.30
	LP/mm	343	867	343	867	343	867	343	867	343	867	343	867
	Mag	1.35	9.75	1.78	12.83	2.38	17.19	3.55	25.66	5.33	38.49	7.11	49.74
	DOF	0.040	0.0063	0.040	0.0063	0.040	0.0063	0.040	0.0063	0.040	0.0063	0.040	0.0063
	FOV	1.06 x 1.42	0.28 x 0.38	1.18 x 1.57	0.22 x 0.29	1.14 x 1.53	0.16 x 0.21	0.77 x 1.02	0.11 x 0.14	0.51 x 0.68	0.07 x 0.10	0.38 x 0.51	0.06 x 0.07
	Cam Format	1/4"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"
	1/4"	-	-	-	-	1.05 x 1.40	0.28 x 0.37	1.24 x 1.66	0.21 x 0.28	1.00 x 1.33	0.14 x 0.19	0.80 x 1.07	0.07 x 0.10
	1/2"	-	-	-	-	-	-	1.11 x 1.47	0.28 x 0.37	1.03 x 1.37	0.19 x 0.25	0.67 x 0.89	0.10 x 0.13
	2/3"	-	-	-	-	-	-	-	-	0.92 x 1.22	0.26 x 0.34	1.22 x 1.63	0.17 x 0.23
	1"	-	-	-	-	-	-	-	-	-	-	1.11 x 1.48	0.25 x 0.33
OPTM 20X M Plan Apo 28-21-11-000 WD: 20 mm	NA	0.24	0.39	0.24	0.39	0.24	0.39	0.24	0.39	0.24	0.39	0.24	0.39
	LP/mm	686	1128	686	1128	686	1128	686	1128	686	1128	686	1128
	Mag	2.70	19.50	3.55	25.66	4.76	34.38	7.11	51.32	10.66	76.97	14.21	99.47
	DOF	0.010	0.0037	0.010	0.0037	0.010	0.0037	0.010	0.0037	0.010	0.0037	0.010	0.0037
	FOV	0.48 x 0.64	0.14 x 0.19	0.51 x 0.68	0.11 x 0.14	0.57 x 0.76	0.08 x 0.11	0.38 x 0.51	0.05 x 0.07	0.26 x 0.34	0.04 x 0.05	0.19 x 0.26	0.028 x 0.037
	Cam Format	1/4"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"	1/3"
	1/4"	-	-	-	-	0.48 x 0.63	0.14 x 0.19	0.50 x 0.66	0.10 x 0.14	0.45 x 0.61	0.07 x 0.09	0.33 x 0.44	0.036 x 0.048
	1/2"	-	-	-	-	-	-	0.47 x 0.63	0.14 x 0.19	0.42 x 0.56	0.09 x 0.12	0.33 x 0.44	0.048 x 0.064
	2/3"	-	-	-	-	-	-	-	-	0.38 x 0.51	0.13 x 0.17	0.56 x 0.74	0.09 x 0.12
	1"	-	-	-	-	-	-	-	-	-	-	0.49 x 0.66	0.12 x 0.17
OPTM 50X M Plan Apo 28-21-50-000 WD: 13 mm	NA	0.48	0.48	0.48	0.48	0.48	0.48	0.48	0.48	0.48	0.48	0.48	0.48
	LP/mm	1440	1440	1440	1440	1440	1440	1440	1440	1440	1440	1440	1440
	Mag	6.75	48.75	9.00	65.00	12.06	87.10	18.00	130.00	27.00	195.00	36.00	252.00
	DOF	0.0025	0.0025	0.0025	0.0025	0.0025	0.0025	0.0025	0.0025	0.0025	0.0025	0.0025	0.0025
	FOV	0.22 x 0.30	0.06 x 0.08	0.24 x 0.31	0.04 x 0.06	0.23 x 0.31	0.03 x 0.04	0.15 x 0.20	0.02 x 0.03	0			

Yellow -- Limited vignetting at low magnifications

Grey -- Configuration not recommended

Tube Clamp Mounts

29-50-10-000
UNIVERSAL CLAMP
(top view)



30-65-34-000
FLAT MOUNT
(for Motorized versions)



30-25-76-000
76mm O.D. TUBE CLAMP
(OLYMPUS, NIKON, LEICA)

30-25-82-000
82mm O.D. TUBE CLAMP
(UNITRON)

30-25-83-000
83mm O.D. TUBE CLAMP
(GENERIC)

30-25-84-000
84mm O.D. TUBE CLAMP
(MEIJI)




Right-Angle TV Tubes

29-90-86-000 - 0.5X RA
29-90-87-000 - 0.67X RA
29-90-74-000 - 1.0X RA
29-90-88-000 - 1.5X RA
29-90-89-000 - 2.0X RA

U-Bend TV Tube

29-90-81-000
1.0X U-BEND
TV TUBE



29-90-26-000
0.5X MINI
TV TUBE



29-90-27-000
0.67X MINI
TV TUBE



29-90-90-000
1.0X MINI
TV TUBE


29-90-28-000
1.5X MINI
TV TUBE




Non-Modular

Dual-Motorized Zooms


39-26-20-000
70XL 12mm FOCUS
DC MOTORIZED BOTH




39-26-37-000
70XL 12mm FOCUS
STEPPING MOTOR BOTH
HALL EFFECT SENSOR



39-28-20-000
70XL 3mm COAX. DC
MOTORIZED BOTH



39-28-37-000
70XL 3mm COAX.
STEPPING MOTOR BOTH
HALL EFFECT SENSOR



29-60-83-000
VSI 110 ILLUMINATOR
110 VOLT / 150 WATT



29-60-84-000
VSI 220 ILLUMINATOR
220VOLT / 150 WATT

30-16-40-000
40" FLEX
BUNDLE

30-16-60-000
60" FLEX
BUNDLE



33-04-00-000
LOW PROFILE, R.A.
COAXIAL LAMPHOUSE
(Replaces standard
bottom module lamphouse)



30-16-01-000
10MM ADAPTER



Illumination

396010-810
SINGLE-CHANNEL
PROGRAMMABLE
LED CONTROLLER



39-60-10-000
STRAIGHT
1-WATT COAXIAL LED



NEW



39-60-20-000
FMOS RIGHT-ANGLE
1-WATT COAXIAL LED

27-60-07-000
110/6 VOLT VARIABLE

27-60-17-000
220/6 VOLT VARIABLE



29-60-09-000
10 WATT
HALOGEN LAMP



Lower Function Modules


29-25-15-000
ZOOM 70XL VARIABLE
AUXILIARY LENS




299510-310
BASIC FUNCTION
MODULE



299591-310
INTERNAL FOCUS
MODULE



296515-310
VERTICAL ILLUM.
MODULE



296591-310
VERTICAL ILLUM &
FOCUS MODULE



33-03-03-000
RING LIGHT ADAPTER



30-16-70-000
VERTICAL RING LIGHT



29-60-81-000
HORIZONTAL
RING LIGHT



29-60-31-000
MINI-ANNULAR
RING PROBE



29-20-09-000
0.18X / 0.25X
CONVERTIBLE
AUX LENS



29-20-11-000
0.5X
AUX LENS



29-20-12-000
0.75X
AUX LENS




29-20-13-000
1.5X
AUX LENS



CAMERA

Standard TV Tubes




Non-Inverting Right-Angle TV Tubes

29-90-92-000 - 0.5X NIRA
29-90-93-000 - 0.67X NIRA
29-90-94-000 - 1.0X NIRA
29-90-95-000 - 1.5X NIRA
29-90-96-000 - 2.0X NIRA



Bayonet TV Tubes

29-77-04-000 - 1.0X F-MOUNT BAYONET
29-77-02-000 - 1.0X 1/2" ENG BAYONET
29-77-03-000 - 1.0X 2/3" ENG BAYONET
29-77-01-000 - 1.0X 1/2" SONY BAYONET

Upper Zoom Modules



Motorization Control & Power Supply



NEW
90°

NEW
x2

29-90-67-000
DUAL MAGNIFICATION
MODULE

Can also integrate FMOS
Upper Fixed Tubes. Refer to
Optem FMOS Brochure

29-90-66-000
LONG WORKING
RIGHT-ANGLE
MODULE

299067-310
COAX
MODULE

299066-304
FMOS / ZOOM 70
LOWER MODULE
COUPLER

299066-306
ZOOM 125
LOWER MODULE
COUPLER

Integrates with Optem Zoom 125C
Lower Function Modules. Refer to Zoom 125C Brochure

299066-302
M25x0.75
ADAPTER

299066-303
OPTEM / MITU.
OBJECTIVES
COUPLER

299066-305
NIKON / OLY
OBJECTIVES
COUPLER

296084-310
VERT. ILLUM.
OBJ. DUAL FOCUS


296085-310
VERT. ILLUMINATOR
OBJECTIVE MODULE

299084-310
OBJECTIVE DUAL
FOCUS MODULE

299085-310
OBJECTIVE ADAPTER
MODULE

29-20-31-000
2X CONVERSION LENS

Auxiliary Lenses



MITUTOYO

Optem® LWD Objectives
28-20-44-000 - 5X HIGH-RES
28-20-45-000 - 10X HIGH-RES
28-20-46-000 - 20X HIGH-RES

Optem® LWD Objectives

28-21-02-000 - 2X M PLAN APO
28-21-05-000 - 5X M PLAN APO
28-21-10-000 - 10X M PLAN APO
28-21-11-000 - 20X M PLAN APO
28-21-50-000 - 50X M PLAN APO

29-20-32-000
RMS ADAPTER

NIKON/OLYMPUS

30-16-50-000
OBJECTIVE RING LIGHT

Substage/Ringlight Illumination Zoom 70XL Optical Performance Matrix

		0.375x TV Tube		0.5x TV Tube		0.67x TV Tube		1.0x TV Tube		1.5x TV Tube		2.0x TV Tube	
		Low Mag	High Mag	Low Mag	High Mag	Low Mag	High Mag	Low Mag	High Mag	Low Mag	High Mag	Low Mag	High Mag
0.18x Aux Lens 29-20-09-000 WD: 468 mm	NA	0.0043	0.015	0.0043	0.015	0.0043	0.015	0.0043	0.015	0.0043	0.015	0.0043	0.015
	LP/mm	13	45	13	45	13	45	13	45	13	45	13	45
	Mag	0.05	0.36	0.07	0.47	0.09	0.63	0.14	0.94	0.20	1.40	0.27	1.90
	DOF	30.67	2.52	30.67	2.52	30.67	2.52	30.67	2.52	30.67	2.52	30.67	2.52
	FOV	1/4"	54.52 x 72.69	7.67 x 10.22	40.89 x 54.52	5.87 x 7.83	30.51 x 40.69	4.38 x 5.84	20.44 x 27.26	2.94 x 3.91	13.63 x 18.17	1.97 x 2.63	10.22 x 13.63
0.25x Aux Lens 29-20-09-000 WD: 310 mm	NA	0.0060	0.020	0.0060	0.020	0.0060	0.020	0.0060	0.020	0.0060	0.020	0.0060	0.020
	LP/mm	18	60	18	60	18	60	18	60	18	60	18	60
	Mag	0.07	0.50	0.09	0.66	0.13	0.88	0.19	1.30	0.28	2.00	0.38	2.60
	DOF	15.75	1.42	15.75	1.42	15.75	1.42	15.75	1.42	15.75	1.42	15.75	1.42
	FOV	1/4"	39.25 x 52.34	5.52 x 7.36	29.44 x 39.25	4.18 x 5.58	21.97 x 29.29	3.14 x 4.18	14.72 x 19.63	2.12 x 2.83	9.81 x 13.08	1.38 x 1.84	7.36 x 9.81
0.5x Aux Lens 29-20-38-000 WD: 178 mm	NA	0.0120	0.040	0.0120	0.040	0.0120	0.040	0.0120	0.040	0.0120	0.040	0.0120	0.040
	LP/mm	36	120	36	120	36	120	36	120	36	120	36	120
	Mag	0.14	0.98	0.19	1.33	0.25	1.75	0.38	2.66	0.56	3.92	0.75	5.25
	DOF	3.94	0.35	3.94	0.35	3.94	0.35	3.94	0.35	3.94	0.35	3.94	0.35
	FOV	1/4"	19.63 x 26.17	2.82 x 3.76	14.72 x 19.63	2.08 x 2.77	10.99 x 14.65	1.58 x 2.10	7.36 x 9.81	1.04 x 1.38	4.91 x 6.54	0.70 x 0.94	3.68 x 4.91
0.75x Aux Lens 29-20-39-000 WD: 114 mm	NA	0.0180	0.060	0.0180	0.060	0.0180	0.060	0.0180	0.060	0.0180	0.060	0.0180	0.060
	LP/mm	54	180	54	180	54	180	54	180	54	180	54	180
	Mag	0.21	1.47	0.28	1.96	0.38	2.66	0.56	3.92	0.85	5.95	1.10	7.70
	DOF	1.75	0.16	1.75	0.16	1.75	0.16	1.75	0.16	1.75	0.16	1.75	0.16
	FOV	1/4"	13.08 x 17.45	1.88 x 2.50	9.81 x 13.08	1.41 x 1.88	7.32 x 9.76	1.04 x 1.38	4.91 x 6.54	0.70 x 0.94	3.27 x 4.36	0.46 x 0.62	2.45 x 3.27
1.0x Config. IO LENS REQ'D WD: 89 mm	NA	0.0240	0.080	0.0240	0.080	0.0240	0.080	0.0240	0.080	0.0240	0.080	0.0240	0.080
	LP/mm	72	240	72	240	72	240	72	240	72	240	72	240
	Mag	0.29	2.03	0.38	2.66	0.50	3.50	0.75	5.25	1.10	7.70	1.50	10.50
	DOF	0.98	0.09	0.98	0.09	0.98	0.09	0.98	0.09	0.98	0.09	0.98	0.09
	FOV	1/4"	9.81 x 13.08	1.36 x 1.81	7.36 x 9.81	1.04 x 1.38	5.49 x 7.32	0.79 x 1.05	3.68 x 4.91	0.53 x 0.70	2.45 x 3.27	0.36 x 0.48	1.84 x 2.45
1.5x Aux Lens 29-20-40-000 WD: 52 mm	NA	0.0360	0.120	0.0360	0.120	0.0360	0.120	0.0360	0.120	0.0360	0.120	0.0360	0.120
	LP/mm	108	360	108	360	108	360	108	360	108	360	108	360
	Mag	0.43	3.00	0.56	3.90	0.75	5.25	1.10	7.70	1.70	11.90	2.30	16.10
	DOF	0.44	0.04	0.44	0.04	0.44	0.04	0.44	0.04	0.44	0.04	0.44	0.04
	FOV	1/4"	6.54 x 8.72	0.92 x 1.23	4.91 x 6.54	0.71 x 0.94	3.66 x 4.88	0.53 x 0.70	2.45 x 3.27	0.36 x 0.48	1.64 x 2.18	0.23 x 0.31	1.23 x 1.64
2.0x Aux Lens 29-20-41-000 WD: 32 mm	NA	0.0480	0.160	0.0480	0.160	0.0480	0.160	0.0480	0.160	0.0480	0.160	0.0480	0.160
	LP/mm	144	480	144	480	144	480	144	480	144	480	144	480
	Mag	0.57	3.99	0.75	5.25	1.00	7.00	1.50	10.50	2.30	16.10	3.00	21.00
	DOF	0.25	0.02	0.25	0.02	0.25	0.02	0.25	0.02	0.25	0.02	0.25	0.02
	FOV	1/4"	4.91 x 6.54	0.69 x 0.92	3.68 x 4.91	0.53 x 0.70	2.75 x 3.66	0.39 x 0.53	1.84 x 2.45	0.26 x 0.35	1.23 x 1.64	0.17 x 0.23	0.92 x 1.23

— Limited vignetting at low magnifications

— Configuration not recommended



서울특별시 금천구 가산동 680 우림라이온스밸리2차 307호

Tel : (02)2027-2078

Web : www.zast.net

Fax : (02)2027-2025

E-mail : info@zast.net



ISO 9001:2000

78 Schuyler Baldwin Drive • Fairport, NY USA 14450-9196

Tel: (585) 223-2370 • Fax: (585) 223-3413 • www.qioptiqimaging.com

Artisan Technology Group is an independent supplier of quality pre-owned equipment

Gold-standard solutions

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

Learn more!

Visit us at [artisan^{tg}.com](https://www.artisantg.com) for more info on price quotes, drivers, technical specifications, manuals, and documentation.

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

We're here to make your life easier. How can we help you today?

(217) 352-9330 | sales@artisan^{tg}.com | [artisan^{tg}.com](https://www.artisantg.com)

