Artisan Technology Group is your source for quality new and certified-used/pre-owned equipment

- FAST SHIPPING AND DELIVERY
- TENS OF THOUSANDS OF IN-STOCK ITEMS
- EQUIPMENT DEMOS
- HUNDREDS OF MANUFACTURERS SUPPORTED
- LEASING/MONTHLY RENTALS
- ITAR CERTIFIED SECURE ASSET SOLUTIONS

SERVICE CENTER REPAIRS

Experienced engineers and technicians on staff at our full-service, in-house repair center

*InstraView*SM REMOTE INSPECTION

Remotely inspect equipment before purchasing with our interactive website at www.instraview.com ↗

WE BUY USED EQUIPMENT

Sell your excess, underutilized, and idle used equipment. We also offer credit for buy-backs and trade-ins. www.artisanng.com/WeBuyEquipment ↗

LOOKING FOR MORE INFORMATION?

Visit us on the web at www.artisanng.com ↗ for more information on price quotations, drivers, technical specifications, manuals, and documentation

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Operating Manual

Please read carefully before use,
and follow all operating and safety
instructions!

BIOHIT

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisanTG.com

Contents

Safety Instructions	4
Components	6
Application and Design	7
Materials	7
Operating Exclusions	8
Operating Limitations	8
Preparing for use	10
Connecting the discharge tube	11
Connecting the filling tube	12
Mounting instrument on bottle	13
Priming the instrument	14
Dispensing	16
Setting the volume	16
Serial dispensing with flexible discharge tube (optional)	18
Connecting drying tube (optional)	20
Cleaning	21
Cleaning/replacing the filling valve	24
Cleaning/replacing the discharge valve	25
Autoclaving	26
Preparation prior to the first autoclaving	26

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Read carefully before use!

Safety Instructions

This instrument may be involved with hazardous materials, operations, and equipment. This manual does not purport to address all of the safety problems associated with its use. It is the responsibility of whomever uses this instrument to consult and establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.

1. **Follow the general instructions for hazard prevention and safety regulations, e.g., wear protective clothing, eye protection and gloves.**
2. Observe the reagent manufacturer's information.
3. Every user must be acquainted with this Operating Manual before operation.
4. Use the instrument only for dispensing liquids, with strict regard to the defined Operating Exclusions and Limitations (► pages 8-9).
If in doubt, contact the manufacturer or supplier.
5. When dispensing flammable media, make provisions to avoid electrostatic charging, e.g., do not dispense into plastic vessels; do not wipe instruments with a dry cloth.

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
Components

- A piston seat
 - B cylinder sleeve
 - C pointer
 - D volume adjustment knob
 - E valve locking ring
 - F* discharge valve
 - G* discharge tube
 - H* closure cap
 - I locking nut
 - J* filling tube
 - K* filling valve
 - L air vent cap
 - M valve block
 - N* glass cylinder
 - O* piston
 - P mounting tool
 - Q drying tube (optional)
 - R seal
- * Parts in contact with the media

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Application and Design

This instrument is designed for dispensing liquids, observing the following limits:

- 15 to 40 °C of instrument and reagent
- vapor pressure up to 500 mbar
- density up to 2.2 g/cm³
- kinematic viscosity up to 500 mm²/s
(dynamic viscosity [mPas] =
kinematic viscosity [mm²/s] x density [g/cm³])

Materials

When the instrument is correctly used, the dispensed liquid comes into contact with only the following chemically resistant materials: Borosilicate glass, Al₂O₃, ETFE, PFA, FEP, PTFE, platinum-iridium; PP (closure cap).

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Operating Exclusions

Never use this instrument for

- liquids attacking ETFE, FEP, PFA, PTFE or Al_2O_3 (e.g., dissolved sodium azide)
- liquids attacking borosilicate glass (e.g., hydrofluoric acid)
- liquids which are decomposed catalytically by platinum-iridium (e.g., H_2O_2)
- fuming acids
- trifluoroacetic acid
- tetrahydrofuran
- explosive liquids (e.g., carbon disulphide)
- suspensions (e.g., of charcoal) as solid particles may clog or damage the instrument.

Operating Limitations

Concentrated hydrochloric acid and concentrated nitric acid, chlorinated and fluorinated hydrocarbons, or chemical combinations which form deposits may make the piston difficult to move or may cause jamming (e.g., concentrated saline solutions, highly crystallizing solutions or highly concentrated alkaline solutions).

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Warning!

If there is a sign of a potential malfunction (e.g., piston difficult to move) never use force. Immediately stop dispensing and follow cleaning instructions (see p. 21) or contact the manufacturer.

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Note:

Compatibility of the instrument for this special application (e. g., trace material analysis) must be checked by the user or contact the manufacturer.

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Preparing instrument for use

Attention!


The adapters supplied with the instrument are made of polypropylene (PP), and can only be applied for media which do not attack PP (► page 35, ETFE-adapters).

Warning!

Wear protective clothing, eye protection and gloves. Follow all Safety Instructions and observe the Operating Exclusions and Limitations (see pages 8, 9).

Warning!

Always wear protective gloves when touching the instrument and the bottle, especially when using dangerous liquids. Carry the mounted instrument only as shown in figure.

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
Connecting the discharge tube

Warning!

Only the slotted-type discharge tube supplied with the instrument allows correct operation. Do not use any other type of discharge tube. Never use damaged or deformed discharge tubes.

1. Push the discharge tube as far as possible into the discharge valve.
2. Firmly finger-tighten the locking nut. Check the discharge tube for a tight fit.

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


Connecting the filling tube

1. For instruments 0-25 ml and 0-50 ml, filling tube can be connected.
2. Cut off filling tube diagonal.
3. Push the filling tube into the valve block as far as possible – with the diagonal end down.

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Mounting instrument on bottle

The instrument can be screwed directly onto a GL 45 screwthread. For other bottles, use the accompanying adapters. Adapters are also available as accessories (see page 35).

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Attention!

The adapters supplied with the instrument are made of polypropylene (PP), and can only be applied for media which do not attack PP (see page 35, ETFE-adapters).

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Priming the instrument

Attention!

For small bottles use a bottle stand to prevent tipping over.

Warning!

Reagent may drip off the discharge tube and closure cap.

- Hold the discharge tube and carefully remove the closure cap.

- Slide the closure cap backwards on the tube support, away from the discharge tube opening.

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


1. To avoid splashes when priming hold discharge tube on inner wall of a suitable recipient.
2. Gently pull the piston up (1 to 2 cm), then press it down rapidly.
3. Repeat this procedure until most of the air bubbles in the glass cylinder and discharge tube have been removed. Pay particular attention to critical areas where bubbles may lodge (see black arrows).

Note:

Before using the instrument for the first time, ensure it is rinsed carefully or discard first few samples dispensed.

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
Dispensing

Warning!

Follow all Safety Instructions (see page 4, 5) and observe the Operating Exclusions and Limitations (see pages 8, 9). The discharge tube must always point away from the user or any other person. Never press down the piston when the closure cap is pushed on. Reagent may drip off the discharge tube and closure cap.

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

Setting the volume

Volumen einstellen

- Hold an appropriate collecting vessel under the discharge tube opening.
- Pull the piston up slowly and evenly to the upper stop to draw reagent into the glass cylinder.
- Push the piston down gently and evenly to dispense the reagent into the vessel.

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Attention!

After use, always leave the piston in the down position.

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Serial dispensing with flexible discharge tube (optional)

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Attention!

Each time you are going to use the tubing, examine it carefully for any damages. To dispense aggressive liquids, e.g. concentrated acids and lyes, you should take safety measures in addition to the normal precautions of wearing protective clothing, eye protection and protective gloves. We recommend that you use a protective screen.

When using the flexible discharge tube, the bottle must be supported using a bottle stand to help prevent the bottle from accidentally tipping over. To help avoid reagent splashing from the tube, always grip the tube firmly by the handle and replace in the holder after use.

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

Operating Exclusions:

The flexible discharge tube is compatible with all liquids suitable for the dispenser.

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Note:

The coil of the tubing can be stretched to a maximum length of 800 mm. The entire coil must lie in regular loops and must not be twisted. When working with the flexible discharge tube, the accuracy and coefficient of variation indicated in the Operating Manual are only obtained for volumes ≥ 2 ml and by smooth and steady operation, gently approaching the upper and lower stops. Rinse the tube carefully! Do not dismantle!

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
Connecting drying tube (optional)

Note:


The screwthreads of the bottle and/or the threaded adapter should be sealed with PTFE tape.

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Cleaning

The instrument must be cleaned in the following situations to assure correct operation:

- if the piston becomes difficult to move
- after dispensing liquids which form deposits
- before changing the reagent
- prior to long term storage
- should liquid have accumulated in the closure cap
- prior to sterilization
- before replacing the valves

Warning!

Follow all Safety Instructions (see pages 4, 5). The glass cylinder, valves, filling and discharge tubes contain reagent. Point the valves and tube openings away from your body. Wear protective clothing, eye protection and appropriate hand protection.

Note:

Never change the pistons of the instruments!

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

1. Push on the closure cap.
2. Empty and unscrew the instrument from the bottle and lift the instrument so that the filling tube is out of the liquid but still in the bottle.
3. Carefully tap the filling tube against the bottle from the inside so that the reagent runs back into the bottle.

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4. Remove the closure cap, and dispense the remaining reagent into the bottle by repeated piston operation.
5. Mount the instrument onto another bottle filled with an appropriate cleaning solution.
6. Flush the instrument.
7. Empty the bottle, flush it and fill it with water. Repeat step 6 and 7.
8. Empty the instrument.

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
9. Unscrew the piston seat (A).
10. Carefully pull the piston (O) out of the glass cylinder (N) by pulling on the piston seat only.
11. Carefully remove deposits at the edge of the glass cylinder with a screw driver. Clean piston and cylinder with a bottle-brush.
12. Carefully reassemble the instrument in the reverse order. Flush instrument with distilled water.

Note:

Insert the piston straight into the cylinder, not at an angle.

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

Cleaning/replacing the filling valve

1. Follow all cleaning instructions (see pages 21-23).
2. Remove the filling tube.
3. Unscrew the filling valve using the mounting tool. Remove the valve with its sealing washer. Make sure the sealing washer does not remain in the valve block.
4. Rinse the filling valve in the cleaning solution, and clean with soft brush. Should valve ball stick, push into filling opening of valve with pointed instrument (e.g., disposable pipette tip), and release valve ball.
5. Screw the cleaned/new filling valve together with its sealing washer into the valve block and tighten it firmly using the mounting tool.

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

Cleaning/replacing the discharge valve

1. Follow all cleaning instructions (see pages 17-19).
2. Unscrew the discharge tube.
3. Remove valve locking ring.
4. Unscrew the discharge valve using the mounting tool. Remove the valve with its sealing washer. Make sure the sealing washer does not remain in the valve block.
5. Rinse the discharge valve in the cleaning solution, and clean with soft brush. Should valve ball stick, push into filling opening of valve with pointed instrument (e.g., disposable pipette tip), and release valve ball.
6. Screw the cleaned/new discharge valve together with its sealing washer into the valve block and tighten it firmly using the mounting tool.
7. Slide on the valve locking ring making sure the two prongs properly align with the notches.

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Autoclaving

This instrument withstands steam sterilization at 121 °C, 2 bar absolute (15 psi) for 20 minutes in accordance with widely used guidelines. It is the user's responsibility to ensure effective autoclaving.

Preparation prior to the first autoclaving

- Unscrew filling valve (K) and discharge valve (F).
- Replace built-in PTFE-rings (R) with the supplied elastic O-rings made of FKM.
- Screw in the filling valve (K) and discharge valve (F) and tighten them firmly using the mounting tool.

Note:

Observe the reduced chemical resistance of FKM!

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

1. Clean the instrument before sterilization (see pages 21-23).
2. Remove the discharge tube (G) with removed closure cap (H), filling tube (J), and drying tube (Q), and put all parts on a towel. **Avoid contact with hot metal surfaces.**
3. Lay the instrument **with piston in down position** on the towel.
4. Autoclave all parts.

Note:

Allow to cool down in the closed autoclave to avoid deformation through too rapid cooling. After every autoclaving, inspect all parts for deformities or damage. If necessary, replace them. Do not re-assemble the instrument until it has cooled down to room temperature (cooling time approx. 2 hours).

Note:

Mounting tool can be steam sterilized at 121 °C.

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Checking the Volume

The accuracy and coefficient of variation of the instrument are determined gravimetrically as follows:

- Set to the nominal volume.
- Dispense distilled H₂O.
- Weigh the dispensed quantity on an analytical balance.
- Calculate the dispensed volume taking the temperature into account.
- Perform at least 10 dispensing and weighing operations.
- Calculate the accuracy (A%) and the coefficient of variation (CV%) by means of the formula used in statistical quality control. The proceeding is described e.g., in DIN EN ISO 8655-6. Observe the operating manual of the balance manufacturer and the corresponding standards.

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Calculations (for the nominal volume)

$$\text{Mean value } \bar{x} = \frac{\sum x_i}{n} \quad \begin{array}{l} x_i = \text{results of weighings} \\ n = \text{number of weighings} \end{array}$$

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$$\mathbf{VK\%} = \frac{100 s}{\bar{V}}$$

Trouble Shooting

Problem	Possible cause	Corrective action
Piston difficult to move	<ul style="list-style-type: none"> - Formation of crystals, dirty 	<ul style="list-style-type: none"> - Stop dispensing immediately. Loosen piston by circular motion, but do not disassemble. Follow all cleaning instructions (► p. 21-23).
Filling not possible	<ul style="list-style-type: none"> - Volume adjusted to minimum setting - Sticking filling valve 	<ul style="list-style-type: none"> - Set to required volume (► page 16). - Clean the filling valve. If necessary replace the valve with sealing washer (► page 24).
Instrument does not dispense liquid	<ul style="list-style-type: none"> - Discharge tube closure cap pushed on - Discharge tube not firmly connected or damaged 	<ul style="list-style-type: none"> - Remove closure cap. - Push on the closure cap. Hold the piston in place, and push the original manufacturer's discharge tube on firmly (► page 11). Replace a deformed or damaged discharge tube (► Accessories and Spare Parts).
Air bubbles in the instrument	<ul style="list-style-type: none"> - Reagent with high vapor pressure has been drawn in too quickly - Seal not inserted, knurled locking nuts not firmly connected 	<ul style="list-style-type: none"> - Slowly draw in reagent. - Double check that seal (R) is in place, and that the knurled locking nuts at the discharge tube are firmly seated and finger tight.

	<ul style="list-style-type: none"> - Air not expelled from the instrument - Filling tube not firmly connected or damaged - Valves dirty, not firmly installed or damaged 	<ul style="list-style-type: none"> - Priming the instrument (► page 14) - Push the filling tube on firmly. If necessary cut off approx. 1 cm of tube at the upper end and re-connect it or replace filling tube. - Cleaning procedure (► page 21-23). Tighten the valves using the mounting tool. If necessary, replace the valves and sealing washers.
<p>Dispensed volume is too low</p>	<ul style="list-style-type: none"> - Discharge tube not firmly connected or damaged - Filling tube not firmly connected or damaged - Valves not firmly installed or damaged - Sticking filling valve 	<ul style="list-style-type: none"> - Push the original manufacturer's tube on firmly (► page 11). - Cleaning procedure (► page 21-23). Push the filling tube on firmly. If this does not suffice, cut off approx. 1 cm of tube at the upper end and re-connect it or replace filling tube (► page 12). - Clean the instrument (► page 21-23). Tighten the valves using the mounting tool, if necessary replace the valves and sealing washers (► page 24) - Clean or replace discharge tube, if necessary.

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* Error limits (A = Accuracy, CV = Coefficient of variation) relative to the nominal capacity (maximum volume) indicated on the instrument, obtained at equal temperature (20 °C) of instrument, ambience and dist. water, and with smooth operation.

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Items supplied:

operating manual, 1 discharge tube, 1 mounting tool,
1 filling tube and adapters of PP of the following sizes:

Nominal vol.	Adapter, PP	Filling tube
0.5 - 5	GL 25, GL 28, GL 32,	250 mm
1 - 10	GL 38, S40	
2.5 - 25	GL 32, GL 38, S 40	500** mm
5 - 50		

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Discharge tube

FEP, complete with tube support, locking nut, and closure cap. Pack of 2.

for nominal volume/size/length		Cat. No.
5, 10 ml	bent, 80 mm	7215 53
25, 50 ml	bent, 120 mm	7215 54
25, 50 ml	bent, 120 mm drawn to fine tip	7215 55
25, 50 ml	bent, 200 mm	7215 56

Filling tubes

with adapter to connect 2 filling tubes, FEP. Pack of 2.

Length x outer diam., mm	Cat. No.
250 x 6.9	7215 59
335 x 6.9	7215 60
390 x 6.9	7215 61

Flexible discharge tube

PTFE, coiled, length 800 mm, with safety handel. Pack of 1.

Nominal volume, ml	Cat. No.
5, 10	7216 33
25, 50	7216 34

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Drying tube

Cat. No. 7215 67

**Valve locking ring
for discharge tube**

Nominal volume 5-10 ml
Nominal volume 25-50 ml

Cat. No. 7215 70

Cat. No. 7215 71

**Filling valve with
sealing washer**

Nominal volume 5-10 ml
Nominal volume 25-50 ml

Cat. No. 7215 72

Cat. No. 7215 73

**Discharge valve with
sealing washer**

Nominal volume 5-10 ml
Nominal volume 25-50 ml

Cat. No. 7215 74

Cat. No. 7215 75

**Seals for filling valve and
discharge valve, pack of 6.**

PTFE
FKM (Fluoro elastomer),
only for sterile applications

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Return for Repair

Attention:

For safety reasons only clean/decontaminated instruments can be checked/repared.

- Therefore: Clean and decontaminate the instrument carefully.
- Complete the “Declaration on the Absence of Health Hazards” and send the instrument to your local Biohit representative and describe
 - the nature of the problem
 - the liquids dispensed.

Return shipment has to be effected at the risk and cost of the sender.

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Warranty

As provided by law, any and all warranties are null and void if the product has been misused, modified or repaired by unauthorized personnel, if the defects are caused by negligence (instruction manual, maintenance) or by normal wear and tear. Use only original manufacturer's accessory/spare parts. For conditions and extent of warranty refer to our General Conditions of Sale.

Subject to technical modification without notice. Errors excepted.

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Dispenser selection chart

Reagent	BIOHIT PROLINE® Prospenser
Acetaldehyde	+
Acetic acid (glacial), 100%	+
Acetic acid, 96%	+
Acetone	+
Acetonitrile	+
Acetylacetone	+
Acrylic acid	+
Acrylonitrile	+
Adipic acid	+
Allyl alcohol	+
Aluminium chloride	+
Amino acids	+
Ammonium chloride	+
Ammonium fluoride	+
Ammonium hydroxide, 30% (Ammonia)	+
Ammonium sulfate	+
n-Amyl acetate	+
Amyl alcohol (Pentanol)	+
Aniline	+
Barium chloride	+
Benzaldehyde	+

Reagent	BIOHIT PROLINE® Prospenser
Chlorobenzene	+
Chlorobutane	+
Chromic acid, 50%	+
Chromosulfuric acid	+
Copper sulfate	+
Cumene (Isopropyl benzene)	+
Cyclohexanone	+
Decane	+
1-Decanol	+
Dibenzylether	+
Dichlorobenzene	+
Diethanolamine	+
Diethylamine	+
1,2 Diethylbenzene	+
Diethylene glycol	+
Dimethyl sulfoxide (DMSO)	+
Dimethylaniline	+
Dimethylformamide (DMF)	+
Diphenyl ether	+
Ethanol	+
Ethanolamine	+
Ethyl acetate	+

Reagent	BIOHIT PROLINE® Prospenser
Monochloroacetic acid	+
Nitric acid, 30%	+
Nitrobenzene	+
Oleic acid	+
Oxalic acid	+
Perchloric acid	+
Petroleum	+
Phenol	+
Phenylethanol	+
Phenylhydrazine	+
Phosphoric acid, 85%	+
Phosphoric acid, 85% + Sulphuric acid, 98%, 1:1	+
Piperidine	+
Potassium chloride	+
Potassium dichromate	+
Potassium hydroxide	+
Potassium permanganate	+
Propionic acid	+
Propylene glycol (Propanediol)	+
Pyridine	+

Benzene (Benzol)	+
Benzoyl chloride	+
Benzyl alcohol	+
Benzylamine	+
Benzylchloride	+
Boric acid, 10%	+
Bromobenzene	+
Bromonaphthalene	+
Butanediol	+
1-Butanol	+
n-Butyl acetate	+
Butyl methyl ether	+
Butylamine	+
Butyric acid	+
Calcium carbonate	+
Calcium chloride	+
Calcium hydroxide	+
Calcium hypochlorite	+
Chloro naphthalene	+
Chloroacetaldehyde, 45%	+
Chloroacetic acid	+
Chloroacetone	+

Ethyl methyl ketone	+
Formaldehyde, 40%	+
Formamide	+
Glycerol	+
Glycol (Ethylene glycol)	+
Glycolic acid, 50%	+
Hexanoic acid	+
Hexanol	+
Hydriodic acid	+
Hydrochloric acid, 20%	+
Isoamyl alcohol	+
Isobutanol (Isobutyl alcohol)	+
Isopropanol (2-Propanol)	+
Isopropyl ether	+
Lactic acid	+
Methanol	+
Methoxybenzene	+
Methyl benzoate	+
Methyl butyl ether	+
Methyl formate	+
Methyl propyl ketone	+
Mineral oil (Engine oil)	+

Pyruvic acid	+
Salicylic acid	+
Salicylaldehyde	+
Scintillation fluid	+
Silver acetate	+
Silver nitrate	+
Sodium acetate	+
Sodium chloride	+
Sodium dichromate	+
Sodium fluoride	+
Sodium hydroxide, 30%	+
Sodium hypochlorite	+
Sulfuric acid, 98%	+
Tartaric acid	+
Tetramethylammonium hydroxide	+
Triethanolamine	+
Triethylene glycol	+
Urea	+
Zinc chloride, 10%	+
Zinc sulfate, 10%	+

The above recommendations reflect testing completed prior to publication. Always follow instructions in the operating manual of the instrument as well as the reagent manufacturer's specifications. In addition to these chemicals, a variety of organic and inorganic saline solutions (e.g., biological buffers), biological detergents and media for cell culture can be dispensed. Should you require information on chemicals not listed, please feel free to contact BIOHIT.

Status as of: 0906/1

The BIOHIT PROLINE® Prospenser is not designed for use with Hydrofluoric acid (HF). Consults the "Operating Exclusions" and "Operating Limitations" sections of the Operating Manual (see pages 9, 10).

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