



Artisan Technology Group is your source for quality new and certified-used/pre-owned equipment

- FAST SHIPPING AND DELIVERY
- TENS OF THOUSANDS OF IN-STOCK ITEMS
- EQUIPMENT DEMOS
- HUNDREDS OF MANUFACTURERS SUPPORTED
- LEASING/MONTHLY RENTALS
- ITAR CERTIFIED SECURE ASSET SOLUTIONS

SERVICE CENTER REPAIRS

Experienced engineers and technicians on staff at our full-service, in-house repair center

*InstraView*SM REMOTE INSPECTION

Remotely inspect equipment before purchasing with our interactive website at www.instraview.com ↗

WE BUY USED EQUIPMENT

Sell your excess, underutilized, and idle used equipment. We also offer credit for buy-backs and trade-ins. www.artisanng.com/WeBuyEquipment ↗

LOOKING FOR MORE INFORMATION?

Visit us on the web at www.artisanng.com ↗ for more information on price quotations, drivers, technical specifications, manuals, and documentation

Contact us: (888) 88-SOURCE | sales@artisanng.com | www.artisanng.com

HiGro®

A high-capacity microwell plate growth system.



The HiGro is a high-capacity, incubating shaker designed to accelerate sample growth in microwell plates. This uniquely engineered instrument combines a small shaking orbital, gas flow and temperature control system in a **user-friendly, compact layout**. Plate-holding cassettes allow for **easy loading** and unloading.

The convenience and control of the HiGro system make it the ideal solution for facilities with high-throughput growth needs as well as for laboratories using applications that are best optimized in a microwell format.

Key Benefits

- Space-saving design
- Optimizes growth in 96- and 384-microwell plate formats
- Easy to load and unload
- Superior control over growth environment

Design Specifications


Dimensions	26" W x 21" D x 22" H (0.66m x 0.53m x .56m) (small footprint allows placement on lab bench)
Capacity	4 growth chambers, each holding 1 cassette with 12 shallow-well or 6 deep-well plates. Total system capacity is 48 shallow-well plates or 24 deep-well plates.
Temperature	Independent control for each chamber. Range is ambient to 55°C. Resolution is 1°C.
Shaking	Rotational orbital is 8.0 mm (ideally suited for 96-well growth) or 4.0 mm (ideally suited for 384-well growth)*. Maximum speed is 650 rpm. Resolution is 1 rpm.
Gas Flow	Oxygen- or air-ready*. Gas flow rate and timing controls located on instrument front panel.

*Depending on Model

Data

HiGro & Floor Shaker Comparison for Plasmid Growth

Figure 1



Figures 1 & 2

Conditions: E.coli strain is DH10B carrying pUC 19; 1x TB media; inoculum: 5 µl glycerol stock culture; selection: Ampicillin; 37°C.


HiGro: 300 µl/well in 96-shallow-well, flat-bottom plate; 550 rpm; Air flow: no delay, 2 sec on, 30 sec off, 6 standard liters per minute; no plate lids.

Figure 1:

Floor Shaker: 1000 µl/well in 96-well Beckman deep-well block, 250 rpm in New Brunswick floor shaker (plates were taped down); aluminum foil placed loosely over each plate; wetted paper towel placed in shaker to maintain ambient humidity.

Cross-Contamination Growth Test with No Plate Lids

Figure 2



Plasmid Yield Time Course from 384-Well Growth*

Figure 3



>1µg plasmid after only 8 hours! (total yield from 90 µl growth volume)

DNA eluted into 50 µl EB

Each band in the ladder represents 25 ng DNA

Figure 3:

Growth Conditions: 90 µl/well Magnificent Media, 384-shallow-well, flat-bottom plate, rounded-square wells; inoculation (1 µl) from O/N culture; high copy plasmid, 3 kb total size; selection: Carbapenem.

HiGro: 37°C; 500 rpm; Air flow: no delay, 0.5 sec on, 0.5 min off, 1.5 standard liters per minute; no plate lids.

Plasmid Prep Conditions: 90 µl (total culture) pipetted into alkaline lysis reagents (100 µl). Neutralization reagents added (100 µl). Pelleted sample. Isolated DNA on commercially available column. Eluted into 50 µl EB.

Gel Conditions: 2 µl & 4 µl of each (50 µl) prep were run on 1% agarose gel. Ladder: two 10-band mass ladders combined (250 ng of each ladder).

Request the latest HiGro Application & Technical Notes.

*Data obtained through customer collaboration

Worldwide Headquarters

Digilab, Inc.
84 October Hill Road
Holliston, MA 01746
USA

Phone: (508) 893-3130
Toll Free: (800) 935-8007
Fax: (508) 893-8011

E-Mail: info@digilabglobal.com

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisanng.com

DIGILAB®



Artisan Technology Group is your source for quality new and certified-used/pre-owned equipment

- FAST SHIPPING AND DELIVERY
- TENS OF THOUSANDS OF IN-STOCK ITEMS
- EQUIPMENT DEMOS
- HUNDREDS OF MANUFACTURERS SUPPORTED
- LEASING/MONTHLY RENTALS
- ITAR CERTIFIED SECURE ASSET SOLUTIONS

SERVICE CENTER REPAIRS

Experienced engineers and technicians on staff at our full-service, in-house repair center

*InstraView*SM REMOTE INSPECTION

Remotely inspect equipment before purchasing with our interactive website at www.instraview.com ↗

WE BUY USED EQUIPMENT

Sell your excess, underutilized, and idle used equipment. We also offer credit for buy-backs and trade-ins. www.artisanng.com/WeBuyEquipment ↗

LOOKING FOR MORE INFORMATION?

Visit us on the web at www.artisanng.com ↗ for more information on price quotations, drivers, technical specifications, manuals, and documentation

Contact us: (888) 88-SOURCE | sales@artisanng.com | www.artisanng.com