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**IKA HBR 4 digital**  
**IKA HB 4 digital**  
**IKA HB 4 basic**

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**Garantie**

You have purchased an original IKA laboratory machine which meets the highest engineering and quality standards.

In accordance with IKA guarantee conditions, the guarantee period is 12 months. For claims under the guarantee please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

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HB4 0596

## Commissioning

Check whether the voltage specified on the type plate matches the mains voltage available. The power socket used must be earthed (protective earth conductor contact). If these conditions are met, the device is ready to operate after plugging in the mains plug. If these procedures are not followed, safe operation cannot be guaranteed and/or the equipment may be damaged.

## Unpacking

Please unpack the equipment carefully and check for any damages. It is important that any damages which may have arisen during transport are ascertained when unpacking. If applicable a fact report must be set immediately (post, rail or forwarder).  
The delivery scope covers: an IKA tempering bath and operating instructions.

## Intended use

Upon continuous operation with high speed and a room temperature of 20°C a warming up of the housing up to 40°C is possible.  
When using an oil bath, only low-viscosity oils must be used, since using high-viscosity oils or solid grease can lead to overheating in localized areas, causing excessive pressure to build up in the bottom of the vessel.  
When using an water bath, we recommend you use water with no mineral content.

## Safety instructions

The safety temperature limit must be set at least 10K under the flashpoint of the tempering fluid used.  
Never operate the heating baths without tempering fluid!  
Filling and emptying when hot can lead to accidents, especially scalds. Please observe the accident prevention regulations applicable for the job and place of work, wear protective goggles and gloves.  
Take care when touching the external jacket, as it can heat up to over 80°C.  
When working with temperature sensitive mediums please consider:

## Contents

| Page | Guarantee  |
|------|--|
| 2    | Safety instructions                              |
| 10   | Intended use                                     |
| 10   | Unpacking  |
| 10   | Commissioning                                    |
| 11   | Starting   |
| 11   | Setting the safety circuit                       |
| 11   | Operating modes                                  |
| 12   | Heating  |
| 12   | Controlling the temperature of the medium        |
| 12   | Stirring   |
| 13   | Error messages                                   |
| 13   | Maintenance and cleaning                         |
| 14   | Technical data                                   |
| 20   | Associated standards and regulations             |
| 20   | Accessories                                      |
| 22   | List of spare parts HB 4 basic                   |
| 23   | Spare parts diagram HB 4 basic                   |
| 24   | List of spare parts HBR 4 digital / HB 4 digital |
| 25   | Spare parts diagram HBR 4 digital / HB 4 digital |

## Starting

The device is switched on with the 2-pole ON/OFF mains switch (A). When switched on the lamp in the rocker switch lights up green.

## Setting the safety circuit

### **IKA HBR4 digital and HB4 digital**

When activating the HBR4 digital and HB4 digital devices an automatic self-test is carried out, all LEDs flash and all segments in **figure 1** are visible on the LCD display. After this, the LCD displays in **figure 2 to figure 5** are visible (standard factory setting).

When setting the safety circuit an upper temperature limit between 50°C and 210°C is set. The safety circuit can only be set in operating mode "A".

◆ Switch on the device.

◆ Within 5 seconds of the display **fig. 6** appearing, keep the mode button (E) pressed down and set the safety temperature desired by simultaneously turning the control knob (B).

◆ Release the mode button (E). For two seconds the safety temperature set and safe OK **fig. 7** is displayed; the safety temperature set is saved.

The safety circuit is automatically checked each time the device is restarted.

### **IKA HB4 basic**

A safety temperature limiter takes precedence over the main controller which works independently from the closed loop control circuit. The safety temperature is set with a small screw driver on the setting axle (K). The end of the slot opposite the flattening is used as a setting mark.

## Operating modes

### **IKA HBR4 digital and HB4 digital**

The device can be operated in three different operating modes (display see **fig. 8 to fig. 10**). To toggle from one operating mode to the other keep the mode button (E) depressed when switching on the device (min. 5 secs). Sequence A-B-C-A-B-C-A- etc.

**Operating mode "A"** This operating mode is the factory setting of the device. When switching on the device the heating and stirring functions are deactivated. The set values last set are saved and are used when the heating or stirring functions are switched on. The set values can be changed.

**Operating mode "B"** When switching on the device the state of the heating and stirring functions are taken over before the device was switched off last as well as the last values set. The set values can be changed.

**Operating mode "C"** When switching on the device the state of the heating and stirring functions are taken over before the device was switched off last as well as the last values set. The set values cannot be changed.

If the heating bath is used as a water bath, the safety temperature limiter is to be set to the symbol **H<sub>2</sub>Omax**.

When the safety temperature is reached the device switches off permanently. Only once the heating bath has cooled down somewhat can it be restarted by pressing in the red safety button (L) for example using a small screwdriver.

## Heating

### IKA HBR4 digital and HB4 digital

With the mode button **(E)** the correcting variable desired (set temperature) can be preselected. The heating bath temperature of the device is kept constant by the closed loop control circuit and is additionally monitored by the safety circuit. In the event of a disturbance in the closed loop control circuit, the heating bath is switched off permanently by the safety circuit. In the event of a fault in the closed loop control circuit or in the safety circuit, the LED "Temp," flashes green and yellow alternately. The LED next to the reset button lights up red. The error is also displayed on the LCD (see chapter Error messages). The heating function cannot be started.

## Controlling the temperature of the medium

### IKA HB4 basic

The desired medium temperature can be set with the control knob **(B)** using the graduation on the front film. The yellow low voltage neon glow lamp **(J)** next to the control knob indicates the energy supply to the heating system of the device.

### IKA HBR4 digital and HB4 digital

With the control knob **(B)** the medium temperature desired can be set between 0 and 200°C, not however higher than the safety temperature set. The value set can be read off on the digital display **(C)** (fig. 11). The heating function is switched on and off by pressing the sensor key "Temp./on/off" **(G)**. If this function is activated the green LED next to the lettering "Temp" lights up. The heating bath is heated up to the temperature set. On the LCD display the set temperature, actual temperature and the PV symbol are displayed. The set temperature and actual temperature relate to the medium.

## Stirring

### IKA HBR4 digital

With the mode button **(E)** the correcting variable desired (set speed) can be preselected. With the control knob **(B)** the desired speed of the stirrer drive can be set between 0 and 800 rpm in steps of 50 rpm. The function is switched on and off by pressing the sensor key "Motor/on/off" **(F)**. If this function is activated the green LED lights up next to the lettering "Mot." The motor runs gently up to the speed set. On the LCD display the set speed, actual speed and the PV symbol are displayed **(fig. 15)**. If both functions (heating and stirring) are switched off, the LCD display always indicates the set temperature set. When the heating function is switched on it is given priority on the LCD display. If the stirring function is started with the sensor key or called up using the mode key **(E)**, the system switches over to speed display for five seconds. With the speed button **(E)** the LCD display can be changed from temperature to speed and vice versa at any time. After five seconds the system switches back to the predominant mode.

With the mode button **(E)** the correcting variable desired (set speed) can be preselected. With the control knob **(B)** the desired speed of the stirrer drive can be set between 0 and 800 rpm in steps of 50 rpm. The function is switched on and off by pressing the sensor key "Motor/on/off" **(F)**. If this function is activated the green LED lights up next to the lettering "Mot." The motor runs gently up to the speed set. On the LCD display the set speed, actual speed and the PV symbol are displayed **(fig. 15)**. If both functions (heating and stirring) are switched off, the LCD display always indicates the set temperature set. When the heating function is switched on it is given priority on the LCD display. If the stirring function is started with the sensor key or called up using the mode key **(E)**, the system switches over to speed display for five seconds. With the speed button **(E)** the LCD display can be changed from temperature to speed and vice versa at any time. After five seconds the system switches back to the predominant mode.

|       |  |   |
|-------|--|---|
| No.   | Error  | Remedy  |
| Er 3  | Internal device temperature above 76°C. Allow the device to cool down.   | Mains switch OFF/ON   |
| Er 4  | Motor speed deviates $\pm 300 \text{ rpm}$ . Reduce load.  | Mains switch OFF/ON   |
| Er 4  | Stirring function turned on from the speed. Set speed greater than $350 \text{ rpm}$ .   | Mains switch OFF/ON   |
| Er 4  | Stirring function turned on. Actual speed = $0 \text{ rpm}$ . Increase the speed or turn the stirring function.                                | Mains switch OFF/ON   |
| Er 20 | Safety relay does not close. Safety circuit is set too low. check mains voltage (mains voltage must be larger than the nominal voltage - 10%). | Mains switch OFF/ON   |
| Er 21 | Safety relay does not open.  | RESET key/service   |
| Er 22 | Disturbance before the last time the device was switched Off or safety temperature is lower than the temperature of the medium.                | RESET button  |
| Er 23 | Adjustable safety circuit faulty.  | Mains switch OFF/ON / service   |
| Er 24 | Safety temperature set exceeded.   | Set safety temperature higher. Mains switch OFF/ON                                  |
| Er 26 | System dry, no medium in the bath container.   | Allow the device to cool down, then refill with medium. Mains switch OFF/ON         |
| Er 27 | Calibration error.   | Mains switch OFF/ON. When switching on, a calibration is automatically carried out. |
| Er 28 | Sensor rupture.  | Mains switch OFF/ON   |
| Er 29 | Short circuit in the safety sensor or controller and safety sensor mixed up.   | Service   |

The following error messages can be displayed (fig. 17) and have the following causes described:

### Error messages

|       |  |         |
|-------|--|---------|
| Er 30 | Short circuit in the controller sensor | Service |
|-------|--|---------|

### Maintenance and cleaning

The safety circuit must be checked over at least once a year by the user. To do this, fill the temperature control bath with 1 l of water and in operating mode „A“ set the safety temperature to 100°C. Then set the target temperature to 80°C and start the heating function by pressing the Temp ON/OFF key. After the target temperature has been reached, switch the moderating bath off and back on again with the power switch and set the safety temperature to 70°C. The temperature of the medium will then be 10°C above the safety temperature, and the safety circuit will be activated. Error message Er 22 or Er 24 must be displayed.

The aparate is maintenance-free. It is subject only to the natural wear and tear of components and their statistical failure rate.

**When ordering spare parts, please give the manufacturing number shown on the type plate, the machine type and the name of the spare part.**

Please send in equipment for repair only after it has been cleaned and is free from any materials which may constitute a health hazard. Use only cleansing agents which have been approved by IKA to clean IKA devices. To remove use:

- Dyes
- isopropyl alcohol
- Construction materials
- water containing tenside / isopropyl alcohol
- Cosmetics
- water containing tenside / isopropyl alcohol
- Foodstuffs
- water containing tenside
- Fuels
- water containing tenside

For materials which are not listed, please request information from IKA. Wear the proper protective gloves during cleaning of the devices. Electrical devices may not be placed in the cleansing agent for the purpose of cleaning. Before using another than the recommended method for cleaning or decontamination, the user must ascertain with the manufacturer that this method does not destroy the instrument.

**Technical data**

|                                | HB4 digital                                  | HB4 digital                                  | HB4 basic                                    |
|--------------------------------|--|--|--|
| Nominal voltage                | VAC 230 ±10%                                 | VAC 230 ±10%                                 | VAC 230 ±10%                                 |
| or                             | VAC 115 ±10%                                 | VAC 115 ±10%                                 | VAC 115 ±10%                                 |
| Frequency                      | Hz 50/60                                     | Hz 50/60                                     | Hz 50/60                                     |
| <b>Heating function</b>        |  |  |  |
| Heating output                 | W 1000                                       | W 1000                                       | W 1000                                       |
| Temperature range (medium)     | °C Room temperature ..200                    | °C Room temperature ..200                    | °C Room temperature ..225                    |
| Temperature range (medium)     | °C Set and actual temperat. LCD              | °C Set and actual temperat. LCD              | °C Set temperat. graduation                  |
| Setting accuracy               | K 1  | K 1  | K ±5   |
| Controller deviation           | K ±0,4                                       | K ±2   | K ±5   |
| (3 l water / 90 °C)            |  |  |  |
| Absolute deviation / mean      | K ±1   | K ±2   |  |
| (3 l water / 90 °C)            |  |  |  |
| Controller deviation           | K ±0,8                                       | K ±1   |  |
| (3 l Marlotherm oil / 150 °C)  |  |  |  |
| Absolute deviation / mean      | K ±2   | K ±10  |  |
| (3 l Marlotherm oil / 150 °C)  |  |  |  |
| <b>Stiring function</b>        |  |  |  |
| Motor                          | ballbearing brush-free                       | ballbearing brush-free                       | EC motor                                     |
| Motor output                   | W 5  | W 5  | W 5  |
| Speed range                    | rpm 100 - 800, adjustable in steps of 50 rpm | rpm 100 - 800, adjustable in steps of 50 rpm | rpm 100 - 800, adjustable in steps of 50 rpm |
| Minimum filling height         | cm 3   | cm 3   | cm 3   |
| Overvoltage class              | % II   | % II   | % II   |
| Duty cycle                     | % 100  | % 100  | % 100  |
| Enclosure rating (DIN 40 050)  | °C IP 21                                     | °C IP 21                                     | °C IP 42                                     |
| Ambient temperature            | % +5 .. 40                                   | % +5 .. 40                                   | % +5 .. 40                                   |
| Ambient humidity (rel.)        | % 80   | % 80   | % 80   |
| Safety class                   | litre DIN 12 877 class 2                     | litre DIN 12 877 class 2                     | litre DIN 12 877 class 2                     |
| Total volume - bath container  | litre 5                                      | litre 5                                      | litre 5                                      |
| Material - bath container      | mm 1.4301                                    | mm 1.4301                                    | mm 1.4301                                    |
| Dimensions (diameter x height) | mm 340 x 250                                 | mm 340 x 250                                 | mm 340 x 250                                 |
| Weight                         | kg 4,4                                       | kg 4,0                                       | kg 3,9                                       |

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### Associated standards and regulations

#### Associated EU guidelines

EMV guidelines: 89/336/EC  
Machine guidelines: 73/023/EWG

#### Construction in accordance with the following safety standards

EN 61 010-1 / VDE 411-1  
EN 50 081  
UL 3101-1  
EN 60 555  
CAN/CSA C22.2 (1010-1)  
EN 55 011

### Accessories

#### H 210 Ring set

**H 159** Intermediate base  
**H 160** Lid

#### H 162 Test tube insert

**IKAFLO** stirrer (dia. 7x60 mm) only for HBR4 digital

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


**HB 4 basic**


| Item | Designation                 |
|------|-----------------------------|
| 9    | Pedestal                    |
| 10   | Braided conductor with fuse |
| 14   | Braided conductor with fuse |
| 18   | Heated pot                  |
| 20   | Rear part                   |
| 21   | Exterior jacket             |
| 22   | Spacer sleeve               |
| 23   | Insulation mat              |
| 24   | Metal plate                 |
| 25   | Device plug                 |
| 26   | Cover                       |
| 27   | Cover                       |
| 29   | Rocker switch               |
| 30   | Base                        |
| 31   | Control knob                |
| 33   | Capillary tube controller   |
| 35   | Front plate                 |
| 37   | Retaining plate             |
| 38   | Upper section               |
| 39   | Lamp                        |
| 2001 | Operating unit compl.       |
| 2002 | Capillary tube limiter      |

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(HBR 4 digital)

**List of spare parts**

| Item | Designation               |
|------|---------------------------|
| 8    | Plastite screw 2.5 x 8    |
| 9    | Flat head screw UNC 3/16" |
| 10   | Securing washer 5.0       |
| 11   | Pedestal                  |
| 12   | Heated pot                |
| 13   | Plug                      |
| 14   | Rocker switch             |
| 16   | EC motor                  |
| 17   | Flange                    |
| 18   | Magnet                    |
| 20   | Spacer sleeve             |
| 21   | Intermediate base         |
| 22   | Control knob              |
| 23   | Rear part                 |
| 24   | Cover                     |
| 27   | Front plate               |
| 28   | BLP operating unit        |
| 29   | BLP control system        |
| 30   | BLP output section        |
| 31   | Exterior jacket           |
| 32   | Base                      |
| 35   | Flat band conductor       |
| 36   | Spacer ring               |
| 47   | Edge protection strip     |
| 2001 | Axial fan compl.          |
| 2002 | Operating unit compl.     |
| 3001 | Fuse 7 A (slow blowing)   |
| 3002 | Fuse                      |

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

**IKÄ**® Works, Inc.

LABORATORY TECHNOLOGY  
ANALYZING TECHNOLOGY  
PROCESSING EQUIPMENT

North America

**IKÄ**® Works, Inc.

2635 North Chase Pkwy SE  
Williamington, NC 28405-7419  
USA  
TEL. +1 800 733-3037  
TEL. +1 910 452-7059  
FAX +1 910 452-7693  
E-mail: usa@ika.net

**IKÄ**® Works, (Asia) Sdn Bhd

LABORATORY TECHNOLOGY  
ANALYZING TECHNOLOGY  
PROCESSING EQUIPMENT

Asia - Australia

**IKÄ**® Works (Asia) Sdn Bhd

Lot 2 Jalan Indah 1/2  
Taman Industri Rawang Indah  
48000 Rawang  
Selangor, Malaysia  
TEL. +60 3 6093-3322  
FAX +60 3 6093-3940  
E-mail: ika@tm.net.my

**IKÄ**® Japan Y.K.

LABORATORY TECHNOLOGY  
ANALYZING TECHNOLOGY  
PROCESSING EQUIPMENT

Japan

**IKÄ**® Japan Y.K.

293-1 Kobayashi-cho  
Yamato Koryama Shi  
639-1026 Japan  
TEL. +81 74358-4611  
FAX +81 74358-4612  
E-mail: japan@ika.de

**IKÄ**® Works do Brasil Ltda.

LABORATORY TECHNOLOGY  
ANALYZING TECHNOLOGY  
PROCESSING EQUIPMENT

South America

**IKÄ**® Works do Brasil Ltda.

Estrada do Guerenquê, 491  
Taquara Jacarepagua, RJ  
Rio de Janeiro  
22713-000 Brasil  
TEL. +55 21 2435-9600  
FAX +55 21 2435-9601  
E-mail: brasil@ika.net

**IKÄ**® Works Guangzhou

LABORATORY TECHNOLOGY  
ANALYZING TECHNOLOGY  
PROCESSING EQUIPMENT

China

**IKÄ**® Works Guangzhou

173-175 Friendship Road  
Guangzhou Economic & Technological  
Development Zone  
Guangzhou 510730, P.R.CHINA  
TEL. +86 20 8222-6772  
FAX +86 20 8222-6776  
E-mail: sales@ikagz.com.cn

## Artisan Technology Group is your source for quality new and certified-used/pre-owned equipment

- FAST SHIPPING AND DELIVERY
- TENS OF THOUSANDS OF IN-STOCK ITEMS
- EQUIPMENT DEMOS
- HUNDREDS OF MANUFACTURERS SUPPORTED
- LEASING/MONTHLY RENTALS
- ITAR CERTIFIED SECURE ASSET SOLUTIONS

### SERVICE CENTER REPAIRS

Experienced engineers and technicians on staff at our full-service, in-house repair center

### *InstraView*<sup>SM</sup> REMOTE INSPECTION

Remotely inspect equipment before purchasing with our interactive website at [www.instraview.com](http://www.instraview.com) ↗

### WE BUY USED EQUIPMENT

Sell your excess, underutilized, and idle used equipment. We also offer credit for buy-backs and trade-ins. [www.artisanng.com/WeBuyEquipment](http://www.artisanng.com/WeBuyEquipment) ↗

### LOOKING FOR MORE INFORMATION?

Visit us on the web at [www.artisanng.com](http://www.artisanng.com) ↗ for more information on price quotations, drivers, technical specifications, manuals, and documentation

**Contact us:** (888) 88-SOURCE | [sales@artisanng.com](mailto:sales@artisanng.com) | [www.artisanng.com](http://www.artisanng.com)