Artisan Scientific

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Contents

	Page
EG- Declaration of conformity	9
Warranty	9
Safety instructions	9
Correct use	11
Unpacking	11
Assembly	12
Commissioning	12
Eliminating malfunctions	13
Maintenance and cleaning	13
Accessories	14
Technical data	14
List of spare parts RV 05 basic	15
Spare parts diagram RV 05 basic	23

EG-Declaration of conformity

We declare under our sole responsibility that this product corresponds to the regulations 89/336/EG and 73/23/EG and conforms with the standards or standardized documents DIN EN IEC 61 010-1 and DIN EN IEC 61 326-1.

Warranty

You have purchased an original **IKA** laboratory machine which meets the highest engineering and quality standards.

In accordance with **IKA** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our works,

enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

If you require servicing, return the appliance in its original packaging. Storage packaging is not sufficient. Please also use suitable transport packaging.

Safety instructions

In order to be able to use the appliance properly and safely, every user must first read the operating instructions and observe the safety instructions contained therein. Take care of these operating instructions and keep them in a place where they can be accessed by everyone.

Only trained staff who are familiar with the appliance and who are authorised to work in this area may operate this appliance, especially under a vacuum, using flammable liquids or explosive or noxious substances.

Please follow the relevant safety instructions and guidelines, and occupational health and safety regulations for use in the laboratory. When working under a vacuum in particular!

When working with the device, the user must select and wear his personal protective equipment according to the mixing hazard category of the medium.

When using media such as toxic liquids, gases, mist, vapours or dusts, or, in some cases, biological or microbiological substances, the user runs the risk of coming into contact with or inhaling the medium.

Please observe the information on the hazardous substance categories for handling solvents and toxic or biological substances.

Always handle distillation fluids in the proper manner.

Adhere to the relevant regulations when working with explosive substances and their compounds.

The appliance must not be operated in explosive atmospheres or under water.

The device must only be used under an all side-closed exhaust, or a comparable protective device.

Ensure that the appliance is on a flat, stable, non-slip surface. Align the feet according to the appliance instructions (see "Assembly" section). Always check before starting up that the cantilever is fitted to the stand in accordance with the assembly instructions for drive **RV 05 basic**. Depending on the structure of the device, there may be a risk to the stationary stability. You should therefore make certain the device is safe and stable each time before placing it in service.

Glass parts should be checked over carefully for damage (for example cracks) each time before placing the unit in service. Damaged glass parts must not be used, as they can break during operation. Glasses must be fitted without tension. When using glass test equipment, it is important to note that bursting can occur because of tension resulting from incorrect assembly, because of external mechanical impact or local temperature peaks. Evacuated glass vessels must not be heated on one side; the evaporating flask must rotate during the heating phase.

Adapt the quantity and the type of distilland to the size of the distillation equipment. The cooler must work properly. Monitor the coolant flow rate at the cooler outlet (see information: "Starting up").

Tasks with the **RV 05 basic** must only be performed when operation is monitored.

Operation with excess pressure is not permitted!

The glass equipment must always be ventilated when working under normal pressure (e.g. open outlet at cooler) in order to prevent a pressure build-up due to extremely volatile gases or an unknown pressure build-up. Please note that dangerous concentrations of gases, vapours or particulate matter can escape through the outlet at the cooler. Take appropriate action to avoid this risk, for example, downstream cold traps, gas wash bottles or an effective extraction system.

The appliance is designed for operation under a vacuum of up to 10 mbar. When carrying out vacuum distillation, uncondensed vapours must be condensed out or safely dissipated. When carrying out vacuum distillation, equipment must be evacuated prior to heating (see "Starting up" section). The equipment must only be aired again after cooling. If there is a risk that the distillation residue could disintegrate in the presence of oxygen, only inert gas must be admitted for stress relief.

Prevent peroxide from forming!

Liquids which are prone to forming organic peroxides must be examined for peroxides prior to distillation and evaporation. Any peroxides detected must be removed. Many organic compounds are prone to the formation of peroxides e.g. dekalin, diethyl ether, dioxane, tetrahydrofuran, as well as unsaturated hydrocarbons, such as tetralin, diene, cumene and aldehydes, ketones and solutions of these substances.

ATTENTION! Organic peroxides can build up in the distillation and waste steam residue of these liquids and solutions and disintegrate explosively! Liquids which are prone to forming organic peroxides must be kept away from light, especially UV radiation.

ATTENTION! Avoid delayed boiling!

Never heat the evaporating flask in the heating bath without switching on the rotary drive! If there are signs that the contents of the flask are beginning to disintegrate in the course of a chemical reaction or distillation (for example, sudden foaming or outgases) switch the heating off immediately and lift the evaporating flask out of the heating bath using the stand. Evacuate the danger zone and warn those in the surrounding area!

ATTENTION! Rotating parts are a source of danger!

Never rotate the evaporating flask on the raised stand. Always lower the stand for rotation. Ensure that parts of the body, hair, items of clothing and jewellery cannot get caught in rotating parts.

There may be dangerous electrostatic activity between the medium and the drive unit.

Check that the filled distillation flask is balanced when setting the speed. Imbalances can cause uncontrolled resonance behaviour in the appliance/equipment which could damage or destroy glass equipment! Set the drive to a speed which ensures that none of the medium is catapulted out by the rotating evaporating flask in the heating bath.

Risk of burns! The motor and glass parts may be hot during operation and for a long time afterwards. Always wear your personal protective equipment!

If the power fails or the motor stops during operation, you must remove the evaporating flask from the heat source immediately by raising the stand.

The machine may only be opened by trained specialists - even during repairs. The machine is to be unplugged from the mains before opening. Live parts inside the machine may still be live for some time after unplugging from the mains.

Correct use

The **IKA RV 05 basic** is designed for evating solutions and for rapid distillation of liquids. The medium in the operating evaporator flask is circulated by rotation. The surface the liquid is increased significantly by wetting the inner surface of the flask with the liquids. This results in a rapid exchange of heat through a large surface.

Operating in a vacuum allows for gentle evaporation.

Unpacking

Please unpack the equipment carefully and check for any damages. It is important that any damages which may have arisen during transport are ascertained when unpacking. If applicable a fact report must be set immediately (post, rail or forwarder).

The delivery scope covers:

One **IKA RV 05 basic**, one extension, one fastening screw, one hexagon socket screw, one ring spanner and one set of operating instructions.

Assembly

Stand:

The stand should be pre-mounted according to directions in the operating instructions. It should be positioned so that the arrow mark (on the base cross) is pointing forward at the lab table. The heat bath should be placed on the right next to the foot standing in front.

Drive:

First, insert the extension into the bore hole on the back of the device and screw it securely in place with the cylinder head screws enclosed.


Using an R 271 cross sleeve, fasten the drive onto the stand. To do this, clamp the cross sleeve in place on the stand tube about 20 cm above the operating knob of the stand.

It should be easy to tip driver counter clockwise, move it into the cross sleeve with the extension and clamp it in place. The drive must be aligned so that it is sitting vertically above the operating knob of the stand.

Placing the vapour tube in the drive:

Turn the ring (R) to the left unit it will go no further. Push the vapour tube in and then turn the ring (R) to the right until it will go no further.

The disassembly follows the direction in the opposite order.

Glass sets:

Drawings with lists are included with the different glass sets. The glass sets should be mounted on the drive according to these drawings.

Commissioning

Check whether the voltage specified on the type plate matches the mains voltage available. The power socket used must be earthed (protective earth conductor contact). If these conditions are met, the device is ready to operate after plugging in the mains plug. If these procedures are not followed, safe operation cannot be guaranteed and/or the equipment may be damaged. Observe the ambient conditions (temperature, humidity, etc.) listed under "Technical Data".

When working with the evaporator, it has proven useful to follow certain orders while placing the unit in service. These orders are designed to largely prevent foaming, which is highly undesirable, or to at least keep it under control. If the evaporation process is running satisfactorily, a vacuum controller can cover subsequent monitoring.

The IKAVAC VC2 vacuum controller allows for recover rates of up to 99%. Working principle: Vacuum control by assigning a defined pressure value.

Recommended start-up sequence

- Set the equipment to the target pressure using the vacuum control before filling the vacuum flask to minimise solvent losses. (It is of course also possible to fill the vacuum flask before applying the vacuum.)
- The distilland is now sucked into the evaporating flask via the feed pipe and the flask is immersed in the heating bath. This is a more environmentally friendly way of filling the flask.
- Switch on the rotary drive and gradually increase the speed. This prevents temperature stress on the glass when you start the heating process and set the required temperature at the heating bath.

Note:

- Please note that the heating process takes some time.
- Only fine tune the pressure/temperature parameters once the bath temperature has been reached.
- Please bear in mind that the cooler load must not exceed approx. 60%. This corresponds to condensation on approx. 2/3 of the cooling coil.
- Adjust the parameters gradually when using unknown solvents or mixtures.
- When using mixtures of substances, first be aware of the low-boiling solvents and work at a low bath temperature initially.

Eliminating malfunctions

The drive is equipped with an overload protection device that will turn the device off in the event of a long-term overload (motor overheating). If the device is turned off during operation and the orange glow discharge light on the front of the device lights up, the drive is overheating.

If the evaporator flask is immersed in the bath, it should be removed from the heat bath immediately by raising the stand. The drive should be turned off with the On/Off switch and must be left to cool off for some time.

Important: The device must be turned off to cool down. It will not turn back on automatically when the motor has cooled off. The device will not start to run again until it has been turned on again. If the cooling off phase was too short, the orange glow discharge light coming on when the device is turned on again will indicate that the motor is still too hot.

Maintenance and cleaning

The **IKA RV 05 basic** is maintenance-free. It is subject only to the natural wear and tear of components and their statistical failure rate.

When ordering spare parts, please give the manufacturing number shown on the type plate, the machine type and the name of the spare part.

Please send in equipment for repair only after it has been cleaned and is free from any materials which may constitute a health hazard. Please tell us if you have used any materials which are dangerous or which may constitute a health hazard.

Use only cleansing agents which have been approved by **IKA** to clean **IKA** devices.

To remove use:

Dyes	isopropyl alcohol
Construction materials	water containing tenside / isopropyl alcohol
Cosmetics	water containing tenside / isopropyl alcohol
Foodstuffs	water containing tenside
Fuels	water containing tenside

For materials which are not listed, please request information from **IKA**.

Wear the proper protective gloves during cleaning of the devices. Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.

Before using another than the recommended method for cleaning or decontamination, the user must ascertain with the manufacturer that this method does not destroy the instrument.

Accessories

IKAVAC VC 2	Vacuum controller
VC 1.1	Water-jet pump
VC 1.3	Vacuum solenoid valve
VC 2.4	Pump controller
MZ-2C	Diaphragm pump
HB 4 basic	IKA Heating bath
RV 05.5 CN	Set of glassware, vertical, 300 ml (1100 cm ² cooling surface)
RV 05.6 CN	IKA Heating bath
RV 06.1	Set of glassware, diagonal (1200 cm ² cooling surface)
RV 06.2	Set of glassware, vertical (1200 cm ² cooling surface)
RV 06.4	Evaporating flask 1 l
RV 06.5	Evaporating flask 2 l
RV 06.6	Evaporating flask 100 ml
RV 06.7	Receiving flask 1 l
RV 06.11	Rotating vapour tube, ϕ 21,6 mm
RV 06.13	Seal (for RV 06.11), material FKM
RV 06.15	Seal (for RV 06.11), material PTFE

Technical data

Design voltage	VAC	220 - 240+10%
or	VAC	100 - 115+10%
Design frequency	Hz	50/60
Input power	W	133
Output power:	W	65
Speed range 50 Hz	rpm	46 ... 260
60 Hz	rpm	55 ... 312
		(infinitely adjustable)
Speed setting		turn dial on front side
Speed display		scale 1 - 10
		speed table on the side
Drive		speed controlled
		asynchronous motor
Overload protection		Self-holding temperature
		switch in the motor winding
		(switching temperature +135 °C).
		malfunction display with
		orange glow discharge
		light
Perm. duration of operation	%	100
Perm. ambient temperature	°C	+5 to +40
Perm. relative humidity	%	80
Protection type to DIN 60 529		IP 21
Protection class		I
Overvoltage category		II
Contamination level		2
Perm. inclination of the device	°	+ 90
Operation at a terrestrial altitude	m	max. 2000 above sea level
Fuses on mains plate	A	4 (Id.Nr. 25 851 00)
Dimensions (W x D x H) (without attachment)	mm	130x200x260
Weight	kg	4,5

List of spare parts RV 05 basic

Item	Designation
1	Connection line
2	Anti-kink protector
3	Cable clip
5	Angel bracket
8	PCB-power
11	Plastite-oval head screw
12	Conveyor plate
13	Rocker switch
14	Glow lamp
20	Grooved ball bearing
21	Snap ring
22	Shaft seal
23	Close-tolerance grooved pin
24	Application
25	Lock ring
26	Oval head screw M2,5x5
27	Plate spring
28	Bearing tube
29	Threaded pin M6x12
30	Worm wheel
31	V - ring SV-35
32	Screwed flange
33	O - ring
34	Snap ring
35	O - ring
36	Ball bearing-spacer ring
37	Motor support
38	Plate spring
39	A. C. motor
40	Cone pulley
41	Friction wheel
42	Operating knob
43	Scale sign
44	Front plate
45	Snap ring
46	Stopper rod
47	Housing bottom part
48	Annular spring
49	Screw cap
52	Grooved ball bearing
53	Worm shaft
54	Snap ring
55	Driving plate
56	Coupling half
57	Pressure spring
58	Condensor
61	Cylinder head screw M5x12
62	Extension
64	Housing upper partI
65	Speed sign
66	Connection clamp
70	Seal RV 06.13
71	Vapour tube
72	Ring spanner
5001	Fan cowl
5002	Fan wing

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IKA® Works, Inc.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

North America

IKA® Works, Inc.

2635 North Chase Pkwy SE
Wilmington, NC 28405-7419
USA
TEL. +1 800 733-3037
TEL. +1 910 452-7059
FAX +1 910 452-7693
E-mail: usa@ika.net

IKA® Works, (Asia) Sdn Bhd

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

Asia - Australia

IKA® Works (Asia) Sdn Bhd

Lot 2 Jalan Indah 1/2
Taman Industri Rawang Indah
48000 Rawang
Selangor, Malaysia
TEL. +60 3 6093-3322
FAX +60 3 6093-3940
E-mail: ika@tm.net.my

IKA® Japan Y.K.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

Japan

IKA® Japan Y.K.

293-1 Kobayashi-cho
Yamato Koriyama Shi
639-1026 Japan
TEL. +81 74358-4611
FAX +81 74358-4612
E-mail: japan@ika.de

IKA® Works Guangzhou

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

China

IKA® Works Guangzhou

173-175 Friendship Road
Guangzhou Economic & Technological
Development Zone
Guangzhou 510730, P.R.CHINA
TEL. +86 20 8222-6772
FAX +86 20 8222-6776
E-mail: sales@ikagz.com.cn

Artisan Scientific

QUALITY INSTRUMENTATION ... GUARANTEED

Looking for more information?

Visit us on the web at <http://www.artisan-scientific.com> for more information:

- Price Quotations
- Drivers
- Technical Specifications, Manuals and Documentation

Artisan Scientific is Your Source for Quality New and Certified-Used/Pre-owned Equipment

- Tens of Thousands of In-Stock Items
- Hundreds of Manufacturers Supported
- Fast Shipping and Delivery
- Leasing / Monthly Rentals
- Equipment Demos
- Consignment

Service Center Repairs

Experienced Engineers and Technicians on staff in our State-of-the-art Full-Service In-House Service Center Facility

InstraView™ Remote Inspection

Remotely inspect equipment before purchasing with our Innovative InstraView™ website at <http://www.instraview.com>

We buy used equipment! We also offer credit for Buy-Backs and Trade-Ins

Sell your excess, underutilized, and idle used equipment. Contact one of our Customer Service Representatives today!

Talk to a live person: 888-88-SOURCE (888-887-6872) | Contact us by email: sales@artisan-scientific.com | Visit our website: <http://www.artisan-scientific.com>