



## Artisan Technology Group is your source for quality new and certified-used/pre-owned equipment

- FAST SHIPPING AND DELIVERY
- TENS OF THOUSANDS OF IN-STOCK ITEMS
- EQUIPMENT DEMOS
- HUNDREDS OF MANUFACTURERS SUPPORTED
- LEASING/MONTHLY RENTALS
- ITAR CERTIFIED SECURE ASSET SOLUTIONS

### SERVICE CENTER REPAIRS

Experienced engineers and technicians on staff at our full-service, in-house repair center

### *InstraView*<sup>SM</sup> REMOTE INSPECTION

Remotely inspect equipment before purchasing with our interactive website at [www.instraview.com](http://www.instraview.com) ↗

### WE BUY USED EQUIPMENT

Sell your excess, underutilized, and idle used equipment. We also offer credit for buy-backs and trade-ins. [www.artisanng.com/WeBuyEquipment](http://www.artisanng.com/WeBuyEquipment) ↗

### LOOKING FOR MORE INFORMATION?

Visit us on the web at [www.artisanng.com](http://www.artisanng.com) ↗ for more information on price quotations, drivers, technical specifications, manuals, and documentation

**Contact us:** (888) 88-SOURCE | [sales@artisanng.com](mailto:sales@artisanng.com) | [www.artisanng.com](http://www.artisanng.com)

**SPECIFICATIONS**

**STANDARD ACCESSORIES**

Optical system	Hg-Xe lamp (150W)
Light source:	
Measurement	220 ~ 420 nm
Measurement range:	0 ~ ±1000 mdeg
Spectral bandwidth:	20 nm
Wavelength accuracy:	±3 nm
Photometric system	PEM modulation method (single beam)
Photometric method:	0.03 mdeg (at 291 nm)
Noise level:	0.05 mdeg/hr (at 291 nm)
Gain:	x1, x10
Response:	FAST, STD, SLOW
Flow through cell	Front-loading cassette cell
Flow cell path length:	25 mm
Optional flow cell:	1, 5 and 10 mm path cell and high-pressure cell for SFE/SFC
Signal input/output	Recorder output: CD output: 10 mV/full scale UV output: 10 mV/full scale (full scale), S (short) Integrator output: 1V/deg
Others:	Marker input, marker output, leak output 1V/AU and 1V/2AU
Flow through cell	Front-loading cassette cell
Flow cell path length:	25 mm
Optional flow cell:	1, 5 and 10 mm path cell and high-pressure cell for SFE/SFC
Program function	Wave length, measurement range, response, auto zero
Spectrum scanning:	220 ~ 420 nm (both of CD and UV)
Self-diagnosis	Memory: 10 files
Self-diagnosis:	Memory (RAM, ROM), lamp and DC power supply
Others	AC 100 ~ 240V ±10% 50/60Hz Power consumption: 400VA Dimension: 300 (W) x 470 (D) x 230 (H) mm Weight: Approx. 21 kg

- Recorder cable 2 pc
- Marker cable 1 pc
- Syringe, 5 mL 1 pc
- Syringe needle 1 pc
- Compression screw, 1/16" 2 pc
- Ferul, 1/16" 2 pc
- Tube, SST, 1/16" O.D. x 0.8 mm I.D. x 300 mm 1 pc
- Screwdriver, Phillips 1 pc
- Wrench, 1/4" x 5/16" 1 pc
- Spare fuses 2 pc
- AC power cable 1 pc



**JASCO INTERNATIONAL CO., LTD.**

4-21, Senhin-cho 2-chome, Hachioji, Tokyo 193-0835, Japan  
 Tel: +81-42-666-1322 Fax: +81-42-665-6512 Internet: <http://www.jascoinc.co.jp/english/index.html>  
 Australia, China, Hong Kong, India, Indonesia, Iran, Korea, Malaysia, New Zealand, Pakistan, Philippines, Russia, Singapore, South Africa, Taiwan, Thailand  
 8649 Commerce Drive, Easton, Maryland 21601-9903, U.S.A  
 Tel: +1-800-333-5272 Tel: +1-410-822-1220 Fax: +1-410-822-7526 Internet: <http://www.jascoinc.com>  
 Canada, Costa Rica, Mexico, Puerto Rico, Argentina, Brazil, Chile, Colombia, Paraguay, Peru, Uruguay  
**JASCO EUROPE s.r.l.**  
 Via Confalonieri 25, 23894 Cremella (Lc), Italy  
 Tel: +39-039-956439 Fax: +39-039-958642 www.jasco-europe.com  
**JASCO Deutschland** www.jasco.de, **JASCO UK** www.jasco.co.uk, **JASCO France** www.jascofrance.fr,  
**JASCO Benelux** www.jasco.nl, **JASCO Spain** www.jasco-spain.com, **JASCO Scandinavia** www.jascoscandinavia.se  
 Austria, Finland, Greece, Hungary, Poland, Portugal, Romania, Switzerland, Algeria, Cyprus,  
 Egypt, Israel, Jordan, Kuwait, Lebanon, Morocco, Saudi Arabia, Syria, Tunisia, Turkey, U.A.E.

**ORDERING INFORMATION**

6838-J002B	Model CD-2095	Circular Dichroism Detector
6837-J002B	Model OR-2090	Chiral Detector
6824-J002B	Model UV-2070	UV/VIS Detector (~ 900 nm)
6825-J002B	Model UV-2075	UV/VIS Detector (~ 600 nm)
6826-J002B	Model UV-2077	Multi-wavelength UV/VIS Detector
6832-J002B	Model MD-2015	Diode Array Detector (~ 900 nm)
6831-J002B	Model MD-2010	Diode Array Detector (~ 600 nm)
6829-J003B	Model FP-2020	Fluorescence Detector
6833-J002B	Model RI-2031	Refractive Index detector
6836-J002B	Model CL-2027	Chemiluminescence Detector
6818-J002B	Model PU-2080	HPLC Pump
6823-J003B	Model PU-2089	Quaternary HPLC Pump
6820-J002B	Model PU-2085	Semi-micro HPLC Pump
6821-J002B	Model PU-2086	Semi-preparative HPLC Pump
6822-J002B	Model PU-2087	Preparative HPLC Pump
6819-J002B	Model PU-20801	Bio-inert HPLC Pump
6823-J004B	Model PU-20891	Bio-inert Quaternary HPLC Pump
6828-J008B	Model AS-2059	Autosampler
6827-J005B	Model AS-2055	Autosampler
6827-J006B	Model AS-2057	Autosampler
6835-J002B	Model CO-2065	Column Oven



Serving the global marketplace  
 with Analytical Instrumentation  
**JASCO Corporation**  
 2967-5, Wakaba-cho, Hachioji, Tokyo 192-0271, Japan  
 http://www.jasco.co.jp  
 LLC09-0607 Printed in Japan

● Specifications are subject to change without notice.

**CD-2095**

**Circular Dichroism Chiral Detector**



Chiral separation has become one of the most important applications in HPLC. The increasing demand for chiral chromatography in various fields, such as drug analysis, drug discovery, biochemical analysis, natural product analysis, organic synthesis, etc. have resulted in a parallel demand for better chiral separation and detection technology. Recent advances in chiral column technology are boosting the use of chiral chromatography in more and more laboratories, however, advances in chiral detection technology have only recently met the demands required in terms of sensitivity and selectivity.

The new CD-2095 has been developed as an innovative CD (Circular Dichroism) detector, using the same technology applied in conventional CD spectropolarimeters, but optimized to meet HPLC requirements. The CD-2095 detector can simultaneously determine both CD and UV absorption of the sample in the same cell for direct determination of optical purity. CD is based on the differential absorption of right and left circularly polarized light due to the existence of a chiral chromophore and is a much more specific technique than OR (Optical Rotation) or Polarimetry. In the past, large and expensive CD spectropolarimeters have often been modified for use as a CD detector in an HPLC system.

The CD-2095 can easily be interfaced with the JASCO HPLC system using the JASCO chromatography data systems.

- Features**
- Unmatched sensitivity and stability
  - Simultaneous CD and UV detection
  - Direct Optical Purity determination
  - Compact, easy to use package
  - Artifact free design
  - Universal HPLC compatibility
  - Spectral scanning capability



**Highly sensitive detection**  
 Chromatograms of Flavonone (0.1 µg injected) for the CD-2095, the OR-2090 OR detector and the UV-2070/2075 UV/VIS detectors are shown in Fig. 1. A 0.1 µg injection still gave a CD chromatogram showing good signal-to-noise ratio (S/N). The UV detector also offered a good chromatogram. However, the

OR detector failed to give any peak. From S/N calculation, the CD-2095 gave approximately 200 times higher sensitivity than the OR detector. Also it is remarkable that both UV and OR detectors gave significant solvent peaks which were not observable in the data from the CD-2095 as a result of the dual beam signal from the difference in left and right absorption. The absorbance change without CD would be canceled.




Fig. 1 Optical Resolution of Flavonone

**Conditions (Flavonone)**

Column: CHIRALCEL OD  
 Eluent: n-Hexane / IPA (90 / 10)  
 Flow rate: 1.0 mL/min  
 Temperature: 25 °C  
 INJ. VOL.: 10 µL  
 CD-1595  
 Scale: 1 mV = 1 mdeg  
 Response: Standard  
 Polarity: +  
 OR-990  
 Scale: 1 mV = 1 mdeg  
 Response: Standard  
 Polarity: +  
 UV 250 nm  
 S/N = 69.5

Upper: CD chromatogram at 310 nm by CD-2095  
 Middle: OR chromatogram by OR-2090  
 Lower: UV chromatogram by UV-2070/2075





Fig. 2 Flavonone CD Spectra



Schematic diagrams of operation principles for both Optical Rotation and Circular Dichroism techniques are shown in Fig. 3. An OR detector measures an angle of rotation when linearly polarized light passes through a flow cell containing optically active compounds. This is due to difference in light propagation velocities, i.e. refractive indices, of right- and left-handed circularly polarized lights. Since the sample is in a liquid stream, multiple reflections off the cell wall and other effects due to gradient changes occur within the flow cell. These effects cause depolarization and gradient artifacts resulting in a reduction of sensitivity and false peaks.





Fig. 3 Optical Rotation and Circular Dichroism Phenomena

A CD detector can differentiate between enantiomers by measuring the difference in absorbance of right and left-handed circularly polarized light. This principle of detection gives much higher intrinsic stability and high sensitivity. Unlike single beam measurements made by OR based detectors the CD measurement of differential absorption is performed within 20 micro seconds, resulting in virtual dual beam detection. A CD detector can thus generate chiral signals with higher sensitivity and stability and is not effected by gradient elution common in drug analyses.

Schematic diagrams of optical system of the OR-2090 and the CD-2095 are detailed in Fig. 4. In the OR-2090, the light from the source is polarized and modulated by a Faraday cell and is then passed through the flow cell. The optical rotation of the sample is measured by the analyzer which determines the change in the intensity of the output light. The change in the intensity is measured by the photodiode and output via an amplifier. In the CD-2095, the light from the source is dispersed by the diffraction grating in order to select an optimum measurement wavelength. The light is directed to a polarizer, and then modulated alternately between right- and left-handed circularly polarized light. The light is then passed through the flow cell and the intensity of each beam is detected by a photodiode. The two opposite polarized light signals are separated by the electronics and the difference in absorption is determined and expressed as the CD signal. The UV absorption is also measured simultaneously and collected as a second data channel.

**CD Spectral collection**  
 Both CD spectra of Flavonone collected with the CD-2095 and a conventional CD spectropolarimeter JASCO Model J-810 are shown in Fig. 2. The major peaks at 245, 310 and 370 nm are in good agreement with each other.

Fig. 4 Optical system of OR-2090 and CD-2095





## Artisan Technology Group is your source for quality new and certified-used/pre-owned equipment

- FAST SHIPPING AND DELIVERY
- TENS OF THOUSANDS OF IN-STOCK ITEMS
- EQUIPMENT DEMOS
- HUNDREDS OF MANUFACTURERS SUPPORTED
- LEASING/MONTHLY RENTALS
- ITAR CERTIFIED SECURE ASSET SOLUTIONS

### SERVICE CENTER REPAIRS

Experienced engineers and technicians on staff at our full-service, in-house repair center

### *InstraView*<sup>SM</sup> REMOTE INSPECTION

Remotely inspect equipment before purchasing with our interactive website at [www.instraview.com](http://www.instraview.com) ↗

### WE BUY USED EQUIPMENT

Sell your excess, underutilized, and idle used equipment. We also offer credit for buy-backs and trade-ins. [www.artisanng.com/WeBuyEquipment](http://www.artisanng.com/WeBuyEquipment) ↗

### LOOKING FOR MORE INFORMATION?

Visit us on the web at [www.artisanng.com](http://www.artisanng.com) ↗ for more information on price quotations, drivers, technical specifications, manuals, and documentation

**Contact us:** (888) 88-SOURCE | [sales@artisanng.com](mailto:sales@artisanng.com) | [www.artisanng.com](http://www.artisanng.com)