



Artisan Technology Group is your source for quality new and certified-used/pre-owned equipment

- FAST SHIPPING AND DELIVERY
- TENS OF THOUSANDS OF IN-STOCK ITEMS
- EQUIPMENT DEMOS
- HUNDREDS OF MANUFACTURERS SUPPORTED
- LEASING/MONTHLY RENTALS
- ITAR CERTIFIED SECURE ASSET SOLUTIONS

SERVICE CENTER REPAIRS

Experienced engineers and technicians on staff at our full-service, in-house repair center

*InstraView*SM REMOTE INSPECTION

Remotely inspect equipment before purchasing with our interactive website at www.instraview.com ↗

WE BUY USED EQUIPMENT

Sell your excess, underutilized, and idle used equipment. We also offer credit for buy-backs and trade-ins. www.artisanng.com/WeBuyEquipment ↗

LOOKING FOR MORE INFORMATION?

Visit us on the web at www.artisanng.com ↗ for more information on price quotations, drivers, technical specifications, manuals, and documentation

Contact us: (888) 88-SOURCE | sales@artisanng.com | www.artisanng.com

T400 SERIES

Precision *In Vitro* & Extracorporeal Tubing Flow Measurement TS410 Flowmeter Module


Whether you require tubing flow measurements in the animal lab or in the R & D engineering suite, Transonic Systems TS410 Tubing Flowmeter Module gives you precision volume flow measurements in tubing with user friendly features applicable to laboratory bench settings, *in vitro* and extracorporeal use.

Compact single-channel TS410 Flowmeter modules plug into 400-Series Instrumentation Consoles (T402/T403) to mix or match measurement capabilities with other Transonic Systems laboratory modules and provide:

- True Volume Flow Measurement in Tubing
- Gold Standard Transit-time Technology
- Measurement Capability for Diverse Fluids: Blood, Saline, Water, Cell Culture, Physiological Buffer Solutions, Blood Analogs such as Glycerine/Water Solutions
- Programmable Alarm Conditions:
 - High / Low Flow Threshold Alarm
 - Bubble / "Received Signal" Alarm
- On-site Recalibration of Flowsensors for changed or new parameters (or conditions)
- 400-Series Instrumentation Console compatibility: T402, T403.

Two tubing flowsensor styles and multiple sizes from 1.2 mm ID to 25.4 mm ID offer versatile functionality to fit every flow circuit application including:

- Circulatory support device development (VAD Performance & Artificial Heart Development)
- Isolated heart and perfused organ studies
- Mock circulatory models
- Non-contact sterile flow measurement of biologicals and bioreactors
- Flow phantoms and any experimental, or non-human clinical or process application where flow measurement in tubing is needed.




TS410 Tubing Flow Module



ME-PXL Clamp-on Sensor



ME-PXN Inline Sensor



ME-PXN inline Sensors



T400 SERIES

TS410 Tubing Module Specifications

GENERAL FEATURES

Weight/ Size 5.125" high x 4" wide x 9.062" deep; 2.3 lbs.
Module is 2 console slots (20HP) wide.

Console Compatibility 400-Series Flowmeter Consoles (T402, T403):
96 pin DIN connector to proprietary backplane.
Thumb-screw to lock into console on front panel

Power: Derives input power from 400-Series Consoles

AUTOMATIC FLOWMETER ADJUSTMENTS

Sensor size identification and corresponding flow output ranges.
Volume flow calibration of applied sensor; serial number displayed of active Flowsensor.

SET/UP STATUS & PROGRAM PARAMETERS

STATUS MODE (White button labels)
Status message displayed on LCD.

Sensor Status: sensor type & calibration
Meter Status: Active flowmeter settings & alarm status
Alarm Mute: Audible alarm On/Off

PROGRAM MODE: (Blue button labels) - EXIT; SCROLL; SELECT.

SENSOR CONTROLS Select factory preprogrammed calibration options; Adjust Flowsensor gain to change calibration on-site

METER CONTROLS

1/4 Flow Scale: increases flow gain by factor of 4 for low flow measurements

Calibrate Scale: sets output to 0 and 1 Volt to calibrate external recording devices with scale factor flow.

Invert flow: inverts polarity of analog outputs & flow display

Alarms Menu: 3 level program to select, set thresholds, and activate Alarms for "Low Flow", "High Flow" and "Received Signal" Interruption

FILTER PROPERTIES: 0.1, 10, 40 Hz: 2nd order Butterworth, with a third passive pole at 160 Hz 160 Hz: 3rd order Butterworth


ZERO FLOW ADJUST: momentary recessed push button to zero flow reading at stopped flow.

DIGITAL DISPLAY: 4-Digit (14 segment) LED

Displays Volume flow / Probe data / Error Messages

Bar Indicator Light: Displays received signal for continuous monitoring of sensor signal quality

LCD DISPLAY: One line 16-character alpha numeric LCD displays program parameters, sensor and meter status, alarm settings.
Default displays sensor serial number.



Zero Offset Adjustment:

Recessed "one press" button adjusts flowmeter to read zero when flow is stopped.

FLOW OUTPUT: Front panel mounted BNC output connector & rear panel terminal block:

Pulsatile/Average Volume Flow

Filtering controlled by front panel selectable filters

Voltage range: -5 to +5 volts

Output resistance: 500 Ohm

Full Range for Flow: -5 to +5 V

(bidirectional flows, with range of 5 x scale factor flow)

SENSOR CONNECTOR Front panel 16-pin connector

Accepts research Inline and Clamp-on Flowsensors and extension cables with male ODU16-pin connectors.

AUTOMATIC DIGITAL SENSOR IDENTIFICATION & CALIBRATION:

TS410 circuitry reads connector operational data (size, scale & calibration) programmed in the sensor's EPROM.

ULTRASONIC FREQUENCY RANGE

600 KHz to 14.4 MHz; Sensor size dependent

FLOWSENSOR COMPATIBILITY

Accepts Transonic "Precision" Tubing Flowsensors with 16-pin ODU connectors

ME-Prefix: Extracorporeal Sensors for 400-Series Flowmeters

PXL-Series: Precision XL-Series Clamp-on Sterile Tubing Flowsensors with 4 transducers

PXN-Series: Precision XN-Series Inline Flowsensors with 4 transducers

SIGNAL OUTPUTS: 8 accessible signals via 400-series Flowmeter

Console's odd numbered back-panel terminal block connections:

Pulsatile Volume Flow; Mean Volume Flow

Received Signal Amplitude (2); Phase (4)

SYNCHRONIZATION: Module rear panel switch for selecting synchronization mode

Self-triggering Mode: "SYNC IN" to "SYNC OUT" jumper on each module:

Sequential Triggering Mode: "SYNC IN" crossed to "SYNC OUT" between modules.



Artisan Technology Group is your source for quality new and certified-used/pre-owned equipment

- FAST SHIPPING AND DELIVERY
- TENS OF THOUSANDS OF IN-STOCK ITEMS
- EQUIPMENT DEMOS
- HUNDREDS OF MANUFACTURERS SUPPORTED
- LEASING/MONTHLY RENTALS
- ITAR CERTIFIED SECURE ASSET SOLUTIONS

SERVICE CENTER REPAIRS

Experienced engineers and technicians on staff at our full-service, in-house repair center

*InstraView*SM REMOTE INSPECTION

Remotely inspect equipment before purchasing with our interactive website at www.instraview.com ↗

WE BUY USED EQUIPMENT

Sell your excess, underutilized, and idle used equipment. We also offer credit for buy-backs and trade-ins. www.artisanng.com/WeBuyEquipment ↗

LOOKING FOR MORE INFORMATION?

Visit us on the web at www.artisanng.com ↗ for more information on price quotations, drivers, technical specifications, manuals, and documentation

Contact us: (888) 88-SOURCE | sales@artisanng.com | www.artisanng.com