## Artisan Technology Group is your source for quality new and certified-used/pre-owned equipment

- FAST SHIPPING AND DELIVERY
- TENS OF THOUSANDS OF IN-STOCK ITEMS
- EQUIPMENT DEMOS
- HUNDREDS OF MANUFACTURERS SUPPORTED
- LEASING/MONTHLY RENTALS
- ITAR CERTIFIED SECURE ASSET SOLUTIONS

### SERVICE CENTER REPAIRS

Experienced engineers and technicians on staff at our full-service, in-house repair center

### *InstraView*<sup>SM</sup> REMOTE INSPECTION

Remotely inspect equipment before purchasing with our interactive website at [www.instraview.com](http://www.instraview.com) ↗

### WE BUY USED EQUIPMENT

Sell your excess, underutilized, and idle used equipment. We also offer credit for buy-backs and trade-ins. [www.artisanng.com/WeBuyEquipment](http://www.artisanng.com/WeBuyEquipment) ↗

### LOOKING FOR MORE INFORMATION?


Visit us on the web at [www.artisanng.com](http://www.artisanng.com) ↗ for more information on price quotations, drivers, technical specifications, manuals, and documentation

**Contact us:** (888) 88-SOURCE | [sales@artisanng.com](mailto:sales@artisanng.com) | [www.artisanng.com](http://www.artisanng.com)

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#### DESCRIPTION

Adjustable volume multichannel micropipette for use with detachable polypropylene tips, for an easy and continuous work with standard 96-well microplates, 9 mm spacing.

The sample is drawn into the polypropylene tip without any contact with the plunger assembly. The integrated tip ejector is located on the handpiece allowing the automatic release of the tips.

Simultaneous delivery of 8 or 12 samples from 5 µl to 200 µl by means of two pipette sizes.

Each instrument is individually calibrated and comes along with its control certificate.

#### HOLDING THE MULTICHANNEL MICROPIPETTE

Recommended way of holding the pipette is shown on fig. 1: place fingerrest on the phalanx of your forefinger with the ejector knob turned towards yourself.

The multichannel casing revolves on the handpiece and allows the selection of the best working position.

#### VOLUME SETTING (fig. 2)

Set desired volume by rotating graduated knob No 31. Graduation mark must be even with vernier ring, vertical line on zero position. By moving vertical line clockwise along vernier graduation, volume will increase by 1 µl per division.

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
Fig. 1 Operation

## TIPS AND FILTERS

For consistent results and proper operation of the ejector, use only our original polypropylene tips, or their equivalent.

These are supplied on 96 pieces racks, with the 9 mm spacing of the barrels of the multichannel micropipette.

Place the instrument on a row of tips, press firmly to ensure a solid fit of the tips on each channel to avoid leakage.

**High-sitting tips** (collar longer than regular tips) can be used after easy replacement of ejector head No 3 with alternative HST-Ejector head available from your distributor.

## PIPETTING OPERATION (fig. 1)

- Phase 1** Depress plunger knob to the intermediate stop ending calibrated stroke before immersing tips into one of the various reservoirs.
- Phase 2** Immerse tips approximately 5 mm and *smoothly* release plunger knob, thus drawing up simultaneously liquid into the 8 or 12 tips. Withdraw tips from sample solution, sliding them along wall of the reservoir.
- Phase 3** Place tips on the wells of the microplate and slowly depress plunger to first stop, thus emptying tips; after approximately one second, depress knob to lower discharge position: the excess air volume produced by depressing the plunger beyond the calibrated stroke blows out possible fluid particles which may cling to the tips. Remove the multichannel pipette from the wells of the microplate, with plunger knob in its lowest position. Release knob slowly. To eject used tips, press ejector knob E. Apply new tips for every liquid; change tips if droplets are retained after blow out. For better accuracy, fill new tip once before sampling.


Fig. 2 Setting of the volume

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## MAINTENANCE

To check tightness of the multichannel pipette, draw liquid into the tips. Droplet forming at the extremity of one of the tips means this channel is leaking. O-ring of the corresponding barrel has to be changed.

Fluid getting accidentally into the multichannel pipette may clog one of the barrels or cause the plungers to drag.

In order to clean the plungers or to change the O-rings, remove the two screws No 4.

Cleaning the plungers: depress the plunger knob No 31: the plungers will protrude from the casing. Clean them, grease slightly and wipe with a clean cloth.

Changing the O-rings: to have access to the barrels No 6, pull out the two ejector shafts No 12 off the ejector head No 3 by pressing together the fingers of the ejector shaft: the ejector head will come out.

Remove the two screws No 5 together with spacers No 18. Take out the leaking barrel and change the two O-rings No 9 and 2 on the pipette 851-200 and the O-ring No 2 on the pipette 851-050.

Depress the plunger knob No 31 when tightening the screws No 4 and 5.

## CALIBRATION

Each Acura® has been calibrated permanently at the factory and individually controlled. The Q.C. certificate accompanying the instrument includes control data, serial number and operator identification. Calibration is carried out at a constant temperature (+/- 0,5°C) comprised between 20 and 25°C with Aqua bidest, according to DIN 12650, part 6.

Under normal conditions, a pipette does not need to be recalibrated. After replacing a part, which could influence the measure, the calibration should be checked, preferably by gravimetry according to the above mentioned norm; an analytical balance is needed. Ten measurements should be performed after checking the seal tightness of the instrument.

**In-lab calibration (fig. 3):** Whenever required, for instance when pipetting viscous or non-aqueous liquids, calibration adjustments can be performed by laboratory personnel. Reach calibration screw located inside plunger button by removing comfort pad a). Unlock side allen screw b), then delicately turn calibration screw c) clockwise to decrease volume, and vice versa. Re-tighten locking screw b). Check resulting volumes after each fine tuning. Cover push-button with comfort pad again.

## WARRANTY

Your multichannel micropipette is guaranteed for a period of two years against any factory defect. Return it to the dealer from whom you purchased it if the instrument presents a defect which you cannot eliminate by applying the ordinary maintenance procedure as described. The identification is facilitated by its serial number.

**Decontaminate instrument prior to returning it.**

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