



Artisan Technology Group is your source for quality new and certified-used/pre-owned equipment

- FAST SHIPPING AND DELIVERY
- TENS OF THOUSANDS OF IN-STOCK ITEMS
- EQUIPMENT DEMOS
- HUNDREDS OF MANUFACTURERS SUPPORTED
- LEASING/MONTHLY RENTALS
- ITAR CERTIFIED SECURE ASSET SOLUTIONS

SERVICE CENTER REPAIRS

Experienced engineers and technicians on staff at our full-service, in-house repair center

*InstraView*SM REMOTE INSPECTION

Remotely inspect equipment before purchasing with our interactive website at www.instraview.com ↗

WE BUY USED EQUIPMENT

Sell your excess, underutilized, and idle used equipment. We also offer credit for buy-backs and trade-ins. www.artisanng.com/WeBuyEquipment ↗

LOOKING FOR MORE INFORMATION?

Visit us on the web at www.artisanng.com ↗ for more information on price quotations, drivers, technical specifications, manuals, and documentation

Contact us: (888) 88-SOURCE | sales@artisanng.com | www.artisanng.com



YSI 2300 STAT PLUS™

Proven Technology for Glucose and Lactate



YSI glucose and lactate analyzers have been recognized standards for the measurement of blood glucose and/or lactate for more than 25 years. The YSI 2300 is known for its exceptional speed, performance, and ease of use. It is a laboratory glucose and lactate analyzer designed for flexibility.

Measure glucose and lactate simultaneously in whole blood or plasma. Alternatively configure the analyzer to measure only glucose or only lactate. Purchase only the reagents you need.

Fast, accurate results in whole blood, plasma, serum, or cerebrospinal fluid

In addition to **speed** (one-minute results), **automatic calibration**, and **small sample size** (25 μ L), the YSI 2300 STAT PLUS provides:

The Standard for Clinical Diagnostics

- STAT lab or central lab
- Diabetes evaluation
- Diabetes/endocrinology research
- Handheld monitor development
- Human performance evaluation
- Exercise/sports physiology
- Cardiac/pulmonary rehabilitation
- Neonatology/pediatrics
- General metabolic research
- Veterinary medicine and research
- Pharmaceutical quality assurance

- **Simple, safe operation:** Menu-driven software guides users through each step. Diagnostics provide performance analysis and troubleshooting. Samples are aspirated from a variety of tube sizes; no needles required.
- **Whole-blood capability:** No need to spin samples. And for lactate samples, lysing agent may be used in the electrode chamber for total blood lactate measurements.
- **Low operating cost:** Immobilized enzyme sensors last 2-3 weeks, minimizing the need for costly reagents and supplies.
- **Batch runs:** Optional 24-position turntable provides batch operation capability for plasma or serum. Each batch and sample is assigned an identification number with date and time.
- **Interference rejection:** Unaffected by many of the physical and chemical interferences that plague other methods.
- **Broad measurement range:** Virtually eliminates the need for sample dilution, measuring glucose to 900 mg/dL (50 mmol/L) and lactate to 30 mmol/L.

www.YSIlifesciences.com



YSI Immobilized Enzyme Biosensor Technology


The key to fast, accurate measurement

An enzyme specific for the substrate of interest—glucose oxidase or L-lactate oxidase—is immobilized between two membrane layers, polycarbonate and cellulose acetate. The substrate is oxidized as it enters the enzyme layer, producing hydrogen peroxide, which passes through cellulose acetate to a platinum electrode, where the hydrogen peroxide is oxidized. The resulting current is proportional to the concentration of the substrate.

YSI membranes contain three layers. The first layer, porous polycarbonate, limits the diffusion of the substrate into the second (enzyme) layer preventing the reaction from becoming enzyme-limited. The third


layer, cellulose acetate, permits only small molecules (such as hydrogen peroxide) to reach the electrode, eliminating many electrochemically-active compounds that could interfere with the measurement.

Immobilized Enzyme Reaction



D-glucose is oxidized, producing hydrogen peroxide and glucono-lactone.

Anode Reaction



Hydrogen peroxide (H_2O_2) is oxidized at the platinum anode, producing electrons. The electron flow is proportional to the H_2O_2 concentration and to the concentration of the substrate.

Call your local
representative or
YSI technical support!

800 659-8895

+1 937 767-7241
support@ysi.com

www.YSIlifesciences.com

YSI Incorporated
1700/1725 Brannum Lane
Yellow Springs, Ohio 45387

CE partner4U BV
Esdoornlaan 13
3951DB Maarn
The Netherlands
+31 (0) 6 516 536 26

YSI UK Ltd.
01252 819844
Fax 01252 819544

YSI (HK) Limited
852 2891 8154
Fax 852 2834 0034

Printed in USA on recycled paper.
0508 B53GL



Instrument Specifications

Sample size: 25 μL (aspirated volume)

Response time (from test tube holder)

Normal mode: Displayed and printed result in 65 seconds. Sample to sample interval is 100 seconds.

Screen mode: Displayed and printed result in 45 seconds. Sample to sample interval is 70 seconds.

Measurement range

Normal mode: Glucose: to 900 mg/dL (9000 mg/L, 50.0 mmol/L)

Lactate: to 30.0 mmol/L (267 mg/dL, 2670 mg/L)

Screen mode: Glucose: to 500 mg/dL (5000 mg/L, 27.8 mmol/L)

Lactate: to 15.0 mmol/L (134 mg/dL, 1335 mg/L)

Precision

Glucose: Whichever is larger, $\pm 2\%$ of reading or 2.5 mg/dL, (25 mg/L or 0.2 mmol/L)

Lactate: Whichever is larger, $\pm 2\%$ of reading or 0.1 mmol/L, (1 mg/dL or 10 mg/L)

Calibration point

Glucose: 180 mg/dL (1800 mg/L, 10.0 mmol/L)

Lactate: 5.00 mmol/L (45 mg/dL, 445 mg/L)

Membrane working life (typical)

Glucose: 21 days

Lactate: 14 days

Operating temperature: 15 to 35° C

Humidity: 10 to 90% noncondensing

Power: 110-120 VAC or 220-240 VAC, 50-60 Hz, 50 Watts nominal

Size: 25.4 x 35.6 x 35.6 cm, 11.4 kg

10 x 14 x 14 inches, 25 lbs

Regulatory Compliance: CE, FDA Registered Class II (862.1345), CSA

CE Device Category: IVD

Note: YSI makes no performance claims for sampling whole blood using the YSI 2710 Turntable.
STAT PLUS is a trademark of YSI Incorporated.



Artisan Technology Group is your source for quality new and certified-used/pre-owned equipment

- FAST SHIPPING AND DELIVERY
- TENS OF THOUSANDS OF IN-STOCK ITEMS
- EQUIPMENT DEMOS
- HUNDREDS OF MANUFACTURERS SUPPORTED
- LEASING/MONTHLY RENTALS
- ITAR CERTIFIED SECURE ASSET SOLUTIONS

SERVICE CENTER REPAIRS

Experienced engineers and technicians on staff at our full-service, in-house repair center

*InstraView*SM REMOTE INSPECTION

Remotely inspect equipment before purchasing with our interactive website at www.instraview.com ↗

WE BUY USED EQUIPMENT

Sell your excess, underutilized, and idle used equipment. We also offer credit for buy-backs and trade-ins. www.artisanng.com/WeBuyEquipment ↗

LOOKING FOR MORE INFORMATION?

Visit us on the web at www.artisanng.com ↗ for more information on price quotations, drivers, technical specifications, manuals, and documentation

Contact us: (888) 88-SOURCE | sales@artisanng.com | www.artisanng.com