

ARTISAN®
TECHNOLOGY GROUP

Your **definitive** source for quality pre-owned equipment.

Artisan Technology Group
(217) 352-9330 | sales@artisantg.com | artisantg.com

Full-service, independent repair center
with experienced engineers and technicians on staff.

We buy your excess, underutilized, and idle equipment
along with credit for buybacks and trade-ins.

Custom engineering
so your equipment works exactly as you specify.

- Critical and expedited services
- In stock / Ready-to-ship
- Leasing / Rentals / Demos
- ITAR-certified secure asset solutions

Expert team | Trust guarantee | 100% satisfaction

All trademarks, brand names, and brands appearing herein are the property of their respective owners.

Find the *Argonaut Technologies Quest 210* at our website: ***Click HERE***

Blank RV's

ARGONAUT TECHNOLOGIES

Blank Reaction Vessels are used to occupy unused reaction vessel positions on the Quest 210. This eliminates wasting reaction vessels when less than 10 reactions are perfomed per reaction vessel bank.

Instructions. Before installing the Blank Reaction Vessels, clean the Quest 210 using the cleaning procedure described in the User's manual. Orient the Blank Reaction Vessel so the flat side is facing down. Insert the Blank Reaction Vessel between the aluminum heating plate and the spring loaded insulators on the Reactor Unit. Seat the Blank Reaction Vessel against the lower manifold. Lower the upper manifold and lock it into place. When the synthesis is complete, raise the upper manifold and remove the Blank Reaction Vessels for reuse (Note: Do not use the reaction vessel extraction pliers to remove the Blank Reaction Vessels from the Reactor Unit. Use of these pliers will damage the Blank Reaction Vessels). When using the Blank Reaction Vessels, solvent delivery and draining will be faster. If used with the Automated Solvent Wash Module, solvent delivery times must be adjusted accordingly to ensure the reaction vessels do not overfill. If a Blank Reaction Vessel leaks, replace the two Teflon® inserts inside the open side of the vessel.

Part Number: 900146 - 5 mL,
900147 - 10 mL
Quantity: 5

Replacement Parts:
900185 Teflon inserts (40/pkg)

Quest Accessories

P/N 102555

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Artisan Technology Group is an independent supplier of quality pre-owned equipment

Gold-standard solutions

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

Learn more!

Visit us at **artisantg.com** for more info on price quotes, drivers, technical specifications, manuals, and documentation.

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

We're here to make your life easier. How can we help you today?

(217) 352-9330 | sales@artisantg.com | artisantg.com

ARTISAN® TECHNOLOGY GROUP