ARTISAN®
TECHNOLOGY GROUP

Your **definitive** source for quality pre-owned equipment.

Artisan Technology Group
(217) 352-9330 | sales@artisantg.com | artisantg.com

Full-service, independent repair center
with experienced engineers and technicians on staff.

We buy your excess, underutilized, and idle equipment
along with credit for buybacks and trade-ins.

Custom engineering
so your equipment works exactly as you specify.

- Critical and expedited services
- In stock / Ready-to-ship
- Leasing / Rentals / Demos
- ITAR-certified secure asset solutions

Expert team | Trust guarantee | 100% satisfaction

All trademarks, brand names, and brands appearing herein are the property of their respective owners.

Find the *Cole-Parmer 12122-36* at our website: *Click HERE*

polystat®

Cole-Parmer®

Cooling/Heating Recirculating Chillers

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Table of Contents

Quick Starts

Preface **i**
Compliance i
Warranty i
Unpacking i

Section 1 **Safety** **1-1**
Safety Warnings 1-1

Section 2 **General Information** **2-1**
Description 2-1
Specifications 2-1
Wetted Parts 2-4

Section 3 **Installation** **3-1**
Ambient Conditions 3-1
Ventilation 3-1
Electrical Requirements 3-2
Plumbing Requirements 3-3
Fluid Considerations 3-4
Approved Fluids 3-5
Filling Requirements 3-7
Draining 3-8
Shipping/Storage 3-8

Section 4 **Operation** **4-1**
Controller 4-1
Setup 4-2
Start Up 4-2
Status Display 4-3
Stand By Mode 4-3
Changing the Setpoint 4-4
Menu Displays 4-5
Menu 4-5
Menu Tree 4-6
Settings 4-7
System 4-12
High Temperature Cutout 4-15
Stopping the Chiller 4-16
Power Down 4-16
Shut Down 4-16
Restarting 4-16

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Section 5 Preventive Maintenance 5-1
Cleaning 5-1
Condenser 5-1
Fluid Maintenance 5-1
Sight Tube 5-1
Hoses 5-1
Testing the Safety Features 5-2
Decommissioning/Disposal 5-2

Section 6 Troubleshooting 6-1
Error Displays 6-1
Check List 6-3

Appendix AC Serial Communications Protocol A-1

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Safety Precautions:

The unit is designed for indoor use only.

Never place unit in a location where excessive heat, moisture, inadequate ventilation, or corrosive materials are present.


Never use flammable or corrosive fluids with this unit.

Never connect process fluid lines to your facility water supply or to any pressurized liquid source.

Before using any fluid or performing maintenance where contact with the fluid is likely refer to the manufacturer's MSDS for handling precautions.

What you need to get started:

- An adjustable wrench
- Appropriate hose or plumbing
- Appropriate size clamps or connection type

1 Ensure that all shipping plugs are removed before installation.


The process fluid connections are located on the rear of the unit and are labeled (PROCESS OUTLET) and (PROCESS INLET).
Connect the PROCESS OUTLET to the fluid inlet on your application. Connect the PROCESS INLET to the fluid outlet on your application. Ensure all connections are secure.

Keep the distance between the unit and the instrument being cooled as short as possible. To minimize back pressure, ensure tubing is as straight as possible. If diameter reductions are required, make them at the inlet and outlet of your application, not at the chiller.

2 Ensure the reservoir overflow located on the rear of the unit is connected to a suitable drain.

⚠ CAUTION **To prevent damage to the unit's plumbing, use a 19 mm backing wrench when removing/installing the external connections. ▲**

Table 1 - Acceptable Fluids:
Use of any fluid not listed below will void the manufacturer's warranty.
5°C to 80°C — Distilled Water or Deionized Water (up to 3 MΩ-cm)
-10°C to 80°C — 50/50 Water with Glycol


3 Lift up on the fill panel on the front of the unit to remove it.

Start with 2.8 liters of fluid.

Insert the supplied funnel into the hose and *slowly* add fluid. **NOTE** Due to the unit's plumbing configuration, adding fluid too quickly may result in overfilling - the fill level indication lags the actual amount added.

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
Avoid overfilling, fluids expand when heated.

4 Verify the appropriate voltage. For units supplied with a line cord, insert female end of power cord into chiller and then insert male end of power cord into power outlet. (The line cord is located under the shipping crate's lid. Do not discard the lid until the cord is located.)

5 Place the circuit protector to the on (I) position.

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Start Up

Units should be left in an upright position at room temperature (~25°C) for 24 hours before starting. This will ensure the lubrication oil has drained back into the compressor.

Before starting the unit, double check all USB (optional), electrical and plumbing connections.

Do not run the unit until fluid is added to the unit. Have extra fluid on hand. If the unit will not start refer to the manual.

- Place the circuit protector located on the front of the unit to the **I** position.


- Press , the Start Display will appear.
- Ensure the start symbol has a highlight box around it, if not use the arrow keys to navigate to the symbol.
- Press . The unit will start and the start symbol will turn into a stop symbol (■).

NOTE It will take 30 seconds for the compressor to start.


After start, check all external plumbing connections for leaks.

Start Display

The **SP1** and **Menu** portions on the top of the display are used to view and/or change the controller's settings. They are explained in detail in the manual.

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Shut Dowm

Ensure the stop symbol has a box around it, if not use the arrow keys to navigate to the symbol.

Press . The unit will stop and the stop symbol will turn into a start symbol (▶).

Press . The thermostat's screen will go blank.

Place the circuit protector on the front of the unit to the **O** position.

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Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Preface

Compliance

Refer to the Declaration of Conformity in the back of this manual.

Warranty

Cole-Parmer warrants to the direct purchaser that the Cole-Parmer Polystat Cooling/Heating Recirculating Chillers, distributed by Cole-Parmer, will be free from defects in material or workmanship for a period of two years from the date of delivery.

Cole-Parmer will repair or replace the product or provide credit, as its sole option, upon prompt notification and compliance with its instructions.

The Distributor warrants to Customer that upon prompt notification and compliance with Distributor's instructions, that the Distributor will repair or replace, at Distributor's sole option, any Product which is defective in material or workmanship.

Distributor expressly disclaims all other warranties, whether expressed, implied or statutory, including the warranties of merchantability, and fitness for a particular purpose. Distributor's sole responsibility and the Customer's exclusive remedy for any claim arising out of the purchase of any Product is repair or replacement, as described above. In no event shall Distributor's liability exceed the purchase price pain therefor; nor shall Distributor be liable for any claims, losses or damage of any third party or for lost profits or any special, indirect, incidental, consequential, or exemplary damages, howsoever arising, even if Distributor has been advised of the possibility of such damages.

Unpacking

Retain all cartons and packing material until the chiller is operated and found to be in good condition. If the chiller shows external or internal damage contact the transportation company and file a damage claim. Under ICC regulations, this is your responsibility.

The chiller does not have handles. Take into account the chiller's weight, 66 pounds (30 kilograms), when unpacking and transporting. We recommend two people lift the chiller from the bottom. ▲

Chillers should be left in an upright position at room temperature for 24 hours before starting. This will ensure the lubrication oil has drained back into the compressor. ▲

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Section 1 Safety

Safety Warnings

Make sure you read and understand all instructions and safety precautions listed in this manual before installing or operating your chiller. If you have any questions concerning the operation of your chiller or the information in this manual, please contact us.

DANGER indicates an imminently hazardous situation which, if not avoided, *will* result in death or serious injury.

WARNING indicates a potentially hazardous situation which, if not avoided, *could* result in death or serious injury.

CAUTION indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury. It is also be used to alert against unsafe practices.

The lightning flash with arrow symbol, within an equilateral triangle, is intended to alert the user to the presence of non-insulated "dangerous voltage" within the chiller's enclosure. The voltage magnitude is significant enough to constitute a risk of electrical shock.

This label indicates the presence of hot surfaces.

This label indicates read the manual.

NOTE The chiller's equipment design incorporates a complete sheet metal enclosure for personnel protection from mechanical and electrical hazards. ▲

Observe all warning labels. ▲

Never remove warning labels. ▲

The chiller construction provides protection against the risk of electrical shock by grounding appropriate metal parts. The protection will not function unless the power cord is connected to a properly grounded outlet. It is the user's responsibility to assure a proper ground connection is provided. ▲

Operate the chiller using only the supplied line cords, never operate equipment with damaged cords. ▲

Ensure all communication and electrical connections are made prior to starting the chiller. ▲

The chiller's power cord is the electrical disconnecting device, it must be easily accessible at all times. ▲

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Always turn the chiller off and disconnect the supply voltage from its power source before moving performing any service or maintenance procedures. Ensure chiller is at a safe temperature before handling. ▲

Never place the chiller in a location or atmosphere where excessive heat, moisture, or corrosive materials are present. ▲

Chillers should be left in an upright position at room temperature (~25°C) for 24 hours before starting. This will ensure the lubrication oil has drained back into the compressor. ▲

The chiller is not designed to be floor mounted. ▲

Other than water, before using any fluid, or when performing maintenance where contact with the fluid is likely, refer to the manufacturer's MSDS and EC Safety Data sheet for handling precautions. ▲

Ensure, that no toxic gases can be generated by the fluid. Flammable gases can build up over the fluid during usage. ▲

Never operate the chiller without fluid in the reservoir. ▲

Never use corrosive or flammable fluids with this chiller. Use of these fluids will void the manufacturer's warranty. ▲

Ensure the fluid is at a safe temperature (20°C to 55°C) before handling or draining. ▲

Drain the chiller before it is transported and/or stored in near or below freezing temperatures, see Draining in Section 3. ▲

Ensure the tubing you select will meet your maximum temperature and pressure requirements. ▲

Never operate damaged or leaking equipment. ▲

Never operate the chiller or add fluid to the reservoir with panels removed. ▲

Transport the chiller with care. Sudden jolts or drops can damage its components. ▲

Do not clean the chiller with solvents, only use a soft cloth and water. ▲

Refer service and repairs to a qualified technician. ▲

Performance of installation, operation, or maintenance procedures other than those described in this manual may result in a hazardous situation and will void the manufacturer's warranty. ▲


Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Section 2 General Information

Description

The Cole Parmer Polystat CR Cooling/Heating Recirculating Chillers are designed to provide a continuous supply of fluid at a constant temperature and flow rate. The chiller consists of an air-cooled refrigeration system, heat exchanger, recirculating pump, polyethylene reservoir, and a microprocessor controller.

Specifications


	CR 250	CR 500
Procedd Fluid Temperature Range	-10°C to +80°C +14°F to +176°F	-10°C to +80°C +14°F to +176°F
Temperature Stability	±0.1°C	±0.1°C
Cooling Capacity[1]		
Heater Capacity at 20°C 50Hz 60Hz	2.0 Kilowatts 1.2 Kilowatts	2.0 Kilowatts 1.2 Kilowatts
Nominal Weight kg lb	30 66	30 66
Pumping Type[2]	Force / Force-Suction	Force / Force-Suction
Maximum Flow Rate lpm gpm	17 / 21 4.5 / 5.5	17 / 21 4.5 / 5.5
Maximum Pressure mbar psi	300 / 805 4.4 / 11.7	300 / 805 4.4 / 11.7
Reservoir Volume liters gallons	2.8 0.7	2.8 0.7

1. Specifications obtained at sea level using water (above 5°C) and 50/50 EG/Water (<5°C) as the recirculating fluid at a 20°C process setpoint, 20°C ambient condition, at nominal operating voltage. Other fluids, process temperatures, ambient temperatures, altitude or operating voltage will affect performance.
2. Nominal values ±10%. Depends on chiller's pump type, force only or force/suction. Derate 100V/50Hz chillers ~15% for maximum pressure and maximum flow.

• Cole Parmer reserves the right to change specifications without notice.


Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Dimensions
(inches/centimeters)


Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Rear View

Side View

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Wetted Materials

Viton
EPDM
Ryton
Ultem
Vectra
Stainless Steel 316
Stainless Steel 304
Ryton

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Section 3 Installation

Ambient Conditions

Ambient Temperature Range*	10°C to 40°C (50°F to 104°F)
Maximum Relative Humidity (Non Condensing)*	80% at 31°C (88°F)
Operating Altitude*	Sea Level to 2000 meters (6560 feet)
Overvoltage Category	II
Pollution Degree	2
Degree of Protection	IP 20

* Maximum relative humidity 80% for temperatures up to 31°C decreasing linearly to 50% relative humidity at 40°C. Derate ambient by 1°C/1000m and capacity by 2%/1000m.

The chiller is designed for continuous operation and for indoor use. Also, it is not designed to be floor-mounted.

⚠ CAUTION

Never place the chiller in a location where excessive heat, moisture, inadequate ventilation, or corrosive materials are present. ▲

⚠ CAUTION

Chillers should be left in an upright position at room temperature (~25°C) for 24 hours before starting. This will ensure the lubrication oil has drained back into the compressor. ▲

Ventilation

The chiller requires clean air for proper operation. Air enters from the rear of the chiller and exits through the sides and front.

The chiller is design to allow blocking of any two sides (except rear) without affecting rated load The ambient temperature is reduced by 3°C if two sides are blocked.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Electrical Requirements

DANGER

The chiller construction provides protection against the risk of electrical shock by grounding appropriate metal parts. The protection will not function unless the power cord is connected to a properly grounded outlet. It is the user's responsibility to assure a proper ground connection is provided. ▲

The chiller is intended for use on a dedicated outlet. All chillers are equipped with a C14 main power inlet and a 20 Amp circuit protector.

NOTE If the circuit protector activates allow the chiller to cool before resetting the protector. Restart the chiller. Contact us if it activates again. ▲

CAUTION

If the chiller's power cord is used as the disconnecting device, it must be easily accessible at all times. ▲

CAUTION

Operate the chiller using only the supplied line cords, never operate equipment with damaged cords. ▲

Refer to the nameplate on the rear, upper-left-hand corner of the chiller for specific electrical requirements. Voltage deviations of ± 10% are permissible. The outlet must be rated as suitable for the total power consumption of the chiller.

The following power options are available:

Volts[1]/Hertz/Phase	Amps[2]	Total Wattage	Plug Type
115/60/1	12	1345	N5-15
100/50-60/1	12	1150	N5-15
220- 230/50/1	12	2395	Country Specific
220/60/1	12	2395	Country Specific

1. Volts ± 10%
2. Maximum amp draw

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Plumbing Requirements


⚠ CAUTION

Ensure that all shipping plugs are removed before installation.

Never connect the process fluid lines to your facility water supply or any pressurized liquid source. ▲

The process fluid connections are located on the rear of the chiller and are labeled (PROCESS OUTLET) and (PROCESS INLET). Connect the PROCESS OUTLET to the fluid inlet on your application. Connect the PROCESS INLET to the fluid outlet on your application. Ensure all connections are secure.

Keep the distance between the chiller and the instrument being cooled as short as possible. To minimize back pressure, ensure tubing is as straight as possible. If diameter reductions are required, make them at the inlet and outlet of your application, not at the chiller.

⚠ CAUTION

To prevent damage to the chiller's plumbing, use a 19 mm backing wrench when removing/installing the external connections. ▲

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Fluid Considerations

DANGER

Only use the approved fluids listed on the next page. Never use corrosive or flammable fluids with this chiller. ▲

CAUTION

Handle and dispose all liquids, other than water, in accordance with the fluid manufacturers specification and/or the MSDS for the fluid used. ▲

CAUTION

Always adjust the chiller's software to the fluid used, see Section 4. ▲

CAUTION

Water/ethylene glycol mixtures require top-offs with pure water, otherwise the percentage of glycol will increase resulting in high viscosity and poor performance. ▲

Cole Parmer takes no responsibility for damages caused by the selection of an unsuitable bath fluid.

Unsuitable bath fluids are fluids which:

- are very highly viscous (much higher than 30 mPas at the respective working temperature)
- have corrosive characteristics or
- tend to break down at high temperatures

For fluid selection consider application requirements, operating temperature range, material compatibility, safety concerns, and environmental issues.

Chlorine

Short term usage of tap water may not cause any adverse affects on the chiller or your application, but in the long term problems may arise. To help alleviate these problems Cole Parmer recommends the use of chlorine.

The duration of time that chlorine remains in solution depends on factors such as water temperature, pH and availability of direct sunlight. We recommend maintaining chlorine levels at proper levels using chlorine test strips, generally 1 to 5 ppm is adequate.

For best results, maintain the pH of the fluid between 6.5 and 7.5. Additional chlorine should not be added without first determining the concentration ratio that already exists in the fluid supply. Corrosion and degradation of the circulation components can result from concentration ratios that are too high. Contact us for additional information.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Approved Fluids

Distilled Water or Deionized Water (up to 3 MΩ-cm)
Normal tap water leads to calcareous deposits necessitating frequent chiller decalcification, see table on next page.

Calcium tends to deposit itself on the heating element. The heating capacity is reduced and service life shortened.

50/50 Water with Laboratory Grade Ethylene Glycol
Below 5°C water has to be mixed with a glycol. The amount of glycol added should cover a temperature range 5°C lower than the operating temperature of the particular application. This will prevent the water/glycol from gelling (freezing) near the evaporating coil.

50% ethylene glycol is the maximum recommendation. Excess glycol deteriorates the temperature accuracy due to its high viscosity.

All heat transfer fluids are supplied with an EC Safety Data Sheet.

Ensure, when selecting the heat transfer fluid, that no toxic gases can be generated. Inflammable gases can build up over the fluid during usage. ▲

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Water Quality and Standards

Process Fluid	Permissible (PPM)	Desirable (PPM)
Microbiologicals		
(algae, bacteria, fungi)	0	0
Inorganic Chemicals		
Calcium	<25	<0.6
Chloride	<25	<10
Copper	<1.3	<1.0
0.020 ppm if fluid in contact with aluminum		
Iron	<0.3	<0.1
Lead	<0.015	0
Magnesium	<12	<0.1
Manganese	<0.05	<0.03
Nitrates\Nitrites	<10 as N	0
Potassium	<20	<0.3
Silicate	<25	<1.0
Sodium	<20	<0.3
Sulfate	<25	<1
Hardness	<17	<0.05
Total Dissolved Solids	<50	<10
Other Parameters		
pH	6.5-8.5	7-8
Resistivity	0.01*	0.05-0.1*

* MΩ-cm (compensated to 25°C)

Unfavorably high total ionized solids (TIS) can accelerate the rate of galvanic corrosion. These contaminants can function as electrolytes which increase the potential for galvanic cell corrosion and lead to localized corrosion such as pitting. Eventually, the pitting will become so extensive that refrigerant will leak into the water reservoir.

As an example, raw water in the United States averages 171 ppm (of NaCl). The recommended level for use in a water system is between 0.5 to 5.0 ppm (of NaCl).

Recommendation: Initially fill the tank with distilled or deionized water. Do not use untreated tap water as the total ionized solids level may be too high. This will reduce the electrolytic potential of the water and prevent or reduce the galvanic corrosion observed.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Filling Requirements

CAUTION


Before using any fluid refer to the manufacturer's MSDS and EC safety data sheets for handling precautions. ▲

NOTE Ensure the reservoir overflow located on the rear of the chiller is connected to a suitable drain. ▲

Start with 2.8 liters of fluid.

Lift up on the fill panel on the front of the chiller to remove it.

Insert the supplied funnel into the hose and *slowly* add fluid. **NOTE** Due to the chiller's plumbing configuration, adding fluid too quickly may result in overfilling - the fill level indication lags the actual amount added. ▲

When pouring fluid into the funnel ensure the fluid contacts the side of the funnel first, A. Pouring the fluid into the center of the funnel, B, will cause air pockets and substantially slow down the filling time.

Using the fill lines, verify the desired fluid level.

Since the reservoir capacity may be small compared to your application and air may need to be purged from the lines, have extra cooling fluid on hand to keep the system topped off when external circulation is started.

CAUTION

Avoid overfilling, fluids expand when heated. ▲

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Draining

CAUTION


Before draining any fluid refer to the manufacturer's MSDS and EC safety data sheets for handling precautions. ▲

Ensure the fluid is at a safe handling temperature (less than 55°C). Wear protective clothing and gloves. ▲

Lift up on the fill panel on the front of the chiller to remove it.

Pull the fill hose off the clips securing it to the chiller.

Keep your thumb on the end of the hose until ready to drain the fluid into a suitable container.

Shipment/Storage


Transporting and/or storing the chiller in near or below freezing temperatures requires draining. Do not store the chiller below -20°C, and <80% relative humidity.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Section 4 Operation


Controller

The Cole Parmer Polystat CR Series of chillers have a digital display and easy-to-use touch pad, five programmable setpoint temperatures, acoustic and optical alarms and adjustable high temperature protection.

Use this button to place the chiller in and out of standby.

Use these navigation arrows to move through the controller displays and to adjust values.

Use this button to start/stop the chiller. This button is also used to make and save changes on the controller's display screen.

Use this button to cancel any changes and to return the controller to its previous display. Canceling a change can only be made before the change is saved. In some cases, it is also used to save changes.

NOTE Holding this button for five seconds resets the display contrast to the default level and also brings up the language menu to change, if needed, the displayed language. See **Settings - Display Options** in this Section. ▲

°C

Use the adjacent dial for adjusting the High Temperature Cutout. Details are explained near the end of this Section.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Setup

CAUTION

Chillers should be left in an upright position at room temperature (~25°C) for 24 hours before starting. This will ensure the lubrication oil has drained back into the compressor. ▲

CAUTION


Before starting the chiller, double check all USB (optional), electrical and plumbing connections. ▲

Start Up

Do not run the chiller until fluid is added. Have extra fluid on hand. If the chiller will not start refer to Section 5 Troubleshooting.

I
O

- Place the circuit protector located on the front of the chiller to the **I** position.
- Press , the Start Display will appear.
- Ensure the start symbol has a highlight box around it, if not use the arrow keys to navigate to the symbol.

Start Display

- Press . The chiller will start and the start symbol will turn into a stop symbol (■).

NOTE It will take 30 seconds for the compressor to start. ▲


NOTE After start, check external plumbing connections for leaks. ▲

The **SP1** and **Menu** portions on the top of the display are used to view and/or change the controller's settings. They are explained in detail later in this Section.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Status Display

If desired, press Esc to toggle between the Start/Status Displays.

Stand By Mode

Press , the display will go blank and the chiller will be in the stand by mode.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com


Changing the Setpoint

NOTE You cannot adjust the setpoint closer than 0.1°C to either of the fluid's system limits, see Fluids Type in this Section, or beyond the chiller's temperature range. ▲

NOTE The setpoint can be changed with the chiller running or not. ▲


The Setpoint is the desired fluid temperature. The controller can store up to five setpoints, **SP1** through **SP5**. The procedure for changing the stored setpoint values is discussed later in this Section.

Use the navigation arrows and move to the **SP1** window and then press ◯ to highlight it as shown below.

Use the up and down navigation arrows to bring up the desired setpoint and then press ◯.

The display on the Setpoint Value Window will now indicate the corresponding setpoint's stored value.

If desired, you can change the displayed setpoint value by using the navigation arrows to highlight the Setpoint Value Window and then pressing ◯. The right-most digit will have a cursor beneath it.

Use the left and right arrows to move the cursor to the desired digit and then use the up and down arrows to change the value. Once all the desired changes are made, press ◯ to save the change.

NOTE Using this procedure also changes the setpoint's stored value. ▲

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com


Menu Displays

The controller uses menus to view/change the chiller's settings.


NOTE The chiller does not need to be running to view/change these settings. ▲

For all Menu displays, once ○ is pressed to change a display, you can press Esc to return to the previous screen.


1. Use the arrow buttons to highlight **Menu** and the controller brings up the Main Menu Display.


2. Use the up and down arrow to highlight the desired setting and then press ○ to bring up additional submenus.

See page 4-7.


See page 4-12.


Menu

The **Menu** line, at the bottom of all the submenu displays, is another way to return the controller back to the Main Menu Display.

1. From any submenu display, use the down arrow button to highlight **Menu**.

2. Press ○ to return to the Main Menu Display.


Pressing Esc from the Menu line returns you tho the previous screen.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Menu Tree


Since the controller can only display five lines of text at a time, keep pressing the down arrow to highlight additional options.


Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Settings - Application Settings is used to view/adjust the controller's five Setpoints and Real Temperature Adjustments (RTA) enable/disable the alarms, change the fluid type, set the pump speed, configure the interfaces (optional), set the clock, turn the timer on or off, and turn auto restart and energy savings on or off.


1. With **Application Settings** highlighted press ○ to view:

Scroll down for additional options:


2. With **Setpoints** highlighted, press ○ to display the list. Use the up/down arrows to highlight the desired **SP**. **NOTE** Use the down arrow to display **SP5**. ▲

3. Press ○.

The Setpoint and RTA are changed using the same procedure. With the desired setpoint highlighted press ↵ to display the submenu.

If the temperatures on the Start/Status Displays do not accurately reflect the actual temperature in the reservoir, an RTA can be applied. The RTA can be set ±10°C (±18°F).

As an example, if the chiller's temperature is stabilized and displaying 20°C but a calibrated reference thermometer reads 20.5°C, the RTA should be set to -0.5°C. After you enter a RTA value allow chiller to stabilize before verifying the temperature. **NOTE** If display accuracy is required, we recommend repeating this procedure at various setpoint temperatures and on a regular basis. ▲

NOTE You cannot adjust the setpoint closer than 0.1°C to either of the fluid's system limits, see Fluids Type in this Section. ▲

4. With the desired line highlighted press ○.

The right-most digit will have a cursor beneath it. Use the left and right arrows to move the cursor to the desired digit and then use the up and down arrows to change the value. Once all the desired changes are made, press ○ to save the change, or Esc to cancel it.


35.0°C

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com


Alarms is used to view/adjust the high and low temperature alarm limits and to enable/disable the audible alarms.

1. With **Alarms** highlighted, press to display:

.

2. With **Temperature Alarms** highlighted, press to display:

3. Highlight the desired limit and press . Follow the same procedure used to change a setpoint.


If the Fault temperature is exceeded the chiller will shut down and, if enabled, the audible alarm will sound. If the Warn temperature is exceeded the chiller will continue to run and, if enabled, the audible alarm will sound. In both cases a message will be displayed.

High Fault cannot be set below **High Warn**.
High Warn cannot be set below **Low Warn**.
Low Fault cannot be set above **High Warn**.

Press Esc, or use the **Menu** window, to save and return to the previous display.

1. With **Audible Alarms** highlighted, press to display the alarms.

Highlight the desired alarm and press to toggle between enable and disable mode.


If **Faults** is enabled the alarm will sound when a fault occurs. If **Warnings** is enabled the alarm will sound when a warning occurs.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Fluids Type is used to identify the type of fluid used. The controller uses the fluid type to automatically set certain operating parameters.

1. With **Fluid Type** highlighted, press ○ to display the list of acceptable fluids.

Highlight the desired fluid and press ○ to select it.


Fluid system limits:	High °C/°F	Low °C/°F
Water	+83/+181	+2/+36
Glycol-Water	+83/+181	-13/+9

2. With the desired fluid selected press Esc , or use the **Menu** window, to save and return to the previous display.

Pump Adjustment is used to review/set the desired pump speed. **NOTE** Force only pumps (FP1) have only two speeds, **High** and **Low**. ▲


1. With **Pump Adjustment** highlighted, press ○ to display the speeds.

Highlight the desired speed and press ○ to select it.


Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Set Clock is used to set the controller's time (hr : min : sec) and date (year - month - day).

On/Off Timer is used to enable and set the controller's timer.

1. With **On/Off Timer** highlighted, press to display the on (**I**) and off (**O**) time as well as the enable box.

After setting the on and off times select **Enable** to activate the timer.

Auto Restart is used to enable the auto restart feature. When enabled, the chiller will automatically restart after a power failure or power interruption condition.

1. With **Auto Restart** highlighted, press to toggle between enable and disable.


Energy Saving is used to enable the energy savings mode. The Energy Saving mode is primarily designed for applications running under a stable load. Enabling the mode saves energy by reducing the chiller's heater power and cooling requirements. This can result in substantial energy savings over the life of the chiller. The default setting is **ON**.

1. With **Energy Saving** highlighted, press to toggle between enable and disable.

Interfaces is used to enable/configure the optional serial communications feature.

1. With **Interfaces** highlighted, press to display the list of parameters.

Highlight the desired parameter and press to view the available options.

Available options:

Serial Type	USB, Off
Baud	19200, 9600, 4800, 2400, 1200, 600, 300
Parity	None, Odd, Even
Data Bits	8
Stop Bits	1, 2

See the Appendix for additional information.


.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Settings - Display Options is used to view/adjust the controller's Temperature Units, the Temperature Resolution, the displayed Language, and the Display Contrast and Display Delay.


1. With **Temp. Unit** highlighted press . Use the up/down arrows to highlight the desired temperature scale.

Press .


2. With **Temp. Resolution** highlighted press .

Use the up/down arrows to highlight the desired resolution. Press .


3. With **Language** highlighted press . Use the up/down arrows to highlight the desired language. (Scroll down for additional languages.)

Press .

4. With **Display Contrast** highlighted press .

Press again and use the up/down arrows keys to change the contrast. With the desired contrast showing, press again.


NOTE Holding Esc for five seconds resets the display contrast to the default level and also brings up the language menu to change, if needed, the displayed language. ▲

5. With **Display Delay** highlighted press to enable/disable it.

Use the up/down arrows to highlight the time and press again.

Use the up/down arrows to change the value. Once the desired delay is displayed press .


With **Display Delay** enabled and the Start Display showing, if no arrows are pressed the Start Display will change to the Status Display after the delay expires. See pages 4-2 and 4-3.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com


System Messages is used to view any Warning or Fault messages.

1. With **Messages** highlighted, press ○ to display the options.


System Run Time is used to view the chiller (**Unit**) and pump operating hours.

1. With **Run Time** highlighted, press ○ to display the times.

System Configuration is used to view the chiller's configuration.

1. With **Configuration** highlighted, press ○ to display the settings.

Pump options:

FP1 Force Pump

FP2 Force/Suction Pump


FP1U Force Pump & USB Option

FP2U Force/Suction Pump & USB Option


Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

System - Password/Reset is used only by a qualified technician. Changing the password enables controller reset options, the temperature sensor calibration procedure and displays PID values.


1. With **Password/Reset** highlighted, press to display:

2. Press and change the number to **1**.


3. Press to display:

Scroll down to display **PID Tuning**.

NOTE The controller resets to the **User** mode whenever the chiller is turned off. The controller also resets to the **User** mode whenever the Start/Status Display is displayed continuously for 10 minutes. ▲

1. If desired, highlight **Reset** and press to display:


Highlight the desired reset option and press .


A confirmation message will appear, press again.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com


1. To calibrate the temperature sensor highlight **Calibration** and press ◯ to display:


2. Press ◯ to display:


3. Highlight the **SP** temperature box and enter either the desired high or low setpoint value.


4. Press ◯ and enter the desired calibration temperature. Allow the fluid temperature to stabilize. Once stable, enter the temperature read off a reference thermometer.

5. Highlight **Calibrate** and the press ◯ to complete the procedure.

6. Repeat for the other calibration point.

Once the calibration is complete you can store it into the controller's memory by selecting **Save User Cal** and pressing ◯ .


You can later restore the same calibration by highlighting **Restore User Cal** and pressing ◯ .

Another option is to restore the factory calibration values by highlighting **Restore Factory Cal** and pressing ◯ .


The **Save Factory Cal** option is designed to be used only by a qualified technician.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com


1. With **PID Tuning** highlighted, press  to display:

2. Highlight the desired PID and press  to display:


3 . If required, press  to change the value.

Factory values are:
P = 05.0
I = 0.10
D = 0.00

High Temperature Cutout

To protect your application, the adjustable High Temperature Cutout (HTC) ensures the heater will not exceed temperatures that can cause serious damage. A temperature sensor is located in the reservoir. A High Temperature fault occurs when the temperature of the sensor exceeds the set temperature limit.

In the event of a fault the chiller will shut down and the controller will display a fault message, see Section 6. The cause of the fault must be identified and corrected before the chiller can be restarted.

The cutout is not preset and must be adjusted. To set the cutout, with the chiller off use a flat head screwdriver to turn the red dial fully clockwise. Start the chiller. Adjust the setpoint for a few degrees higher than the highest desired fluid temperature and allow the chiller to stabilize at the temperature setpoint. Then slowly turn the dial counterclockwise until the chiller shuts down and the message appears. Press the enter key to clear the message.

Before you can restart the chiller it has to cool down a few degrees. To restart press the black ring surrounding the red dial - and then press the enter key again. If Auto Restart is enabled the chiller will restart, if disabled use the Start Up procedure.


NOTE: We recommend periodically rechecking operation or if the chiller is moved. ▲

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Stopping the Chiller

Ensure the stop symbol is highlighted, if not use the arrow keys to navigate to the symbol.

Press  . The chiller will stop and the stop symbol will turn into a start symbol (▶).

Power Down

Press  . The controller's screen will go blank .

Shut Down

Place the circuit protector on the front of the chiller to the **O** position.

⚠ CAUTION

Using any other means to shut the chiller down can reduce the life of the compressor. ▲

⚠ CAUTION

Always turn the chiller off and disconnect it from its supply voltage before moving. ▲

⚠ CAUTION

The circuit protector located on the front of the chiller is not intended to act as a disconnecting means. ▲

Restarting


NOTE When quickly restarting, the compressor may take up to 10 minutes before it starts to operate. Limit chiller cycling to 8 times/hour. ▲

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Section 4 Operation


Controller

The Cole Parmer Polystat CR Series of chillers have a digital display and easy-to-use touch pad, five programmable setpoint temperatures, acoustic and optical alarms and adjustable high temperature protection.

Use this button to place the chiller in and out of standby.

Use these navigation arrows to move through the controller displays and to adjust values.

Use this button to start/stop the chiller. This button is also used to make and save changes on the controller's display screen.

Use this button to cancel any changes and to return the controller to its previous display. Canceling a change can only be made before the change is saved. In some cases, it is also used to save changes.

NOTE Holding this button for five seconds resets the display contrast to the default level and also brings up the language menu to change, if needed, the displayed language. See **Settings - Display Options** in this Section. ▲

°C

Use the adjacent dial for adjusting the High Temperature Cutout. Details are explained near the end of this Section.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Setup

CAUTION

Chillers should be left in an upright position at room temperature (~25°C) for 24 hours before starting. This will ensure the lubrication oil has drained back into the compressor. ▲

CAUTION


Before starting the chiller, double check all USB (optional), electrical and plumbing connections. ▲

Start Up

Do not run the chiller until fluid is added. Have extra fluid on hand. If the chiller will not start refer to Section 5 Troubleshooting.

I
O

- Place the circuit protector located on the front of the chiller to the **I** position.
- Press , the Start Display will appear.
- Ensure the start symbol has a highlight box around it, if not use the arrow keys to navigate to the symbol.

Start Display

- Press . The chiller will start and the start symbol will turn into a stop symbol (■).

NOTE It will take 30 seconds for the compressor to start. ▲


NOTE After start, check external plumbing connections for leaks. ▲

The **SP1** and **Menu** portions on the top of the display are used to view and/or change the controller's settings. They are explained in detail later in this Section.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Status Display

If desired, press Esc to toggle between the Start/Status Displays.

Stand By Mode

Press , the display will go blank and the chiller will be in the stand by mode.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com


Changing the Setpoint

NOTE You cannot adjust the setpoint closer than 0.1°C to either of the fluid's system limits, see Fluids Type in this Section, or beyond the chiller's temperature range. ▲

NOTE The setpoint can be changed with the chiller running or not. ▲


The Setpoint is the desired fluid temperature. The controller can store up to five setpoints, **SP1** through **SP5**. The procedure for changing the stored setpoint values is discussed later in this Section.

Use the navigation arrows and move to the **SP1** window and then press ◯ to highlight it as shown below.

Use the up and down navigation arrows to bring up the desired setpoint and then press ◯.

The display on the Setpoint Value Window will now indicate the corresponding setpoint's stored value.

If desired, you can change the displayed setpoint value by using the navigation arrows to highlight the Setpoint Value Window and then pressing ◯. The right-most digit will have a cursor beneath it.

Use the left and right arrows to move the cursor to the desired digit and then use the up and down arrows to change the value. Once all the desired changes are made, press ◯ to save the change.

NOTE Using this procedure also changes the setpoint's stored value. ▲

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com


Menu Displays


The controller uses menus to view/change the chiller's settings.

NOTE The chiller does not need to be running to view/change these settings. ▲


For all Menu displays, once ◯ is pressed to change a display, you can press Esc to return to the previous screen.

1. Use the arrow buttons to highlight **Menu** and the controller brings up the Main Menu Display.


2. Use the up and down arrow to highlight the desired setting and then press ◯ to bring up additional submenus.

See page 4-7.


See page 4-12.

Menu


The **Menu** line, at the bottom of all the submenu displays, is another way to return the controller back to the Main Menu Display.

1. From any submenu display, use the down arrow button to highlight **Menu**.

Pressing Esc from the Menu line returns you tho the previous screen.


2. Press ◯ to return to the Main Menu Display.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Menu Tree


Since the controller can only display five lines of text at a time, keep pressing the down arrow to highlight additional options.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Settings - Application Settings is used to view/adjust the controller's five Setpoints and Real Temperature Adjustments (RTA) enable/disable the alarms, change the fluid type, set the pump speed, configure the interfaces (optional), set the clock, turn the timer on or off, and turn auto restart and energy savings on or off.


1. With **Application Settings** highlighted press ◯ to view:


Scroll down for additional options:

2. With **Setpoints** highlighted, press ◯ to display the list. Use the up/down arrows to highlight the desired **SP**. **NOTE** Use the down arrow to display **SP5**. ▲

3. Press ◯.

If this temperature on the Start/Status Displays does not accurately reflect the actual temperature in the reservoir, an RTA can be applied. The RTA can be set ±10°C (±18°F).

As an example, if the chiller's temperature is stabilized and displaying 20°C but a calibrated reference thermometer reads 20.5°C, the RTA should be set to -0.5°C. After you enter a RTA value allow chiller to stabilize before verifying the temperature. **NOTE** If display accuracy is required, we recommend repeating this procedure at various setpoint temperatures and on a regular basis. ▲

NOTE You cannot adjust the setpoint closer than 0.1°C to either of the fluid's system limits, see Fluids Type in this Section. ▲

4. With the desired line highlighted press ◯.

The right-most digit will have a cursor beneath it. Use the left and right arrows to move the cursor to the desired digit and then use the up and down arrows to change the value. Once all the desired changes are made, press ◯ to save the change, or Esc to cancel it.


Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com


Alarms is used to view/adjust the high and low temperature alarm limits and to enable/disable the audible alarms.

1. With **Alarms** highlighted, press to display:

.

2. With **Temperature Alarms** highlighted, press to display:

3. Highlight the desired limit and press . Follow the same procedure used to change a setpoint.


If the Fault temperature is exceeded the chiller will shut down and, if enabled, the audible alarm will sound. If the Warn temperature is exceeded the chiller will continue to run and, if enabled, the audible alarm will sound. In both cases a message will be displayed.

High Fault cannot be set below **High Warn**.
High Warn cannot be set below **Low Warn**.
Low Fault cannot be set above **High Warn**.

Press Esc, or use the **Menu** window, to save and return to the previous display.

1. With **Audible Alarms** highlighted, press to display the alarms.

Highlight the desired alarm and press to toggle between enable and disable mode.


If **Faults** is enabled the alarm will sound when a fault occurs. If **Warnings** is enabled the alarm will sound when a warning occurs.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Fluids Type is used to identify the type of fluid used. The controller uses the fluid type to automatically set certain operating parameters.

1. With **Fluid Type** highlighted, press ○ to display the list of acceptable fluids.

Highlight the desired fluid and press ○ to select it.


Fluid system limits:	High °C/°F	Low °C/°F
Water	+83/+181	+2/+36
Glycol-Water	+83/+181	-13/+9

2. With the desired fluid selected press Esc , or use the **Menu** window, to save and return to the previous display.

Pump Adjustment is used to review/set the desired pump speed. **NOTE** Force only pumps (FP1) have only two speeds, **High** and **Low**. ▲


1. With **Pump Adjustment** highlighted, press ○ to display the speeds.

Highlight the desired speed and press ○ to select it.


Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Set Clock is used to set the controller's time (hr : min : sec) and date (year - month - day).

On/Off Timer is used to enable and set the controller's timer.

1. With **On/Off Timer** highlighted, press to display the on (**I**) and off (**O**) time as well as the enable box.

After setting the on and off times select **Enable** to activate the timer.

Auto Restart is used to enable the auto restart feature. When enabled, the chiller will automatically restart after a power failure or power interruption condition.


1. With **Auto Restart** highlighted, press to toggle between enable and disable.

Energy Saving is used to enable the energy savings mode. The Energy Saving mode is primarily designed for applications running under a stable load. Enabling the mode saves energy by reducing the chiller's heater power and cooling requirements. This can result in substantial energy savings over the life of the chiller. The default setting is **ON**.

1. With **Energy Saving** highlighted, press to toggle between enable and disable.

Interfaces is used to enable/configure the optional serial communications feature.

1. With **Interfaces** highlighted, press to display the list of parameters.
Highlight the desired parameter and press to view the available options.

Available options:

Serial Type	USB, Off
Baud	19200, 9600, 4800, 2400, 1200, 600, 300
Parity	None, Odd, Even
Data Bits	8
Stop Bits	1, 2

See the Appendix for additional information.


.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Settings - Display Options is used to view/adjust the controller's Temperature Units, the Temperature Resolution, the displayed Language, and the Display Contrast and Display Delay.


1. With **Temp. Unit** highlighted press . Use the up/down arrows to highlight the desired temperature scale.

Press .


2. With **Temp. Resolution** highlighted press .

Use the up/down arrows to highlight the desired resolution. Press .


3. With **Language** highlighted press . Use the up/down arrows to highlight the desired language. (Scroll down for additional languages.)

Press .

4. With **Display Contrast** highlighted press .

Press again and use the up/down arrows keys to change the contrast. With the desired contrast showing, press again.


NOTE Holding Esc for five seconds resets the display contrast to the default level and also brings up the language menu to change, if needed, the displayed language. ▲

5. With **Display Delay** highlighted press to enable/disable it.

Use the up/down arrows to highlight the time and press again.

Use the up/down arrows to change the value. Once the desired delay is displayed press .


With **Display Delay** enabled and the Start Display showing, if no arrows are pressed the Start Display will change to the Status Display after the delay expires.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com


System Messages is used to view any Warning or Fault messages.

1. With **Messages** highlighted, press ○ to display the options.


System Run Time is used to view the chiller (**Unit**) and pump operating hours.

1. With **Run Time** highlighted, press ○ to display the times.

System Configuration is used to view the chiller's configuration.

1. With **Configuration** highlighted, press ○ to display the settings.

Pump options:

FP1 Force Pump

FP2 Force/Suction Pump


FP1U Force Pump & USB Option

FP2U Force/Suction Pump & USB Option


Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

System - Password/Reset is used only by a qualified technician. Changing the password enables controller reset options, the temperature sensor calibration procedure and displays PID values.


1. With **Password/Reset** highlighted, press to display:

2. Press and change the number to **1**.


3. Press to display:

Scroll down to display **PID Tuning**.

NOTE The controller resets to the **User** mode whenever the chiller is turned off. The controller also resets to the **User** mode whenever the Start/Status Display is displayed continuously for 10 minutes. ▲

1. If desired, highlight **Reset** and press to display:


Highlight the desired reset option and press .


A confirmation message will appear, press again.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com


1. To calibrate the temperature sensor highlight **Calibration** and press ◯ to display:


2. Press ◯ to display:


3. Highlight the **SP** temperature box and enter either the desired high or low setpoint value.


4. Press ◯ and enter the desired calibration temperature .

5. Highlight **Calibrate** and the press ◯ to complete the procedure.

6. Repeat for the other calibration point.

Once the calibration is complete you can store it into the controller's memory by selecting **Save User Cal** and pressing ◯ .


You can later restore the same calibration by highlighting **Restore User Cal** and pressing ◯ .

Another option is to restore the factory calibration values by highlighting **Restore Factory Cal** and pressing ◯ .


The **Save Factory Cal** option is designed to be used only by a qualified technician.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com


1. With **PID Tuning** highlighted, press  to display:

2. Highlight the desired PID and press  to display:


3 . If required, press  to change the value.

Factory values are:
P = 05.0
I = 0.10
D = 0.00

High Temperature Cutout

To protect your application, the adjustable High Temperature Cutout (HTC) ensures the heater will not exceed temperatures that can cause serious damage. A temperature sensor is located in the reservoir. A High Temperature fault occurs when the temperature of the sensor exceeds the set temperature limit.

In the event of a fault the chiller will shut down and the controller will display a fault message, see Section 6. The cause of the fault must be identified and corrected before the chiller can be restarted.

The cutout is not preset and must be adjusted. To set the cutout, with the chiller off use a flat head screwdriver to turn the red dial fully clockwise. Start the chiller. Adjust the setpoint for a few degrees higher than the highest desired fluid temperature and allow the chiller to stabilize at the temperature setpoint. Then slowly turn the dial counterclockwise until the chiller shuts down and the message appears. Press the enter key to clear the message.

Before you can restart the chiller it has to cool down a few degrees. To restart press the black ring surrounding the red dial - and then press the enter key again. If Auto Restart is enabled the chiller will restart, if disabled use the Start Up procedure.


NOTE: We recommend periodically rechecking operation or if the chiller is moved. ▲

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Stopping the Chiller

Ensure the stop symbol is highlighted, if not use the arrow keys to navigate to the symbol.

Press  . The chiller will stop and the stop symbol will turn into a start symbol (▶).

Power Down

Press  . The controller's screen will go blank .

Shut Down

Place the circuit protector on the front of the chiller to the **O** position.

⚠ CAUTION

Using any other means to shut the chiller down can reduce the life of the compressor. ▲

⚠ CAUTION

Always turn the chiller off and disconnect it from its supply voltage before moving. ▲

⚠ CAUTION

The circuit protector located on the front of the chiller is not intended to act as a disconnecting means. ▲

Restarting

NOTE When quickly restarting, the compressor may take up to 10 minutes before it starts to operate. Limit chiller cycling to 8 times/hour. ▲

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Section 5 Preventive Maintenance

⚠ WARNING **Laboratory Grade Ethylene glycol (EG) is poisonous and flammable. Before performing any preventive maintenance refer to the manufacturer's most current MSDS for handling precautions. ▲**

⚠ CAUTION **Disconnect the power cord prior to performing any maintenance. Ensure chiller is at a safe temperature before handling. ▲**

⚠ CAUTION **Handle the chiller with care. Sudden jolts or drops can damage its components. ▲**

There are no user serviceable components within the equipment cabinet. Only Cole Parmer should provide any required replacement parts listed below.

Cleaning

Clean the chiller's surface with a soft cloth and warm water only.

Condenser

Clean the condenser by running a vacuum over the grid on the lower rear section of the chiller.

⚠ CAUTION **Failure to clean the condenser will cause a loss of cooling capacity and lead to premature failure of the cooling system. ▲**

Clean/replace the optional condenser filter.

Fluid Maintenance

Check the fluid level, concentration and pH on a regular basis, see Section 3. Change the fluid if it is discolored.

Sight Tube

As required, clean the sight tube using a small soft brush.

Hoses

Inspect and tighten the chiller's external hoses and clamps daily.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Testing the Safety Features

°C

High temperature protection
Using a flat-head screw driver, set a cut-off temperature that is lower than the desired setpoint temperature.

Switch on the circulator and ensure the chiller shuts down at the set cut-off temperature. Press the black ring surrounding the dial to reset the HTC.

If the chiller does not shut down have it checked by a qualified technician.

Reset the safety to the desired temperature.

Decommissioning/ Disposal

Decommissioning prepares equipment for safe and secure transportation.

Laboratory Grade Ethylene glycol (EG) is poisonous and flammable. Before disposing refer to the manufacturer's most current MSDS for handling precautions. ▲

Decommissioning must be performed only by qualified dealer using certified equipment. All prevailing regulations must be followed. ▲

Consider decommissioning the chiller when:

- It fails to maintain desired specifications
- It no longer meets safety standards
- It is beyond repair for its age and worth

Refrigerant (R134A) and oil (Polyol Ester) must be recovered from equipment before disposal.

NOTE Keep in mind any impact your application may have had on the chiller. ▲

Direct questions about chiller decommissioning or disposal to our Sales, Service and Customer Support.

Handling and disposal should be done in accordance with the manufacturers specification and/or the MSDS for the material used. ▲

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Section 6 Troubleshooting

Error Displays

Error messages are cleared by pressing the enter key (). Once the cause of the error message is identified and corrected, to restart the chiller - on chillers with an HTC error press the black ring surrounding the red dial - and then for all chillers press the enter key again. If **Auto start** is enabled the chiller will restart, if disabled use the Start Up procedure.

FAULT:
HTC
HIGH TEMP. CUTOUT-
PRESS ENTER

Error Message (Typical)

Fault Displays

The heating element, pump and, if applicable, refrigeration shut down with a fault. A fault also sounds the alarm, if enabled.

Message	Action
High Fixed Temp.	• chiller's nonadjustable high temperature protection limit exceeded • check fluid selection • check environmental conditions
High Temperature	• adjustable high temperature protection limit exceeded • check limit setting • check fluid selection • ensure chiller has adequate ventilation
High Temperature Refrigeration	• check voltage supply • the refrigeration may need servicing
HPC High Press. Cutout	• check for obstructions to air flow • the refrigeration may need servicing
HTC High Temp. Cutout	• high temperature protection limit exceeded • check limit setting • allow chiller to cool down • reset HTC by pressing black ring
Low Fixed Temp.	• chiller's nonadjustable low temperature protection limit exceeded • check fluid selection

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Low Temperature	• adjustable high temperature protection limit exceeded • check limit setting • check fluid selection
Motor Fault	• it can take over 10 minutes for the motor temperature to get low enough before the chiller can be restarted
MOL Motor Overload	• allow chiller to cool down
Open RTD1 Internal	• open internal temperature sensor
Shorted RTD1 Internal	• shorted internal temperature sensor

Warning Displays

The chiller will continue to run with a warning. A warning also sounds the alarm, if enabled.

Bad Calibration	• redo calibration
High Temperature	• adjustable high temperature protection limit exceeded • check limit setting • check fluid selection
Low Temperature	• adjustable low temperature protection limit exceeded • check limit setting • check fluid selection

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Checklist

Chiller will not start

Check the controller for error displays, see Error Displays in this section.

Ensure the circuit protector is in the on (I) position.

Ensure supply voltage is connected and matches the chiller's nameplate rating ±10%.

No display on controller

Recycle the circuit protector on the front of the chiller.

Chiller will not circulate process fluid

Check the reservoir level. Fill, if necessary.

Check the application for restrictions in the cooling lines.

The pump motor overloaded. The pump's internal overtemperature overcurrent device will shut off the pump causing the flow to stop. This can be caused by low fluid, debris in system, operating chiller in a high ambient temperature condition or excessively confined space. Allow time for the motor to cool down.

Make sure supply voltage matches the chiller's nameplate rating ±10%.

Inadequate temperature control

Verify the setpoint.

Make sure the condenser is free of dust and debris.

Check the fluid concentration.

Ensure chiller installation complies with the site requirements in Section 3.

Make sure supply voltage matches chiller nameplate rating ±10%.

If the temperature continues to rise, make sure your application's heat load does not exceed the rated specifications.

Check for high thermal gradients (e.g., the application load is being turned on and off or rapidly changing).

Chiller shuts down

Ensure  button wasn't accidently pressed.

Ensure the circuit protector is in the on (I) position.

Check the controller for error codes.

Ensure supply voltage is connected and matches the chiller's nameplate rating ±10%.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

USB Driver Not Recognized

If your operating system does not automatically recognize the optional driver log on to:

http://www.ftdichip.com/FTDrivers.htm

for instructions.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Appendix AC Serial Communications Protocol

Serial communication is accomplished through the USB port on the chiller. If your operating system does not automatically recognize the optional driver log on to: http://www.ftdichip.com/FTDrivers.htm for instructions.

Note This appendix assumes you have a basic understanding of communications protocols. Information on the NC, STANDARD and NAMUR protocols is available upon request. ▲

Note Keypad operation is still available with serial communications enabled. ▲

All commands must be entered in the exact format shown in the tables on the following pages. The tables show all commands available, their format and responses. Chiller responses are either the requested data or an error message. The chiller response must be received before the host sends the next command.

The host sends a command embedded in a single communications packet, then waits for the chiller's response. If the command is not understood, the chiller responds with an error command. Otherwise, the chiller responds with the requested data.

Commands are not case sensitive. Upper or lower case letters may be used. Commands are listed in the Commands Table, error responses are given in the Errors Table, and symbols are shown in the Key Table.

Key	
Symbol	Meaming
[B]	A binary value 0 or 1 (0 = Off, FALSE or Disable(d); 1 = On, TRUE or Enable(d)).
[CR]	Carriage return – used as the termination character.
[U]	Text representing the units associated with a value.
[V]	A value that can be requested in a read command or sent as part of a set command.
[V_{MAX}]	Maximum allowed value. Part of error message when set value is too high.
[V_{MIN}]	Minimum allowed value. Part of error message when set value is too low.

Value: Read commands return analog [V] or bit [B] values or settings, while set commands send analog or bit settings. Read commands return values with the same displayed precision. Set command messages missing the space character between the command and the setting will be rejected, as the user's intent is unclear.

Units: A read command returning an analog [V] value or setting, will include the units [U] associated with that value or setting. A set command sending an analog value will not include the units. The units returned by the complementary read command are assumed.

Termination character: A carriage return [CR] is used to terminate command and response messages. (Typically the "Enter" key on the keyboard.)

Note The inter-character timeout (time between transmitted characters) is set to 30 seconds. Exceeding the timeout will clear the receiver buffer and require the message to re retransmitted.

Note Special characters (backspace, delete, insert, etc.) are not recognized by the protocol and generate error responses.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Commands Table:

Commands *All messages from master and slave are terminated with a carriage return [CR]*		
Command Description *Notes*	Master Sends	Sample Slave Response *(echo off)* *Alternate units*
Read Temperature *Internal*	RT	*[V]*C *F K*
Read Displayed Setpoint	RS	*[V]*C *F K*
Read Setpoint X (X = 1 to 5)	RSX	*[V]*C *F K*
Read High Temperature Fault	RHTF	*[V]*C *F K*
Read High Temperature Warn	RHTW	*[V]*C *F K*
Read Low Temperature Fault	RLTF	*[V]*C *F K*
Read Low Temperature Warn	RLTW	*[V]*C *F K*
Read Proportional Cool Band Setting	RPC	*[V]*%
Read Integral Cool Band Setting	RIC	*[V]*Repeats per minute
Read Derivative Cool Band Setting	RDC	*[V]*Minutes
Read Temperature Precision	RTP	*[V]*
Read Temperature Units	RTU	*[V]* *C,F,K*
Read Unit On	RO	*[B]*
Read Auto Restart Enabled	RAR	*[B]*
Read Energy Saving Mode	REN	*[B]*
Read Time	RCK	*hh:mm:ss*
Read Date	RDT	*mm/dd/yyyy or dd/mm/yyyy*
Read Date Format	RDF	*mm/dd/yyyy or dd/mm/yyyy*
Read Firmware Version	RVER	*[V]*
Read Firmware Checksum	RSUM	*[V]*
Read Unit Fault Status	RUFS	*[V_1, V_2 , V_3, V_4, V_5] See page 5*

Commands	*All messages from master and slave are terminated with a carriage return [CR]*	
Command Description *Notes*	Master Sends	Sample Slave Response *(echo off)* *Alternate units*
Set Displayed Setpoint	SS *[V]*	OK
Set Setpoint X (X = 1 to 5)	SSX *[V]*	OK
Set High Temperature Fault	SHTF *[V]*	OK
Set High Temperature Warning	SHTW *[V]*	OK
Set Low Temperature Fault	SLTF *[V]*	OK
Set Low Temperature Warning	SLTW *[V]*	OK
Set Proportional Cool Band Setting	SPC *[V]*	OK
Set Integral Cool Band Setting	SIC *[V]*	OK
Set Derivative Cool Band Setting	SDC *[V]*	OK
Set Temperature Resolution	STR *[V]*	OK
Set Temperature Units	STU *[V] C,F,K*	OK
Set Unit On Status	SO *[B]*	OK
Set Auto Restart Enabled	SAR *[B]*	OK
Set Energy Saving Mode	SEN *[V]*	OK
Set Pump Speed	SPS *[V] L,M,H*	OK

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Errors Table:

Errors	
Error Description *Notes*	Slave Responds
Not defined, not implemented or incorrectly formatted	? Unsupported command
Extra characters…	? Format error
Set value too high	? Maximum allowed is *[V_{MAX}]*
Set value too low	? Minimum allowed is *[V_{MIN}]*
Argument to binary set command not 0 or 1	? Value must be 0 or 1
Set command attempted while in read only mode	? Mode is read only
Set command failed (e.g. SO 1 with low level)	? Failed

RUFS Read Unit Fault Status

This command returns 5 values. These are decimal representations of hexadecimal values. Each individual bit of the value represents a different warning, fault or status.

decimal	hex	B7	B6	B5	B4	B3	B2	B1	B0
1	1	0	0	0	0	0	0	0	1
2	2	0	0	0	0	0	0	1	0
4	4	0	0	0	0	0	1	0	0
8	8	0	0	0	0	1	0	0	0
16	10	0	0	0	1	0	0	0	0
32	20	0	0	1	0	0	0	0	0
64	40	0	1	0	0	0	0	0	0
128	80	1	0	0	0	0	0	0	0

Value	Description of bits
V_1	B0 – b5 unused B6 rtd1 shorted B7 rtd1 open
V_2	B0 HTC fault B1 high RA temperature fault B2 – b7 unused
V_3	B0 low level warn B1 lo temperature warn B2 high temperature warn B3 low level fault B4 lo temperature fault B5 high temperature fault B6 low temperature fixed fault B7 high temperature fixed fault
V_4	B0 PWM heat duty cycle > 0 B1 compressor On/Off B2 Pump On status B3 Unit On status B4 Unit Stopping B5 Unit fault status B6 unused B7 Beeper On status
V_5	B0 Pump speed fault B1 MOL fault B2 HPC fault B3 Cool Icon On steady (unit is cooling at max capacity) B4 Cool Icon flashing (unit is B5 Heat Icon On steady B6 Heat Icon flashing B7 External sensor controlling

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Refer to Key table on page 1 for explanation of symbols and their meanings.

Examples:

Read Temperature:

Host

R	T		**CR**
Command			[CR]

Controller:

2	0	.	0	C	**CR**
[V]				[U]	[CR]

Set Setpoint:

Host

S	S		2	0	**CR**
Command			[V]		[CR]

Controller:

O	K	**CR**
Command Accepted		[CR]

Read Temperature 2:

Host:	R	T	2	**CR**		
Controller:	2	0	.	0	C	[CR]

Set Setpoint to -22°C when minimum allowed is -20°C: Minimum allowed is *[V_{MIN}]*

Host:	S	S		-	2	2	**CR**																		
Controller:	?		M	i	n	i	m	u	m		a	l	l	o	w	e	d		i	s		-	2	0	**CR**

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

DECLARATION OF CONFORMITY

Manufacturer: Cole Parmer

Address: 625 East Bunker Court
Vernon Hills, IL 60061

CE

Year of inception 2013

Products: Heated and Refrigerated Liquid Circulators.

We declare that the following products conform to the Directives and Standards listed below:

Cole Parmer Polystat models CR250 & CR500.

All rated :

1. 100 Volts, 50 Hz & 100 Volts, 60 Hz or
2. 115 Volts, 60 Hz or
3. 220 Volts, 60 Hz or
4. 220 - 230 Volts, 50 Hz

Equipment Class: Measurement, control and laboratory

Directives and Standards:

2004/108/EC – Electromagnetic Compatibility (EMC Directive):

EN 61326-1: 2006 – Electrical equipment for measurement, control, and laboratory use – EMC requirements, EMC Class A

2006/95/EC – Low Voltage Directive (LVD):

EN 61010-1: 2010 – Safety requirements for electrical equipment for measurement, control, and laboratory use – Part 1: General requirements.

EN 61010-1: 2004 – Safety requirements for electrical equipment for measurement, control, and laboratory use – Part 1: General requirements.

EN 61010-2-010: 2003 – Safety requirements for electrical equipment for measurement, control, and laboratory use – Part 2-010: Particular requirements for laboratory equipment for the heating of materials.

2011/65/EU - Restriction of the Use of Certain Hazardous Substances In Electrical and Electronic Equipment (ROHSD).

EN 50581: 2012 - Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances.

2012/19/EU - Waste from Electrical and Electronic Equipment (WEEED).

Manufacturer's Authorized Representative: Date:

Robin Wiley

10 July 2014

Robin Wiley Compliance Engineering

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Artisan Technology Group is an independent supplier of quality pre-owned equipment

Gold-standard solutions

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

Learn more!

Visit us at **artisantg.com** for more info on price quotes, drivers, technical specifications, manuals, and documentation.

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

We're here to make your life easier. How can we help you today?

(217) 352-9330 | sales@artisantg.com | artisantg.com

ARTISAN® TECHNOLOGY GROUP