ARTISAN®
TECHNOLOGY GROUP

Your **definitive** source for quality pre-owned equipment.

Artisan Technology Group
(217) 352-9330 | sales@artisantg.com | artisantg.com

Full-service, independent repair center
with experienced engineers and technicians on staff.

We buy your excess, underutilized, and idle equipment
along with credit for buybacks and trade-ins.

Custom engineering
so your equipment works exactly as you specify.

- Critical and expedited services
- In stock / Ready-to-ship
- Leasing / Rentals / Demos
- ITAR-certified secure asset solutions

Expert team | Trust guarantee | 100% satisfaction

All trademarks, brand names, and brands appearing herein are the property of their respective owners.

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SAFETY PRECAUTIONS

DANGER: High voltages exist and are accessible in the Console Drive. Use extreme caution when servicing internal components.

WARNING: PRODUCT USE LIMITATION

This product is not designed for, nor intended for use in patient connected applications; including, but not limited to, medical and dental use, and accordingly has not been submitted for FDA approval.

Variety of Pump Heads Accepted

Mount up to 1 (600 rpm) or 2 (200 rpm) MASTERFLEX® L/S® pump heads and all MASTERFLEX-compatible pump heads.

Setup and Drive Operation

1. Mount pump head and load tubing. (See pump head manual.)
2. Adjust flow rate with the 1-turn potentiometer speed control.

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
SAFETY PRECAUTIONS
SETUP AND DRIVE OPERATION

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MOTOR BRUSH CHECK/REPLACEMENT

NOTE: Brushes should be checked every 6 months or 2000 operating hours if erratic operation occurs.

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1. Place the POWER SWITCH in the off position.
2. Disconnect the AC power input line cord from the AC receptacle.
3. Remove the screws from each side of the housing and lift off the housing.
4. Carefully unscrew each brush cap. Withdraw the brush, and examine it for wear.

NOTE: Replace both brushes, if either brush is less than 7.6 mm (0.300 in) long from base to point.

5. Insert brushes and install brush cap.
6. Install housing and secure with the screws on each side.
7. Connect the AC power input line cord to the AC receptacle.

Fuses

NOTE: These drives have current limiting which cause drive to stop under overload conditions. If drive does not start return to dealer for servicing.

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MOTOR BRUSH CHECK/REPLACEMENT CLEANING

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
Power Output: 37 W (0.05 hp)
Maximum Current:
- **115 V Units:** 1.5A, shorted output conditions
- **230 V Units:** 0.9A, shorted output conditions

Speed Regulation:
- **Line:** ±2%
- **Load:** ±3%
- **Temperature Drift:** ±10%

Enclosure Rating: IP22 per IEC529
Humidity (non-condensing): 10% to 90%
Operating Temperature: 0˚ to 40˚C
Storage Temperature: –45˚ to 60˚C
Chemical Resistance: Exposed material is painted steel housing
Altitude: Less than 2000 m
Display: Green LED
Dimensions (L × W × H): 220 mm × 180 mm × 135 mm (8 11/16 in × 7 1/8 in × 5 5/16 in)
Weight: 4.1 kg (9 pounds)
Compliance: 115V: UL508, CSA C22.2, No. 14-M91
230V (For CE Mark):
EN61010-1/A2: 1995 (EU Low Voltage Directive) and
EN61326-1/A1: 1998 (EU EMC Directive)
Pollution Degree: Pollution Degree 2 per IEC 664 (Indoor usage—lab, office)
Installation Category: Installation Category II per IEC 664 (Local level—appliances, Portable Equipment, etc.)

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EU Declaration of Conformity

Name of Apparatus: MASTERFLEX® Economy Analog Console L/S® Pump Drives
Model Numbers: 7554-85 & -95, 77910-25
Description of Apparatus: Variable speed, peristaltic pump motor drive. Used with pump head and peristaltic tubing to pump fluids.
7554-85 (200 rpm) & -95 (600 rpm);
77910-25; (Bundled 7554-95)

Barnant Company declares that the above models are in conformity to the following harmonized standards and directives

Applicable Directives	Applicable Specifications	Manufacturer's Report Number
73/23/EEC 93/68/EEC	EN61010-1/A2:1995	TR9651
89/336/EEC 92/31/EEC 93/68/EEC	EN61326-1/A1:1998	TR9652

The last two digits of the year in which the current configuration of the above models were assessed per the Low Voltage Directive is :00.

Manufacturer: Barnant Company Division
Cole-Parmer Instrument Company
28W092 Commercial Avenue
Barrington, IL 60010-2392
USA
Tel: 847-381-7050

Manufacturer's Signature: James W. Doll, 17 July, 2000
James W. Doll — Date
Vice President, Engineering

MASTERFLEX® L/S®

WARRANTY

Use only MASTERFLEX precision tubing with MASTERFLEX pumps to ensure optimum performance. Use of other tubing may void applicable warranties.

The Manufacturer warrants this product to be free from significant deviations from published specifications. If repair or adjustment is necessary within the warranty period, the problem will be corrected at no charge if it is not due to misuse or abuse on your part, as determined by the Manufacturer. Repair costs outside the warranty period, or those resulting from product misuse or abuse, may be invoiced to you.

The warranty period for this product is noted on the Warranty Card.

PRODUCT RETURN

To limit charges and delays, contact the seller or Manufacturer for authorization and shipping instructions before returning the product, either within or outside of the warranty period. When returning the product, please state the reason for the return. For your protection, pack the product carefully and insure it against possible damage or loss. Any damages resulting from improper packaging are your responsibility.

TECHNICAL ASSISTANCE

If you have any questions about the use of this product, contact the Manufacturer or authorized seller.

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ISO 9001 REGISTERED

Printed in U.S.A.
030801

Artisan Technology Group is an independent supplier of quality pre-owned equipment

Gold-standard solutions

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

Learn more!

Visit us at **artisantg.com** for more info on price quotes, drivers, technical specifications, manuals, and documentation.

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

We're here to make your life easier. How can we help you today?

(217) 352-9330 | sales@artisantg.com | artisantg.com

ARTISAN® TECHNOLOGY GROUP

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Trademarks bearing the ® symbol in this publication are registered in the U.S. and in other countries.
Barrington, Illinois U.S.A. 60010-2392