ARTISAN®
TECHNOLOGY GROUP

Your **definitive** source for quality pre-owned equipment.

Artisan Technology Group
(217) 352-9330 | sales@artisantg.com | artisantg.com

Full-service, independent repair center
with experienced engineers and technicians on staff.

We buy your excess, underutilized, and idle equipment
along with credit for buybacks and trade-ins.

Custom engineering
so your equipment works exactly as you specify.

- Critical and expedited services
- In stock / Ready-to-ship
- Leasing / Rentals / Demos
- ITAR-certified secure asset solutions

Expert team | Trust guarantee | 100% satisfaction


All trademarks, brand names, and brands appearing herein are the property of their respective owners.

Visit our website - *Click HERE*

DOC023.97.80196

Sample Heater

02/2012, Edition 3

User Manual
Manuel d'utilisation

English 3
Français 20

Specifications

Specifications subject to change without notice.


Application	For water samples. Do not use the application for concentrated acids, bases or any type of volatile sample.
Output temperature precision	±1.5 ºC (2.7 ºF) under reference conditions [1]
Display resolution	1 ºC (1 ºF)
Set point range	20 to 50 °C (68 to 122 °F) Refer to operating charts for more information.
Flow rate range	50-300 mL/min; Sample Heater internal volume is 39 mL
Power requirements for single phase (with grounded neutral) only	120 Vac, 50/60 Hz, 11A 230 Vac, 50/60 Hz, 5.5A
Installation category	II
Inlet pressure	689.5 kPa (100 psi) max. internal pressure. This rating is for the Sample Heater and does not affect the Analyzer requirements.
Sample fittings	Two ¼" OD Tubing Compression fittings
Weight	Net—4.98 kg (11 lb) Shipping—6.35 kg (14 lb)
Wetted parts	316 Stainless Steel
Ambient operating temperature	0 to 50 °C (32 to 122 °F)
Pollution Degree	2
Protection Class	I
Maximum Altitude	2000 m (6560 ft)
Enclosure Rating	NEMA 4X, IP 65
Compliance Information	Listed for use in general locations to UL and CSA Standards by ETL. CE certified.
Storage Temperature	-20 to 70 ºC (-4 to 158 ºF) with fluids drained.
Humidity	Maximum relative humidity 80% for temperatures up to 31 ºC decreasing linearly to 50% R.H. at 40 ºC.

[1] Within proportional control range (see charts) and input temperature ±5 ºC, flow rate ±10%, 115 VAC, 60 Hz ±2 Hz and ambient temperature ± ºC. Instrument tuned at operating conditions (Autotune 3). Cal offset adjusted to yield the set temperature at the point of use.

Sample heater operating range


Figure 1 shows in graph form the sample heater operating range in terms of sample output temperature and minutes. Figure 2 shows in graph form the sample heater operating range in terms of sample output temperature and line AC voltage.

Figure 1 Temperature-minutes

1 Sample output temperature, ºC	**2** Minutes

Figure 2 Temperature-line AC voltage

1 Sample output temperature, ºC	**2** Line voltage, VAC

**milliliters per minute

General information

In no event will the manufacturer be liable for direct, indirect, special, incidental or consequential damages resulting from any defect or omission in this manual. The manufacturer reserves the right to make changes in this manual and the products it describes at any time, without notice or obligation. Revised editions are found on the manufacturer's website.

Safety information

NOTICE
The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

Please read this entire manual before unpacking, setting up or operating this equipment. Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

Make sure that the protection provided by this equipment is not impaired. Do not use or install this equipment in any manner other than that specified in this manual.

Use of hazard information

⚠ DANGER
Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.

⚠ WARNING
Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.

⚠ CAUTION
Indicates a potentially hazardous situation that may result in minor or moderate injury.

NOTICE
Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

Precautionary labels

	This is the safety alert symbol. Obey all safety messages that follow this symbol to avoid potential injury. If on the instrument, refer to the instruction manual for operation or safety information.
	This symbol indicates that a risk of electrical shock and/or electrocution exists.
	Electrical equipment marked with this symbol may not be disposed of in European public disposal systems after 12 August of 2005. In conformity with European local and national regulations (EU Directive 2002/98/EC), European electrical equipment users must now return old or end-of-life equipment to the Producer for disposal at no charge to the user.

Certification

Canadian Radio Interference-Causing Equipment Regulation, IECS-003, Class A:

Supporting test records reside with the manufacturer.

This Class A digital apparatus meets all requirements of the Canadian Interference-Causing Equipment Regulations.

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FCC Part 15, Class "A" Limits

Supporting test records reside with the manufacturer. The device complies with Part 15 of the FCC Rules. Operation is subject to the following conditions:

1. The equipment may not cause harmful interference.
2. The equipment must accept any interference received, including interference that may cause undesired operation.

Changes or modifications to this equipment not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference, in which case the user will be required to correct the interference at their expense. The following techniques can be used to reduce interference problems:

1. Disconnect the equipment from its power source to verify that it is or is not the source of the interference.
2. If the equipment is connected to the same outlet as the device experiencing interference, connect the equipment to a different outlet.
3. Move the equipment away from the device receiving the interference.
4. Reposition the receiving antenna for the device receiving the interference.
5. Try combinations of the above.

Product overview

⚠ DANGER

The Sample Heater Model 48685 is not designed for use with flammable or explosive samples. If a sample solution other than water is used in this product, test the sample/product compatibility to assure user safety and correct product performance.

Figure 3 shows the sample heater. The sample heater has a thermostatically controlled outlet temperature and is designed for use with water analysis instruments that require an elevated and controlled sample temperature.

Figure 3 Sample heater

The heater contains a heater module, temperature controller, temperature sensor, thermal cut-off switch, circuit breaker and power-connection terminal. The housing consists of a gasket-sealed plastic enclosure.


A single piece of stainless steel tubing cast in an aluminum block and two 700 watt electric heater elements make up the heater module. The stainless steel tube has two stainless steel adapter fittings brazed on each end.

The temperature controller is a Watlow EZ-ZONE PM Express with microprocessor-based auto-tuning control.

Figure 3 shows the normal operating display, the control panel and the controller keys, and describes the controller key functions.

The lower window displays the set-point temperature and the upper window displays the process temperature. The operator enters the setpoint and the controller heats the sample to the setpoint. The 1 indicator light is on when the heating element turns on. A thermocouple attached to the outlet tube of the heater module makes up the temperature sensor. The sensor indicates the sample temperature at the outlet. A thermostat mounted on the heater block protects the product and user from a thermal run-away fault (Wire the instrument on page 11).

Figure 4 Sample heater display and control panel

1 Infinity key--clears and silences alarms. This key is factory locked by default. When the menu is unlocked, this key will clear and silence alarms, back up one level, or if held for 2-3 seconds return to the Operation Menu.	**5** Lower display--shows setpoint temperature
2 Upper display--shows process temperature	**6** Up and Down arrow keys--Adjusts the setpoint in the lower display. In other unlocked pages, changes the upper display to a higher or lower value, or changes a parameter selection.
3 Output 1 indicator--shows output activity	**7** Advance key--moves through parameter prompts
4 Temperature unit indicator (Celsius)	**8** A/M key--toggle the control mode. Locked by default. When the Operation Menu is unlocked, this key will toggle the control mode, and set mode to Manual when the A/M key is pressed for 3 seconds.

Product components

Make sure that all components have been received. If any items are missing or damaged, contact the manufacturer or a sales representative immediately.

Installation

Mount the instrument

⚠DANGER

Burn, electrical, fire and chemical hazards. Only qualified personnel should conduct the installation tasks in this section of the manual.

Figure 5 and Figure 6 show the sample heater dimensions. The sample heater has four mounting feet. Install the instrument vertically on a wall with the sample connections on the bottom. Refer to the appropriate manual for more information.

Figure 5 Sample heater (front view)



Figure 6 Sample heater bottom and side view

1 .875" (22.23 mm) O.D. for ½" NPT conduit fitting	**3** Sample inlet--0.25 (6.35 mm) O.D. tubing stainless steel compression fitting
2 Sample outlet--0.25 (6.35 mm) O.D. tubing stainless steel compression fitting	

Sample connections

⚠ DANGER

The Sample Heater Model 48685 is not designed for use with samples that are flammable or explosive in nature. If any sample solution other than water is used in this product, test the sample/product compatibility to assure user safety and proper product performance.

⚠ CAUTION

The instrument is factory tested to withstand 100 psi. Pressures greater than 100 psi may result in operator injury and instrument damage.

Sample connections are stainless steel compression type fittings sized for ¼" OD stainless steel or rigid wall plastic tubing. No tubing is supplied with the sample heater. Sufficient tubing is normally supplied with analyzers to allow for sample heater installation.

The heater tubing is ¼" OD stainless steel cast in position in the heater block. The threaded end pieces are brazed to the tubing. Twisted tubing may cause irreversible damage to the instrument. Keep the tubing from becoming twisted and prevent the ends from bending when tightening the connecting tubing. Make sure to plumb the heater so that the heater tubing is free from stress and vibration.

If a leak develops between the brazed-on end fitting and the tubing adapter, hold the end fittings while tightening the tubing adapter fitting. Use a 5/16" wrench on the flats machined on the fitting.

To make the sample connections:

1. Put the tubing in the fitting assembly and push it firmly against the internal shoulder.
2. Loosen the compression nut and make sure the tubing passes through the ferrules.

3. Finger tighten the nut. Use a 9/16" wrench to hold the adapter fitting and a second 9/16" wrench to tighten the nut 1¼ turns.
 Note: *Never tighten the tubing nut without holding the hex portion of the tubing adapter.*

Test the sample connections

Make sure that all the tubing connections are correct. Switch the sample on and gradually increase flow through the instrument. Make sure that no leaks exist. If leaks exist, tighten the connections in small increments. If a leak continues, loosen the fitting. Make sure that the ferrule is correctly installed and the fitting contains no foreign material. Tighten the fitting again.

Connect to AC power

⚠ DANGER

Potential Electrocution Hazard. Only qualified personnel should conduct the tasks described in this section of the manual.

⚠ DANGER

To meet UL, CSA and other applicable instrument safety standards, an external power disconnect switch must be installed. Install the device in a location and position that gives easy access to the disconnect device and its operation.

⚠ DANGER

Connect equipment in accordance with local, state, and national electrical codes.

Table 1 Power configuration and approved use by model number

Sample heater model number	Power connections	Approved use
48685-00 (48685-60 on reference serial number label)	115 VAC, 50/60 Hz line power ***Note:*** *The 115 VAC model cannot be modified for 230 VAC operation.*	United States, Canada
48685-02 (48685-62 on reference serial number label)	230 VAC, 50/60 Hz line power for single phase (with grounded neutral) only	United States, Canada, European Community and other countries with 230 V single phase

Wire the instrument

The instrument is designed to be hard-wired. A ½" NPT conduit opening ($^{7}/_{8}$" diameter) is provided for electrical connection. Use the correct sealing hardware to keep the NEMA 4 rating.
Note: *Do not apply power to the heater until sample flow has been established and sample is in the heater tubing. Sample must be flowing through the heater at all times when the heater is on.*

Requirements for the power cord installation

To use a power cord for the connection to AC mains, make sure that the:

- Cord is less than 3 m (9 ft) in length
- Cord wire gauge is less than 16 AWG
- Outer power cord insulation is suitable for end use location
- Cord has a 3-circuit cord plug that has a ground connection that is necessary for a supply connection
- Cord wires have the correct size of spade tongue terminals on the connection end. The terminals must have insulated barrels.
- Cord has a cord cable gland that has the correct rating for the installation location. The cable gland must attach the cord and seal to the enclosure to keep environmental ratings within specifications.
- Cord plug serves as the power supply disconnect. Make sure that it is close to the instrument and easily accessible by the operator.

- Cord is wired according to Color codes for the power cord on page 12.

Color codes for the power cord

Table 2 shows the color codes for the power cord.


Table 2 Color codes for the power cord

Conductor Type	North American	IEC
Line conductor (Line 1 (Hot))	Black	Brown
Neutral conductor (Line 2 (Neutral))	White	Blue
Ground conductor	Green	Green with yellow stripe

Terminal block connections


Connect the ground line to the heater block terminal that shows the ground symbol (Figure 7). Connect the hot and neutral power supply lines to TB1 (Figure 7). The terminal block accepts up to 12 gauge wire. The manufacturer recommends 16–12 AWG wire based on installation requirements. The wiring insulation must be rated to 90 °C minimum.
Figure 8 shows the interior of the sample heater with connections.

Figure 7 Power connections at the J3 terminal block

1 Line 1 (Hot)	**3** Ground
2 Line 2 (Neutral)	**4** Top of instrument

Figure 8 Sample heater interior with connections

1 Power connections (top to bottom: Ground, Neutral, Hot)	**3** Circuit breaker for 115 V (230 V has two switches)
2 Thermostat reset switch	

Operation

Sample heater start up

Before startup, find the required sample temperature and sample flow rate of the analyzer. Start the sample flow through the heater, then apply power to the instrument. Adjust the sample conditioning system on the analyzer and make sure there are no leaks.

Instrument power

Apply power at the external (user-supplied) switch and the display lights up. If the display does not light, qualified service personnel must check the circuit breaker and the thermostat reset switch inside the heater enclosure. The area that contains the circuit breaker is not user-serviceable.

Set the temperature

1. Push the UP or DOWN arrows to set the desired temperature. The lower display shows the temperature setpoint (e.g., 50 °C).

2. Use the thermometer included with the Sample Heater to measure the temperature of the sample in the analyzer sample cell. Put the thermometer in the cell and wait for the temperature to stabilize. The difference between the temperature in the cell and the temperature shown on the display is due to loss of heat between the cell and the heater. For operator convenience, enter a calibration offset to adjust for the temperature loss (Enter the calibration offset on page 14). If necessary, insulate the sample lines to decrease the amount of heat loss.
3. The factory sets the PID parameters for the heater. The factory-set values are appropriate for most applications. If it is necessary to change the PID parameters, refer to Autotune on page 14.

Enter the calibration offset

The calibration offset is the difference between the required analyzer set-point temperature and the temperature in the sample cell. For example, if the required temperature is set to 50 ºC and the measured temperature of the sample cell in the analyzer is 47 ºC, the calibration offset is -3 ºC.

Wide variations in the temperature of the incoming sample or sample flow rate may cause a difference in the indicated process temperature. If temperature changes occur in the analyzer sample cell, make sure that the temperature setting on the sample heater and the offset agree with the required values.

If variations in the sample-intake temperature or flow cause temperature instability, change the offset to adjust. Refer to Sample Heater Specifications on page 3 to find the capability for the operating conditions.

To enter a calibration offset:

1. Push the ADVANCE key until **CAL** shows in the lower display. Push the DOWN arrow to decrease the value of the CAL offset by 3 ºC.
 Note: *If no CAL offset value is entered, the CAL offset value will be 0.*
2. Push the ADVANCE key.
 The sample heater returns to normal operation. The process temperature and the setpoint indicator should read the same temperature value.

Autotune

The Watlow EZ-ZONE PM ExpressMicroprocessor-based Auto-tuning Control has PID (proportional band, resent/integral and rate derivative) control capability. The factory sets the PID parameters for the heater. These values suit most applications. If it is necessary to change these parameters, follow the procedure below. This procedure allows the controller to automatically tune the parameters.

Note: *For successful auto-tuning, the setpoint temperature should be set at least 10 ºC above the inlet temperature.*


1. Push the UP and DOWN arrows and hold them for 3 seconds to access the ***Loc*** menu.
2. Push the UP arrow and change the upper display to 5.
3. Push the INFINITY key. The display returns to the Operation menu.
4. Push the ADVANCED key and scroll to the ***AUto*** menu.
5. Push the UP arrow and highlight ***YES***.
 The upper display flashes ***tune*** and then ***YES*** while the sample heater does the auto-tune procedure. The temperature must cross the auto-tune setpoint five times to complete the procedure. The controller installs the appropriate PID parameters and saves them in nonvolatile memory. The sample heater will auto-tune itself around 90% of the set temperature value.

Controller setup menu



The operational setup of the PID controller depends on specifics of the use-case. The manufacturer has set the controller to operate with the manufacturer's sample heater.

The tables show the values in the upper and lower displays and give information about the values. The values shown are default values. If necessary, use the UP and DOWN arrows to change a value and then move to the next step.


1. Push and hold the up and down arrow keys for approximately 3 seconds to go to the Lockout menu. Make the upper display show **5** to allow read and write access.

5	Operations menu and setup menu with full R/W access
	Lockout menu


2. Push **ADVANCE**.

	Sensor type set to Thermocouple
	Sensor type


3. Push **ADVANCE**.

	J-type thermocouple as input measuring device
	Linearzation



4. Push **ADVANCE**.

	Temperature input and display units in whole numbers
	Decimal



5. Push **ADVANCE**.

	Temperature input and display units in Celsius
	Display units



6. Push **ADVANCE**.

	The lower range set point limit is 20 ºC
	Range Low



7. Push **ADVANCE**.

	The upper range set point limit is 50 ºC
	Range High



8. Push **ADVANCE**.

	Output set to heating
	Function of Output 1



9. Push **ADVANCE**.

	PID proportional-integral derivative
	Heat Algorithm



10. Push **ADVANCE**.

	Alarm set to off
	Alarm type



11. Push **ADVANCE**.

	Ramp action set to off
	Ramp action



12. Push **ADVANCE**.

	Set maximum value output 1 range
	Power scale high output 1


13. Push **ADVANCE**.

	Active process value
	Upper or Left display


14. Push **ADVANCE**.

	Active Set Point
	Lower or Right display

15. Push **ADVANCE**.

	Set to 1
	Zone Address–Standard Bus Communication

16. Push **ADVANCE** and then the DOWN arrow until the upper display shows **2**.

2	A/M button disabled, set point R/W only.
	Lockout Menu


17. Push INFINITY.
The sample heater returns to normal operations.

Sample heater operation menu

The upper display window shows the process temperature. The lower display window shows the temperature setpoint. The 1 indicator lamp is on when the heater is on. No other display indicators show.


If necessary, use the up and down arrow keys to change a value shown in the display, then move to the next step.

1. To go to the Lockout menu, push and hold the up and down arrow keys for approximately 3 seconds. Change the upper display to 5.



5	Operations menu and setup menu with full R/W access
	Lockout menu

2. Push INFINITY.
The sample heater enters the operation menu.



3. Push **ADVANCE**.

	Make sure that this is set to **no**. (To initially tune the heater, select **YES**. Refer to Autotune on page 14.)
	Autotune


4. Push **ADVANCE**.

	Make sure that this is set to **Auto**
	Control Mode Active


5. Push **ADVANCE**.

	Set to 9
	Heat Proportional Band



6. Push **ADVANCE**.

	Set to 109
	Time Integral



7. Push **ADVANCE**.

	Set to 18
	Time Derivative

8. Push **ADVANCE**.

	Set to 7
	Time Base Output 1

9. Push **ADVANCE**.

	Set to -3
	Calibration Offset

The calibration offset is shown in the upper display window. The value entered here corrects for any temperature offsets that exist in the sampling application. Change this value to make the process display agree with the assessment of the process temperature. Refer to Enter the calibration offset on page 14.

10. Push the INFINITY key to go to the Operation menu.
11. Push and hold the UP and DOWN arrow keys for approximately 3 seconds to go to the Lockout menu.
12. Change the upper display to **2**, then push INFINITY.
 The sample heater returns to normal operation.

Maintenance

Instrument maintenance and cleaning

The manufacturer designed the instrument to be maintenance-free except for exterior cleaning. If the instrument exterior becomes dirty, clean the exterior surfaces with a damp cloth.

Troubleshooting

⚠ WARNING
Multiple hazards. Do not disassemble the instrument for maintenance or service. If the internal components must be cleaned or repaired, contact the manufacturer.

Common problems

Problem	Possible solutions
The PID display does not light	• Make sure the external power switch is in the On position • Check the wiring connections • Make sure the circuit breaker is in the On position • Make sure the thermostat cutoff is not tripped (Wire the instrument on page 11). If the cutoff has tripped, the middle button will snap down when lightly pressed.
The PID display is on but the process temperature does not drive to the set temperature	• Make sure the sample flow rate and temperature are within the operational range • Check all electrical connections • Make sure the temperature controller setup is correct • Make sure the temperature of the sample temperature is not higher than the set temperature of the heater

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© Hach Company/Hach Lange GmbH, 2011–2012. All rights reserved. Printed in U.S.A.

Artisan Technology Group is an independent supplier of quality pre-owned equipment

Gold-standard solutions

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

Learn more!

Visit us at **artisantg.com** for more info on price quotes, drivers, technical specifications, manuals, and documentation.

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

We're here to make your life easier. How can we help you today?

(217) 352-9330 | sales@artisantg.com | artisantg.com

ARTISAN® TECHNOLOGY GROUP

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(800) 227-4224 (U.S.A. only)