ARTISAN®
TECHNOLOGY GROUP

Your **definitive** source for quality pre-owned equipment.

Artisan Technology Group
(217) 352-9330 | sales@artisantg.com | artisantg.com

Full-service, independent repair center
with experienced engineers and technicians on staff.

We buy your excess, underutilized, and idle equipment
along with credit for buybacks and trade-ins.

Custom engineering
so your equipment works exactly as you specify.

- Critical and expedited services
- In stock / Ready-to-ship
- Leasing / Rentals / Demos
- ITAR-certified secure asset solutions

Expert team | Trust guarantee | 100% satisfaction


All trademarks, brand names, and brands appearing herein are the property of their respective owners.

Find the *Hach H160* at our website: *Click HERE*

DOC022.97.80058

H160, H170

10/2013, Edition 1

User Manual
Manuel d'utilisation

English 3
Français 21

Table of contents

Specifications on page 3
General information on page 4
Installation on page 7
User interface and navigation on page 8
Startup on page 10
Operation on page 11
Maintenance on page 14
Troubleshooting on page 17

Specifications

Specifications are subject to change without notice.

Specification	Details
Dimensions (W x D x H)	9 x 20 x 5 cm (3.5 x 8 x 2 in.)
Enclosure	Water resistant
Weight	1300 g (3.0 lb.)
Power requirements	4-ANSI 15 A or IEC-LR6 (AA Alkaline)
Power consumptions	1 W with backlight on
Battery life	Up to 200 hours continuous operation, 10-hour advance of low battery warning
Operating temperature	5 to 40 °C (41 to 104 °F)
Storage temperature	–20 to +40 °C (4 to 140 °F)
Operating humidity	Relative humidity: 50% maximum at 25 ºC (77 ºF), non-condensing
Input connectors	8-pin ISFET, BNC with phono jack temperature, 12-pin conductivity (H170 only)
Languages	English
Temperature correction	Off, automatic and manual (parameter dependent)
Measurement	Continuous measurement
Protection rating	IP67
Certifications	CE
pH	
Calibration	Up to five points: 1.68, 4.01, 6.86, 7.00, 9.18, 10.01, 12.45
Accuracy	±0.01 pH
Resolution	0.01 pH
Range	0.00 to 14.00
mV	
Calibration	None
Accuracy	0.1 mV
Resolution	Autoranging, 0.1 and 1
Range	Autoranging, ±199.9 mV to ±1999 mV
Temperature	
Calibration	None
Accuracy	±0.5 °C

Specification	Details
Resolution	0.1 °C (0.1 °F)
Range	–0 to 100 °C (32 to 212°F)
Conductivity	
Calibration	Up to five points
Accuracy	±1% ppt (–2 to +35 °C)
Resolution	0.01 iS, 0.1 iS ,1 iS, 0.01 mS, 0.1 mS
Range	Autoranging: 0.00 to 19.99 iS, 20.0 to 199.9 iS, 200 to 1999 iS, 2.00 to 19.99 mS, 20.0 to 199.9 mS
TDS	
Calibration	Up to five points
Accuracy	±1% full scale or ±1 digit
Resolution	0.01 ppm, 0.1 ppm, 1 ppm, 0.01 ppt, 0.1 ppt, 1 ptt, 0.1 mg/L, 1mg/L ,0.01 gal/L, 0.1 gal/L
Range	Autoranging: 0.00 to 9.99 ppm, 10.0 to 99.9 pm, 100 to 999 ppm,1.00 to 9.99 ppt, 10.0 to 99.9 ppt, 100 to 200 ppt mg/L: 0.00 to 199.9 mg/L, 200 to 1999 mg/L, 2.00 to 19.99 gal/L, 20.0 to 50 gal/L
Salinity	
Calibration	None (derived from conductivity)
Accuracy	±0.1 ppt, –2 to +35 °C (28.4 to 95 °F)
Resolution	0.1 ppt 1%
Range	0 to 42 ppt

General information

In no event will the manufacturer be liable for direct, indirect, special, incidental or consequential damages resulting from any defect or omission in this manual. The manufacturer reserves the right to make changes in this manual and the products it describes at any time, without notice or obligation. Revised editions are found on the manufacturer's website.

Safety information

NOTICE

The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

Please read this entire manual before unpacking, setting up or operating this equipment. Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

Make sure that the protection provided by this equipment is not impaired. Do not use or install this equipment in any manner other than that specified in this manual.

Use of hazard information

⚠ DANGER

Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.

⚠ WARNING

Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.

⚠ CAUTION

Indicates a potentially hazardous situation that may result in minor or moderate injury.

NOTICE

Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

Precautionary labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed. A symbol on the instrument is referenced in the manual with a precautionary statement.

	This is the safety alert symbol. Obey all safety messages that follow this symbol to avoid potential injury. If on the instrument, refer to the instruction manual for operation or safety information.
	This symbol indicates that a risk of electrical shock and/or electrocution exists.
	Electrical equipment marked with this symbol may not be disposed of in European public disposal systems after 12 August of 2005. In conformity with European local and national regulations (EU Directive 2002/96/EC), European electrical equipment users must now return old or end-of-life equipment to the Producer for disposal at no charge to the user. ***Note:*** *For return for recycling, please contact the equipment producer or supplier for instructions on how to return end-of-life equipment, producer-supplied electrical accessories, and all auxiliary items for proper disposal.*

Certification

Canadian Radio Interference-Causing Equipment Regulation, IECS-003, Class A:

Supporting test records reside with the manufacturer.

This Class A digital apparatus meets all requirements of the Canadian Interference-Causing Equipment Regulations.

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FCC Part 15, Class "A" Limits

Supporting test records reside with the manufacturer. The device complies with Part 15 of the FCC Rules. Operation is subject to the following conditions:

1. The equipment may not cause harmful interference.
2. The equipment must accept any interference received, including interference that may cause undesired operation.

Changes or modifications to this equipment not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference, in which case the user will be required to correct the interference at their expense. The following techniques can be used to reduce interference problems:

1. Disconnect the equipment from its power source to verify that it is or is not the source of the interference.
2. If the equipment is connected to the same outlet as the device experiencing interference, connect the equipment to a different outlet.
3. Move the equipment away from the device receiving the interference.
4. Reposition the receiving antenna for the device receiving the interference.
5. Try combinations of the above.

Product overview

NOTICE


Always disconnect power to the meter before an electrode is changed. Only use the meter as instructed in this manual or instrument damage can occur.

The meter is available in two models:

- H160—pH and ORP
- H170—pH, ORP, conductivity, TDS and salinity

The H-series waterproof handheld meters use glass sensor electrodes with BNC connectors and non-glass probes with ISFET (ion sensitive field effect transistor) silicon chip sensors. The meters use one pH electrode (a BNC pH electrode or an ISFET pH probe) at a time. When a meter is set to on, the meter automatically identifies the type of electrode that is attached.

Figure 1 Handheld meter


1	Rubber dust caps	**5**	ISFET pH probe connector (8-pin)
2	Conductivity probe connection (12-pin, H170 only)	**6**	LCD display
3	BNC probe connector	**7**	Power button
4	3.5 mm phono jack for glass pH electrode, ISE, ORP or DO temperature sensors	**8**	Battery cover

Product components

Make sure that all components have been received. Refer to Figure 2. If any items are missing or damaged, contact the manufacturer or a sales representative immediately.

Figure 2 Handheld meter components

1 Waterproof handheld meter (H160 shown)	**2** AA Alkaline batteries (4x)

Installation

⚠ CAUTION

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

Electrical installation

Install the batteries

⚠ CAUTION

Explosion hazard. Incorrect battery installation can cause the release of explosive gases. Be sure that the batteries are of the same approved chemical type and are inserted in the correct orientation. Do not mix new and used batteries.

NOTICE

Do not tighten the screws too much or instrument damage can occur.

NOTICE

Do not complete this operation with probes attached to the meter. Remove all probes from the meter or instrument damage can occur.


The meter uses AA alkaline or rechargeable NiMH batteries. To conserve the battery life, refer to Configure the meter on page 10 to configure automatic shutdown of the meter.

Refer to Figure 3 to install the batteries.

Items to collect:

- Phillips screwdriver
- AA Alkaline batteries (4x)

Figure 3 Battery installation

Probe and connection types

NOTICE

Make sure that all the rubber plugs are fully installed before the handheld meter is put in water or the interior of the 3.5 mm phono jack will fill with water. However, the interior of the handheld meter will stay sealed. If water gets into the 3.5 mm phono jack, immediately invert the meter to let water come out of the phone jack, then let the connector dry.

NOTICE

Always set the power to off before a probe connection is changed or damage to the instrument can occur.

The meter can connect to two pH-sensing technologies: ISFET pH probes and glass electrodes. When the meter is set to on, it automatically recognizes the type of probe attached. The meter only uses one type of pH probe at a time.

Note: *If a probe is not attached to the meter, the selected mode is inactive.*

User interface and navigation

Keypad description

Figure 4 and Table 1 show the function of the keys for each meter.

Figure 4 Keypad descriptions



Table 1 Keypad functions

Key	Action	H160	H170
	ON/OFF: Set the power to on or off.	x	x
pH	pH mode	x	x
mV	mV mode	x	x
	Set the backlight to on. After 2 minutes without a keystroke, the backlight is set to off.	x	x
TDS	TDS mode		x
CON	Conductivity mode		x
SAL	Salinity mode		x
HOLD	Hold mode: Operates only in measurement modes. All readings freeze when hold mode is set to on.	x	x
CAL	Calibration mode	x	x
	Scroll through values, setup screens and options.	x	x
SETUP	Setup mode	x	x
ENTER	Select an option, setting or value.	x	x

Display description

Figure 5 shows measurement modes and values, data connection and storage options, battery status, temperature values, stabilization lock and connection statuses shown on the display.

Figure 5 Display

1 Setup mode	**7** Measurement mode
2 Temperature and data values	**8** Measurement units
3 Measured value	**9** Automatic Temperature Compensation (ATC) value (pH, conductivity, TDS or salinity)
4 Stabilization lock	**10** Data storage memory location
5 Hold indicator	**11** Temperature unit
6 ISFET probe indicator	**12** Low battery indicator

Startup

Set the power to on

Push the power button to set the power to on or off. Make sure that the power supply (AC power or battery power) is correctly installed.

Configure the meter

Complete the operations in order:

1. Set the power to on.
2. Push **SETUP** to put the meter in setup mode.
3. Use the arrow keys to select an option, then push **ENTER**:

Option	Function	Description
READY	Stabilization lock	Sets the stabilization lock to on or off. • The "READY" icon shows when an endpoint occurs and the stabilization lock is set to on. • After stabilization, the display locks the value and the meter ignores slight changes in the measurement of pH, conductivity or TDS. • Sets the stabilization lock to off during titrations or slight change detection. The display automatically unlocks after a significant measurement change is found.
ºC ºF	Temperature display	Sets the temperature unit to ºC or ºF.
Snd	Sound options	Sets the sound alerts to on or off. Three sounds tell the user about different functions: • **One beep**: Stabilization occurred with the stabilization lock on. • **Two beeps**: An error occurred. The error code shows on the display. Refer to Error codes on page 18. • **Three beeps**: Measurement stability in calibration mode, regardless of the stabilization lock setting.
oFF	Automatic shutdown	Sets shutdown parameters from 1 minute to 000 minutes (continuous power). The meter beeps 1 minute before shutdown. Make sure to power cycle for the automatic shutdown to occur. Automatic shutdown disables during: data logging, Bluetooth™ transfers and when connected to the USB port.
CLr ALL CAL	Clear calibration points	Erases all calibration points in storage. Calibrate the meter again.[1]

[1] This does not erase calibration data for an ISFET probe.

4. Push **ENTER** to keep the changes and go back to setup mode.
 Note: *To exit and not keep the changes, push* ***SETUP****.*

Operation

Override an error

The user can override some errors (e.g., poor electrode condition) until a replacement electrode is ordered and received.
Note: *It is not possible to override all errors. The manufacturer does not recommend that the user override error codes. Contact technical support to override error codes.*

To cancel or override an error:

1. Find the cause of the error. Refer to Error codes on page 18. If possible, push **SETUP** or **pH, COND, TDS** or **SAL** to cancel the calibration or other error.
2. To override the error, do not complete step 1. Instead, push **ENTER**.
3. Use the arrow keys until "YES" shows.
4. Push **ENTER**.

Set the pH operations

Set the calibration alarm, the pH resolution, the calibration points and the calibration methods. Complete the operations in order.

1. Set the power to on.
2. Push **SETUP**, then **pH**.
3. Use the arrow keys to select an option, then push **ENTER**:

Operation	Function
pH resolution	Sets the pH resolution to 0.1 or 0.01 pH. The stabilization lock (the "READY" icon) is faster at a 0.1 pH resolution.
Set the 2 or 3-point pH calibration	Shows the current number of calibration points. Use the arrow keys to set the calibration points.
pH calibration method	• Sets the pH buffer. The default is a 2-point calibration (pH 7.00 then pH 4.00). • Push **ENTER** again to set the first calibration buffer. The display shows the next buffer. Do this procedure until all buffers are set. Previously-selected buffers do not show.

4. Push **ENTER** to keep the changes and go back to setup mode.
 Note: *To exit and not keep the changes, push* ***SETUP****.*

Set the manual temperature compensation

Glass pH electrodes are automatically temperature compensated if a temperature probe is sensed in the 3.5 mm phono jack. If no probe is sensed, the meter defaults to the manual temperature compensation at 25.0 ºC. Manual temperature compensation is always in ºC. Refer to Set the manual temperature compensation for glass electrodes.

pH buffer values at various temperatures							
25 ºC	0 ºC	5 ºC	10 ºC	20 ºC	30 ºC	40 ºC	50 ºC
1.68	1.67	1.67	1.67	1.67	1.68	1.69	1.71
4.01	4.00	4.00	4.00	4.00	4.01	4.03	4.06
6.86	6.98	6.95	6.92	6.87	6.85	6.84	6.83
7.00	7.12	7.09	7.06	7.01	6.99	6.97	6.97
9.18	9.46	9.40	9.33	9.23	9.14	9.07	9.02
10.01	10.32	10.25	10.18	10.06	9.97	9.98	9.83
12.45	13.42	13.21	13.00	12.63	12.29	12.04	11.70

1. Push **SETUP**, then **pH** to put into pH setup mode.
2. Use the arrow keys until the manual temperature shows "- - -.- ºC" on the display.
3. Push **ENTER** to see the manual temperature value.
4. Use the arrow keys to set the manual temperature compensation.
5. Push **ENTER** to keep the changes and go back to setup mode.
 Note: *To exit and not keep the changes, push* ***SETUP****.*

Calibrate for pH measurement

The meter accepts 1,-2-or 3-point calibrations. The default is a 2-point calibration.

1. Set the power to on.
2. Connect the probe to the meter.
3. Push **pH** to put the meter in pH mode.
4. Put the probe in the buffer solution.
 Note: *To erase all previous calibration points, slopes and offsets, push* ***CAL****, then* ***SELECT****. To update a previous calibration with a single-point calibration, ignore this step.*
5. Push **CAL**, then the "CAL" icon flashes on the display.
 Note: *Stop a calibration when the pH value is in flash mode. Push* ***pH*** *or* ***SETUP*** *to go back to the measurement mode.*
6. If the buffer is correct, push **ENTER**, then the pH buffer flashes until the calibration is complete.

7. If the pH buffer is not correct, use the arrow keys to select another buffer.
8. Use the arrow keys to see the buffers: 1.68, 4.01, 6.86, 7.00, 9.18, 10.01 and 12.45. During "CAL 2" or "CAL 3," the meter ignores the values used in previous calibrations.
9. The meter beeps three times when a calibration is complete. "CAL 3" and the next pH buffer value (default = 4.01) flashes on the display.
10. Rinse the probe in deionized water. Put the probe in the second buffer solution.
11. If OK, push **ENTER**, or use the arrow keys to change the value. Push **pH** to accept a 1-point calibration and to go back to the pH measurement mode. If **ENTER** is selected, the buffer value flashes until calibration is complete.
12. If the meter is set for a 2-point calibration, the meter beeps three times when the calibration is complete. The display shows the slope, the mV offset at pH 7.00 and the number of calibration points. Review the calibration data, then push **pH** or **STORE** to start reading the pH value of another solution.
13. If the meter is setup for a 3-point calibration, do steps 1 through 12 until all the calibration points are set.

Set the conductivity and TDS options

Set the options for conductivity and TDS: conductivity calibration points, conductivity temperature compensation, conductivity temperature, cell constant, TDS conversion factor, TDS units and salinity units in ppm or as a percentage. Complete the operations in order:

1. Set the power to on.
2. Push **SETUP**, then **COND** or **TDS** to put the meter in conductivity or TDS setup modes.
3. Use the arrow keys to select an option, then push **ENTER**:

Operation	Function
Conductivity calibration point	Keeps up to five calibration points. • Use the arrow keys to select each option, then push **ENTER**. • The first active calibration point shows. Use the arrow keys to see other in-use values. ***Note:*** *"- - - -" shows on the display if there are no stored conductivity calibration points.*
Conductivity temperature[1]	• Shows the conductivity normalization value. • Use the display settings between 20 to 25 °C (default = 25 °C).
Conductivity probe cell constant	The cell constant is automatically found and shown on the display. Sets the cell constants to: • K = 0.5 for general procedures • K = 1 for low conductivity procedures (e.g., drinking water) • K = 10 for salinity or brackish water • If no conductivity probe is attached, the display shows "000."
TDS conversion factor	• TDS values are related to conductivity. Calibrate with a TDS standard or with conductivity standards, then set the meter with a conversion factor. • Sets the conversion factor to 0.40 to 1.00 (default = 0.50).

Operation	Function
TDS units	Use the arrow keys to set the TDS units to mg/L (gal/L) or ppm (ppt).
Salinity units in ppm (ppt) or %	Use the arrow keys to set the salinity units to ppm (ppt) or as a percentage.

[1] Conductivity varies a lot with temperature. The conductivity automatic temperature compensation (ATC) adjusts conductivity measurements to calculate the conductivity changes in the reading that are caused by temperature. The readings are referenced to or "normalized" at a standard temperature. The ATC gives the equivalent conductivity or TDS of a solution normalized at 20 or 25 °C.

4. Push **ENTER** to keep the changes and go back to setup mode.
 Note: *To exit and not keep the changes, push* ***SETUP****.*

Calibrate for conductivity and TDS measurement

Calibrate the meter with up to five calibration points with conductivity or TDS standards. If the calibration point is 20% larger than a calibration point in storage, both conductivity points are kept. If the calibration point is 20% smaller than a calibration point in storage, only the new calibration point is kept. If there are five calibration points in storage, the new calibration data replaces values near the same range.

By default, the meter uses the temperature coefficient of 2% per °C. The user can set the coefficient in setup mode. Review calibration points in setup mode.

Note: *Salinity measurements are identified from conductivity.*

1. Set the power to on.
2. Push **COND** or **TDS**.
3. Put the probe in the first standard solution. The meter measures the standard.
4. Use the arrow keys to select the value of the conductivity or TDS standard.
5. Push **ENTER**. The main display shows, then the "WAit" icon flashes.
 Note: *Push* ***TDS*** *or* ***COND*** *to stop calibration. The meter stops the current calibration and keeps the previous calibration data.*
6. Push **CAL**. The main value flashes.
7. To erase all previous calibration points, slopes and offset values, push **CAL**, then **DELETE**. Ignore this step to update a previous calibration with a single point calibration.
8. Push **COND** or **TDS** to read the conductivity or TDS value of any solution.
9. Do steps 2 through 9 again to complete multiple calibration points (maximum of 5).
10. When the calibration is complete, the meter shows the slope as compared to the nominal sensitivity for a probe of that cell constant and the total number of conductivity and TDS calibration points in the meter.
 Note: *The number at the bottom right corner of the display is the calibration points that are in the meter. The value above the slope is the slope relative to the nominal sensitivity.*

Maintenance

Replace the batteries

The screen reads "bAt" when the batteries are too low to give a reliable measurement. The low battery icon shows when there is approximately 25 hours of battery power. Measurement errors are possible when the batteries are low. Refer to Install the batteries on page 7 to replace the batteries.

Clean the instrument

Clean the exterior of the instrument with a moist cloth and a mild soap solution and then wipe the instrument dry.

ISFET pH probe maintenance

⚠ CAUTION

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

⚠ CAUTION

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

The expected life of an ISFET probe is approximately 18 months. The reference electrode has a KCl gel that is diluted over time. The reference is sealed and is non-refillable. Replace the probe when it becomes difficult to calibrate.

- Keep the probe dry with the protective shroud on when not in use.
- Clean oil, fat, food particles, starch, protein or other materials from the probe tip after use.
- Never use sharp metal objects (e.g., a needle, a pin, etc.) to clean the sensor surface.
- Remove the probe from environments with static electricity. Electrostatic discharge (ESD) can permanently damage the probe.
- Remove the probe from environments that will damage the sensor, such as hydrofluoric acid or abrasive samples.
- Remove the probe from environments that will damage the epoxy materials used in the probe tip (e.g., keep the probe away from acetone, toluene, methylene chloride, xylene and other strong organic solvents).
- Do not use the probe in temperatures more than 60 ºC (140 ºF). Thermal cycling can decrease the life of the probe.
- For semi-solids use, carefully twist the probe to make sufficient contact with the sample to the sensor.
- When semi-solids are tested, make sure that solid objects (i.e., bone or gristle) do not scratch the sensor.
- Cool samples to room temperature to maximize probe life.
- Always use new buffers and new rinse solutions.

Prepare the probe for storage

NOTICE

Do not use the probe for long-term pH measurement applications.

Note: *No electrode storage solution is necessary.*

A new probe usually has visible reference gel at the tip of the probe, as well as in the interior of the rubber dust cap. To remove the gel, carefully clean with a soft-bristled toothbrush and mild soapy water (a few drops of dish soap in a warm cup of water). The gel can continue to show for two to five days. Do not put the rubber dust cap on the probe until all of the reference gel is removed. Do not complete the steps that follow until all reference gel is removed from the probe and the interior and exterior of the rubber dust cap.

Prepare the probe for storage

1. Put the new probe (or the probe that was in extended storage) in pH 7 buffer for at least 5 minutes.
2. Stir the probe in pH 7 buffer solution to dislodge air bubbles.
3. Use fresh deionized water to rinse the probe.
4. Dry the probe with a lint-free cloth.

5. Put the rubber dust cap back on.
6. Keep the probe in dry storage when not in use.

Calibrate the ISFET probe

Refer to Prepare the probe for storage on page 15 before calibration.

NOTICE

Do not use the probe for applications that cycle between hot and room temperature samples.

Note: *Remove other active or non-active measurement devices during pH or conductivity measurement. Other devices, even AC power, can cause interference.*

Note: *Only do calibrations away from the sun. Direct sunlight can cause unstable readings or difficulty in calibration.*

1. Connect the probe to the meter.
2. Set the meter power to on.
3. Clean the probe with new deionized water and dry with a lint-free cloth.
4. Put the probe in the pH 7 buffer.
5. Clean the probe with deionized water and dry with a lint-free cloth.
6. Put the probe in the second buffer (pH 4 or pH 10).
7. Read the results.
8. If the result is not correct, the probe is not correctly hydrated.
 Soak the probe for another 5 minutes in pH 7 buffer, then do the calibration again.

Clean the ISFET probe

NOTICE

Do not use sharp metal objects (a needle, a pin, etc.) to clean the sensor. This can scratch the sensor and cause permanent damage to the probe.

For use in dairy, cheese or meat applications, soak the probe in Pepsin Cleaning Solution for 15 minutes before the probe is cleaned.

Regularly clean a non-glass probe:

1. Remove the rubber dust cap from the probe, then rinse with new deionized water.
2. Use a soft-bristle toothbrush and a mild detergent (a few drops of dish soap in a warm cup of water) to carefully clean the probe.
3. Rinse with new deionized water to remove all debris from the sensor surface.
4. Dry the probe with a lint-free cloth.
5. Calibrate the probe again. Refer to Calibrate the ISFET probe on page 16, then Prepare the probe for storage on page 15.

Repair the ISFET probes

NOTICE

Do not use the probe to find out if the buffer is above 60 °C (140 °F). If the probe is suddenly put into very hot liquid, the probe can be permanently damaged.

Monitor how long the probes are in dry storage. If the probes are in storage for an extended period of time, the KCl gel at the reference junction can crystallize.

1. Heat pH 7 buffer to approximately 45 to 60 °C (115 to 140 °F).
2. Soak the probe for 2 minutes.
3. Put the probe in room temperature pH 7.00 buffer and let cool.

BNC electrode maintenance

Prevent unstable readings:

- Keep the probe in an electrode storage solution.
- Start measurements with at least a 2-point calibration. Update often with 1, 2 or 3-point calibrations.
- Use new buffers and new deionized rinse solution.
- Use buffer solutions with pH values no larger than 3 pH units apart. Ideally, the buffers bracket the anticipated pH values of the unmeasured samples.
- Use deionized water to rinse residual buffer and sample solutions from the probe after calibration and measurement.
- Calibrate at the same temperature as the sample solution. Although the meter has an ATC, get the best results when the calibration buffers and the sample are the same temperature.
- Keep the connectors clean and dry. Dirty or damp connectors can cause unstable readings.

Troubleshooting

Problem	Possible cause	Solution
No display	Automatic shutdown set the power to off.	Set the power to on again.
	There is no power.	Replace the batteries.
Unstable reading	The probe is dirty.	Clean the probe.
	The probe or the meter connectors are dirty.	Clean the probe contacts on the cable connector. Clean the meter with methanol on a cotton swab. Let dry completely. Connect the probe to the meter again.
	No flow in the reference junction	Clean the warm buffer.
	ISFET probe is not properly hydrated.	Soak the probe for at least 5 minutes in a pH 7.00 buffer.
	Interference from other devices	Remove other devices from the solution.
	Interference from direct sunlight	Use protection for the probe from direct sunlight.
	The probe is in a very low ionic strength solution.	A stable reading is not possible.
	The pH or temperature of the solution changes.	A stable reading is not possible until pH and temperature are constant.

Problem	Possible cause	Solution
The pH or temperature of the solution changes.	The probe sensor surfaces are dirty or the probe needs to be conditioned again.	Clean and condition the probe again.
	No flow in the reference junction	Clean the warm buffer.
	The buffers are contaminated or expired.	Calibrate with new buffers again.
	Interference from other devices in solution	Remove other devices from the solution.
	The battery is low.	Replace the batteries if the battery icon shows low battery power.
	Interference from direct sunlight. Out of calibration.	Use protection for the probe from direct sunlight.
	The probe is too old.	Replace the probe.
The meter continually shows 14.00 or 0.00 with an electrode attached. The "ISFET" icon does not show when an ISFET probe is attached.	The ISFET probe is not sensed by the meter or the probe.	Set the meter power to off. If the ISFET icon does not show, replace the ISFET probe.
	Out of calibration.	Calibrate the ISFET probe on page 16
	The probe is not in the solution.	Put the probe in liquid. Carefully shake the probe to make sure that air bubbles are not caught on the sensor surface.
	The probe is dirty.	Clean the probe.
	The probe is damaged.	Replace the probe. Contact technical support.

Error codes

Table 2 shows what codes can occur for various reasons. Error codes can show instrument malfunction or user error.

Table 2 Error code descriptions

Error code	Description	Solution
E03	Glass electrode offset error.	Clean the probe. If the error continues, contact technical support to replace the probe.
E04	Glass electrode slope error.	Clean the probe. If the error continues, contact technical support to replace the probe.
E06	ISFET pH probe slope error.	Clean the probe. Soak in 40 ºC (113 ºF) pH 7.00 buffer for 2 minutes. Try calibration again. If the error continues, contact technical support to replace the probe.
E08	Too long to calibrate. Signal is not stable during calibration.	Clean the probe. Disconnect stirrers and other AC power sources. Make sure that calibration solution temperature is constant. If the error continues, contact technical support to replace the probe.

Table 2 Error code descriptions (continued)

Error code	Description	Solution
E13	The ISFET pH probe temperature sensor is out of range.	Replace the probe. Contact technical support for probe replacement information.
E20	The conductivity probe temperature sensor is out of range.	Replace the probe. Contact technical support for probe replacement information.

Replacement parts and accessories

⚠ WARNING

Personal injury hazard. Use of non-approved parts may cause personal injury, damage to the instrument or equipment malfunction. The replacement parts in this section are approved by the manufacturer.

Note: *Product and Article numbers may vary for some selling regions. Contact the appropriate distributor or refer to the company website for contact information.*

Replacement parts

Description	Quantity	Item no.
AA Alkaline batteries	4	1938004

Accessories

Description	Quantity	Item no.
ISFET pH probes		
General purpose round tip probe	1	PH77-SS
Micro pH probe	1	PH17-SS
Piercing tip micro pH probe	1	PH37-SS
pH probe for NMR tubes	1	PH47-SS
Heavy duty piercing tip micro pH probe	1	PH57-SS
Conductivity probes		
General purpose conductivity probe, epoxy, K = 0.5	4 graphite contacts	CDW97-KP5
LIS conductivity probe, glass, K = 1	4 platinum contacts	CDW97-K1
Conductivity probe for salinity and brackish water, epoxy, K = 10	4 contacts	CDW97-K10
Glass pH probes		
Gel-filled combination pH, BNC and 3.5 mm phone	1	5193511
Pt Series combination pH, BNC and 3.5 mm phone	1	5191011
Refillable combination pH, BNC and 3.5 mm phone	1	5194011

Consumables

Description	Quantity	Item no.
Probe cleaning solution for removal of protein deposits on pH probes	500 mL	2964349
Rinse solution, non-ionic surfactant	500 mL	2964449
Buffer solution, red, pH 4.01	500 mL	2283449
Buffer solution, yellow, pH 7.00	500 mL	2283549
Buffer solution, blue, pH 10.01	500 mL	2283649
Conductivity standard solution, 100 µS/cm	500 mL	2971849
Conductivity standard solution, 1000 µS/cm	500 mL	1440049
Conductivity standard solution, 10,000 µS/cm	500 mL	2972249
Conductivity standard solution of 35 ppt salinity, 53 mS/cm	500 mL	2714349
Conductivity standard solution, 100 µS/cm	100 mL	2971842
Conductivity standard solution, 1000 µS/cm	100 mL	1440042
Conductivity standard solution, 10,000 µS/cm	100 mL	2972242

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© Hach Company/Hach Lange GmbH, 2013.
All rights reserved. Printed in U.S.A.

Artisan Technology Group is an independent supplier of quality pre-owned equipment

Gold-standard solutions

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

Learn more!

Visit us at **artisantg.com** for more info on price quotes, drivers, technical specifications, manuals, and documentation.

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

We're here to make your life easier. How can we help you today?

(217) 352-9330 | sales@artisantg.com | artisantg.com

ARTISAN® TECHNOLOGY GROUP

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(800) 227-4224 (U.S.A. only)