ARTISAN®
TECHNOLOGY GROUP

Your **definitive** source for quality pre-owned equipment.

Artisan Technology Group
(217) 352-9330 | sales@artisantg.com | artisantg.com

Full-service, independent repair center
with experienced engineers and technicians on staff.

We buy your excess, underutilized, and idle equipment
along with credit for buybacks and trade-ins.

Custom engineering
so your equipment works exactly as you specify.

- Critical and expedited services
- In stock / Ready-to-ship
- Leasing / Rentals / Demos
- ITAR-certified secure asset solutions

Expert team | Trust guarantee | 100% satisfaction

All trademarks, brand names, and brands appearing herein are the property of their respective owners.

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2669600

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Operating Instructions

IKA COLOR SQUID

The IKA Color Squid is a magnetic stirrer suitable for mixing thinly fluid media up to 800 ml by using a magnetic bar up to a length of 30 mm. The speed can be adjusted between 0 and 1500 1/min.

Unpacking: Please unpack the instrument carefully and look out for any damage. It is important to recognize a transportation damage already when unpacking. If this is the case, an immediate factual report is required (mail, rail or forwarding agency).

Safety Warnings: The instrument may not be operated in rooms with explosion hazards. Furthermore, it has to be placed upon a fire-proof and/or not burnable surface. **Attention: Magnetism!** Effects of the magnetic field have to be taken into account (e.g. data carriers, cardiac pavemakers...). When using PTFE-coated magnetic bars, the following has to be noted: Chemical reactions of PTFE occur in contact with molten or dissolved alkaline and alkaline earth metals, as well as with fine-particled powders of metals of the 2. and 3. group of the periodical systems at temperatures above 300–400°C. Only elementary metals do attack it; halogen hydrocarbons have a reversibly swelling effect.

(Source: Römpps Chemie-Lexikon and „Ullmann" Vol. 19)

Attention: Surface should not be in contact with aggressive medium! Do not place under water!

Plug the instrument into a socket with a voltage of 100 ... 240 V ± 10% and a frequency of 50/60 Hz.

2669600

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Starting-Up:

Put a magnetic bar into a vessel permeable for magnetic field lines, and fill the vessel with the medium to be processed. Please make sure that the rotary knob for speed adjustment is at the left-hand stop, and then put the vessel on the glass hood. Now you can set the requested speed on the rotary knob.

Technical Data:

Max. stirring quantity: 800 ml
speed range: 0 ... 1500 1/min
voltage: 100 ... 240 V AC ± 10 %
frequency: 50/60 Hz
power: 4 watts
degree of contamination: 2
overvoltage category: II
protection class DIN 40050: IP 54
weight: 520 g
dimensions ∅ x H: 130 x 48 mm
preventive measure class: 0

– No additional protective device – Only for instruments involving no danger if control fails – The operating condition of the instruments has to be checked in regular intervals – For cleaning use only water with a tenside containing detergent or isopropyl alcohol (if very soiled).

Guarantee:

You have acquired an original IKA-laboratory instrument that meets highest requirements in techno-logy and quality. According to the IKA-terms of warranty, the guarantee period is 24 months. In the event of a claim under guarantee, please contact your local dealer. You can also send the instrument directly to our Works, enclosing the invoice and stating your complaints. Freight charges are at your expense.

IKA®

Reg.-No.4343

IKA®-WERKE GMBH & CO.KG
Janke & Kunkel-Str. 10 · D-79219 STAUFEN
TEL. 07633/831-0 · FAX 07633/831-98
E-mail: sales@ika.de Internet: http://www.ika.net

2669600

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

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Artisan Technology Group is an independent supplier of quality pre-owned equipment

Gold-standard solutions

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

Learn more!

Visit us at **artisantg.com** for more info on price quotes, drivers, technical specifications, manuals, and documentation.

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

We're here to make your life easier. How can we help you today?

(217) 352-9330 | sales@artisantg.com | artisantg.com

ARTISAN® TECHNOLOGY GROUP