ARTISAN®
TECHNOLOGY GROUP

Your **definitive** source for quality pre-owned equipment.

Artisan Technology Group
(217) 352-9330 | sales@artisantg.com | artisantg.com

Full-service, independent repair center
with experienced engineers and technicians on staff.

We buy your excess, underutilized, and idle equipment
along with credit for buybacks and trade-ins.

Custom engineering
so your equipment works exactly as you specify.


- Critical and expedited services
- In stock / Ready-to-ship
- Leasing / Rentals / Demos
- ITAR-certified secure asset solutions

Expert team | Trust guarantee | 100% satisfaction

All trademarks, brand names, and brands appearing herein are the property of their respective owners.

Find the *IKA HKD T 06 D* at our website: ***Click HERE***

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Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Table of Contents

	Page
Description	6
Drive	6
Trough	6
Kneading Blades	6
Bearing of Kneading Blades	7
Sealings	7
Vacuum	7
Starting-Up *Electrical Function*	7
The Kneading Process	7
Measuring the Temperature in the Product	8
Maintenance and Cleaning	8
Technical Data	8
List of Parts	9
Accessories	9
Guarantee	9

Description

The IKA-high-efficiency laboratory kneader HKD - T 0,6 is suitable for processing non-flowable, highly-viscous media. The uniform mixing is based upon an intensive processing by the wide-bladed kneading elements. The kneading medium is moved inside the trough horizontally as well as vertically. Thus, even smallest agglomerations are separated and distributed equally. The kneading elements intermesh deeply and narrowly. The drive is by a geared engine that is secured against overload by a temperature measuring sensor in the motor winding.
The kneader consists of a cast-aluminium base, a flangeconnected geared engine, the kneading gear unit as well as the kneading trough with cover and side part. The item numbers in the text following the individual part descriptions are repeated in the list of parts as well as in the diagrams.

Drive

The torque of appr. 40 Nm released by the geared engine (39) of 0,18kW, 36 min^{-1}, is transfered via a coupling sleeve (40) to the kneading gear unit. The kneading gear unit reduces the speed (2:1) of the rear kneading blade (118). The front kneading blade (117) is driven directly.
The shafts (112) and (113) of the kneading gear unit are supported in ball bearings. The gear unit has a permanent lubrication (46) and is sealed with radial sealing rings (33) and (34).
The right-hand side of the kneading gear unit is closed by a stainless steel bearing housing (70). The side wall of this bearing housing is at the same time the unremovable left side wall of the trough. When looking into the large bore holes of the bearing housing, with the trough removed, one recognizes the plug couplings into wich the kneading blades are plugged.

Trough

If the trough (116) is filled to the top, the filling volume if appr. 600ml. The actual kneading area (range of the rotating kneading blades) is appr. 300ml.
The bottom of the trough (can be heated and/or cooled) is double-walled. The direction of flow is optional. supply and discharge lines are fitted to the connections - R1/8" - of the trough wall. Here the permissible maximums are 210°C and 3bar.

Kneading Blades

The kneading blades are burnished. The rear kneading blade (118) has a longer shaft end on the plug-coupling side. If the two kneading blades are inadvertently mixed up when fixing, the right trough side cannot be fitted tightly.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Bearing and Kneading Blades

At the left-hand side the kneading blades are fitted in plug-coupling side. At the right-hand side the kneading blades are supported in ball bearings in the side wall. The axial clearence is adjusted by the adjusting screws (14). When removing the right side wall, the kneadingblades are generally left in the plug-coupling. They can be drawn off together with the trough.
When re-installing them (before installing the trough), the kneading blades must be placed in the direction of the slot in the plug coupling (view to the inside wall of the trough).

Sealings

All rubberlike sealings are made of Viton A.

PTFE - sealing strips (110) are in the grooves of the trough side walls. New PTFE - sealing strips (110) are simply pressed into the grooves, the ends of wich may overlap but should not projekt into the trough.
When tightening the three clamping screws (114), the PTFE - sealing strips get flatz also in the place they overlap.

The left pivots of the kneading blades are sealed by O-rings (37) and shaft sealing rings (35) towards the inside of the trough. The shaft sealing rings (35) are installed in such a way that they are vacuum-tight.

For exchanging the shaft sealing ring (35), the sealing bushing (87) must be draw n out of the bearing housing (70) after removing the small threaded pin (12). Then the shaft sealing ring can be ejected over the washer (88).The sealing bushings (87) are sealed outside with O-rings (38).

The right pinions of the kneading blades are sealed by O-rings (37). These O-rings sit directly in the removable side wall (115). For exchanging the O-rings (37), the protective covers (126) and (127) have to be unscrewed and the ball bearing (32) drawn out. After removing the supporting rings (120), the O-rings (37) are accessible.

Vacuum

The instrument is equipped for vacuum operation. The vacuum is connected at the R1/8" connection of the trough cover
Obtainable final pressure: appr. 50 mbar.

The viewing glass cover requiredfor vacuum operation is enclosed loosely. It is simply put on the O-ring in the trough cover. For the sake of safety, the opening in the trough cover has been provided with a guard grille.

When using the optionally available pestle, no vacuum operation is possible.

Starting-Up

Electrical Function The kneader is supplied with a standard three-phase motor of 3x400V. The voltage indicated on the identification plate hasto conform witch the mains voltage of the user.
If two voltages (e.g. 230/400V) are indicated on the identification plate, the instrument was switched to the higher voltage at works. On the buttom plate inside the instrument there is a circuit diagram.
With three-phase motors, attention has be paid to the sense of rotation of the motor. If the sense of rotation does not conform with the arrow of the sense of rotation on the motor, two phases in the plug have to be exchanged.
In line with the motor relay there is a push button (55) activated by the cover (119). It prevents the kneader from being started-up while it is still open. This switch also opens the relay current circuit when the trough is tilted. This safety device may not be deactived or circumvented.

The Kneading Process

During the kneading process, the front kneading blade partly wipes the kneading medium off the rear kneading blade (Duplex-principle). A perfect kneading function can only be optained if sufficient material has been added for filling the voids forming around the front kneading blade. For this purpose appr. 100ml
kneading material are required. The deeply intermeshing kneading blades transport the
kneading material from the side walls to the center of the trough.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

A further improvement of the kneading process can be obtained by using a double-walled pestle that is available on request. This pestle is adapted to the blade course and is used instead of the cover, wich implies that no vacuum operation is possible.
After the kneading process, flowable products can be discharged by tilting the trough. The arrangement of the kneading trough allows for a receiving tank to be placed below.

Non-flowable products have to be removed after the trough has been dismantled. After the three star screws have been unscrewed, the right side part (115) can be drawn off. Then the trough can be removed from the adjusting pins

The operating elements are designed in such a way that during high-temperatures protective gloves ensure operations

Measuring the Temperature in the Product

In the rear kneading blade shaft there is a deep boring reaching from the front side right to the middle of the trough. Into this boring the PT100 thermocouple (available as an accessory part) can be plugged in after the locking screw (5) has been unscrewed from the protectiv cover (126). The thermocouple is then connected to a suitable temperature measuring device (e.g. DTM 11 IKATRON digital temperature measuring devive).

Maintenance and Cleaning

The motor with spur-gear system and the kneading gear unit are maitenance-free. In case of repair, the instrument may only be openend by an expert.
When ordering spare parts, please always quote type of instrument, series number and ordering number of the spare part.
Please submit only instruments for repair that have been cleaned and are free from substances dangerous to health. The instrument surface can be cleaned with water with a tenside-containing detergent or isopropyl alcohol (e.g. if very soiled).

Technical Data

Type of kneading blade	Duplex
Trough	
Useful capacity min/max	100 / 300ml
total capacity	600ml
attainable vacuum	50mbar
trough saddle heatable up to	210 °C
permissible pressure in trough saddle	3bar
Drive	
Motor principle / type of power connection	asynchronous / 3
motor power (input / output)	320 / 180 watts
motor protection	thermal contact
nominal torque	40Nm
speed of front kneading blade	36 1/min
speed of rear kneading blade	18 1/min
safety device	cover contact
mains voltage and frequency	3x400V 50Hz *or* 3x230V 60Hz
Dimensions (W x T x H)	660 x 250 x 380mm
Gewicht	27kg
Zulässige Umgebungstemperatur	5 - 40 °C
Zulässige Feuchte	80%
Schutzart nach DIN 40 050	IP 54

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

List of Parts

Item	Designation	Item	Designation
5	locking screw	67	operating plate
6	cheese head screw	70	bearing housing
12	threaded pin	71	ring bearing
14	threaded pin	72	handle
16	hexagonal nut	75	guide bushing
20	sealing ring	76	shifter rod
22	safety washer	77	pressure spring
23	guard ring	78	toothed wheel
24	cylindrical pin	79	toothed wheel
25	cylindrical pin	81	protective hose
27	tapered pin	84	connection cable
31	grooved ball bearing	85	circuit diagram
32	grooved ball bearing	86	stop sleeve bush
33	radial shaft sealing ring	87	sealing bushing
34	radial shaft sealing ring	88	contact washer
35	radial shaft sealing ring	89	bearing foil
36	O-ring	103	glue-in-plate
37	O-ring	104	circuit diagram
39	three-phase engine	110	sealing strip
40	coupling bushing 12"	111	base plate
41	sealing ring	112	front gear shaft
42	grease, Q-paste NB 50	113	rear gear shaft
45	rubber foot	114	clamping screw
46	pinion grease Centoplex HO	115	side part
48	arrow plate	116	trough
54	small motor relay	117	front kneading blade
55	inst. limit switch	118	rear kneading blade
56	push botton „ON"	119	cover
57	push botton „OFF"	120	intermediate ring
59	stuffing-box stud	126	protective cover
63	screwed pipe joint	127	protective cover

Accessories


HKD 06.1	pestle	LT 5.1	hose, Viton
DTM 11	IKATRON digitales temp. meas.device	VC 2	IKAVAC vacuum controller
PT 100.7	temp. meas.sensor	VC 2.4	pump control
LT 5 digital	IKA circulating thermostat	MZ-2C	chemical diaphragm pump


Guarantee

You have purchased an original IKA laboratery machine which meets the highest engineering and quality standards.


In accordance with IKA guarantee conditions, the guarantee period is 12 months. For claims under the guarantee please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com


Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Schnitt C - D

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

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Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Supplement to operating instructions

HKD-T 06 D

The IKA high performance laboratory kneader meets European standards and has the CE mark.
Some kneader functions have been modified in order to ensure compliance with CE mark certification standards. The accompanying Instruction Manual does not include these modifications and does not therefore fully describe the current version of the IKA high-performance laboratory kneader.

This supplement to the instruction manual describes the functions which have been modified.

2. Installation

2.1. Electrical functions

The kneader is fitted as standard with a 3 x 400V three-phase AC motor. Before connecting up you must ensure that your mains power supply voltage matches that shown on the motor rating plate.
The rating plate shows two voltages (e.g. 230V/ 400V); the motor is switched to the higher voltage at our works prior to delivery. A circuit diagram is provided inside the unit on the base plate.
On AC drives check the direction of rotation of the motor. The motor must rotate in the direction of the arrow shown on the motor casing. If this is not the case, invert the two phase wires in the plug.
Warning: This operation is for an authorised, qualified electrician only.

The following controls are located on the unit base:


Illuminated pushbutton	(56)	**"ON – I"**
Pushbutton	(57)	**"OFF – 0"**
Mushroom-head button	(58)	**"EMERGENCY STOP"**

The motor protection consists of two limits switches (55 and 170) actuated by the kneader trough cover (119). These limit switches trigger and prevent the motor from starting when the kneader trough is opened or tipped up. This safety device is necessary because the torque inside the trough during a standard kneading cycle is high.


Do not disable or bypass this safety device.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com


400V 3/N/PE


230V 3/PE

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

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IKA® Works, Inc.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

North America

IKA® Works, Inc.
2635 North Chase Pkwy SE
Wilmington, NC 28405-7419
USA
TEL. +1 800 733-3037
TEL. +1 910 452-7059
FAX +1 910 452-7693
E-mail: usa@ika.net

IKA® Works (Asia) Sdn Bhd

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

Asia - Australia

IKA® Works (Asia) Sdn Bhd
No. 17 & 19, Jalan PJU 3/50
Sunway Damansara Technology Park
47810 Petaling Jaya
Selangor, Malaysia
TEL. +60 3 7804-3322
FAX +60 3 7804-8940
E-mail: sales@ika.com.my

IKA® Japan Y.K.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

Japan

IKA® Japan Y.K.
293-1 Kobayashi-cho
Yamato Koriyama Shi
639-1026 Japan
TEL. +81 74358-4611
FAX +81 74358-4612
E-mail: japan@ika.de

IKA® Works do Brasil Ltda.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

South America

IKA® Works do Brasil Ltda.
Estrada do Guerenguê, 491
Taquara Jacarepagua, RJ
Rio de Janeiro
22713-000 Brasil
TEL. +55 21 2435-9600
FAX +55 21 2435-9601
E-mail: brasil@ika.net

IKA® Works Guangzhou

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

China

IKA® Works Guangzhou
173-175 Friendship Road
Guangzhou Economic & Technological
Development Zone
Guangzhou 510730, P.R.CHINA
TEL. +86 20 8222-6772
FAX +86 20 8222-6776
E-mail: sales@ikagz.com.cn

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Artisan Technology Group is an independent supplier of quality pre-owned equipment

Gold-standard solutions

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

Learn more!

Visit us at **artisantg.com** for more info on price quotes, drivers, technical specifications, manuals, and documentation.

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

We're here to make your life easier. How can we help you today?

(217) 352-9330 | sales@artisantg.com | artisantg.com

ARTISAN® TECHNOLOGY GROUP