ARTISAN®
TECHNOLOGY GROUP

Your **definitive** source for quality pre-owned equipment.

Artisan Technology Group
(217) 352-9330 | sales@artisantg.com | artisantg.com

Full-service, independent repair center
with experienced engineers and technicians on staff.

We buy your excess, underutilized, and idle equipment
along with credit for buybacks and trade-ins.

Custom engineering
so your equipment works exactly as you specify.

- Critical and expedited services
- In stock / Ready-to-ship
- Leasing / Rentals / Demos
- ITAR-certified secure asset solutions

Expert team | Trust guarantee | 100% satisfaction

All trademarks, brand names, and brands appearing herein are the property of their respective owners.

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CE-DECLARATION OF CONFORMITY ***GB***

We declare under our sole responsibility that this product corrosponds to the regulations 89/336/EG; 92/31/EG; 94/9/EG and 73/23/EG and conforms with the standards or standardized documents EN 61010-1 (IEC/CEI 61010-1); EN 61326-1 (IEC 61326-1) EN 50021 and EN 60079-10.

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Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Contens

	Page
Safety instructions	6
Correct use	6
Unpacking	6
Commissioning	6
Maintenance and cleaning	7
Guarantee	7
Technical data	8

Safety instructions

This instrument is designed for use in Explosion Hazard Zone 2 areas. Do not use in Explosion Hazard Zone 1 areas.

If you are working in explosion hazard environments, do not use plastic beakers. Plastic beakers can create static electrical charges and can ignite explosive substances. The power supply plug must be connected up outside the explosion hazard zone. The housing must be earthed. To earth the instrument, use a earthing cable with a commercial 6.3 quick-connect terminal and plug this into the socket on the rear of the instrument (see figure).

Damage to the protective membrane sheet will reduce the IP and explosion hazard zone rating of the instrument.

Hazard warning – Magnetism. The instrument creates a magnetic field. Items which are sensitive to magnetism (e.g. data storage devices, pacemakers...). must be kept away from the instrument.

Correct use

The **IKA® lab disc** is a magnetic stirrer designed for mixing up to 800 ml of thinly fluid media with magnetic bars up to 30 mm long. Stirring speed can be adjusted in the range from 15 to 1500 rpm.

Unpacking

Please unpack the device carefully and check it for damage. It is important that any transport damage is detected when the device is unpacked. If necessary, an inventory of the damage should be made immediately (postal service, railway, haulage company).
As supplied, the package contains:
The **IKA® lab disc**, a plug-in power supply unit and an instruction manual.

Commissioning

Plug the instrument into a power supply socket. (For the power supply specifications, see the section "Technical Data").

Place the magnet bar and the medium to be mixed in a container which is transparent to magnetic fields. Set the stirrer speed to minimum. To do this, switch on the instrument by pressing the ⏻ button and then press the ⊖ button until the flashing LED between the minus and plus buttons stops flashing and lights permanently. You can now place the container on the mark in the centre of the set-up surface.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Now press the ⊖ button and the ⊕ button to set the speed you require. The LED between the buttons will flash during setting.

To select the reverse mode and invert the direction of stirring at regular intervals, press the button. In the reverse mode the stirring direction will be inverted every 30 seconds; there will be pause of 5 seconds between each changeover. The LED above the button will flash when the reverse mode is active.

The speed setting and reverse mode selection will remain unchanged when the instrument is unplugged.
Caution: At some speeds, irregularities in the shape of the stirring bar and of the bottom of the container can cause the bar to move irregularly during stirring. Increase or decrease the speed until the stirring bar begins to rotate evenly.

Maintenance and cleaning

The **IKA® lab disc** is maintenance-free. It is subjekt only to the natural wear and tear of components and their statistical failure rate. Use only cleansing agents which have been approved by IKA to clean IKA devices.
To remove use:

Dyes	isopropyl alcohol
Constructions materials	water containing tenside / isopropyl alcohol
Cosmetics	water containing tenside / isopropyl alcohol
Foodstuffs	water containing tenside
Fuels	water containing tenside

For materials which are not listed, please request information from IKA. Wear the proper protective gloves during cleaning of the devices. Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.
Before using another then the recommended method for cleaning or decontamination, the user must ascertain with the manufacturer that this method does not destroy the instrument.

Guarantee

You have purchased an original IKA laboratory machine which meets the highest engineering and quality standards.

In accordance with IKA guarantee conditions, the guarantee period is 24 months. For claims under the guarantee please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Technical data

Plug-in power supply unit:		
Rated voltage:	**V AC**	100 - 240
Rated frequency:	**Hz**	50 - 60
Output voltage:	**VDC**	12
Absorbed power:	**W**	5
Output power:	**W**	min. 3
Magnetic stirrer:		coil drive 4-pole
Speed range:	**1/min**	15 - 1500
Stirred quantity (max.):	**ml**	500 (at 1500 rpm with 25mm stirring bar)
Operating time:	**%**	100
Protection class:		II
Pollution degree:		2
Overvoltage category:		II
EN 60529 protection class:		IP 65
Ambient humidity (relative)	**%**	80
Ambient temperature:	**°C**	5 to 40 (in the operation)
Transport temperature:	**°C**	-25 to +70
Storage temperature:	**°C**	-10 to +70
Material: Housing		PA GF TPU
Front membrane sheet		Polyester
Weight:	**kg**	0,3
Installed weight (max.):	**kg**	2,0
Set-up surface:	**mm**	ø 90
Dimensions: (W x D x H)	**mm**	114 x 161 x 12
Explosion hazard class:		Ⓔx II 3G EEx nC IIB T6

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

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IKA®-Works, Inc.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

North America

IKA® Works, Inc.
2635 North Chase Pkwy SE
Wilmington, NC 28405-7419
USA
TEL. +1 800 733-307
TEL. +1 910 452-7059
FAX +1 910 452-7693
E-mail: usa@ika.net

IKA®-Works, (Asia) Sdn Bhd

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

Asia - Australia

IKA® Works (Asia) Sdn Bhd
No. 17 & 19, Jalan PJU 3/50
Sunway Damansara Technology Park
47810 Petaling Jaya
Selangor, Malaysia
TEL. +60 3 7804-3322
FAX +60 3 7804-8940
E-mail: sales@ika.com.my

IKA® Japan Y.K.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

Japan

IKA® Japan Y.K.
293-1 Kobayashi-cho
Yamato Koriyama Shi
639-1026 Japan
TEL. +81 74358-4611
FAX +81 74358-4612
E-mail: japan@ika.de

IKA® Works do Brasil Ltda.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

South America

IKA® Works do Brasil Ltda.
Estrado do Guerenguê, 491
Taquara Jacarepagua, RJ
Rio de Janeiro
22713-000 Brasil
TEL. +55 21 2435-09600
FAX +55 21 2435-09601
E-mail: ika@ika.com.br

IKA® Works Guangzhou

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

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Artisan Technology Group is an independent supplier of quality pre-owned equipment

Gold-standard solutions

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

Learn more!

Visit us at **artisantg.com** for more info on price quotes, drivers, technical specifications, manuals, and documentation.

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

We're here to make your life easier. How can we help you today?

(217) 352-9330 | sales@artisantg.com | artisantg.com

ARTISAN® TECHNOLOGY GROUP