ARTISAN®
TECHNOLOGY GROUP

Your **definitive** source for quality pre-owned equipment.

Artisan Technology Group
(217) 352-9330 | sales@artisantg.com | artisantg.com

Full-service, independent repair center
with experienced engineers and technicians on staff.

We buy your excess, underutilized, and idle equipment
along with credit for buybacks and trade-ins.

Custom engineering
so your equipment works exactly as you specify.

- Critical and expedited services
- In stock / Ready-to-ship
- Leasing / Rentals / Demos
- ITAR-certified secure asset solutions

Expert team | Trust guarantee | 100% satisfaction

All trademarks, brand names, and brands appearing herein are the property of their respective owners.

Find the *IKA MS 3* at our website: ***Click HERE***

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

CE-DECLARATION OF CONFORMITY *EN*

We declare under our sole responsibility that this product corresponds to the regulations 73/23/EEC and 89/336/EEC and conforms with the standards or standardized documents DIN EN IEC 61010-1 and DIN EN IEC 61326-1.

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Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Contents

EN

	Page
CE-Declaration of conformity	2
Safety instructions	9
Unpack	10
Correct use	10
Commissioning	10
Accessories	12
Maintenance	13
Technical data	14
Warranty	14

Safety instructions

- **Read the operating instructions in full before starting up and follow the safety instructions.**
- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- The voltage stated on the nameplate must correspond to the mains voltage.
- The device must only be operated with the original plug-in power supply unit.
- Wear your personal protective equipment in accordance with the hazard category of the medium to be processed. Otherwise there is a risk of:
 - splashing liquids
 - projectile parts
 - body parts, hair, clothing and jewellery getting caught.
- Set up the appliance in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- The feet of the appliance must be clean and undamaged.
- Position the knob (C) at the left stop before starting up. Gradually increase the speed.
- Keep hand/arm vibration to a minimum (Directive 2002/44/EC) when operated manually:
 - Keep the sample container as vertical as possible.
 - Only exert the minimum pressure necessary for the mixing process.
 - Use attachments for multiple samples when there is a large number of samples.
- If the shaking movement of the appliance increases (resonance), reduce the speed or pass through the critical phase as quickly as possible.
- Reduce the speed if
 - the medium splashes out of the vessel because the speed is too high
 - the appliance is not running smoothly
 - the appliance begins to move around because of dynamic forces.
- Firmly secure the accessories and vessels in place, otherwise shaking vessels could be damaged or projected out.
- Check the appliance and accessories beforehand for damage each time you use them. Do not use damaged components.
- Position one single shaking vessel in the centre and several shaking vessels so that they are evenly spread out.
- Sharp edged vessels cause wear debris on the attachments.
- Beware of the risk of
 - flammable materials
 - glass breakage as a result of mechanical shaking power
- Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.
- **Do not** operate the appliance in explosive atmospheres, with hazardous substances or under water.
- Safe operation is only guaranteed with the accessories described in the "Accessories" chapter.
- Always disconnect the plug before fitting accessories.
- The appliance does not start up again automatically following a cut in the power supply.
- Protect the appliance and accessories from bumps and impacts.
- The appliance may heat up when in use.
- The appliance may only be opened by experts.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Unpack

- **Unpack**
 - Please unpack the device carefully
 - In the case of any damage a fact report must be set immediately (post, rail or forwarder)
- **Delivery scope**

MS 3 basic	**MS 3 digital**
- Shaker **MS 3 basic**	- Shaker **MS 3 digital**
- Power supply	- Power supply
- Standard attachment MS 3.1	- Standard attachment MS 3.1
- Universal attachment MS 3.3	- Universal attachment MS 3.3
- One-hand insert MS 1.21	- One-hand insert MS 1.21
- Operating instructions	- Microtiter attachment MS 3.4
	- Test tube insert MS 1.32
	- Operating instructions


Correct use

- **Use**
 - ▸ For mixing liquids
 - Touch mode for single test tubes
 - Continuous mode for single or several vessels
- **Range of use**
 - Laboratories
 - Schools
 - Pharmacies

This device is suitable for use in all areas except:
- Residential areas
- Areas that are connected directly to a low-voltage supply network that also supplies residential areas.

The protection for the user is not ensured no more, if the equipment with accessories is operated, which are supplied or are not recommended not by the manufacturer or if the equipment in not correct use is operated against the manufacturer default.





Commissioning

	Setting	Action	Display
		 Position the rotating knob at the left stop before starting up	
		 If these conditions are met, the unit is ready for service when the mains plug has been plugged in	
1	**Switching on**	 Mode A Operating mode "Touch operation" Speed upto 3000 rpm *(= Factory setting)*	 
	Change to continuous operation	 Mode A Operating mode "Continuous operation" Speed upto 1300 rpm	

Operating modes

	Modus A *(with speed limiter)*	Modus B *(without speed limiter)*
Touch mode *(with standard-attachment MS 3.1)*	max. 3000 rpm	max. 3000 rpm
Continuous mode *(with all attachments)*	max. 1300 rpm	max. 3000 rpm

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

	Setting	Action		Display
2	**Setting timer / Change timer setting** *Not available on MS 3 BASIC* **Change to minutes mode**		Seconds mode The value previously set will be displayed *(Factory setting=000)*	
			Countdown running: An audible signal will be given when the countdown is complete	
			Switch off the device	
			Switching on the device: Press and hold button "1" and also press the "Power" button	
			Minutes mode	
			Countdown running: An audible signal will be given when the countdown is complete	
3	**Speed adjustment**			
4	**Change mode A → B**		Switch off the device	
			Switching on the device: Press and hold button "Start/Stop" and also press the "Power" button	
			Mode B Operation modes "Touch" and "Continuous operation" Speed upto 3000 rpm	 
	B → A		Switch off the device and switch on again, see **1**	

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Accessories


Possible uses and permitted speed ranges of attachments:

Designation	Description	Id.- No.	Touch mode	Continuous mode	Speed range rpm	Mode
 MS 3.1 Standard attachment	For test tubes and small vessels upto ø 50 mm	3426300	X - -	- X X	0/200-3000 0/200-1300 0/200-3000	A and B A B
 MS 3.3 Universal attachment	For different rubber foam inserts	3426600	-	X	0/200-1300	A
 MS 3.4 Microtiter attachment	For one microtiter plate	3426400	-	X	0/200-1300	A
 MS 1.21 One-hand attachment	For inserting into the universal attachment	L001540	-	X X	0/200-1300 0/200-3000	A B
 MS 1.31 Test tube insert	- For inserting into the universal attachment - For 14 test tubes ø 10 mm	L001840	-	X	0/200-1300	A
 MS 1.32 Test tube insert	- For inserting into the universal attachment - For 6 test tubes ø 12 mm	L001850	-	X	0/200-1300	A
 MS 1.33 Test tube insert	-For inserting into the universal attachment - For 4 test tubes ø 16 mm	L001860	-	X	0/200-1300	A
 MS 1.34 Test tube insert	- For inserting into the universal attachment - You can make holes as you like	L001830	-	X	0/200-1300	A


Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Changing attachments


MS 3.1

MS 3.3
MS 3.4

Using inserts

Maintenance

The appliance is maintenance-free.

Cleaning

Only use cleansing agents which have been recommended by **IKA**

Use to remove:

Dyes	isopropyl alcohol
Construction materials	water containing tenside/ isopropyl alcohol
Cosmetics	water containing tenside/ isopropyl alcohol
Foodstuffs	water containing tenside
Fuels	water containing tenside

- Do not allow moisture to get into the appliance when cleaning
- Wear protective gloves during cleaning the devices.
- Before using another than the recommended method for cleaning or decontamination, the user must ascertain with **IKA** that this method does not destroy the instrument.

Spare parts order

When ordering spare parts, please give:

- Machine type
- Manufacturing number, see type plate
- Item and designation of the spare part, see **www.ika.de**, spare parts diagram and spare parts list

Repair

In case of repair the device has to be cleaned and free from any materials which may constitute a health hazard.

If you require servicing, return the appliance in its original packaging. Storage packaging is not sufficient. Please also use suitable transport packaging.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Technical data

Power supply		
Input	V	100 - 240
	A	0,8
	Hz	50-60
Output	Vdc	24
		24W.LPS
		(Limited power source)
	Hz	50-60
Protection class		2 (double insulated) ⧈
Shaker		
Operating voltage	Vdc	24
	mA	800
Power consumption, normal operation	W	20
Power consumption, standby operation	W	2
Motor-output power	W	8
Drive		EC - motor
Speed range	rpm	0/200 – 3000
		infinitely adjustable
Speed adjustment		rotating knob in front
		of the appliance
Speed display		scale 0 – 3000
Agitation stroke	mm	4,5
Shaking motion		horizontal, circular
Perm duration of operation	%	100
Timer mode seconds		1 - 999 sec
mode minutes		1 - 999 min
Timer display		digital
Perm ambient temperature	°C	+5 bis +40
Perm. relative humidity	%	80
Protection type acc. to DIN EN 60529		IP 21
Contamination level		2
Overvoltage category		II
Operation at a terrestrial altitude	m	max. 2000
Dimensions total W x D x H without attachment	mm	148 x 205 x 63
Weight without supported load	kg	2,9
Max. supported load incl. attachment	kg	0,5

Subject to technical changes!

Warranty

In accordance with **IKA** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

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IKA® Works, Inc.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

North America

IKA® Works, Inc.
2635 North Chase Pkwy SE
Wilmington, NC 28405-7419
USA
TEL. +1 800 733-3037
TEL. +1 910 452-7059
FAX +1 910 452-7693
E-mail: usa@ika.net

IKA® Works, (Asia) Sdn Bhd

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

Asia - Australia

IKA® Works (Asia) Sdn Bhd
Lot 2 Jalan Indah 1/2
Taman Industri Rawang Indah
48000 Rawang
Selangor, Malaysia
TEL. +60 3 6093-3322
FAX +60 3 6093-3940
E-mail: ika@tm.net.my

IKA® Japan K.K.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

Japan

IKA® Japan K.K.
293-1 Kobayashi-cho
Yamato Koriyama Shi
639-1026 Japan
TEL. +81 74358-4611
FAX +81 74358-4612
E-mail: japan@ika.de

IKA® Works Guangzhou

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

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Artisan Technology Group is an independent supplier of quality pre-owned equipment

Gold-standard solutions

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

Learn more!

Visit us at **artisantg.com** for more info on price quotes, drivers, technical specifications, manuals, and documentation.

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

We're here to make your life easier. How can we help you today?

(217) 352-9330 | sales@artisantg.com | artisantg.com

ARTISAN® TECHNOLOGY GROUP