ARTISAN®
TECHNOLOGY GROUP

Your **definitive** source for quality pre-owned equipment.

Artisan Technology Group
(217) 352-9330 | sales@artisantg.com | artisantg.com

Full-service, independent repair center
with experienced engineers and technicians on staff.

We buy your excess, underutilized, and idle equipment
along with credit for buybacks and trade-ins.

Custom engineering
so your equipment works exactly as you specify.

- Critical and expedited services
- In stock / Ready-to-ship
- Leasing / Rentals / Demos
- ITAR-certified secure asset solutions

Expert team | Trust guarantee | 100% satisfaction

All trademarks, brand names, and brands appearing herein are the property of their respective owners.

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
PET042004

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Contens

1 Guarantee 27
2 Quality assurance 28
3 Safety instructions 28
4 Unpacking and scope of delivery 28
5 Operating principle 29
6 Correct use 29
6.1 Use 29
6.2 Care and Maintenance 30
7 Description and material specification of IKA®-*PET* pipette 30
8 Operating the IKA®-*PET* pipette 30
8.1 Volume setting 30
8.2 Pipetting or Dosing a Test Volume in a Weighing Vessel using the IKA®-*PET* pipette 31
8.3 Steps in the Testing Process for Air Displacement Pipettes / Checking the Measurement Value 33
8.4 Fitting the Pipette Tips 33
8.5 Recalibrating the IKA®-*PET* Pipette 33
9 Assembling and Disassembling the Pipette 35
10 Cleaning / Maintenance 37
11 Sterilization 37
12 Accessories 37
12.1 IKA-Pipette Tips 37
12.2 Pipette Stands IKA®-*PET* Rack 1 or 3 38
13 Sources of error, Error Detection and Error Elimination 39
13.1 Air Displacement Piston Stroke Pipettes 39
13.2 Error Detection and Error Elimination 39
14 Special Notes on Checking and Calibrating the Pipette 41
14.1 Test room 41
14.2 Minimum requirement of measuring instruments 41
14.2.1 Balance resolution / accuracy 41
14.2.2 Timer 42
14.2.3 Thermometer 42
14.2.4 Hygrometer 42
14.2.5 Barometer 42
14.3 Weighing procedure 42
14.4 Test Conditions During the Weighing Procedure 42
14.5 Test liquid 42
14.6 Vaporization 42
14.7 Test volume 43
14.8 Temperature differences 43
14.9 Leakage test 43
14.10 Duration of the Test Cycle 43
14.11 Number of measurements per test 43
14.12 Margins of error 44
14.12.1 Margins of error of IKA®-*PET* pipettes 44
14.13 Evaluation 45
14.13.1 Calculating the Systematic Measurement Deviation in % and µl (accuracy) 45
14.13.2 Calculating the Systematic Measurement Deviation in % and µl (coefficient of variation 45
14.13.3 Correction Factor Z for Distilled Water 46
14.14 Readjusting for Liquids with Densities different to that of Twice-Destilled, Degassed Water 47
14.15 Measurement data login 47
15 Used Symbols and Units 48
16 Spare parts ordering 48

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

DECLARATION OF CONFORMITY ***EN***

We hereby certify with sole responsibility that this product (the air displacement piston stroke pipette) conforms to the requirements stipulated in standard EN ISO 8655 ff and DIN 12600.

The exact conformity certification process is described in the standard. By affixing the conformity symbol, we, the manufacturer, certify that the piston stroke air displacement devices described in these operating instructions meet the requirements of the relevant standards.

1 Guarantee

You have purchased an original IKA laboratory machine which meets the highest engineering and quality standards.

In accordance with IKA sales and delivery conditions, the product is guaranteed against material or manufacturing defects for a period of 24 months. Please contact your dealer for any guarantee claims. You may also send the decontamination machine direct to our works, enclosing the delivery invoice and giving reasons for the claim.
You will be liable for freight costs.
The guarantee does not extend to wearing parts, nor to defects that are the result of improper use, insufficient care and maintenance or failure to follow the instructions in this operating manual.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

2 Quality assurance

IKA is a manufacturer of certified volumetric measuring devices and as such its quality assurance system and product quality are monitored by the German Weights and Measures Inspectorate. This ensures that all the requirements regarding inspection equipment monitoring and their realignment with the national measurement standards, as well as the skills of the staff are fulfilled. The stated margins of error only apply to use of original IKA accessories which are suitable for the respective measuring process (e.g. IKA pipette tips).

We recommend that the IKA pipettes be checked every 3 to 12 months. However, the frequency can be chosen to suit your individual requirements. They should be checked more frequently if the frequency of usage increases or if exposed to caustic media. For regular checks in accordance with ISO 9000 and as required by the GLP guidelines, you can also use the IKA calibration service or the calibration software. The device sent to us for calibration will be returned to you complete with its test report within a few days. You can obtain more detailed information about this process from your dealer.

3 Safety instructions

To ensure safe and problem-free dosing of liquids using this pipette, please read the operating instructions carefully.

In order to be able to use the appliance properly and safely, every user must first read the operating instructions and observe the safety instructions contained therein. Take care of these operating instructions and keep them in a place where they can be accessed by everyone.

Follow the safety precautions and health and safety guidelines specified for your type of work.
Covers or parts that are removed from the device must later be refitted to ensure safe operation. This will prevent the infiltration of foreign objects, liquids and other contaminants. Check the nominal volume after replacing parts that could affect the nominal volume.
When working with the pipette, users must wear the suitable personal protective equipment specified for the hazard class of the substance to be pipetted. Defective or inadequate personal protective equipment can expose the user to the risk of personal injury from caustic or contaminated media. Contaminated pipettes must not be sent (passed on) to third parties. Even though the pipette has been developed in accordance with the latest ergonomic findings -as is generally known in the field of liquid handling - complaints may arise when the pipette is in permanent uniform motion.

4 Unpacking and scope of delivery

Please unpack the equipment carefully and check for any damages. It is important that any damages which may have arisen during transport are ascertained when unpacking. If applicable a fact report must be set immediately (post, rail or forwarder).
The delivery scope covers:


- (1x) IKA®-*PET* pipette (with specified nominal volume)
- (1x) Sample pack of pipette tips (3 units)
- (1x) Test log
- (1x) Calibration and assembly tool
- (1x) Tube of silicone grease
- (1x) Operating instructions
- (1x) IKA pipette calibration software

PET042004

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

5 Operating principle

A plastic or glass tip is fitted to the piston stroke pipette. With the piston at the lower stop, the tip is immersed in the liquid that is to be measured and dosed. Retracting the piston draws in the liquid. Pressing down on the piston expels the volume of the liquid to be dosed. Piston stroke air displacement pipettes have an additional air volume that can be used to expel the last drops of liquid (blow-out function).
The advantage of pipetting with an air displacement pipette is that no wetting of the devices occurs as the liquid only enters the tip. The tip is only used once. This ensures completely residue-free work, which is especially important for laboratory work carried out in a sterile or contamination-free environment.

6 Correct use

An air displacement piston stroke pipette is used for collecting or dispensing defined volumes of liquids (liquid handling).
It is manually operated and mechanically controlled. The single-channel piston stroke pipette has a lubricated piston and can dose the volume specified by the manufacturer. In addition, all IKA® *PET* pipettes are fitted with a pipette tip ejector and blow-out function. For safe, straightforward work using the IKA® *PET* pipette, the following points must be observed.

6.1 Use

- The pipettes may only be used when a pipette tip is fitted.
- Before operating, you should always check for leakages (see section 14.9).
- Do not store the pipette with a filled pipette tip.
- Moisten the pipette tip before use
- Keep the pipette upright when immersing it in a liquid.
- Press down on the piston button smoothly and not too quickly, then slowly release it.
- Wipe off the pipette tip on the side of the vessel.
- There should be approximately 3 mm of test liquid in the weighing vessel.
- Avoid sudden piston movements and non-uniform rhythm and time intervals when pipetting.
- Avoid temperature differences between pipettes and pipette tips, as well as the liquid used. Temperature differences can cause mistakes in the volumes dispensed. This can also happen when using liquids with a high vapor pressure.
- When using infected, radioactive, poisonous, or other

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

solutions which present a health hazard, observe the relevant safety regulations.


- When using organic solvents and caustic chemicals, check their compatibility with the pipette tips (polypropylene-PP) and parts of the pipette housing.
- When using solutions that differ greatly from water in terms of their physical behavior (i.e. in terms of density, vapor pressure and wetting behavior, for example: glycerine), check the dispensed volume according to section 8.5 Recalibration and/or 14.14 Readjustment.

6.2 Care and Maintenance

- Do not allow any liquids to enter the pipette.
- Do not clean the pipette with acetone or caustic solutions.
- Use only original spare parts and accessories (pipette tips).
- When not in use, store the pipette in an upright position.
- If your appliance requires servicing, return it in its original packaging. Storage packaging is not sufficient. Also use appropriate transport packaging.

7 Description and material specification of IKA®-*PET* pipette

Name	Material
a.) Piston button	>PP<
b.) Moulding ring	>PP<
c.) Tip ejector	>PP<
d.) Body (hand grip)	>PP<
e.) Tip ejector sleeve	>PP<
f.) Suction tube (pipette cone)	>PP<


8 Operating the IKA®-*PET* pipette

8.1 Volume setting

With the IKA®-*PET* pipette, it is not possible to vary the nominal volume. This type of pipette is a fixed volume air displacement piston stroke pipette.
If the pipette volume is incorrect, you can, of course, readjust it using the tool provided (see section 8.5).

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

8.2 Pipetting or Dosing a Test Volume in a Weighing Vessel using the IKA®-*PET* Pipette


1	2	3	4	5
Hold the pipette almost upright. Press the piston button down to the first stop.	Vertically immerse the tip a maximum of 2 to 3 mm in the liquid*	Allow the piston button to slide back steadily until all of the liquid is drawn into the tip.	Wait approximately 1 to 2 seconds.	Pull the tip out of the liquid vertically

* The depth to which you immerse the tip can have a decisive effect on experimental results. If the tip is immersed too deep, droplets will form on the outside of the tip and will then be expelled together with the sample. If the tip is not immersed deep enough, this will cause an eddy effect, preventing the pipette from drawing in the desired volume.


Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

6


Lay the tip at a 10° to 45° angle on the inner wall of the weighing vessel.

7


Press the piston button down slowly and steadily to the first stop and wait at least 1 second. Then press the piston button past the first stop (blow-out function) to expel the remaining liquid.

8


9


Hold the piston button down while you extract the pipette tip from the vessel along its inner wall Now release the piston button.

10

Press the tip ejector button to eject the tip.
Always discard used tips. The residual liquid could damage the pipette if stored improperly

11


Use a pipette stand for storing your pipette

Additional note to 10

Caution:
Used tips contain liquid residue. Always follow the necessary safety precautions and regulations when disposing of the disposable parts. To avoid contact with contaminated tips, hold the pipette over a waste container and press the tip ejector push-button.


Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

8.3 Steps in the Testing Process for Air Displacement Pipettes / Checking the Measurement Value

8.4 Fitting the Pipette Tips

The IKA®-*PET* pipette will only function when a pipette tip is attached. The liquid is drawn into the tip.
Insert the suction tube (pipette cone) firmly into the tip by turning slightly to ensure a good seal between the pipette tip and the cone of the piston stroke pipette (see fig.).
The tip immersion depth in the liquid should be in the range of 2 to 3 mm. Use IKA filter tips to prevent sample migration and contamination of the pipette and user as far as possible. When pipetting liquids with wetting properties that differ from those of water, please note the recommendation in section 14.14 The sterile IKA filter tips available are free of detectable nucleic acids (RNA, DNA) and RNAses, DNAses and proteases.


8.5 Recalibrating the IKA®-*PET* Pipette

If your pipette was disassembled for repair or autoclaving or if the margins of error of your pipette are no longer correct, you will need to recalibrate the pipette. Before adjusting the pipette, you should check that all the pipette parts work properly. Before readjusting a pipette, please see the points in sections 6.1, 13 and 14.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

After correctly assembling the pipette and adhering to the testing and test room conditions specified in the standard and highlighted in these operating instructions, you can set the pipette piston stroke using the tool provided in the kit contents. First, pull out the piston button to remove it.

Turn the adjustment screw counter-clockwise to increase the stroke. Turn the screw clockwise, to decrease the stroke and volume to be adjusted.
The following table shows how much the stroke or volume changes per turn for a given pipette type.
As a guide you can use the markings on the moulding ring. The "Preliminary measurement" option in the IKA calibration software menu also offers specific help regarding the necessary direction of rotation of the adjustment screw and thus the correct stroke setting.

Pipette model	Stroke adjustment at a rotation **[mm]**	Volume adjustment at a rotation approx. **[µl]**
IKA®-PET 25µl–1ch	0.75	2.35
IKA®-PET 50µl–1ch	0.75	5.30
IKA®-PET 100µl–1ch	0.75	8.50


After you have set the pipette stroke, put the piston button back on the pipette and carry out the testing and evaluation as described in sections 8.2, 8.3 and 14.
If the margins of error for your pipette are still incorrect, you will have to change the stroke again using the adjustment screw and retest the pipette.
High quality pipettes which undergo regular maintenance do not have to be readjusted very often. Nevertheless, the accuracy (margins of error) of all pipettes should be tested regularly for safety reasons (recommended every 3-12 months).

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

9 Assembling and Disassembling the Pipette

The following figures show how to dismantle the pipette using the tool provided. The pipette is assembled as shown in the reverse order.
Only reassemble appliances that have dried out properly. Faulty parts (e.g. O-rings, piston etc.) must be replaced.
After cleaning the pipette, the piston must be slightly lubricated again (use the silicone grease provided).


Removing the piston/sealing unit (see section 17):

PET042004

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Removing the adjustment screw and stroke limiter (see section 17):

PET042004

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

10 Cleaning / Maintenance


As with any precision instrument, the outside of the pipette should be cleaned daily. In particular, care should be taken to ensure that the tip cone is clean and undamaged. We recommend you regularly clean and decontaminate all pipette parts. Only use IKA-approved cleaning agents when cleaning your pipette.
Use water containing detergent for cleaning
(soap solution) or 75% ethanol or 60% isopropanol. Then rinse the parts with distilled water and if necessary dry off with a soft, lint-free cloth, or allow them to air dry.
Wear protective gloves when cleaning the pipette. Heavy soiling due to in-seepage of a medium can only be removed after disassembling (see section 9).

11 Sterilization

For sterilization, only the bottom half of the pipette's tip ejector sleeve can be vapor-autoclaved (121°C, 1bar, 20min) (see section 9 for information on how to remove). The autoclaved parts must be dried afterwards at room temperature. Only reassemble the pipette parts after they have completely dried and cooled.

CAUTION: decontamination is the responsibility of the user. Piston stroke devices that have not been decontaminated must not be sent or passed on to third parties.
After dismantling the pipette you must check and, if necessary, recalibrate it to ensure accuracy.

Note: sterilization by autoclaving is often not enough. For example, this does not kill RNAses.

12 Accessories

Please only use our recommended original accessories. You should be aware that the precision, accuracy and life span of our devices can be adversely affected if you do not use recommended spare parts and disposables. We accept no responsibility for damage caused by the failure to use recommended original spare parts and material. Such damage is not covered by guarantee.


12.1 IKA-Pipette Tips:

IKA pipette tips are affordable, high-quality, universal, precision pipette tips designed to meet the most exacting requirements of research and routine laboratories.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Table 1: IKA-recommended Pipette Tips with Universal Fit :

Order.-Nr.	Name	Pack	Steril / Color
33 307 00	IKA-TIP 200µl - bag	1000	no / natural
33 307 01	IKA-TIP 200µl – box	2x1000	no / natural
33 307 02	IKA-TIP 200µl – 96x-rack	10x96	no / natural
33 307 03	IKA-TIP 200µl – 96x-rack	10x96	yes / natural
33 308 00	IKA-TIP Filter 1-100µl - bag	1000	no / natural
33 308 01	IKA-TIP Filter 1-100µl - bag	1000	yes / natural
33 309 00	IKA-TIP Filter 1-200µl - bag	1000	no / natural
33 309 01	IKA-TIP Filter 1-200µl – bag	1000	yes / natural
33 310 00	IKA-TIP Filter 100µl – 96x-rack	10x96	yes / natural
33 311 00	IKA-TIP Filter 200µl – 96x-rack	10x96	yes / natural


The conformity-certified pipette tips recommended by IKA can be autoclaved at 121°C for 15 minutes at 1 bar.
The margins of error (see section 13) always apply to the entire system, consisting of piston stroke pipette and a pipette tip. The supplier's declaration of conformity or certification does not apply if pipette tips other than those recommended by IKA are used. In this case, the test must be carried out adhering to the margins of error (using the alternate tips) in accordance with the test procedures specified in ISO 8655-6.
Plastic tips for piston stroke air displacement pipettes are designed for disposal after one single use. They must not be cleaned for re-use. Disposable tips lose their measuring characteristics after a single use.
The snap-on disposable tips are made of uncolored polypropylene. The tips are available with or without filters.

12.2 Pipette Stands IKA®-*Pet* Rack 1 or 3

To prevent corrosion, contamination and damage to your pipette, you should not place it on the work bench. It could come into contact with chemicals on the bench or get damaged by falling off.

For routine work, use IKA® *Pet* rack 1 or IKA® *Pet* rack 3 pipette stands

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

13 Sources of error, Error Detection and Error Eliminationg

13.1 Air Displacement Piston Stroke Pipettes

If problems occur, please read the following information in Table 2. It will help you to identify and eliminate the errors. The sources of error for air displacement piston stroke pipettes are also described in detail in standard EN ISO 8655-2 Piston-Operated Volumetric Apparatus in Table B.1 (influencing values, resulting measurement deviations and counter measures).

13.2 Error Detection and Error Elimination

Table 2: List of causes of errors and their elimination

Error	Causes of erreo	Error elimination
Droplets remain on the pipette tip	• Uneven wetting of the plastic walls	• Put on a new pipette tip
	• Pipette tips with poor wetting properties were used	• Use original IKA tips
Dripping tip; the dosed volume is incorrect	• Tip loose	• Press firmly on the tip
	• Poorly fitting pipette tip used	• Use original IKA tips
	• Tip ejector incoorectly assembled or damaged	• Correctly assemble the tip ejector or change it
	• Piston sealing system leaks	• Remove the piston and the seal from the bottom part of the pipette (see section 9). Change the seal and, if necessary, the piston.
	• A liquid with a vapor pressure different from that of water was pipetted	• In a case like this, we recommend pipetting with a positive displacement pipette
	• A liquid with a density different from that of water was pipetted	• Put on a new pipette tip
	• The viskosity and/or flow behavior and wetting characteristics of the liquid different from those of water	• Use original IKA tips
	• Measuring apparatus not clean enough	• Clean pipette again

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Error	Causes of erreo	Error elimination
Push button sticks, does not run smoothly	• Pipette incorrectly assembled • Piston dirty • Seal dirty • Solvent vapors have penetrated	• Reassemble • Clean the piston and lubricate it slightly (see section Assembling / Disassembling) • Remove the piston and seal from the bottom part of the pipette. (see section 9). Rinse bottom part toroughly with warm water. Rinse again with distilled water and allow the dry. Replace the seal if necessary. • Remove the piston from the bottom of the pipette (see section 9) and allow the bottom to vent. Clean the piston and lubricate it slightly.
Pipette has no suction	• The pipette is incorrectly assembled • Bottom part is not firmly screwed onto sealing system	• Reassemble • Tighten bottom part
Pipette clogged, not enough liquid drawn in	• Liquid has penetrated the pipette suction tube and dried there	• Remove the piston from the bottom part of the pipette (see section 9). Rinse bottom with warm water. Rinse again with distilled water and allow to dry.
Irregular stroke action	• Push button pressed and released too quickly	• Check your own pipetting technique
Temperature differences	• Pipette, tip and test liquid have different temperatures	• The pipette to be tested and the tip should be left in the test room for approximately 2 hours before the measurement

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

14 Special Notes on Checking and Calibrating the Pipette

Note: the pipette calibration software developed by IKA allows you to automatically record the weight value and the test room condition when using balances, thermometers, hygrometers and barometers with an RS232 interface. The software also allows you to carry out evaluations according to DIN EN ISO 8655, record manufacturer-independent master data, monitor the pipette at scheduled times and much more.

14.1 Test room

Testing must be conducted in a draught-free room at constant conditions. The test room must have a relative humidity RH of over 50% and a constant temperature of between 20°C and 25°C (±0.5°C). Before testing, the device to be tested and the test liquid (water) must be stored for a sufficient amount of time in the test room (at least 2 hours) in order to adjust to the room conditions.

14.2 Minimum requirement of measuring instruments:

14.2.1 Balance resolution/accuracy

The laboratory or analytical balance or equivalent weighing appliance required for calibration should at least fulfil the requirements specified in Table 3. The scale division value of the balances used for testing should correspond to the selected volume of the device to be tested.

Table 3: Minimum requirement of balances

Selected volume of the device to be tested V	Display resolution [mg]	Repeatable precision and linearity [mg]	Standard measurement uncertainty [µl]
1 µl **V** 10 µl	0,001	0,002	0,002
10 µl **V** 100 µl	0,01	0,02	0,02
100 µl **V** 1000 µl	0,1	0,2	0,2
1 ml **V** 10 ml	0,1	0,2	0,2
10 ml **V** 200 ml	1	2	3

If the standard measurement uncertainty of the balance is known (from the balance calibration certificate, for example), it can be used instead of the repeatable precision and linearity. The standard measurement uncertainty should not be greater than 2 to 3 times the resolution. This ensures that the spread within the series of measurements for estimating systematic and chance measurement deviation as per DIN EN ISO 8655 is recorded with adequate precision.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

The manufacturer's specifications are within these standard limiting values and can, therefore, require a higher degree of balance precision. For test volumes < 50µl electrostatic charges must be prevented from building up on all components that come into direct or indirect contact with the pipette. We recommend that you use vessels and humidity traps with glass covers to prevent static charge as far as possible.
The photometric measuring procedure can also be used for volumes <10µl if necessary.

14.2.2 Timer

The timer can have a maximum standard measurement uncertainty of 1 second (1).

14.2.3 Thermometer

The thermometer for ambient temperature and test liquid temperature can have a maximum standard measurement uncertainty of 0.2°C (1)..

14.2.4 Hygrometer

The hygrometer for relative ambient humidity can have a maximum standard measurement uncertainty of 10%(1).

14.2.5 Barometer

The barometer can have a maximum standard measurement uncertainty of 0.5kPa(1).

NOTE(1): all measurement uncertainties are stated for a coverage factor k=1.

14.3 Weighing procedure

IKA pipettes are tested for their ability to expel (Ex) test liquids into the weighing vessel.

14.4 Test Conditions During the Weighing Procedure

Before beginning and after concluding the weighing procedure, the temperatures of the test liquid in the storage container must be documented at an accuracy of 0.2°C. The air pressure in the test room must be recorded at an accuracy of 1kPa and the relative humidity RH at an accuracy of 10%.

14.5 Test liquid

Quality 3 distilled or de-ionized water as per ISO 3696, which is degassed or in equilibrium with the air. The test liquid must be at the same temperature as the test room.

14.6 Vaporization

For small volumes in particular (under 50 µl), errors caused by test liquid vaporization during weighing must be taken into consideration. The shape of the weighing vessel and the duration of the test cycle are also important. In order to keep vaporization errors to a minimum, the following additional points can be observed when testing volumes under 50 µl.

- Balance with vaporization trap
- Dose the test liquid into a capillary tube (**note:** this is not customary practice).
- There should be approximately 3 mm of test liquid in the weighing vessel.
- When testing devices with small volumes in particular, we recommend a weighing vessel height to diameter ratio of at least 3:1. Otherwise use weighing vessels with covers.

In addition to the above, errors due to vaporization can be determined through experimentation in a series of tests and then compensated for through calculations. The compensation uncertainty should be added to the measurement uncertainty.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

14.7 Test volume

For devices with fixed volumes, the test volume is the nominal volume

14.8 Temperature differences

Before testing the device to be tested, the
pipette tip and the test liquid must be left for an adequate period of time (at least 2 hours) in the test room to allow them to adjust to the room conditions. It is essential that the pipettes and the measuring area are not exposed to direct sunlight or other influences.

14.9 Leakage test

Before working with the pipette, following maintenance or repair work and before calibration, you should check that it does not contain any leaks as follows:

- Attach a pipette tip of the correct size to the shaft
- Draw in the max. volume of distilled water (see section 14.5) from a vessel
- Hold the pipette upright and wait 20 seconds
- If a water droplet appears on the end of the tip, the pipette seal is faulty. (see section 13)
- If no water droplet forms, immerse the tip in the vessel filled with water.
- The liquid level in the tip must not change. If the liquid level drops, the pipette has a faulty seal (see section 13)
- To facilitate seal testing, use a colored liquid.

14.10 Duration of the Test Cycle

The duration of the test cycle must be kept to a minimum. It should not last longer than 60 seconds if possible. It is important for the duration to be the same both during each cycle and, if possible, from cycle to cycle, so that a reliable calculated compensation of the error due to vaporization can be conducted during the series of measurements.

14.11 Number of measurements per test

When using the gravimetric test procedure, ten measurements must be conducted per test volume. The systematic and chance measurement deviations are calculated from these measurements as per DIN ISO 8655-6.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

14.12 Margins of error

14.12.1 Margins of error of IKA®-*PET* pipettes

Table 3: Margins of error of IKA®-PET Pipettes

Pipette model	Nominal volume	Accuracy of systematic measuring deviation e_s in [%]	Precision of change measuring deviation CV in [%]
IKA®-*PET* 25µl - 1ch	25 µl	± 1,0 %	0,3 %
IKA®-*PET* 50µl - 1ch	50 µl	± 0,7 %	0,3 %
IKA®-*PET* 100µl - 1ch	100 µl	± 0,6 %	0,2 %
Calculation of eS in [%]	see Calculation section 14.13	$\mathbf{e_s}=100 \bullet \frac{(V - V_S)}{V_S}$	
Calculation of CV in [%]	see Calculation section 14.13		$\mathbf{CV}=100 \bullet \frac{\sqrt{\frac{\sum_{i=1}^{n} (V_i - V)^2}{n - 1}}}{\overline{V}}$

The margins of error can only be guaranteed when using original IKA pipette tips. The test conditions must be as specified in EN ISO 8655-6 Piston-Operated Volumetric Apparatus, tests must use precision balances with vaporization protection tested by a standards and calibration certification authority.
The basis for the adjustment (Ex) are 10 measured values; degassed, twice-distilled water, test room constant at 20-25°C±0.5°C; with pre-wetted pipette tip, input into the weighing vessel (vessel wall); at a volume of 10µl adjustment (IN) i.e. output from the weighing vessel (cause: vaporization)

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

14.13 Evaluation

(For symbols and unit see section 15)

14.13.1 Calculating the Systematic Measurement Deviation in % and µl (accuracy)

The systematic measurement deviation e_S of the piston stroke device is to be calculated in µl as in equation **(1)** and in % as in equation **(2)**

$$\mathbf{e_S} = \overline{\mathbf{V}} - \mathbf{V_S} \quad \text{in [µl]} \qquad \mathbf{(1)}$$

$$\mathbf{e_S} = 100 \bullet \frac{(V - V_S)}{V_S} \quad \text{in [%]} \qquad \mathbf{(2)}$$

with

- the target volume $\mathbf{V_S}$ with fixed volume pipettes the test volume Vs and the nominal volume V_O are identical
- the correction factor **Z** in $\frac{\mu l}{mg}$ see DIN EN ISO 8655-6 Table A.1
- the single volume $\mathbf{V_i = m_i \bullet Z}$ in [µl] **(3)**
- Number of measurements **n** (at least 10)
- the mean value of the single volume $\overline{\mathbf{V}} = \frac{1}{10} \bullet \sum_{i=1}^{n} \mathbf{V_i}$ in [µl] **(4)**

14.13.2 Calculating the Systematic Measurement Deviation in % and µl (coefficient of variation)

The chance measurement deviation s is to be calculated as a repeatable standard deviation in µl as in equation **(5)** or in % as a coefficient of variation CV as in equation **(6)** to n.

$$\mathbf{s} = \sqrt{\frac{\sum_{i=1}^{n} (V_i - \overline{V})^2}{n - 1}} \quad \text{in [µl]} \qquad \mathbf{(5)}$$

$$CV = 100 \bullet \frac{\sqrt{\frac{\sum_{i=1}^{n} (V_i - \overline{V})^2}{n - 1}}}{\overline{V}} \quad \text{in [%]} \qquad \mathbf{(6)}$$

with

- the correction factor **Z** in $\left[\frac{\mu l}{mg}\right]$ see DIN EN ISO 8655-6 Table A.1
- the single volume $\mathbf{V_i = m_i \bullet Z}$ in [µl] **(3)**
- Number of measurements **n** (at least 10)
- the mean value of the single volume $\overline{\mathbf{V}} = \frac{1}{10} \bullet \sum_{i=1}^{n} V_i$ in [µl] **(4)**

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

14.13.3 Correction Factor Z for Distilled Water

Values for the correction factor Z in [µl/mg] as per Table A.1 of standard DIN EN 8655-6, dependent on the test temperature and air pressure

Temperature in [°C]	Air pressure in [hPa]							
	800	850	900	900	950	1000	1013	1050
15,0	1,0017	1,0018	1,0019	1,0019	1,0019	1,0020	1,0020	1,0020
15,5	1,0018	1,0019	1,0019	1,0019	1,0020	1,0020	1,0020	1,0021
16,0	1,0019	1,0020	1,0020	1,0020	1,0021	1,0021	1,0021	1,0022
16,5	1,0020	1,0020	1,0021	1,0021	1,0021	1,0022	1,0022	1,0022
17,0	1,0021	1,0021	1,0022	1,0022	1,0022	1,0023	1,0023	1,0023
17,5	1,0022	1,0022	1,0023	1,0023	1,0023	1,0024	1,0024	1,0024
18,0	1,0022	1,0023	1,0023	1,0023	1,0024	1,0025	1,0025	1,0025
18,5	1,0023	1,0024	1,0024	1,0024	1,0025	1,0025	1,0026	1,0026
19,0	1,0024	1,0025	1,0025	1,0025	1,0026	1,0026	1,0027	1,0027
19,5	1,0025	1,0026	1,0026	1,0026	1,0027	1,0027	1,0028	1,0028
20,0	1,0026	1,0027	1,0027	1,0027	1,0028	1,0028	1,0029	1,0029
20,5	1,0027	1,0028	1,0028	1,0028	1,0029	1,0029	1,0030	1,0030
21,0	1,0028	1,0029	1,0029	1,0029	1,0030	1,0031	1,0031	1,0031
21,5	1,0030	1,0030	1,0031	1,0031	1,0031	1,0032	1,0032	1,0032
22,0	1,0031	1,0031	1,0032	1,0032	1,0032	1,0033	1,0033	1,0033
22,5	1,0032	1,0032	1,0033	1,0033	1,0033	1,0034	1,0034	1,0034
23,0	1,0033	1,0033	1,0034	1,0034	1,0034	1,0035	1,0035	1,0036
23,5	1,0034	1,0035	1,0035	1,0035	1,0036	1,0036	1,0036	1,0037
24,0	1,0035	1,0036	1,0036	1,0036	1,0037	1,0037	1,0038	1,0038
24,5	1,0037	1,0037	1,0038	1,0038	1,0038	1,0039	1,0039	1,0039
25,0	1,0038	1,0038	1,0039	1,0039	1,0039	1,0040	1,0040	1,0040
25,5	1,0039	1,0040	1,0040	1,0040	1,0041	1,0041	1,0041	1,0042
26,0	1,0040	1,0041	1,0041	1,0041	1,0042	1,0042	1,0043	1,0043
26,5	1,0042	1,0042	1,0043	1,0043	1,0043	1,0044	1,0044	1,0044
27,0	1,0043	1,0044	1,0044	1,0044	1,0045	1,0045	1,0045	1,0046
27,5	1,0045	1,0045	1,0046	1,0046	1,0046	1,0047	1,0047	1,0047
28,0	1,0046	1,0046	1,0047	1,0047	1,0047	1,0048	1,0048	1,0048
28,5	1,0047	1,0048	1,0048	1,0048	1,0049	1,0049	1,0050	1,0050
29,0	1,0049	1,0049	1,0050	1,0050	1,0050	1,0051	1,0051	1,0051
29,5	1,0050	1,0051	1,0051	1,0051	1,0052	1,0052	1,0052	1,0053
30,0	1,0052	1,0052	1,0053	1,0053	1,0053	1,0054	1,0054	1,0054

PET042004

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

14.14 Readjusting for Liquids with Densities different to that of Twice-Destilled, Degassed Water

It is possible to adjust the pipettes described above for a volume of a liquid with a different density to water in such a way that the volume displayed on the pipette corresponds to the pipetted volume. The adjustment process is comparable to the process described in sections 7.1, 9.5, 14 and 15. The only difference is that the mean value of the weightings is converted in microlitres according to the formula:

$$\text{pipette volume} = \frac{\text{mean value of the weights measured}}{\text{density of liquid weighed}}$$

The value calculated in this way is the actual value to which the pipette volume is set.

CAUTION:
The value set to this volume must be checked gravimetrically (DIN ISO 8655-6). The pipette set in this way only provides corresponding dosing volumes for this liquid. The setting may have to be revised and checked again.

After the pipette has been adjusted to another test liquid, the measuring appliance must be clearly marked, for example, using a label. The user must clearly state the name of the liquid to which the nominal volume applies on the outside of the pipette.

CAUTION:
The manufacturer's warranty expires following adjustment!

14.15 Measuremet data loggin

For electronic (digital) measurement data logging and management, we recommend the pipette calibration software developed by IKA in accordance with the latest ergonomic and technical findings.
The user-friendly software automatically records the measurements determined on a balance and performs the evaluation for both fixed volumes and variable pipettes according to DIN EN ISO 8655. The software also allows you to manage all the pipettes used in your laboratory with their respective margins of error and individual test cycles. The software will notify you automatically whenever a given pipette is due for calibration (according to your scheduled preferences). All the individual values from a measuring series are saved and assigned to a database for whichever pipette they were generated. This means that you can trace back measurements years later in accordance with ISO and GLP guidelines.

- Testing as per DIN EN ISO 8655
- Automatic measurement transfer for balances, hygrometers, thermometers and barometers using a RS 232 interface
- Automatic correction of the ambient conditions of air pressure and temperature
- Master data for more than 580 pipettes by various manufacturers is already saved
- Manufacturer-independent inventory of individual pipettes
- Calibration schedule monitoring
- Measured results saved in a database
- Report as per GLP guidelines with specifications of all the individual measured values, mean value, accuracy and precision.
- Can be run on all Microsoft Win32 operating systems (Win95/98/ME/NT/2000/XP etc.)

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

15 Used Symbols and Units

Physical size	Symbols	Characters
Individual mass	m_i	mg
Correction factor	Z	µl/mg
Individual volumes	V_i	µl
Mean value of individual volumes	$\overline{V}$	µl
Target volume (test volume)	V_S	µl
Systematic measurement deviation	e_S	µl
Systematic measurement deviation	e_S	%
Change measuring deviation Standard repeatable deviation	s	µl
Change measuring deviation Coefficient of variation	CV	%

16 Spare parts ordering

When ordering spare parts, please specify the manufacturer's number (on the device), the device type, the position number and the description of the spare part. Please only send devices in for repair that have been cleaned and are free of materials which might present health hazards. Please inform us if you have worked with hazardous or dangerous materials.

CAUTION: before returning a pipette, please ensure that it is completely free of chemical, biological and radioactive materials!

To maintenance, calibration or repair of pipettes you send these on:


GENEO BioProducts GmbH
Sillemstrasse 60A
D-20257 Hamburg
Germany
Tel.: +49 040 4019 7717
Fax: +49 040 4019 7719
Email: info@geneo.de
Internet:www.geneo.de

oder

MÜLLER **ratiotec** GmbH
Hildastrasse 17
D-77948 Friesenheim b. Lahr/Schwarzw.
Germany
Tel.: +49 07821 63704
Fax: +49 07821 67076
Email: info@mueller-ratiotec.de
Internet:www.mueller-ratiotec.de

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Assembly sequence
Service drawing

Assembly step	Pos.	Name	Material
1	1	Gehäuse / housing	PP
3	2	Spitzenabwurf / tips dropping	PP
21	3	Spitzenabwurfhülse / tips dropping bushing	PP
10	4	Zierring / ornamental ring	PP
11	5	Kappe / cap	PP
16	6	Abschlusskappe / end cap	PP
17	7	Kolbenführung / piston guide	POM
12	8	Saugrohr / suction pipe	PP
14	9	Dichtungshülse / sealing sleeve	1.4571
20	11	Kolben kmpl. / piston	1.4571/coatet
6	12	Justierschraube / adjusting screw	POM
8	13	Schubstange / push rod	1.4571
15	15	Gewindehülse / threaded sheath	POM
18	16	Gewindebuchse / threaded bushing	POM
8	17	Druckscheibe / thrust washer	1.4571
4	18	Sicherungsscheibe	DIN 6799-2.3/1.4571
2	19	Druckfeder / compression spring	1.4310
9	20	Druckfeder / compression spring	1.4310
19	21	Druckfeder / compression spring	1.4310
15	25	O-Ring	EPDM
13	26	O-Ring	EPDM

PET042004

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

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IKA® Works, Inc.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

North America

IKA® Works, Inc.
2635 North Chase Pkwy SE
Wilmington, NC 28405-7419
USA
TEL. +1 800 733-3037
TEL. +1 910 452-7059
FAX +1 910 452-7693
E-mail: usa@ika.net

IKA® Works, (Asia) Sdn Bhd

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

Asia - Australia

IKA® Works (Asia) Sdn Bhd
No. 17 & 19, Jalan PJU 3/50
Sunway Damansara Technology Park
47810 Petaling Jaya
Selangor, Malaysia
TEL. +60 3 7804-3322
FAX +60 3 7804-8940
E-mail: sales@ika.com.my

IKA® Japan Y.K.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

Japan

IKA® Japan Y.K.
293-1 Kobayashi-cho
Yamato Koriyama Shi
639-1026 Japan
TEL. +81 74358-4611
FAX +81 74358-4612
E-mail: japan@ika.de

IKA® Works do Brasil Ltda.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

South America

IKA® Works do Brasil Ltda.
Estrada do Guerenguê, 491
Taquara Jacarepagua, RJ
Rio de Janeiro
22713-000 Brasil
TEL. +55 21 2435-9600
FAX +55 21 2435-9601
E-mail: brasil@ika.net

IKA® Works Guangzhou

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

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Artisan Technology Group is an independent supplier of quality pre-owned equipment

Gold-standard solutions

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

Learn more!

Visit us at **artisantg.com** for more info on price quotes, drivers, technical specifications, manuals, and documentation.

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

We're here to make your life easier. How can we help you today?

(217) 352-9330 | sales@artisantg.com | artisantg.com

ARTISAN® TECHNOLOGY GROUP