ARTISAN®
TECHNOLOGY GROUP

Your **definitive** source for quality pre-owned equipment.

**Artisan Technology Group**
(217) 352-9330 | sales@artisantg.com | artisantg.com

**Full-service, independent repair center**
with experienced engineers and technicians on staff.

**We buy your excess, underutilized, and idle equipment**
along with credit for buybacks and trade-ins.

**Custom engineering**
so your equipment works exactly as you specify.


- Critical and expedited services
- In stock / Ready-to-ship
- Leasing / Rentals / Demos
- ITAR-certified secure asset solutions

**Expert team | Trust guarantee | 100% satisfaction**

All trademarks, brand names, and brands appearing herein are the property of their respective owners.

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***CE-DECLARATION OF CONFORMITY EN***

*We declare under our sole responsibility that this product corrosponds to the regulations 89/336/EEC; 92/31/EEC und 73/23/EEC and conforms with the standards or standardized documents EN 61010-1 (IEC/CEI 61010-1) and EN 61326-1 (IEC 61326-1).*

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## Warranty

You have purchased an laboratory machine which meets the highest engineering and quality standards.
In accordance with guarantee conditions, the guarantee period is 24 months. For claims under the guarantee please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.
The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

If you require servicing, return the appliance in its original packaging. Storage packaging is not sufficient. Please also use suitable transport packaging.

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# Contents

| | Page |
|---|---|
| Warranty | 2 |
| Safety instructions | 7 |
| Correct use | 7 |
| Unpacking | 8 |
| Commissioning | 8 |
| Function Stirring | 8 |
| Function Heating | 8 |
| Maintenance and cleaning | 9 |
| Associated standards and regulations | 9 |
| Accessories | 9 |
| Technical data | 9 |
| List of spare parts RO / RT power | 10 |
| Spare parts diagram RO / RT power | 15 |

# Safety instructions

In order to be able to use the appliance properly and safely, every user must first read the operating instructions and observe the safety instructions contained therein. Take care of these operating instructions and keep them in a place where they can be accessed by everyone.

Only staff who have been trained accordingly, know the appliance and are authorised to carry out work in this field should use this appliance.

The machine may only be opened by trained specialists - even during repairs. The machine is to be unplugged from the mains before opening.

When using PTFE-coated magnetic bars, please pay attention to the following: *Chemical reactions vof PTFE occur in contact with molten or dissolved alkaline and alkaline earth metals, as well as with fine-particled powders of metals of the 2. and 3. group of the periodical system at temperatures above 300 - 400°C. Only elementary fluorine, chlorine trifluoride and alkaline metals do attrack PTFE, halogen hydrocarbons have a reversibly swelling effect.*

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**Attention: Magnetism!** Effects of the magnetic field have to be taken into account (e.g. data carriers, cardiac, cardiac pacemakers...).

If the power cable is replaced, it must be by one of the same quality.

Please put the supplied transparent foil on the supporting surface. The foil prevents that the glas vessels don't slip to the side.

**RT:** Caution: Never heat up any medium with ignition point higher than 180°C.

The instrument may not be operated in rooms with explosion hazard.

Furthermore, it has to be placed upon a fire-proof and/or not burnable surface.

**RT:** Please make sure that the mains cable does not contact the heating plate!

**RT:** Exercise caution when touching the housing parts and the heating plate!

**Risk of burns!**

The heating plate can reach temperaturesin excess of 125°C.

**RT:** Please consider possible endangeringdue to release of toxid or flamable gaser caused by heated medium.

# Correct use

IKAMAG RO 5 power, IKAMAG RO 10 power and IKAMAG RO 15 power are laboratory devices. They are used in laboratories in the chemical industry, schools and pharmacies to name a few. The devices are suitable for simultaneous mixing of substances which are placed on mounting plates. The built-in stirrer makes it possible to mix 5, 10 or 15 stirring positions at the same time and with equal intensity. The mixing intensity depends on the motor speed and the size of the magnetic stirring bar. The model variants IKAMAG RT 5 power, IKAMAG RT 10 power and IKAMAG RT 15 power are also suitable for controlling temperature.

## Unpacking


Please unpack the equipment carefully and check for any damages. It is important that any damages which may have arisen during transport are ascertained when unpacking. If applicable a fact report must be set immediately (post, rail or forwarder).
The delivery scope covers: an IKAMAG RO 5, 10 or 15 power and operating instructions bzw.an IKAMAG RT 5, 10 or 15 power and operating instructions.

## Commisioning

Check whether the voltage specified on the type plate matches the mains voltage available. The power socket used must be earthed (protective earth conductor contact). If these conditions are met, the device is ready to operate after plugging in the mains plug. If these procedures are not followed, safe operation cannot be guaranteed and/or the equipment may be damaged.
Observe the ambient conditions (temperature, humidity, etc.) listed under Technical Data.

## Function Stirring


The device is switched ON/OFF by means of the appliance switch **(A)** located on the right. The LED for "Power" **(B)** becomes green when the switch is activated.
The speed of the motor is set on the control knob "Mot" **(C)** with the scale from 1 to 10 **(D)**, which roughly corresponds to a range from 0 to 1200 rpm. The maximum speed of 1200 rpm is obtained when the control knob is turned all the way to the right in no-load operation.
If the viscosity of the medium increases during processing, then the speed drops slightly. If the viscosity decreases, e.g. due to the addition of low-

viscosity media, then the speed rises slightly. In either case the control knob can be adjusted to compensate for this speed change, provided the required torque does not exceed the maximum torque.
Note that fluctuations in the mains voltage within the permissible tolerances will also result in slight variations in speed.

## Function Heating

The device is switched ON/OFF by means of the appliance switch **(A)** located on the right. The LED for "Power" **(B)** becomes green when the switch is activated.
The temperature of the heater plate can be adjusted on the „Temp" turn dial **(E)** with the corresponding scale from 1 to 10 **(F)**.
If energy is being directed to the heater plate, the „Temp" LED **(G)** lights up yellow. The scale values 1 through 10 correspond to temperatures on the surface from 0°C trough 120°C.

## Maintenance and cleaning

The IKAMAG RO 5 power, RO 10 power or RO 15 power and IKAMAG RT 5 power, RT 10 power or RT 15 power are maintenance-free. It is subject only to the natural wear and tear of components and their statistical failure rate.

**When ordering spare parts, please give the manufacturing number shown on the type plate, the machine type and the name of the spare part.**

Please send in equipment for repair only after it has been cleaned and is free from any materials which may constitute a health hazard. Use only cleansing agents which have been approved by IKA to clean IKA devices. To remove use:

| | |
|---|---|
| Dyes | isopropyl alcohol |
| Construction materials | water containing tenside / isopropyl alcohol |
| Cosmetics | water containing tenside / isopropyl alcohol |
| Foodstuffs | water containing tenside |
| Fuels | water containing tenside |

For materials which are not listed, please request information from IKA. Wear the proper protective gloves during cleaning of the devices. Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.

## Associated standards and regulations

**Construction in accordance with the foolwing safety standards**
EN 61 010-1 UL 3101-1 CAN/CSA C22.2 (1010-1)

**Construction in accordance with the foolwing EMC standards**
EN 50 082-1 EN 55 014 EN 60 555-2; -3

**Associated EU guidlines**
EMC-guidlines: 89/336/EWG
Machine guidelines: 73/023/EWG

## Accessories

**IKAFLON** stirrers, PTFE, 8 mm dia length 30mm
**RSE** Stirrer remover

**Recommended IKA accessories**

**EH 4 basic** suspension thermostat
**RO 5.1** bath vessel
**RO 10.1** bath vessel
**RO 15.1** bath vessel

## Technical data

| | | | |
|---|---|---|---|
| Design voltage: | **VAC** | 230±10% | |
| or | **VAC** | 115±10% | |
| Design frequency: | **Hz** | 50/60 | |
| Motor: | | ball-bearing brushless EC motor | |
| Speed range: | **rpm** | 0 to 1100 infinitely variable | |
| Deviation between the stirrer points | **%** | 5 | |
| Perm. ambient temperature: | **°C** | +5 to +40 | |
| Perm. relative humidity: | **%** | 80 | |
| Perm. ON duration: | **%** | 100 | |
| Interference suppression: | | to VDE 0871 | |
| Protection to DIN 40 050: | | IP 42 | |
| Mounting surface | **mm** | 120 x 440 | (RO 5) |
| | **mm** | 180 x 440 | (RO 10) |
| | **mm** | 270 x 440 | (RO 15) |
| **RT**: Max surface temperature | **°C** | 120 | |
| **RT**: Temperature constancy in the medium | **K** | ± 2 depending on the container | |

| | | | |
|---|---|---|---|
| Max. stirring volume (water) | | 200 ml per stirring position in 250 ml Erlenmeyer flask, for other vessels | |
| | | 2ltr | (RO 5) |
| | | 4ltr | (RO 10) |
| | | 6ltr | (RO 15) |
| Stirring points | | 5 | (RO 5) |
| | | 10 | (RO 10) |
| | | 15 | (RO 15) |
| Stirring point distance | **mm** | 90 | (RO 5) |
| | **mm** | 90x90 | (RO 10) |
| | **mm** | 90x90 | (RO 15) |
| Max. power consumption: | **W** | 10 | (RO 5) |
| | **W** | 20 | (RO 10) |
| | **W** | 30 | (RO 15) |
| Max. power consumption: | **W** | 185 | (RT 5) |
| | **W** | 395 | (RT 10) |
| | **W** | 600 | (RT 15) |
| Dimensions: (B x T x H ) | **mm** | 122 x 552 x 65 | (RO 5) |
| | **mm** | 182 x 552 x 65 | (RO 10) |
| | **mm** | 272 x 552 x 65 | (RO 15) |
| | **mm** | 138 x 552 x 65 | (RT 5) |
| | **mm** | 198 x 552 x 65 | (RT 10) |
| | **mm** | 288 x 552 x 65 | (RT 15) |
| Weight: | **kg** | 2,3 | (RO 5) |
| | **kg** | 3,2 | (RO 10) |
| | **kg** | 4,7 | (RO 15) |
| | **kg** | 3,0 | (RT 5) |
| | **kg** | 4,2 | (RT 10) |
| | **kg** | 6,0 | (RT 15) |

## List of spare parts

**RO 5 / 10 / 15 power**
**RT 5 / 10 / 15 power**

| Item | Designation |
|---|---|
| 1 | Spacer |
| 6 | Deep-groove ball bearing |
| 8 | Protective cap |
| 9 | Plug |
| 12 | Clamping bush |
| 13 | Motor retainer |
| 14 | Potentiometer shaft |
| 15 | Rocker switch |
| 16 | BLP-Motor |
| 17 | EC - motor |
| 18 | Front plate |
| 19 | Hood |
| 21 | Bearing pins |
| 22 | Reflector |
| 25 | Miniature fuse |
| 27 | Locking ring |
| 28 | Pedestal |
| 29 | Driving belt |
| 30 | BLP Motor (small) |
| 40 | Cable clamp |
| 41 | Cable binder |
| 42 | Reference plate |
| 52 | Snap cable feed-through |
| 53 | Heater hood |
| 2001 | Magnet holder bearing pedestal (complete) |
| 2002 | Magnet holder motor (complete) |

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## IKA® Works, Inc.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

North America

**IKA**® Works, Inc.
2635 North Chase Pkwy SE
Wilmington, NC 28405-7419
USA
TEL. +1 800 733-3037
TEL. +1 910 452-7059
FAX +1 910 452-7693
E-mail: usa@ika.net

## IKA® Works, (Asia) Sdn Bhd

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

Asia - Australia

**IKA**® Works (Asia) Sdn Bhd
No. 17 & 19, Jalan PJU 3/50
Sunway Damansara Technology Park
47810 Petaling Jaya
Selangor, Malaysia
TEL. +60 3 7804-3322
FAX +60 3 7804-8940
E-mail: sales@ika.com.my

## IKA® Japan Y.K.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

Japan

**IKA**® Japan Y.K.
293-1 Kobayashi-cho
Yamato Koriyama Shi
639-1026 Japan
TEL. +81 74358-4611
FAX +81 74358-4612
E-mail: japan@ika.de

## IKA® Works Guangzhou

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

China

**IKA**® Works Guangzhou
173-175 Friendship Road
Guangzhou Economic & Technological
Development Zone
Guangzhou 510730, P.R.CHINA
TEL. +86 20 8222-6772
FAX +86 20 8222-6776
E-mail: sales@ikagz.com.cn

**Artisan Technology Group** is an independent supplier of quality pre-owned equipment

## Gold-standard solutions

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

## We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

## Learn more!

Visit us at **artisantg.com** for more info on price quotes, drivers, technical specifications, manuals, and documentation.

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

We're here to make your life easier. How can we help you today?

(217) 352-9330 | sales@artisantg.com | artisantg.com

ARTISAN® TECHNOLOGY GROUP