ARTISAN®
TECHNOLOGY GROUP

Your **definitive** source for quality pre-owned equipment.

**Artisan Technology Group**
(217) 352-9330 | sales@artisantg.com | artisantg.com

**Full-service, independent repair center**
with experienced engineers and technicians on staff.

**We buy your excess, underutilized, and idle equipment**
along with credit for buybacks and trade-ins.

**Custom engineering**
so your equipment works exactly as you specify.

- Critical and expedited services
- In stock / Ready-to-ship
- Leasing / Rentals / Demos
- ITAR-certified secure asset solutions

**Expert team | Trust guarantee | 100% satisfaction**

All trademarks, brand names, and brands appearing herein are the property of their respective owners.

Find the *IKA RV 10* at our website: ***Click HERE***

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**CE-DECLARATION OF CONFORMITY** **EN**
We declare under our sole responsibility that this product corresponds to the regulations 89/336/EC and 73/23/EC and conforms with the standards or standardized documents DIN EN IEC 61010-1, -2-010 and DIN EN IEC 61326-1.

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Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

# Contents

EN

| | Page |
|---|---|
| CE-Declaration of conformity | 3 |
| Safety instructions | 12 |
| Correct use | 13 |
| Unpacking | 13 |
| Useful information | 13 |
| Setting up | 14 |
| Drive RV 10 basic/ digital | |
| Heating bath | |
| Glassware | |
| Hose system | |
| Interfaces and outputs | 15 |
| Commissioning | 16 |
| Maintenance and cleaning | 17 |
| Accessories | 17 |
| Error codes | 17 |
| Technical data | 18 |
| Solvent table (excerpt) | 18 |
| Warranty | 18 |
| Spare parts list | 33 |
| Spare parts diagram | 34 |

# Safety instructions

For your protection

- **Read the operating instructions in full before starting up and follow the safety instructions.**
- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations. **When working under a vacuum in particular!**
- Wear your personal protective equipment in accordance with the hazard category of the medium to be processed. Otherwise there is a risk of:
  - splashing liquids,
  - body parts, hair, clothing and jewellery getting caught,
  - injury as a result of glass breakage.
- **CAUTION!** Inhalation of or contact with media such as poisonous liquids, gases, spray mist, vapours, dusts or biological and microbiological materials can be hazardous to user.
- Set up the device in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- Ensure that there is sufficient space above the device as the glass assembly may exceed the height of the device.
- Prior to each use, always check the device, accessories and especially the glass parts for damage. Do not use damaged components.
- Ensure that the glass assembly is tension-free! Danger of cracking as a result of:
  - stress due to incorrect assembly,
  - external mechanical hazards,
  - local temperature peaks.
- Ensure that the stand does not start to move due to vibrations respectively unbelance.
- Beware of hazards due to:
  - flammable materials,
  - combustible media with a low boiling temperature,
  - glass breakage.
- CAUTION! Only process and heat up media that has a flash point higher than the adjusted safe temperature limit of the heating bath that has been set.
  The safe temperature limit of the heating bath must always be set to at least 25 °C lower than the fire point of the media used.
- Do **not** operate the appliance in explosive atmospheres, with hazardous substances or under water.
- Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.
- Tasks with the device must only be performed when operation is monitored.
- Operation with excess pressure is not permitted.
- Do not cover the ventilation slots of the device in order to ensure adequate cooling of the drive.
- There may be electrostatic discharges between the medium and the drive which could pose a direct danger.
- The appliance is not suitable for manual operation.
- Safe operation is only guaranteed with the accessories described in the "Accessories" chapter.
- Refer to the operating instructions for the HB 10 basic/ digital heating bath.
- Refer to the operating instructions for the accessories.
- Use the HB 10.2 cover or the HB 10.1 spray guard.
- Only use the device under an all side-closed exhaust, or a comparable protective device.
- Adapt the quantity and the type of distilland to the size of the distillation equipment. The cooler must work properly. Monitor the coolant flow rate at the cooler outlet.
- The glass equipment must always be ventilated when working under normal pressure (e.g. open outlet at cooler) in order to prevent a pressure build-up.
- Please note that dangerous concentrations of gases, vapours or particulate matter can escape through the outlet at the cooler. Take appropriate action to avoid this risk, for example, downstream cold traps, gas wash bottles or an effective extraction system.
- Evacuated glass vessels must not be heated on one side; the evaporating flask must rotate during the heating phase.
- The glassware is designed for operation under a vacuum of up to 10 mbar. The equipment must be evacuated prior to heating (see chapter "Commissioning"). The equipment must only be aired again after cooling. When carrying out vacuum distillation, uncondensed vapours must be condensed out or safely dissipated. If there is a risk that the distillation residue could disintegrate in the presence of oxygen, only inert gas must be admitted for stress relief.
- **CAUTION!** Avoid peroxide formation. Organic peroxides can accumulate in distillation and exhaust residues and explode while decomposing!
  Keep liquids that tend to form organic peroxides away from light, in particular from UV rays and check them prior to distillation and exhaust for the presence of peroxides. Any existing peroxides must be eliminated. Many organic compounds are prone to the formation of peroxides e.g. dekalin, diethyl ether, dioxane, tetrahydrofuran, as well as unsaturated hydrocarbons, such as tetralin, diene, cumene and aldehydes, ketones and solutions of these substances.

- **DANGER OF BURNING!** The heating bath, tempering medium, evaporator piston and glass assembly can become hot during operation and remain so for a long time afterwards! Let the components cool off before continuing work with the device.
- **ATTENTION!** Avoid delayed boiling! Never heat the evaporating flask in the heating bath without switching on the rotary drive!Sudden foaming or exhaust gases indicate that flask content is beginning to decompose. Switch off heating immediately. Use the lifting mechanism to lift the evaporator piston out of the heating bath. Evacuate the danger zone and warn those in the surrounding area!
- When the device is switched off or the power supply disconnected, the internal safety lift removes the evaporator piston from the heating bath. The maximum load (evaporator piston and contents) is 1.5 kg for a permissible safety lift with no power.
- **CAUTION!** Never operate the device when the evaporator piston is rotating and the lift is raised. Always lower the evaporator piston into the heating bath first before starting the rotation drive. Otherwise hot tempering medium may be sprayed out!
- Set the speed of the drive so no tempering medium is sprayed out as a result of the evaporator piston rotating in the heating bath. If necessary reduce the speed.
- Do not touch rotating parts during operation.
- Imbalance may result in uncontrolled resonance behaviour of the device or assembly. Glass apparatus may be damaged or destroyed. In the event of unbalance or unusual noises, switch off the appliance immediately or reduce the speed.
- The appliance does not start up again automatically following a cut in the power supply.

062008

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

- The device is only disconnected from the power supply network if the device power switch is off or the plug is pulled out.
- The socket for the mains cord must be easily accessible.

For protection of the equipment

- The voltage stated on the type plate must correspond to the mains voltage.
- Socket must be earthed (protective ground contact).
- Removable parts must be refitted to the appliance to prevent the infiltration of foreign objects, liquids etc.
- Protect the appliance and accessories from bumps and impacts.
- The appliance may only be opened by experts.

## Correct use

**• Use**

▸ Together with the accessories recommended by **IKA®** , the device is suitable for:

- quick and gentle distillation of liquids
- evaporation of solutions and suspensions
- crystallisation, synthesis or cleaning of fine chemicals
- drying of powder and granulate material
- recycling of solvents

Mode of operation: Tabletop device

**• Range of use**

- Laboratories
- Pharmacies
- Schools
- Universities

The protection for the user is not ensured no more, if the equipment with accessories is operated, which are supplied or are not recommended not by the manufacturer or if the equipment in not correct use is operated against the manufacturer default.

## Unpacking

**• Unpacking**

- Please unpack the device carefully
- In the case of any damage a fact report must be set immediately (post, rail or forwarder)

**• Delivery scope**

see table

| | Drive RV 10 basic | Drive RV 10 digital | Heating bath HB 10 basic | Heating bath HB 10 digital | Vertical glassware RV 10.1 | Vertical glassware RV 10.10, coated | Diagonal Glassware RV 10.2 | Diagonal Glassware RV 10.20, coated | Woulff bottle | Bracket | Cooler locking device, compl. | Operating instructions |
|---|---|---|---|---|---|---|---|---|---|---|---|---|
| RV 10 basic V | x | | x | | x | | | | x | x | x | x |
| RV 10 basic VC | x | | x | | | x | | | x | x | x | x |
| RV 10 basic D | x | | x | | | | x | | x | x | x | x |
| RV 10 basic DC | x | | x | | | | | x | x | x | x | x |
| RV 10 digital V | | x | | x | x | | | | x | x | x | x |
| RV 10 digital VC | | x | | x | | x | | | x | x | x | x |
| RV 10 digital D | | x | | x | | | x | | x | x | x | x |
| RV 10 digital DC | | x | | x | | | | x | x | x | x | x |

## Useful information

Distillation is a thermal separating process for liquid compounds based on substance-specific, pressure-dependent boiling points through evaporation and subsequent condensation.
The boiling point temperature decreases with decreasing external pressure which means that work is usually done under reduced pressure. In this way the heating bath can be maintained at a constant temperature (e.g. 60 °C). Using the vacuum, the boiling point is set with a steam temperature of approx. 40 °C. The cooling water for the condensation cooler should not be warmer than 20 °C (60-40-20 rule).

A chemical resistant membrane pump with a vacuum controller should be used to create the vacuum. The pump is protected from solvent residue by the addition of a Woulff bottle.
Working with a jet pump to create a vacuum can only be recommended to a limited extent as the solvents may contaminate the environment when using these systems

Speed, temperature, piston size and system pressure all affect the evaporator air capacity.
The optimum capacity of the flow-through condenser is approx. 60%.

This corresponds to condensation on approx. 2/3 of the cooling coil. With larger capacities there is the risk that the uncondensed solvent vapour will be extracted.

The device is equipped with a piston-operated safety mechanism. If the power cuts out, the evaporator piston is automatically lifted out of the heating bath by an integrated gas spring.

**CAUTION!** The maximum load (evaporator piston and contents) is 1.5 kg for a permissible safety lift with no power for glass apparatus with diagonal or vertical condenser.
When using other types of condensers such as dry ice or intensive condensers as well as when using return distillation distributors with slip-on condensers, the max. load can be reduced by the added weight of the glass apparatus. Thus, prior to distillation, check whether the lift goes up without power when laden with the glass and distillation material.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

## Drive RV 10 basic/ digital

**CAUTION!** Loosen transportation lock (fig. 4a)!

☞ Hold the lift with your hand on the height position and remove the thumb screw on the back of the appliance (A).

➢ Once the transportation lock has been removed, the lift moves slowly to its upper end position. The distance is approx. 140 mm.



Fig. 4a

☞ Connect the device to the power supply (B) using the supplied power cord (B).

Adjustable base (Fig. 4b)

Fig. 4b


Fig. 5

Fig. 6

☞ Attach the bracket for the Woulff bottle to the left side of the lift (fig. 5).

☞ Insert the bottle and attach the supplied hose connectors to the bottle (fig. 6).

☞ Remove the clamping device for the angle setting of the rotation drive on the right side of the lift by rotating the knurled screw counter clockwise (by lightly pressing and turning at the same time, the knurled screw comes out farther).

☞ Set the drive at an angle of approx. 30° (fig. 7).


Fig. 7

☞ Then secure the rotation drive from being accidentally turned by tightening the knurled screw in a clockwise direction.

## Heating bath HB 10 basic / digital

**CAUTION!** Refer to the chapter entitled "Commissioning" in the heating bath instruction manual!

☞ Place the heating bath on the stand of the rotation drive and push it into the left position (fig. 8).

Fig. 8

## Glassware

☞ Open the locking device on the drive head by turning it 60° counter clockwise (D), (fig. 9).

☞ Feed the steam pipe in until it stops.

☞ Then lock the locking device by turning it clockwise by 60° (fig. 10).

☞ Check the axial locking device on the steam pipe.



Fig. 9

Fig. 10

☞ Place the RV 06.15 or RV 06.13 condenser seal in the condenser receptacle and fit the glassware to the device according to the assembly instructions (Fig. 11 a and b).

Fig. 11a

Prior to commissioning - Mounting the seal RV 06.15

☞ Insert vapour tube (1).

☞ Insert seal RV 06.15 (2).

☞ Mount the cooler (3a).

☞ Tighten the cap nut hand-screwed. Start the device at 120 rpm for 20 minutes. Then tighten the cap nut firmly (3b).


Fig. 11b

Push-off mechanism

☞ The plastic screw nut (E) helps loosen tight-fitting piston ground-in connections (Fig. 9).

☞ To do so, hold the tight-fitting evaporator piston and turn the plastic screw nut (E) until the evaporator piston neck.

*Note: Prior to commissioning, hand-tighten the plastic screw nut (E) left-aligned.*

Assembling the vertical glassware cooler locking device

☞ Assemble the condenser locking device according to the diagram (fig. 12).

- Mount the plate using the thumb screw.
- Put the support rod on the plate and attach it with the screw nut.
- Attach the rubber protector.
- Fasten the Velcro band to the support rod.

☞ Secure the vertical glassware with the Velcro band.

Fig. 12

## Hose system

☞ Connect the water hoses to the condenser according to the counter-flow principle (fig. 13).

☞ - Install the vacuum connections to the condenser, Woulff bottle, vacuum controller with valve and vacuum pump.

☞ Always connect the vacuum hose to the condenser at the high-test point to minimise solvent losses during suctioning.


Fig.13

☞ Use typical laboratory vacuum hoses with an inside diameter of 8 mm and a wall thickness of 5 mm (see accessories).

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Data transfer

Fig.14

**HB 10 digital**
The heating bath transfers data via IR interfaces. These are located on the left display side of the heating bath or on the right side of the drive unit. Do not place any objects between the two operating units as otherwise the data transfer may be interrupted (fig. 14)!

**RV10 digital**
The RS 232 interface at the back of the device is fitted with a 9-pole SUB-D jack and can be connected to a PC. The pins have serial signals. Using the laboratory device software "labworldsoft", the device can be operated in "Remote" mode.
Note: Please note the system requirements as well as the operating instructions and help section included with the software.

Serial interface RS 232 (V24)

Configuration
- The functions of the interface lines between the device and the automation system are a selection from the signals specified in the EIA standard RS232 C, as per DIN 66 020 Part 1.
- For the electrical characteristics of the interface and the allocation of signal status, standard RS 232 C applies in accordance with DIN 66 259 part 1.
- Transmission procedure: asynchronous character transmission in start-stop mode.
- Type of transmission: full duplex.
- Character format: character representation in accordance with data format in DIN 66 022 for start-stop mode. 1 start bit; 7 character bits; 1 parity bit (even); 1 stop bit.
- Transmission speed: 9600 bit/s.
- Data flow control: none
- Access procedure: data transfer from the device to the computer takes place only at the computer's request.

Command syntax and format

The following applies to the command set:
- Commands are generally sent from the computer (Master) to the device (Slave).
- The stirrer machine sends only at the computer's request. Even fault indications cannot be sent spontaneously from the stirrer machine to the computer (automation system).
- Commands are transmitted in capital letters.
- Commands and parameters including successive parameters are separated by at least one space (Code: hex 0x20).
- Each individual command (incl. parameters and data) and each response are terminated with Blank CR Blank LF (Code: hex 0x20 hex 0x0d hex 0x20 hex 0x0A) and have a maximum length of 80 characters.
- The decimal separator in a number is a dot (Code: hex 0x2E).

The above details correspond as far as possible to the recommendations of the NAMUR working party (NAMUR recommendations for the design of electrical plug connections for analogue and digital signal transmission on individual items of laboratory control equipment, rev. 1.1).

The NAMUR commands and the additional specific **IKA**® commands serve only as low level commands for communication between the rotary evaporator and the PC. With a suitable terminal or communications programme these commands can be transmitted directly to the rotary evaporator. The **IKA**® software package, labworld*soft*, provides a convenient tool for controlling rotary evaporator and collecting data under MS Windows, and includes graphical entry features, for motor speed ramps for example.

The following table summarises the (NAMUR) commands understood by the **IKA**® control equipment.

Abbreviations used

| | | | |
|---|---|---|---|
| X,y | = | | Numbering parameter (integer) |
| m | = | | Variable value, integer |
| n | = | | Variable value, floating-point number |
| X | = | 4 | Speed |
| X | = | 60 | Interval time (1-99 seconds, 1 <= m >=99) |
| X | = | 61 | Timer (1-199 minutes, 1 <= m >=199) |
| X | = | 62 | Direction of lift (m=2-> lift up , m=1-> lift down) |

| NAMUR Commands | Function |
|---|---|
| IN_NAME | Request designation |
| IN_PV_X<br>X = 4 | Read actual value |
| IN_SOFTWARE | Request software Id-number, date and version |
| IN_SP_X<br>X = 4 | Read target value input |
| OUT_SP_X m<br>X = 1,60,61,62 | Set target value to m |
| RESET | Switch to normal operation |
| START_X<br>X = 4,60,61,62 | Switch on appliance (remote) function |
| STATUS | Status output<br>0:Manual operation without interruption<br>1:automatic operation Start (without interruption)<br>ERROR z (z error number see table) |
| STOP_X<br>X = 1,60, 61, 62 | Switch off appliance function. The variable set with OUT_SP_X remains. |

**PC 1.1 Cable (Fig. 15)**
This cable is required to connect the 9-pin connector to a PC.


Fig. 15

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

The unit is ready for service when the mains plug has been plugged in.

On/Off switch (Fig. 16)
☞ Switch on the device on the right side of the device.
➢ Device functions activated.

Fig.16

**Function keys** ***(see page 4)***

**• "Power" key (1)**
☞ Press the "Power" key once
➢ "Stand-by" mode
☞ Press the "Power" key again
➢ Device functions reactivated
"Software version" display changes to "Parameter" display.

*Note: The device must first be switched on at the on/off switch on the right side of the device.*

*Note: The device is not dead or deenergised until the on/off switch has been switched to off (right side of device)!*

**• Lift position key "▲" (2)**
☞ When you press the "▲"key, the lift goes up.
☞ When you release the "▲" key, the lift stops and remains in that position.

**• Lift position key "▼" (3)**
☞ When you press the "▼"key, the lift goes down.
☞ When you release the "▼" key, the lift stops and remains in that position.

**• Rotating knob (4)**
Setting the speed
☞ Turn the rotating knob
➢ Speed setting accuracy: ± 5 rpm,
Speed range 20-270 rpm
Note: When you select the speed > 100 rpm, smooth start is automatically activated.

Start / Stop Rotation drive
☞ Press the rotating knob
➢ Rotation drive starts
☞ Press the rotating knob again
➢ Rotation drive stops

**• Timer key (5)**
☞ Press the "Timer" key.
➢ The timer block appears on the display, the "TIMER" indicator flashes.
☞ Set the target value "Timer" on the knob by turning (1 to 199 minutes).
☞ Save the "Timer" target value by pressing the "Timer" key again.
☞ Starting the "Timer" function:
- Automatic when rotation drive is running,
- Start the rotation drive.

*Note: Once the timer has run out, there is a beep, the rotation finishes and the lift goes up.*

➢ The previously saved value appears on the display.

☞ Deactivate the timer by setting the target value to "0".

**• Int key (6)**
Interval cycle for right-left running (change of direction) of the evaporator piston in seconds.
☞ Press the "Int" key.
➢The interval block appears on the display, the "INT" indicator flashes.
☞ Set the target value "Int" on the knob by turning (1 to 60 minutes).
☞ Starting the "Int" function
- Automatic when rotation drive is running
- Start the rotation drive

*Note: In interval mode the maximum speed is limited to 200 rpm.*
☞ Save the "Int" target value by pressing the "Int" key again.
☞ Deactivate the interval by setting the target value to "0".

**Display** ***(see page 4)***

**• Display "REMOTE" (10)**
(only for RV 10 digital)
Remote control via PC in connection with software "labworld*soft*" is active.

**• Display "888 rpm" (20)**
Current target and actual speed indicator in revolutions per minute [rpm]. When the device is switched on the most recently selected target value "speed" appears.

**• Display ◯ (30)**
Rotation drive active.

**• Display "TIMER" (40)**
Time lapse control active.
Indication of remaining distillation time in minutes [min].
*Note: With the software "labworld*soft*", RV 10 digital can also switch off the heating bath once the timer has run out.*

**• Display "INT" (50)**
Interval mode activated.

**• Display "E01" (20)**
--- not pictured ---
Error codes, see chapter entitled "Error Codes".

**Setting the lower end stop**
**CAUTION!** Depending on the size of the piston, the setting angle of the rotation drive and the position of the heating bath and lift, the evaporator piston can be in contact with the bottom of the heating bath. Risk of glass breakage!
Limit the lower lift position with the variable end stop.

☞ Press the "▼" key until the lift reaches the desired position.
*Note: The evaporator piston should be 2/3 immersed in the heating bath.*
☞ To adjust the stop element, press the centre button on the front of the lift (fig. 17).
☞ Move the stop to the desired position (fig. 18).
☞ Press the "▲" key until the lift has reached the upper end stop.
*Note: The path is limited from 0 – 6 cm.*

Fig.17

Fig.18

**Check that the limit stop is functioning properly:**
☞ Lower the lift by holding down the "▼" key.
➢ The drive stops automatically when the preset desired lower end position has been reached.
☞ Move the drive back to the upper position.

**Filling the evaporator piston**
Prior to filling the evaporator piston, a vacuum controller is used to regulate the glass apparatus to the target pressure.
☞ Now fill the evaporator piston using the backfeed line.
➢ Due to the vacuum present, the solvent is suctioned into the evaporator piston. This enables you to keep solvent loss due to suctioning to a minimum.

☞ You can also fill the evaporator piston manually prior to creating the vacuum. The evaporator piston should not be filled more than half its volume.

*Note: **CAUTION!** The maximum allowed load (evaporator piston and contents) is 3 kg.*

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

**Setting up the heating bath**

**Refer also to the operating instructions for the HB 10 basic /digital heating bath!**

☞ Move the lift to the bottom position and check the position of the heating bath in relation to the evaporator piston. When using larger evaporator pistons (2 or 3 litres) or depending on the angle of the rotation drive, you can move the heating bath 50 mm to the right.

☞ Fill the heating bath with the tempering medium until the evaporator piston is surrounded by tempering medium to 2/3 of its volume.

☞ Switch on the rotation drive and slowly increase the speed.

*Note: Avoid creating waves.*

☞ Switch on the heating bath using the main on/off switch.

*Note: Avoid stress on the glass due to different evaporator piston and heating bath temperatures when lowering the evaporator piston into the heating bath!*

## Maintenance and cleaning

The device is maintenance-free. It is subject only to the natural wear and tear of components and their statistical failure rate.

For cleaning disconnect the main plug.

To clean the drive use only water with a detergent that contains tensides, or use isopropylalcohol for stubborn soiling.

Spare parts order

When ordering spare parts, please give:

- Machine type
- Manufacturing number, see type plate
- Item and designation of the spare part, see spare parts diagram and spare parts list

Repair

In case of repair the device has to be cleaned and free from any materials which may constitute a health hazard.

If you require servicing, return the appliance in its original packaging. Storage packaging is not sufficient. Please also use suitable transport packaging.

## Accessories

| | |
|---|---|
| IKAVAC VC2 | Vacuum controller |
| VC1.1 | Jet pump |
| VC1.3 | Electric vacuum valve |
| VC2.4 | Pump control |
| RV 10.1 | Vertical glassware |
| RV 10.10 | Vertical glassware ,coated |
| RV 10.2 | Diagonal glassware |
| RV 10.20 | Diagonal glassware, coated |
| RV 10.3 | Intensive condenser glassware |
| RV 10.30 | Intensive condenser glassware, coated |
| RV 10.4 | Dry ice condenser glassware |
| RV 10.40 | Dry ice condenser glassware, coated |
| RV 10.5 | Reflux distillation glassware |
| RV 10.50 | Reflux distillation glassware, coated |
| RV 10.6 | Reflux distillation glassware, intensive condenser |
| RV 10.60 | Reflux distillation glassware, intensive condenser, coated |
| HB 10.1 | Spray guard (for heating bath HB10 basic/ digital/ control) |
| HB 10.2 | Cover (for heating bath HB10 basic/ digital/ control) |
| NS 29/32 | Evaporator piston 100 ml |
| NS 29/32 | Evaporator piston 1000 ml |
| NS 29/32 | Evaporator piston 2000 ml |
| NS 29/32 | Evaporator piston 3000 ml |
| NS 29/32 | Powder piston 1000 ml |
| NS 29/32 | Evaporator cylinder 1500 ml |
| NS 29/32 | Distillation spider with 6 sleeves 20 ml |
| NS 29/32 | Distillation spider with 12 sleeves 20 ml |
| NS 29/32 | Distillation spider with 20 sleeves 20 ml |
| NS 29/32 | Distillation spider with 5 pistons 100 ml |
| S35/20 | Receiving flask 1000 ml |
| S35/20 | Receiving flask 1000 ml, coated |
| S35/20 | Receiving flask 2000 ml |
| S35/20 | Receiving flask 2000 ml, coated |
| NS 29/32 | Steam pipe |
| NS 24/29 | Steam pipe |
| NS 24/29 | Evaporator piston 100 ml |
| NS 24/29 | Evaporator piston 1000 ml |
| NS 24/29 | Evaporator piston 2000 ml |
| NS 24/29 | Evaporator piston 3000 ml |
| NS 24/29 | Powder piston 1000 ml |
| NS 24/29 | Evaporator cylinder 1500ml |
| NS 24/29 | Distillation spider with 6 sleeves 20 ml |
| NS 24/29 | Distillation spider with 12 sleeves 20 ml |
| NS 24/29 | Distillation spider with 20 sleeves 20 ml |
| NS 24/29 | Distillation spider with 5 pistons 100 ml |
| RV 06.13 | Seal: FKM material |
| RV 06.15 | Seal; PTFE material |
| PC 1.1 | Cable |

## Error codes

Any malfunctions during operation will be identified by an error message on the display (only for RV 10 digital).
Once an error message has been displayed, the lift moves to the top end position and the device can no longer be operated. The lift can be operated again.

Proceed as follows in such cases:

☞ Switch off device using the main switch (A)
☞ Carry out corrective measures
☞ Restart device

| Error code | Effect | Cause | Correction |
|---|---|---|---|
| **E01** | No or too little rotation | - Cable break<br>- Load too large<br>- Motor blocked | Reduce the volume in the evaporator piston |
| **E02** | No communication with the heating bath | Interface hidden or dirty | Check / clean the interface |

If the actions described fail to resolve the fault or another error code is displayed then take one of the following steps:
- Contact the service department
- Send the device for repair, including a short description of the fault.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

## Technical Data

| | | RV 10 basic | RV 10 digital |
|---|---|---|---|
| Operating voltage range | Vac | 100 - 230 ± 10% | |
| Rated voltage | Vac | 100 - 230, 50/60 Hz | |
| Frequence | Hz | 50 / 60 | |
| Power input without heating bath | W | 75 | |
| Power input in operation "stand by" | W | 3,3 | |
| Speed | rpm | 20-270 | |
| Speed display | | digital | |
| Dimensions of visible display area (W x H) | mm | 48 x 35 | |
| Display | | 7 segment display | |
| Multiple languages | | - | |
| Right and left movement/ interval operation | | yes | |
| Smooth start | | ja | |
| Lift | | automatic | |
| Lifting speed | mm/s | 50 | |
| Stroke | mm | 140 | |
| Setting of lower end stop | | 60 mm, contactless | |
| Head angle adjustable | | 0° - 45° | |
| Timer | | 1-199 Minuten | |
| Interval | | 1-60 sec | |
| Integrated vacuum controller | | - | |
| Interface | | - | RS 232 |
| Ramp programming | | - | with accessory labworld*soft* |
| Distillation processes programmable | | - | with accessory labworld*soft* |
| Remote control | | - | with accessory labworld*soft* |
| Cooling surface | $cm^2$ | 1200 | |
| Fuse | | 2xT1 ; 6A 250V 6,35x19,05 | |
| Perm. On- time | % | 100 | |
| Perm. ambient temperature | °C | 5-40 | |
| Perm. humidity | % | 80 | |
| Protection acc. to DIN EN 60529 | | IP 20 | |
| Protection class | | I | |
| Overvoltage category | | II | |
| Contamination level | | 2 | |
| Weight (no glassware; no heating bath) | kg | 16,8 | 17,5 |
| Dimensions (W x D x H) | mm | 500 x 410 x 430 | |
| Operation at a terrestrial altitude | m | max. 2000 above sea level | |

Subject to technical changes!

## Solvent table (excerpt)

| Solvent | Formula | Pressure for boiling point 40 °C in mbar |
|---|---|---|
| Acetic acid | $C_2H_4O_2$ | 44 |
| Acetone | $C_3H_6O$ | 556 |
| Acetonitrile | C2H3N | 226 |
| N-Amylalcohol, n-Pentanol | $C_5H_{12}O$ | 11 |
| n-Butanol | $C_4H_{10}$ | 25 |
| tert. Butanol, 2-Methyl-2-Propanol | $C_4H_{10}O$ | 130 |
| Butylacetate | $C_6H_{12}O_2$ | 39 |
| Chlorobenzene | $C_6H_5Cl$ | 36 |
| Chloroform | $CHCl_3$ | 474 |
| Cyclohexane | $C_6H_{12}$ | 235 |
| Dichloromethane, Methylenechloride | $CH_2Cl_2$ | atm.press. |
| Diethylether | $C_4H_{10}O$ | atm.press. |
| 1,2,-Dichloroethylene (trans) | $C_2H_2Cl_2$ | 751 |
| Diisopropylether | $C_6H_{14}O$ | 375 |
| Dioxane | $C_4H_8O_2$ | 107 |
| Dimethylformamide (DMF) | $C_3H_7NO$ | 11 |
| Ethanol | $C_2H_6O$ | 175 |
| Ethylacetate | $C_4H_8O_2$ | 240 |
| Ethylmethylketone | $C_4H_8O$ | 243 |
| Heptane | $C_7H_{16}$ | 120 |
| Hexane | $C_6H_{14}$ | 335 |
| Isopropylalcohol | $C_3H_8O$ | 137 |
| Isoamylalcohol, 3-Methyl-1-Butanol | $C_5H_{12}O$ | 14 |
| Methanol | $CH_4O$ | 337 |
| Pentane | $C_5H_{12}$ | atm.press. |
| n-Propylalcohol | $C_3H_8O$ | 67 |
| Pentachloroethane | $C_2HCl_5$ | 13 |
| 1,1,2,2,-Tetrachloroethane | $C_2H_2Cl_4$ | 35 |
| 1,1,1,-Trichloroethane | $C_2H_3Cl_3$ | 300 |
| Tetrachloroethylene | $C_2Cl_4$ | 53 |
| Tetrachloromethane | $CCl_4$ | 271 |
| Tetrahydrofurane (THF) | $C_4H_8O$ | 357 |
| Toluene | $C_7H_8$ | 77 |
| Trichloroethylene | $C_2HCl_3$ | 183 |
| Water | $H_2O$ | 72 |
| Xylene | $C_8H_{10}$ | 25 |

## Warranty

In accordance with **IKA**® warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

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
| Item | Designation |
|---|---|
| 5 | Foot |
| 6 | Metal sheet |
| 14 | Damping foil |
| 15 | Metal sheet |
| 16 | Metal sheet |
| 17 | Metal sheet |
| 19 | Foil |
| 20 | PCB power supply RV 10 |
| 21 | PCB power & logic |
| 24 | Cover *(only RV 10 basic)* |
| 26 | Rocker switch, 2-pin, compl. |
| 27 | Connector plug |
| 29 | Angle plug |
| 30 | Connection cable |
| 31 | Flat strip conductor, 16-pin. |
| 32 | Connection cable |
| 33 | Device fuse |
| 35 | Filter foil |
| 36 | Base |
| 37 | Cover |
| 38 | PCB display |
| 39 | Front foil |
| 40 | Operating knob |
| 41 | Display cover |
| 44 | Bracket |
| 45 | Blind nut M4, 4x |
| 46 | Stepping motor |
| 47 | Column, 2x |
| 48 | Reed sensor |
| 49 | Clamp |
| 53 | Sliding carriage |
| 54 | Bearing |
| 55 | Magnet (3x) |
| 57 | Gear rod |
| 58 | Bearing bushing |
| 59 | Threaded bushing |
| 62 | Detent pin |
| 63 | Detent sleeve |
| 64 | Pressure spring |
| 67 | Bridge |
| 68 | Reed sensor |
| 69 | Metal sheet |
| 70 | Flat strip bracket |
| 71 | Gas pressure spring |
| 72 | Pressure spring |
| 73 | Slide |
| 74 | Knob |
| 75 | Detent plate |
| 76 | Pressure spring |
| 78 | Guiding |
| 79 | Spring nut |
| 86 | Cover |
| 87 | Flange |
| 89 | Rear panel |
| 108 | Retaining ring |
| 117 | Cheese-head screw |
| 122 | Label |
| 125 | Label |
| 129 | Knurled screw |
| 134 | Ferrite core |
| 1954 | RV 06.15 |
| 1955 | Support rod |
| 1956 | Protection hose |
| 1957 | Bracket |
| 1958 | Knurled screw |
| 1959 | Velcro® strip |
| 1960 | Bracket |
| 1961 | Condensate bottle |
| 1964 | Screw plug |
| 1965 | Screw coupling |
| 1966 | Hose connection |
| 1967 | Steam pipe |
| 3004 | Housing - bottom part |
| 3007 | Gehäuse - upper part |
| 3009 | Gear motor |
| 3011 | Bearing tube |
| 3012 | O-ring |
| 3013 | Flat spring |
| 3014 | Insert |
| 3015 | Gear belt disc |
| 3016 | Gear belt disc |
| 3017 | Gear belt |
| 3018 | Pressure ring |
| 3019 | Locking ring |
| 3020 | Connection cable |
| 3021 | Cap nut |
| 3022 | Threaded flange |
| 3023 | O-ring |
| 3025 | Circlip |
| 3034 | Ball bearing |
| 3035 | Retaining ring |
| 3036 | Retaining ring |
| 3037 | V-ring |
| 3038 | Rotary shaft seal |
| 3039 | Grooved ball bearing |
| 3041 | Close tolerance grooved pin |
| 3055 | Caution label |

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




## India

IKA® - Werke
GmbH & Co.KG
Liaison Office India
No. 31 (Old No. 264)
1st Floor, 10th Cross
1st "N" Block, Rajajinagar
560 010 Bangalore
Tel.: +91 80-41157736
Fax: +91 80-41157735
E-Mail: info@ikaindia.com

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**Artisan Technology Group** is an independent supplier of quality pre-owned equipment

## Gold-standard solutions

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

## We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

## Learn more!

Visit us at **artisantg.com** for more info on price quotes, drivers, technical specifications, manuals, and documentation.

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

We're here to make your life easier. How can we help you today?

(217) 352-9330 | sales@artisantg.com | artisantg.com

ARTISAN® TECHNOLOGY GROUP

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IKA® Works Inc.