ARTISAN®
TECHNOLOGY GROUP

Your **definitive** source for quality pre-owned equipment.

Artisan Technology Group
(217) 352-9330 | sales@artisantg.com | artisantg.com

Full-service, independent repair center
with experienced engineers and technicians on staff.

We buy your excess, underutilized, and idle equipment
along with credit for buybacks and trade-ins.

Custom engineering
so your equipment works exactly as you specify.

- Critical and expedited services
- In stock / Ready-to-ship
- Leasing / Rentals / Demos
- ITAR-certified secure asset solutions

Expert team | Trust guarantee | 100% satisfaction

All trademarks, brand names, and brands appearing herein are the property of their respective owners.

Find the *IKA RW 14* at our website: ***Click HERE***

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CE-DECLARATION OF CONFORMITY *EN*

We declare under our sole responsibility that this product corresponds to the regulations 89/336/EG, 98/37/EG and 73/23/EG and conforms with the standards or standardized documents DIN EN IEC 61 010-1; DIN EN IEC 61 326-1; EN 292-1 -2 and EN 60 204-1.

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Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Guarantee

You have purchased an original **IKA** laboratory machine which meets the highest engineering and quality standards.
In accordance with **IKA** guarantee conditions, the guarantee period is 24 months. For claims under the guarantee please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

If you require servicing, return the appliance in its original packaging. Storage packaging is not sufficient. Please also use suitable transport packaging.

Contents

CE-Declaration of conformity	2
Guarantee	12
Safety instructions	13
Correct use	14
Unpacking	15
Useful facts	15
Drive	15
Motor circuit breaker	15
Speed display	16
Commissioning	16
Securing the laboratory stirrer to the stand	
Securing the stirring element	
Securing the stirring shaft protection	
Switching on the machine	17
Output shaft	17
Signals and visual display	17
Maintenance and cleaning	18
Accessories	19
Permitted **IKA** stirrer tools	19
Associated standards and regulations	19
Spare parts list	19
Technical data	20
Spare parts diagram/ wiring plan	111

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Safety instructions

In order to be able to use the appliance properly and safely, every user must first read the operating instructions and observe the safety instructions contained therein. Take care of these operating instructions and keep them in a place where they can be accessed by everyone. Only staff who have been trained accordingly, know the appliance and are authorised to carry out work in this field should use this appliance.

Never operate the equipment with the stirrer tools rotating freely. Ensure that parts of the body, hair or items of clothing cannot be trapped by the rotating parts. Pay attention when setting the speed to any imbalance of the stirrer tools and possible spraying of the medium to be stirred. Please use stirring shaft protective equipment (see fig.1).

The equipment is not suitable for manual operation. Please follow the relevant safety instructions and guidelines, and occupational health and safety regulations for use in the laboratory.

When working with the stirring unit, the user must select and wear his personal protective equipment according to the mixing hazard category. Defective or inappropriate protective equipment can expose the user to the risk of splashing or getting caught in the stirring element.


Fig. 1: points of danger

It is important to note that the surfaces of the motor (cooling fins) and certain parts of the roller bearing can get very warm during operation.

Avoid knocks and impacts on the lower end of the shaft and the chuck gear teeth. Even minor, invisible damage can lead to imbalance and uneven shaft action. Careful handling guarantees safe working and the long life of the equipment.

Imbalance of the output shaft, the chuck and in particular the stirring tools can lead to uncontrolled resonant vibrational behaviour of the equipment and the whole assembly. Glass apparatus and stirrer containers can be damaged or shattered by this. It can cause injury to the operator, as can the rotating stirring tool. If vibration of the equipment is noticed, the speed must be reduced immediately in all cases until no more vibration occurs. In this case exchange the stirring tool for one without imbalance or remedy the cause of the imbalance.

Operation with a freely rotating stirring element without a stirring vessel is dangerous. You must ensure that the stirring element is securely clamped in the chuck.

The machine may only be opened by trained specialists - even during repairs. The machine is to be unplugged from the mains before opening. Live parts inside the machine may still be live for some time after unplugging from the mains.

The high torque developped by the **RW 14 basic** requires particular care in the choice of stand, cross sleeve and anti-rotation element for the stirring bowl.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Note ! Covering or parts that are capable of being removed from the unit without accessory equipment have to be reattached to the unit for safe operation in order to prevent, for example, the ingress of fluids, foreign matter, etc..

The agitated vessels used for stirring have to be secured. Ensure that the whole assembly is stable.

Additional hazards to the user may occur if inflammable materials are used during the stirring operation or, if glass equipment is used, the energy transmitted from the stirring device results in glass breakage.

Furthermore, the stirring unit may only be used to stir those materials or material mixtures that the user knows will not react dangerously to the extra energy produced by the stirring. This also applies to extra energy by means of solar radiation during the stirring procedure.

Make certain that the unit is set at the lowest speed before commissioning; otherwise, the unit will begin running at the speed last set. Stirring at high speeds can cause the medium to splash out of the mixing vessel. There may be dangerous electrostatic activity between the medium and the output shaft.

The appliance starts up again automatically following a cut in the power supply or mechanical failure during a stirring procedure.

Stirring unit **RW 14 basic** may not been used in explosive atmospheres, for mixing dangerous substances or under water.

Accessories may only be assembled once the plug has been disconnected.

Correct use

The **RW 14 basic** laboratory stirrer can be used to stir and mix liquids of low to medium viscosity with various stirring tools. It is designed for use in laboratories.


Fig. 2

In order to use the appliance properly, it must be carefully secured to the stand using the cantilever (see fig. 2) and a cross sleeve. Assemble the appliance on an even, stable base. Ensure that the mounting equipment (e.g. fixed plate) is assembled on a stable base and does not start to move during the stirring procedure.

Use stirring shaft protective equipment (e.g. R 301 with R 301.1 see fig. 6).

The safety of the user cannot be guaranteed if the appliance is operated with accessories that are not supplied or recommended by the manufacturer or if the appliance is operated improperly contrary to the manufacturer's specifications.

Do not touch rotating parts during operation.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Unpack

Please unpack the machine carefully and inspect for damage. It is important that any transit damage should be noted at the time of unpacking. In certain circumstances it may be necessary to investigate immediately (post, rail or freight forwarder).
The guarantee covers: One **RW 14 basic** stirring machine in accordance with the type ordered, a chuck key and operating instructions.

Useful facts

In buying this machine you have chosen a high quality product. It is particularly useful for reproducibly setting the speed of stirring procedures or for producing media in which high viscosity values may occur temporarily. The appliance is manageable and simple to operate, thus ensuring ease of use.

You can set speeds of 100 to 2000 rpm with the appliance. The motor heats up during use. The outer surface of the motor acts as a cooling surface distributing and radiating the heat as evenly as possible.

You can also use high - sided mixing vessels together with **IKA** stands without any problem.

The chuck accommodates stirrer shanks of up to max. 8 mm in diameter. You must ensure that the stirrer shank is inserted into the drill chuck as far as the stop. You must also ensure that the stirrer shank is securely clamped in the chuck.

Drive

The speed of the directly powered **RW 14 basic** laboratory stirrer can be infinitely adjusted by electronic means. The quiet drive (DC motor equipped with ball bearings) is maintenance - free and has an anti - overload system. The constant torque over the whole speed range is particularly suitable for reproducibly stirring and mixing media of medium viscosity and for test series carried out successively and in parallel under the same conditions.
If the max. output torque is exceeded because of a process-related high increase in viscosity, the appliance will switch to overload mode (motor current limit) thus reducing the appliance's output power (red LED flashes). The speed of the stirring shaft will also decrease. Please note that fluctuations in the mains voltage within the permitted tolerances will also cause minor fluctuations in the speed.

Motor circuit breaker

The DC motor, which is equipped with ball bearings, is maintenance-free and is suitable for continuous operation. The motor has an anti-stall and anti - overload system which is operated by an electronic motor current limit. If a fault occurs, the motor is immediately switched off permanently. This will be displayed by means of a red LED on the front of the appliance (see "C" fig. 3). Please check that the appliance is still operable by switching it off and then on again.


Fig. 3

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Speed display

The speed of the **RW 14 basic** stirring unit is set using the knob (D) on the front of the appliance (see fig. 7). The scale on the front panel displays the set speed (The values only apply when the stirring shaft is in normal mode and not in overload mode).

Commissioning

Assemble the test texture on a stable, even, non- slip surface. The **RW 14 basic** laboratory stirrer must be secured to a stable stand (e.g. R 1825, R 1826 or R 1827) with a cross - sleeve (e.g. R 182) (see fig 1). The stirring bowl must always be securely fixed for safety reasons. You must also ensure that the mounting equipment (e.g. fixed plate) is set up in such a way that it is not liable to topple and does not start to move during the stirring procedure.

Accessories must be assembled according to the following assembly instructions:

Securing the laboratory stirrer to the stand

First secure the cross sleeve (K) to the stand leg (I).
Then secure the cantilever (L) of the laboratory stirrer in the open side of the cross sleeve which is facing upwards. Once you have set the position required for the stirring procedure, firmly tighten the two clamping screws (M).


Fig. 4

Securing the stirring element

Slide the stirring element (G) into the chuck (F) as far as the stop as shown and firmly tighten using the chuck key (H) (⊕ chuck locks/clamps; ⊖ chuck opens).


Fig. 5

Securing the stirring shaft protection

You must use stirring shaft protection when working to prevent injury. This is secured to the stand leg (I) using a cross sleeve (K) and a supporting bar (N). The plastic half shells are secured to the supporting bar (N) as shown in fig. 6 using the screws (P) provided (+ tighten screw; ⊖ loosen screw). You can alter the length of the stirring shaft protection using the screw (Q).

Check each time before starting up the appliance and at regular intervals that the laboratory stirrer, the stirring element and the stirring shaft protection are securely clamped. You may only change the position of the laboratory stirrer and the stirring shaft protection when the appliance is not running and plug is disconnected. The chuck key must be removed before switching on.


Fig. 6

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Switching on the machine

Check whether the voltage given on the type plate corresponds to the available mains voltage. The socket used must be earthed (fitted with earth contact). If these conditions have been met, the machine is ready to operate when the mains plug is plugged in. If these conditions are not met, safe operation is not guaranteed and the machine could be damaged.

The appliance will begin to run when you set the on/ off switch (A) to position "I". The appliance is 'ON' when the control light (B) is green.

When you start up the appliance, the drive shaft will begin to run at the last set speed. It is therefore important to check the position of the knob/ the set speed before switching on. Also make sure that the set speed is safe for the selected test texture. If you are not sure, set the speed knob (D) to the lowest speed (left stop), (see fig. 7).


Fig. 7

Only switch the appliance on once the designated working height has been set and the stirring instrument is immersed in the mixing vessel.

Output shaft

The chuck on the off- drive shaft allows stirring elements with shank diameters of 0.5 to 8 mm to be clamped.

The stirring elements must be clamped deeply enough in the chuck. You must ensure that the stirring element shank is securely clamped in the chuck (see fig. 5).

Stirring elements which have not been properly secured may be catapulted out of the test texture and could damage glass structures and cause injury.

Operation with a freely rotating stirring element without a stirring vessel is dangerous and is not permitted.

Signals and visual display

If the max. permitted torque is exceeded or the stirring shaft stalls, the laboratory stirrer will switch off automatically and a red LED "C" (overload display) will illuminate on the front of the appliance (see fig. 7).

Once the cause of the mechanical overload has been removed, please press the ON/ OFF button (A) to start the appliance again.

If the appliance still does not work, even after switching off and switching on again following a lengthy break, please contact our service department. Please describe the fault to us in as much detail as possible. This makes it easier for us to locate the fault and give an initial response.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Maintenance and cleaning

The stirrer is maintenance - free. It is subject only to the natural wear and tear of components and their statistical failure rate.

When ordering spare parts, please give the manufacturing number shown on the type plate, the machine type and the name of the spare part.

Please send in equipment for repair only after it has been cleaned and is free from any materials which may constitute a health hazard. Please tell us if you have used any materials which are dangerous or which may constitute a health hazard.
Only use cleansing agents which have been approved by **IKA** to clean **IKA** appliances.

Use to remove:

Dyes	isopropyl alcohol
Construction materials	water containing tenside/ isopropyl alcohol
Cosmetics	water containing tenside/ isopropyl alcohol
Foodstuffs	water containing tenside
Fuels	water containing tenside

For materials which are not listed, please request information from **IKA**.

Wear the proper protective gloves during cleaning of the devices. Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Accessories

R 1825	Plate stand	R 301	Stirring shaft protection
R 1826	Plate stand	R 301.1	Stand holder
R 1827	Plate stand	RH 3	Tensioner
R 182	Cross sleeve	FK 1	Flexible coupling

Permitted IKA stirrer tools

	Max. speed rpm
R 1300 Dissolver stirrer	2 000
R 1303 Dissolver stirrer	2 000
R 1311 Turbine stirrer	2 000
R 1312 Turbine stirrer	2 000
R 1330 Anchor stirrer	1 000
R 1332 Anchor stirrer	800
R 1352 Centrifugal stirrer	2 000
R 1373 Surface stirrer	1 000
R 1342 Propeller stirrer 4fl	2 000
R 1381 Propeller stirrer 3fl	2 000
R 1382 Propeller stirrer 3fl	2 000
R 1389 Propeller stirrer 3fl	800

Associated standards and regulations

Associated EU guidelines

73/ 23/ EG	Low voltage guideline
89/ 336/ EG	EMC - guideline
98/ 37/ EG	Machine guidelines

Construction in accordance with the following safety standards
DIN EN IEC 61010 – 1, EN 292-1,
EN 292-2, EN 60 204-1,
UL 61010A-1; UL 61010A-2.051; CAN/CSA C22,2 No 1010-1

Construction in accordance with the following EMV standards
DIN EN IEC 61326 + A1 + A2

Spare parts list

1 Fillister head screw
2 Counter sunk screw
4 Counter sunk screw for thermoplastics
5 Threaded pin
10 Housing- upper part incl. extension
11 Housing base
12 BLP - control **RW 14 basic**
13 Connection line
14 Cable strap
15 Flange
16 DC motor
17 Connecting terminal 3-pole
18 Pulse ring
19 Front panel
20 Front plate
21 Rotary control knob
22 Chuck incl key

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Technical Data

Speed range:	rpm	100 – 2000	
Max. torque stirrer shaft:	Ncm	8	
Permitted on - time:	%	100	
Speed adjustment:		Rotary control knob (potentiometer) - front of appliance	
Speed display:		Scale - front of appliance	
Max. speed deviation (at nominal torque):	%	± 5	
Design voltage:	VAC	220-240	100-120
Design frequency:	Hz	50/60	50/60
Input power (motor):	W	28.5	28.5
Output power (motor):	W	17	17
Max. output at stirrer shaft :	W	17	17
Max. power input:	W	45	45
Protection type to DIN 60 529:		IP 20	
Excess voltage category:		II	
Protection class:		I (protective earth)	
Contamination level:		2	
Protection at overloaded:		motor current limit	
Perm ambient temperature:	°C	+5 bis+40	
Perm. ambient humidity (rel.):	%	80	
Operating position:		on stand, chuck facing downwards	
Drive:		DC motor direct drive	
Maximum stirring quantity- water:	ltr	8	
For medium viscosity (MV):	mPas	1000-10000	
Clamping chuck- clamping range:	mm	0.5 - 8	
Hollow shaft internal diameter Ø:	mm	no hollow shaft	
Stand clamp:		bowl clamp	
Housing:		thermoplastic synthetic	
Max. dimension (W x D x H) without extension arm	mm	70 x 176 x 197	
Max. dimension (W x D x H) with extension arm	mm	70 x 336 x 197	
Weight with extension arm and clamping chuck:	kg	2	
Operation at a terrestrial altitude:	m	max. 2000 above sea level	

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

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IKA® Works, Inc.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

North America

IKA® Works, Inc.
2635 North Chase Pkwy SE
Wilmington, NC 28405-7419
USA
TEL. +1 800 733-3037
TEL. +1 910 452-7059
FAX +1 910 452-7693
E-mail: usa@ika.net

IKA® Works, (Asia) Sdn Bhd

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

Asia - Australia

IKA® Works (Asia) Sdn Bhd
No. 17 & 19, Jalan PJU 3/50
Sunway Damansara Technology Park
47810 Petaling Jaya
Selangor, Malaysia
TEL. +60 3 7804-3322
FAX +60 3 7804-8940
E-mail: sales@ika.com.my

IKA® Japan Y.K.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

Japan

IKA® Japan Y.K.
293-1 Kobayashi-cho
Yamato Koriyama Shi
639-1026 Japan
TEL. +81 74358-4611
FAX +81 74358-4612
E-mail: japan@ika.de

IKA® Works Guangzhou

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

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Artisan Technology Group is an independent supplier of quality pre-owned equipment

Gold-standard solutions

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

Learn more!

Visit us at **artisantg.com** for more info on price quotes, drivers, technical specifications, manuals, and documentation.

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

We're here to make your life easier. How can we help you today?

(217) 352-9330 | sales@artisantg.com | artisantg.com

ARTISAN® TECHNOLOGY GROUP