ARTISAN®
TECHNOLOGY GROUP

Your **definitive** source for quality pre-owned equipment.

Artisan Technology Group
(217) 352-9330 | sales@artisantg.com | artisantg.com

Full-service, independent repair center
with experienced engineers and technicians on staff.

We buy your excess, underutilized, and idle equipment
along with credit for buybacks and trade-ins.

Custom engineering
so your equipment works exactly as you specify.

- Critical and expedited services
- In stock / Ready-to-ship
- Leasing / Rentals / Demos
- ITAR-certified secure asset solutions

Expert team | Trust guarantee | 100% satisfaction

All trademarks, brand names, and brands appearing herein are the property of their respective owners.

Find the *IKA 8015701* at our website: ***Click HERE***

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CE-DECLARATION OF CONFORMITY *GB*

We declare under our sole responsibility that this product corresponds to the regulations 89/336EEC and 73/023EEC and conforms with the standards or standardized documents EN 61 010; EN 50 081; EN 50 082; EN 55 011 and EN 60 555.

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Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Contents

	Page
Safety Precautions	15
Proper use	15
Unpacking	16
General information	16
Setup	16
Mounting the Basic Unit onto the Stand	17
Mounting Agitators	17
Setting Up the Display Monitor	18
Mounting the Monitor onto the Stand	18
Interfaces and Outputs	18
Commissioning	20
Carrying out Measurements	20
Evaluation of Measured Values	21
Maintenance and Cleaning	21
Accessories	21
Adapters / Agitators	23
Technical Specifications	24
Associated standards and regulations	24
Guarantee	24
List of spare parts	36
Spare parts diagram	37

Safety Precautions

VK2500196

When assembling equipment: Always insert the jack plug of the power supply unit into the monitor jack first, and then plug the power supply unit into the power outlet.

When disassembling equipment: Always remove the jack plug of the power supply unit from the power outlet and then, if necessary, from the monitor jack.

As the IKA Viscoklick VK 250 and VK 600 are both extremely accurate measuring instruments, careful handling of the equipment represents the basic prerequisite for preserving the unit's precision.

Important: Please be sure that the agitating unit is capable of moving freely in the basic unit after being clamped.

Agitating tools may only be aligned by service personnel. Do not touch any rotating parts (agitating tools). Use an agitator shaft guard where possible.

Be sure to always prevent the agitated vessel from rotating.

Do not touch any of the moveable parts of the measuring unit (agitator and arm, Viscoklick - upper part, cable between the agitator and clamping plate) during measurement.

The clamp sleeve may be clamped with the clamping lever only if the agitator is fitted.

Fit only stirrer attachments onto the basic unit which either have the corresponding flange diameter or which are provided with a suitable adapter.

Connecting cables may only be removed or inserted when the unit is switched off.

Proper Use

The measuring units IKA Viscoklick VK 250 and VK 600 function as torque meters for laboratory agitators. Absolute torque and changes in torque can be measured and displayed during the agitation process.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Unpacking

Unpack the machine carefully and check it thoroughly for damage. It is important that any damage resulting from transit is determined when unpacking. An immediate damage report should be filed if necessary (post office, rail authority or shipping company).
The machine is supplied complete with: one IKA Viscoklick VK 250.1 or VK 600.1 with power supply unit, VK1 display monitor, VK1.2 connecting cable, one double offset screwdriver and operating instructions (stand, cross sleeves and other accessories have to be ordered separately).

❖ **NOTE!** *As the IKA Viscoklick VK 250 and VK 600 are both extremely accurate measuring instruments, careful handling of the equipment represents the basic prerequisite for preserving the unit's precision.*

General Information

The measuring unit transforms the dynamic torque, which acts on the agitating element, into a static torque. This static moment is measured by a force transducer and digitally displayed by means of display electronics.
The measuring unit is suitable for many laboratory agitators (for mounting requirements see the chapter on **Technical Specifications**) so that torque can be measured on already available agitators.
With a measuring range of 0 - 270 Ncm, the IKA Viscoklick VK 250 is sufficient for most laboratory agitators and tasks. The VK 600, with a range of 0 - 600 Ncm, is available for higher torque.
As the instrument is extremely accurate the viscosity of the medium can be determined by the torque measured, which is directly proportional to the viscosity. Viscosity conversion is carried out by means of the torque and the rotational speed, additional parameters include: temperature, form and size of the agitating element as well as the form and size of the vessel.
IKASOFT rheo software calculates, evaluates and displays the viscosity from both the prescribed and the measured parameters.

Setup

cf. Fig. 1

First of all, the two arms (8) supplied with the basic unit have to be attached. Insert the two arms into the bore holes on the back of the unit and, by turning, align them in such a way that both fillister head screws are capable of being screwed into the threaded holes of the arms. Now tighten the fillister head screws by using the tool supplied. After fastening the arms close the bore hole located on the upper surface of the unit by using the black stopper in order to prevent fluids from penetrating the unit.
The basic unit is fastened to a stabile stand by using two cross sleeves. The unit is only in a guaranteed vertical position if the two cross sleeves used are of the same size. Original IKA stands with a leg diameter of 34 mm and IKA R260 cross sleeves should be used in order to ensure accuracy.
We recommend a telescopic stand for easier handling.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Fastening the Basic Unit to the Stand

First attach both cross sleeves (9) to the stand (10) at a distance of 7 cm apart. Then slide the basic unit with the two arms into the cross sleeves. The easiest way to do this is if one of the cross sleeves on the stand leg is slightly moveable. Tighten the cross sleeves alternately when the basic unit has thus been brought into the desired position.
Attachment of the unit by means of the two cross sleeves ensures that the measuring unit stands straight and provides it with the required stability.

Fastening Agitator Attachments

For agitators having a flange diameter of 60 mm and a flange height of at least 10 mm for VK 250 and a flange diameter of 62 mm and a height of 10 mm for VK 600 (cf. page 13).

The agitator does not require an arm in order to be mounted onto the basic unit. It can be removed for easier handling. A stopper is supplied with EUROSTAR and RW16 agitators in order to close the arm bore holes. If the arm is left on the agitator then care should be taken to prevent it from coming into contact with anything.
The mounting ring (brass part) must be free from dirt and fat. If necessary, the ring should be cleaned.
Open the clamping lever (12) on the basic unit until the clamping ring completely rests on the housing. Now insert the agitator with the flange into the opening of the clamping ring.
The base of the agitator may, but does not have to, stand on the basic unit.
Agitators having a flange height of more than 15 mm (e.g. IKA agitator RE 162) are inserted into the clamping ring of the basic unit until stopping. Now align the agitator such that its sides and the basic unit are roughly parallel. Be sure that if an arm is attached to the agitator that it is not in contact with the stand leg or any other object. There must be a space of at least 1 cm between the stand leg and the arm. If this is the case then the agitator can be connected to the basic unit with a sturdy pull of the clamping lever.
Further tightening is not necessary. Over-tightening can lead to unit malfunctioning (bearing strain; the required easy movement is no longer ensured; measuring errors).
The agitator may only be connected to the basic unit via the flange; an attached arm may not be connected to the stand under any circumstances, as measurement is otherwise impossible.
The agitator mains cable (21) must be locked to the clamping plate (6) in such a way that it slightly sags between the clamping plate and the agitator without touching anything else. The clamping plate can be raised by pressing the button (5) found on the vertical section at the back of the basic unit. The mains cable can be inserted into the left or right clamping plate guide, depending on the type of agitator. By releasing the button the mains cable is then held in the corresponding position. The same procedure shall be followed if an interface cable is used.
Assembly of the agitator to the basic unit has been correctly executed if the agitator and the upper part of the basic unit can be moved easily (approximately 1 mm deflection of the upper part of the basic unit). Once assembly has been completed a check should be made to ensure that the agitator, the agitator arm and the cable leading away from the agitator, between the agitator and the clamping plate, are not touching.
An adapter is required for agitators without a cylindrical agitator flange diameter of 60 mm or 62 mm (e.g. IKA RW20 series).

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Display Monitor Setup

The display monitor (14) for reading the torque can be set up anywhere, independent of the basic unit. Connect the power supply unit (20) to the monitor first. In order to do so insert the plug jack into the socket provided on the monitor. Now connect the monitor to the basic unit by inserting and screwing in the round, 8-pin connecting cable plug (13) into the jack (16) located on the right side of the monitor. Now insert the other side of the connecting cable into the jack (11) located on the side of the basic unit, screwing it in as well.

❖ **NOTE!** *The connecting cable may only have a maximum length of 2 m. Proper functioning is not guaranteed if longer cables are used. The connecting cable should only be removed or inserted when the unit is switched off.*

Mounting the Monitor on to the Stand

A cross sleeve and an IKA H36 stay bar are required in order to fasten the monitor onto the stand. If an original IKA stand with a stand leg diameter of 34 mm is used, then an R260 cross sleeve is required. Fasten the stay bar with the cross sleeve to the stand, the plastic cylinder of the stay bar must point upward. The monitor can now be slid onto this plastic cylinder. The mount for this, a plastic clamping piece, is located on the back of the monitor housing.
Be sure that none of the monitor parts touch the moveable components of the basic unit and the agitator.

Interfaces and Outputs

The VK 1 unit has a 15-pin interface for transferring data. The interface jack (17) is located on the side of the monitoring housing and is covered with a protecting cap. The jack should always be covered with the protecting cap when not in use. The interface can be used to transmit analog and digital signals. By using IKA adapter PC5.1 (available as accessory equipment) the 15-pin interface can be split into an analog output (round, 7-pin) and a 9-pin RS232 interface. The analog output can be used, for example, to connect recording instruments in order to document the torque. The RS232 interface is used to transfer measured values to a PC. Software can then be used to further process the data (e.g. graphical representation of torque characteristics).

❖ **NOTE!** *Do not connect two interface cables for extension purposes. The signal lines will cross each other and prevent proper functioning.*
See the chapter on **Technical Specifications** for interface output signals and **Accessories** for available software and extender modules.

- Analog output (**Fig. 8**)

A voltage value for the measured quantity (torque) is available at the analog assigned pins.

(9)	measured quantity (torque)	1 VDC/100Ncm (VK250) 1 VDC/200Ncm (VK600)
(10)	analog GND	

- RS 232 (V24) serial interface (**Fig, 9**)

Configuration of the RS 232 C interface

- The function of the interface cables between the measuring unit and the automation system represent a selection of the signals

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

specified in EIA Standard RS232 C, in accordance with DIN 66 020 Part 1.
- Standard RS 232 C, in accordance with DIN 66 259 Part 1, applies to the electrical properties of the interface cables and assignment of the signal states.
- Transmission method: asynchronous character transmission in start-stop mode.
- Transmission type: full-duplex.
- Character format: character imaging in accordance with data format in DIN 66 022 for start - stop operation. One start bit; seven character bits; one parity bit (even); one stop bit.
- Transfer rate: 9600 bps.
- Data flow control: RTS/CTS hardware handshake
 RTS: (pin 7) LOW (positive voltage): PC may transmit
 RTS: (pin 7) HIGH (negative voltage): PC may not transmit
 CTS: (pin 8) LOW (positive voltage): PC ready-to-receive
 CTS: (pin 8) HIGH (negative voltage): PC not ready-to-receive
- Access method: data transmission from the measuring unit to the computer only takes place after a computer prompt.

Command Syntax and Format

The following applies to the command syntax:
- commands are generally sent from the computer (master) to the measuring unit (slave).
- the measuring unit only transmits at the computer's request. Even error messages cannot be sent spontaneously from the measuring unit to the computer (automation system).
- commands are transmitted in upper case letters.
- commands and parameters as well as successive parameters are separated by at least one blank space (code: hex 0x20).
- each individual command (including parameters and data) and each answer is concluded with CR blank LF (code: hex 0x20 hex 0x0d hex 0x20 hex 0x0A) and has a maximum length of eighty characters.
- the period is the decimal separator in a floating-point digit (code: hex 0x2E).

In general, the above-mentioned corresponds to the recommendations of the NAMUR Working Group. (NAMUR - Recommendations for the design of electrical plug-and-socket connections for analog and digital signal transmission to individual systems in laboratory control and instrumentation technology, Item 1.1).
The NAMUR commands and the additional commands specific to IKA serve only as low level commands for the purpose of communication between the measuring unit and the PC. These commands can be directly transmitted to the measuring unit by means of a suitable terminal or communications program.
The following chart provides an overview of the (NAMUR) commands understood by IKA Viscoklick.
Abbreviations used: X = 5 : torque

NAMUR commands	**Function**
IN_PV_X X = 5	actual value reading
IN_TYPE	laboratory equipment identification request
IN_NAME	designation request
OUT_NAME name	display of name designation (6 characters max., default: IKA_VK)

- PC 5.1 Adapter IKA Control (**Fig. 10**)

An adapter cable is available as an accessory to the Viscoklick series. It splits the analog and serial signals. The analog output signals are connected across a 7-pin jack in accordance with the NAMUR Recommendation and the serial signals across a 9-pin sub-D-female connector (RS 232 C).

(2) torque
(3) analog GND

(2) R x D
(3) T x D
(5) RS 232 GND
(7) RTS
(8) CTS

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

- PC 1.2 Adapter

This adapter is required in order to connect the 9-pin jack to an 8-fold serial interface (25-pin plug).

- PC 2.1 Cable (**Fig. 11**)

This cable is required in order to connect the 9-pin jack to a PC.

- AK 2.1 Cable

This cable is required in order to connect the 7-pin jack to a recording instrument (4 mm banana plug).

- AK 2.2 Cable

This cable is required in order to connect the 15-pin jack to a recording instrument (4 mm banana plug).

Commissioning

Before putting the measuring unit into service make certain that the agitator exhibits rotational accuracy. A balance test is best performed at a lower rotational speed. Large unbalances would affect the measurement result. A slight unbalance can also be caused by the agitator clamping chuck and thus it may be advisable to open the chuck and to twist the agitator element and refit it. Agitator attachments which are strongly unbalanced have to be aligned. As a rule, small agitating elements with a small agitator shaft diameter center themselves, especially when the medium has a high viscosity and the unbalance is slight.

Measuring

Dip the agitating element into the vessel containing the medium to be stirred up to the desired depth and as close to the center as possible. Make certain that the agitated vessel cannot be twisted.
Activate the measuring unit by switching on the ON-OFF button (15). After a self-test (entire display flashes three times), (**Fig. 2**) the display switches to 0.0 Ncm (**Fig. 3**), and the LCD shows 0.0 at the upper right. The smaller display value is used to store measured data. Now turn on the agitator and set the desired rotational speed. After a brief run-up period a particular torque value is set which is then visible on the display. A slight fluctuation of the display value can result from a non-homogeneous medium or the untrue running of a stirrer attachment.
The torque displayed now corresponds to the actual torque acting on the agitator or generated by it. If the change in torque is important (increase or decrease as a result of admixture to the medium) in your stirring procedure, rather than the absolute torque, then the measuring system can be set to output value zero by pressing the offset button (18) at the desired moment.
The small display at the upper right shows the stored measured value which was valid at the moment of resetting, in this case absolute torque.
A decreasing torque is represented by a negative torque display, an increasing torque by a positive torque display. The absolute torque can be determined at any time by means of addition in accordance with the specified sign.
Absolute torque can, for example, be displayed at the end of a mixing operation by renewed pressing of the reset button. The absolute torque value then appears as a small display value, the large display value then switches again to 0.0 Ncm (**Fig. 4 / Fig 5**)
The original, absolute torque measured value, visible as a small display value at the top of the screen is added, in accordance with the specified sign, to the current differential value, which appears as the large display value. This makes comprehensive reading of absolute torque possible at any time and observation of changes in torque (viscosity) with the correct signs.
If the measured value exceeds 199,9 Ncm, the display will change from Ncm to Nm (**Fig. 6**).

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Measured Value Evaluation

The measured value in the run-up period (increase in the agitator torque) and during powering-down (decrease in the agitator rotational speed), as well as all the changes in torque set on the agitator, cannot be used in measurement evaluation. The highly increased values during revving up occur as a result of inertia and acceleration, they are not justified in the viscosity change of the medium. The same applies to the lower values during revving down.
Only those measured values obtained at constant speed or at very slight speed variation (e.g. during admixture of substances into the medium) can be utilized. Ideally, the agitator should maintain constant speed, although the viscosity is subject to change and thus too the torque on the agitator. This is the case for all IKA-EUROSTAR agitators.
No moveable parts of the measuring unit (agitator and arm, upper part of Viscoklick, cable between the agitator and the clamping plate) may be touched during the measurement.
Touching the moveable parts results in measuring errors due to generally high or low peak measured values. This applies particularly to operating functions carried out on the agitator (setting the speed with the rotary button, switching the digital display, etc.) which require contact with it. Peak measured values thus generated cannot be utilized for the evaluation of a series of measurements.

IKA Software (cf. chapter on Accessories) offers the possibility of controlling EUROSTAR power control agitators from a PC, thus eliminating contact with the agitator.

Maintenance and Cleaning

The measuring system requires no maintenance. It is only subject to the natural aging of the components and their statistical failure rate.
When ordering spare parts always provide the serial number given on the rating plate, the machine type and item number, as well as the designation of the spare part.
Please send in only those machines which have been cleaned and are free from materials which present health hazards. Only water and a tenside detergent additive should be used for cleaning or isopropyl alcohol in the event of heavy soiling.

Accessories

PC 5.1 Adapter IKA-Control

The adapter is used to split the 15-pin jack into a 9-pin RS232 digital output and a 7-pin analog output. The analog output can be used, for example, to connect recording instruments. The digital output is used to transmit measured data, e.g. to a PC for further processing.

IKASOFT

Suitable software is required in order to record the torque which is measured and then transmitted via the RS232 interface. IKA Software makes it possible to record torque characteristics and to represent them graphically. If there is an interface on the agitator as well, then this software can be used to record the torque and speed (the agitator speed is measured).

IKASOFT rheo

Several process parameters are necessary in order to determine the viscosity of a medium.
Torque, speed and temperature are the parameters that have to be measured. These parameters are recorded by the software and processed with the other predetermined parameters. Predetermined parameters include: sample vessel,

agitating element and the density of the medium at a particular temperature. The software uses this data to calculate the viscosity and provides graphical representation of the results.

R 301.1 Stand Mounting

If the HWS agitator shaft protector is to be attached in order to avoid injuries, then the R301.1 stand mounting is required because the agitator flange screw-on thread will be covered by the IKA-Viscoklick clamping ring.
While assembling make certain that the agitator shaft protector and the stand mounting do not come into contact with any part of the measuring system.

VK60/01 Adapter ***cf. Fig. 7***

This adapter is required for attaching the RW20 type series (cf. page 13).

Attachment to a VK250.1 measuring system: Place the adapter into the inner sleeve of the basic unit. The pin on the edge of the adapter has to catch in the inner sleeve slot and thereby prevent twisting, otherwise installation is not possible.

Attachment to a VK600.1 measuring system: The VK62/60 adapter is required in addition to the VK60/01 adapter for attachment to a VK600.1. Place the VK60/01 adapter into the VK62/60 adapter, the pin of the VK60/01 adapter must catch in the slot of the VK62/60 adapter. Both adapters are thus inserted into the inner sleeve of the basic unit. The pin on the edge of the adapter has to catch in the inner sleeve slot and thereby prevent twisting, otherwise installation is not possible.

An RW20 agitator can now be attached. Insert the ring for gear shifting (located above the chuck) into the adapter ring as far as it will go. The agitator has to be in alignment with the basic unit; there must be enough space (approximately 10 mm) between the agitator arm and stand bar. Never turn the agitator clockwise in order to adjust it as a gear speed change or neutral position will otherwise occur on the agitator. Assuming a frontal view, attach the agitator somewhat swung-out to the basic unit. Adjust the agitator while the gearing flange is in the adapter by turning clockwise and brace it with the clamping lever. Make certain that the gear shift ring lies evenly on the adapter at the bottom and that the agitator is thus in horizontal alignment with the basic unit. Before switching on the machine make certain that the agitator cannot be easily turned in a counterclockwise direction and that the arm is not on the stand bar.
Set the desired gear speed before assembling the agitator.
Should it be necessary to change gears during a series of measurements, then the agitator in gear range **II** (gear shift all the way to the top) has to be mounted onto the basic unit.
In order to shift to gear range **I** switch off the agitator and twist it clockwise without disengaging the clamp lever. Now pull the agitator upwards and turn it back counterclockwise as far as it will go. Now the agitator is no longer in alignment with the basic unit, but this has no effect on the measurement results. The measurement can be continued after the mains cable is clamped again.

VK62/60 Adapter ***cf. Fig. 7***

This adapter is used to attach agitators having a flange diameter of 60 mm to measuring system VK600.1 (cf. page 13).

Place the adapter ring into the inner sleeve of basic unit VK600.1, the adapter ring groove must be in more or less alignment with the groove on the inner sleeve (cf. Fig. 7 for permission deviation). The agitator used is then mounted onto the basic unit as described in the above and clamped using the clamping lever.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Adapters / Agitators

Agitators	VK250.1 0 - 270 Ncm	VK600.1 0 - 600 Ncm
RW 16 basic	•	VK62/60
EUROSTAR basic	•	VK62/60
EUROSTAR digital	•	VK62/60
EUROSTAR digi-visc	•	VK62/60
EUROSTAR Power basic	•	VK62/60
EUROSTAR Power digi-visc	•	VK62/60
EUROSTAR Power control-visc	•	VK62/60
EUROSTAR Power digi-visc-P4	•	VK62/60
EUROSTAR Power digi-visc-P7	◊	VK62/60
RW 20	VK60/01	VK60/01 + VK62/60
RW 20 DZM	VK60/01	VK60/01 + VK62/60
RW 20 DZM P4	—	•
RW 20.n	VK60/01	VK60/01 + VK62/60
RW 20 DZM.n	VK60/01	VK60/01 + VK62/60
RW 20 DZM.n P4	—	•
RE 162 analog	•	VK62/60
RE 162 P4 analog	—	•
RE 162 P7 analog	—	•
Outside products having a flange diameter of **60 mm** and a flange height of approx. 10 mm	•	VK62/60
Outside products having a flange diameter of **62 mm** and a flange height of approx. 10 mm	—	•

• agitators fit without an adapter
— agitators do not fit or exceed the specified torque measuring range of the measuring system
◊ agitator fits without an adapter, but exceeds the specified torque measuring range of the measuring system

VK2500196

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Technical Specifications

Mains voltage:	**VAC**	230
Frequency:	**Hz**	50 / 60
Operating voltage:	**V**	15
Acceptance power:	**W**	5
Percentage duty cycle:	**%**	100
Protective system:		IP21
Dimensions		
basic unit W x L x H	**mm**	110 x 216 x 110
monitor W x L x H	**mm**	80 x 180 x 79
Weight		
basic unit	**kg**	1.61
monitor	**kg**	0.24
Protection class:		2
Overvoltage category:		II
Contamination level:		2
Permissible ambient temperat.:	**°C**	5 - 40
Humidity: (relative)	**%**	80
Display:		LCD, 3 1/2 digit
Resolution:	**Ncm**	0.1 (0 - 199 Ncm)
		1 (as of 199 Ncm)

Interfaces: ***VK250.1***

analog, output signal: 0V - 5V = -100 Ncm — +400 Ncm (1V = 0 Ncm; Δ1V = Δ100Ncm)

RS232

Specified measuring range:	**Ncm**	0 - 270
Accuracy of measurement:	**Ncm**	± 0.5 (0 - 60 Ncm)
		± 1 (60 - 270 Ncm)
Hub flange Ø:	**mm**	60 ±0.5

Interfaces: ***VK600.1***

analog, output signal: 0V - 5V = -199 Ncm — +800 Ncm (1V = 0 Ncm; Δ1V = Δ200Ncm)

RS232

Specified measuring range:	**Ncm**	0 - 600
Accuracy of measurement:	**Ncm**	± 0.5 (0 - 60 Ncm)
		± 1 (60 - 600 Ncm)
Hub flange Ø:	**mm**	62 ±0.5

The measuring system has overload protection.

The following agitators are suitable without accessories:

All IKA-EUROSTAR IKA RE 162 agitators
All outside products having an agitator flange with a diameter of 60 mm and a minimum height of approximately 10 mm.

Associated standards and regulations

Associated EU guidelines

EMV guidelines: 89/336/EC
Low voltage guidelines: 73/023/EWG

Construction in accordance with the following safety standards:

EN 61 010-1 / VDE 411-1
EN 50 081
EN 50 082
EN 55 011
EN 60 555
UL 3101-1
CAN/CSA C22.2 (1010-1)

Guarantee

You have purchased an original IKA laboratory machine which meets the highest engineering and quality standards.
In accordance with IKA guarantee conditions, the guarantee period is 12 months. For claims under the guarantee please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

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List of spare parts

Item	Designation
	VK250.1 / VK600.1
1	pressure screw
2	PCB plug
4	clamping lever
5	pressure spring
8	clamping plate
11	force sensor
15	PCB-MD measuring amplifier
18	thin-ring four-point contact bearing
20	retaining ring SB 79
21	retaining ring SW 63
22	bearing cover
23	housing - lower part
24	bolt
25	pressure spring
26	cylindrical pin
27	cylindrical pin
28	flat ribbon conductor
2001	housing - upper part, complete

Item	Designation
2003	arm
2014	plug
	VK 1
2	switch
3	PCB logistics
4	clip
5	curved washer
6	rubber base
7	sealing strip
8	PCB plug 8-/15-pin
9	PCB plug
10	membrane keyboard
11	housing - upper part
12	housing - lower part
13	locking mechanism
14	flat ribbon conductor
15	flat ribbon conductor
16	cover

Renewed calibration is required if parts 11 (force sensor) or 15 (PCB-MP measuring amplifier) are replaced. Please send the device to our service department.

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

LABORTECHNIK
JANKE & KUNKEL GMBH
& CO. KG

IKA® WORKS, INC.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

IKA® WORKS, (Asia) Sdn Bhd

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

Europa - Africa

IKA LABORTECHNIK
JANKE & KUNKEL-STR. 10
79 219 STAUFEN
GERMANY
TEL (07633) 831-0
FAX (07633) 831-98

America

IKA Works, Inc.
2635 NORTH CHASE PKWY, SE
WILMINGTON, NC 28405-7499
TEL. 800/733-3037
TEL. 910/452-7059
FAX 910/452-7693

Asia - Australia

IKA Works (Asia)Sdn Bhd
No. 3, Lot PT 3403,
Batu 22 3/4 Lebuh Raya,
Kuala Lumpur Serendah,
48200 Serendah
Selangor, Malaysia
TEL :03-601 3122
FAX :03-601 3126

Artisan Technology Group is an independent supplier of quality pre-owned equipment

Gold-standard solutions

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

Learn more!

Visit us at **artisantg.com** for more info on price quotes, drivers, technical specifications, manuals, and documentation.

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

We're here to make your life easier. How can we help you today?

(217) 352-9330 | sales@artisantg.com | artisantg.com

ARTISAN® TECHNOLOGY GROUP