ARTISAN®
TECHNOLOGY GROUP

Your **definitive** source for quality pre-owned equipment.

Artisan Technology Group
(217) 352-9330 | sales@artisantg.com | artisantg.com

Full-service, independent repair center
with experienced engineers and technicians on staff.

We buy your excess, underutilized, and idle equipment
along with credit for buybacks and trade-ins.

Custom engineering
so your equipment works exactly as you specify.

- Critical and expedited services
- In stock / Ready-to-ship
- Leasing / Rentals / Demos
- ITAR-certified secure asset solutions

Expert team | Trust guarantee | 100% satisfaction

All trademarks, brand names, and brands appearing herein are the property of their respective owners.

Find the *Thermo Spectronic Helios* at our website: ***Click HERE***

Heλios Gamma & Delta User Manual

9499 230 45811
Issue 1a
001201

Copyright © 2000 Unicam Limited, Registration No. 441506
All rights reserved

Thermo Spectronic (Europe, Middle East and Africa)
Mercers Row, Cambridge CB5 8HY, UK
Tel: +44 (0)1223 446600
Fax: +44 (0)1223 446644
E-mail: info@thermospectronic.co.uk
www.thermospectronic.com

Thermo Spectronic (North America and Latin America)
820 Linden Avenue, Rochester, NY 14625, USA
Tel: (800) 654-9955 or (716) 248-4000
Fax: (716) 248-4014
E-mail:info@thermospectronic.com
www.thermospectronic.com

For Indian Sub-continent, Far East and Australasia
Contact Thermo Spectronic (Europe, Middle East and Africa)

Thermo Spectronic

A Thermo Electron business

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Environmental and Electrical Requirements

Your Helios Gamma or Delta spectrophotometer has been designed to operate under the environmental and electrical requirements listed below.

Line Voltages 100-240V ±10% (max)
50-60Hz ±10%

Power input 180VA

Operating Environment

Temperature	5°C to 40°C
Rate of change of temperature	<1°C/hour
Relative humidity	20% to 90% non condensing
Absolute humidity	3g/m^3 to 20g/m^3
Air pressure	68kPa to 106kPa
Air velocity	0.5m/s

Transport and storage requirements

Temperature	-20°C to 70°C
Rate of change of temperature	5°C/hour
Relative humidity	5% to 95% non condensing
Absolute humidity	0.1g/m^3 to 35g/m^3
Air pressure	25kPa to 110kPa
Air velocity	< 30m/s

For Indoor Use only.

Installation Category II.

The following are important for the safe operation of the instrument:

When connecting the instrument into a mains socket, ensure the socket is earthed. Position the instrument sufficiently away from other objects, walls etc, to ensure a free airflow around the instrument and also to allow access to the mains switch on the rear of the product.

The mains switch is marked “0” and “I”. When “0” is depressed the instrument is off. When ”I” is depressed, the instrument is switched on.

The fuse on the internal power supply is rated at 3.15A and is a type T on all products with a Deuterium lamp fitted and a 3.15A type F on Visible only variants.

The following connectors are provided on the rear of the instrument:

Sipper:	Use to connect a Supersipper or Minisipper
RS232:	Use to connect a computer
Parallel:	Use to connect a printer with a parallel interface
Recorder:	Use to connect a chart recorder

Ensure that any equipment such as printers which are connected to the sockets on the rear panel are to relevant IEC safety standards.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

USER MANUAL

CONTENTS

Section **Page**

INSTALLATION
Unpacking, locating. 1
Making connections, Initializing 1

SAFETY
Safety notes 3

GENERAL
User Interface 4
Text entry 6
Data Store 6

HOME PAGE
Live Display 7

METHODS LIBRARY
Methods Library Page 8

FIXED
Fixed Method 10
Fixed Results 12

QUANT
Quant Method 14
Quant Standards 16
Quant Calibration 16
Quant Results 17

RATIO
Ratio Method 18
Ratio Results 20

SCAN
Scan Method 21
Scan Graph 24

RATE
Rate Method 28
Rate Graph 30
Rate Results 33
Parallel Rate 34

SETUP
Overview, Clock 35
Printers, Environment 36
Wavelength Calibration and Optical Initialization ... 38
White Light, Recorder…. 39
Lamps 40

CELL PROGRAMMER
Cell Programmer modes of operation 42

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Removal and Refit 43

SUPERSIPPER

SuperSipper modes of operation 44
Sipper Calibration 46

MINISIPPER

MiniSipper modes of operation 47
MiniSipper Calibration 48

INTERNAL PRINTER

Internal Printer operation 50

MAINTENANCE

Routine Maintenance 51
Removal and Replacement of Lamps 53

FAULT FINDING


Fault Finding Guide 55
Connecting Helios to a PC 57

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

THIS SYSTEM IS DESIGNED TO BE USER INSTALLED -


So, use the following procedures to quickly get your new spectrophotometer system in position and running the way you want it.

1. Read the Safety instructions on page 3 of this manual.


2. Check the component parts of the system against the Despatch Note and Packing List. **Immediately report by telephone (confirmed in writing) any discrepancies**.

3. Find a suitable location. (see below)


4. Connect the supplied power cord to the spectrophotometer.

5. Ensure NO SAMPLE(S) are in place in the sample compartment .
- Turn the spectrophotometer on.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com


6. Local Control display shows this sequence if OK
- ✓appear sequentially as each test is completed.


Approximately

3 minutes

THIS COMPLETES THE INSTALLATION OF A STANDALONE SYSTEM - IF YOU WISH TO INSTALL A PRINTER GO TO POINT 7.

Connecting up a Printer to a Standalone system

7. **ENSURE THAT THE SPECTROPHOTOMETER IS TURNED OFF -**

Now, connect the printer to the Parallel port on the spectrophotometer using the cable supplied.

8. To set up the Printer in the software:

POWER UP THE SYSTEM AND ALLOW TO INITIALISE
- From the HOME page, press MENU
- from the MENU page, press SETUP
- select PRINTER using the cursor keys,
- press ENTER. The PRINTERS page is displayed with the default printer selected.
- press ENTER again to display the Printer Menu and using the cursor keys, select the required printer.
- press ENTER.

THIS COMPLETES THE INSTALLATION OF A STANDALONE SYSTEM

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

SAFETY

- THIS SPECTROPHOTOMETER SYSTEM HAS BEEN DESIGNED TO BE USER INSTALLED.

 READ THIS PAGE CAREFULLY BEFORE INSTALLING AND USING THE INSTRUMENT AND ITS ACCESSORIES.

- The safety statements in this manual comply with the requirements of the HEALTH AND SAFETY AT WORK ACT 1974.

- The instrument and accessories described in this manual are designed to be used by properly trained personnel only. Adjustment, maintenance and repair of exposed equipment must be carried out only by qualified personnel who are aware of the hazards involved. Where indicated in the relevant manual, certain maintenance processes may be carried out by the user, who must be fully aware of, and apply, the following safety precautions.

- For the correct and safe use of this instrument and its accessories it is essential that both operating and servicing personnel follow generally accepted safety procedures in addition to the safety precautions specified in this manual.

- Specific warning and caution statements, where applicable, can be found throughout the manual. Warning and caution statements and/or symbols are marked on the apparatus where necessary.

- The instrument covers and accessories should only be removed by personnel who have been trained to avoid the risk of electric shocks. The mains electricity supply to the instrument must be disconnected at the mains supply connector and at least three minutes allowed for capacitors to discharge.

- Some of the chemicals used in spectrophotometry are corrosive, and/or flammable and samples may be radioactive, toxic or potentially infective. Care should be taken to follow the normal laboratory procedures for handling chemicals.

- The UV radiation from a deuterium lamp can be harmful to the skin and eyes. Always view the lamp through protective glasses/goggles that will absorb the UV radiation and avoid looking directly at the deuterium arc. Do not expose the skin to direct or reflected UV radiation.

- Whenever it is likely that safety protection has been impaired, the instrument and/or accessory must be made inoperative and secured against any unintended operation. The matter should then be referred to the appropriate servicing authority. Safety protection is likely to be impaired if for example, the instrument fails to perform the intended measurements or shows visible damage.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

GENERAL

Introduction


- The range comprises a group of UV-Visible and Visible spectrophotometers which can be controlled independently via an integral keypad and LCD display, or externally from a computer.

- The system is composed of a spectrophotometer with integral keypad, LCD display, Local Control Software and output device.

LCD display contrast can be controlled during initialisation, or from the MENU page using the left and right arrow keys.

- The Local Control software controls all aspects of the systems operation.

User Interface

Key	Operation
Arrow Keys	Highlight menu options, or move track cursor, or move Cell Programmer, depending on page in use.
Change	display contrast with < > from the Menu Page.
Numerical Keys	Enter numbers, minus and decimal point.
Function Keys	Access and perform system functions. Operation depends on screen in use, and is indicated by labels at bottom of screen.
Clear	Clears entry leaving field or parameter ready for new entry, clears pop-up, and clears error messages.
Enter	Enters changes to field or parameter.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Key	Operation
Run	Starts instrument measurement according to current method.
Home	Returns to Home page.
Zero/Base	Performs a zero or baseline as appropriate to application.


REMOVE THE SAMPLE AND ENSURE THAT THE SAMPLE BEAM IS CLEAR OR CONTAINS THE BASELINE SAMPLE RELEVANT TO THE ANALYSIS BEFORE ZEROING THE INSTRUMENT OR PERFORMING A BASELINE SCAN.

Software

The Local Control Software is organised in a tree structure with all functions accessed initially from the MENU Page. Movement between the software pages is controlled by function keys or by highlighting with cursor keys and pressing ENTER. Pressing the Home key will always display the HOME page.

The FIXED, QUANT, RATIO, SCAN, and RATE applications are entirely separate. Only one application can be operational at any one time. Loading another application will overwrite any current data.

Local and Computer Control

REMOTE

- From switch-on the instrument will automatically be under the control of the LOCAL CONTROL software. To enable control from an external computer press REMOTE on the MENU page. Control will go to the computer providing the instrument is idle.
- To return control to the Local Control software press the HOME key. Control will revert to the Local Control software providing that the instrument is idle, and the PC software has relinquished control.


Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

TEXT ENTRY Page

This page is used for entry of units and method identities in all applications. The title displayed at the top of the page depends on which parameter is being edited. On entry to this page the name field is filled with the current value.

METHOD ID Enter up to 11 characters for method description.
UNITS Enter up to 11 characters for Conc units.

- Use the Arrow keys to move the cursor to the required character and press ENTER. Once all the required characters have been entered press ACCEPT. Numbers are entered using the numeric keypad. If you make a mistake ○○○○● will act as a backspace or CLEAR will clear the whole entry.

(←) - Acts as backspace and clears the last character in the line.
ACCEPT - Enters the names and returns to the calling screen.
CANCEL - Abandons the naming operation and returns the calling page. The original entry is not changed.


SAVE / RENAME

This page is used for saving or renaming files from any part of the LOCAL CONTROL software. The page is used to name a file and/or change the ID.

- On coming to the page the Filename field is ready for editing. Use the Arrow keys to move the cursor to the required character and press ENTER. Once all the required characters have been entered press ACCEPT. Numbers are entered using the numeric keypad. If you make a mistake ○○○○● will act as a backspace or CLEAR will clear the whole entry.

Filenames are limited to 8 characters.

- The ID field will contain the characters entered for the method. To edit the ID highlight the ID field and press ENTER, change as necessary then press ACCEPT.

(←) - Acts as backspace and clears the last character in the line.


ACCEPT/SAVE - Accepts the edited field ready for saving. Only reads ACCEPT if edit is in progress.
SAVE - saves file (not present during an edit).

CANCEL - Abandons the naming operation. No files will be changed and the software returns to the calling screen.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

HOME Page

- This page provides a large display of the live reading at the current wavelength.
- To access the HOME page from any other page, press HOME.

UNITS – Changes the concentration units for the live display.

λ nm – Changes the current wavelength.

MODE – Selects Absorbance, %Transmittance, or Concentration mode for the live display.

PRINT - Prints the current reading


MENU – Displays the MENU page which provides access to all other areas of the software.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

METHODS LIBRARY Page

- To access the Methods Library select METHODS on the MENU page and press ENTER.

From the Methods Library, method files can be loaded, renamed, and deleted. Saving files to the library is done from the method pages of the application.

- The METHODS LIBRARY Page lists all the files, including the file type, description, and file name.

FILE TYPE

This describes the type of information contained in the file.

M FIXED	Fixed wavelength method
M QUANT	Quant method, including calibration
M RATIO	Ratio or corrected ratio method
M SCAN	Scan method
M RATE	Rate method


DESCRIPTION

This is the description entered when the file was saved, and corresponds to the Method ID.

FILENAME

This is the DOS compatible filename used by the instrument.


Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

FORMAT LIBRARY - Clears all files from the library. *THIS WILL DELETE ALL THE FILES IN THE LIBRARY.*

MENU – Displays the MENU page.

Using the Methods Library

To perform an operation on a method file first select the file. To do this move the cursor to highlight the required file using the Up/Down Arrow keys.

If the file does not appear on the list the > and < arrow keys will scroll one page at a time down and up the list. There is a slight delay after the key is pressed while the next section of the directory is loaded. Once the required file is in the window, move the cursor to it using the Up/Down arrow keys as above.

With the required file highlighted, press ENTER to display the file pop-up menu. Use the Up/Down arrow keys to highlight the required operation and press ENTER.

LOAD Loads the file from the Methods Library

RENAME Switches to the SAVE/RENAME page where file name and description can be changed.

DELETE Removes the highlighted file from the Methods Library.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

FIXED

- Instrument and analysis parameters are set up on the FIXED page. Move the cursor to the required parameter using the Up/Down Arrow keys. Change the parameter by pressing the ENTER key.

- Once the method has been set up press ZERO/BASE to zero the instrument for the current method (See page 5) and press RUN. The spectrophotometer will perform a measurement and display the result on the FIXED RESULTS page.

FIXED METHOD Page

```
          * FIXED *
MODE                : ABS

ID                  :
λ SELECT            : SINGLE λ
WAVELENGTH(S)       : 550.0 nm
INTEGRATION         : 1 s
DELAY TIME          : 00:00

AUTOPRINT           : OFF
CELL PROG           : AUTO
REF MODE            : OFF
```

MODE : Select from ABS / %T / CONC.

ABS	Selects Absorbance.
%T	Selects % transmittance.
CONC	Selects Concentration mode. (Only available when SINGLE λ is selected)

- CONC allows the result to be automatically multiplied by a factor to give concentration measurements. If selected two further parameters are added to the method, FACTOR and UNITS.

ID : Enter a description using the TEXT ENTRY screen.

λ SELECT : Select from SINGLE λ or MULTI λ

SINGLE λ	This option is used to measure each sample at a single wavelength which is the same for each sample.
MULTI λ	This option allows each sample to be measured at up to 20 wavelengths, which are the same for each sample.

WAVELENGTH(S) :

SINGLE λ	Use the numeric key pad to enter the required wavelength into the pop-up box. Press ENTER when finished.
MULTI λ	Use the up and down arrow keys to move to the wavelength to be entered or edited and press ENTER to display the edit box. Use the numeric keypad to enter the wavelength and press ENTER when finished. The instrument returns to the MULTI λ screen with the next wavelength in the list highlighted. Up to 20 wavelengths may be entered. When the list is finished press the ACCEPT function key to accept the new list or the CANCEL function key to return to the FIXED METHOD page without changing the wavelength list.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

The Deuterium lamp will be switched on if a wavelength less than 325.0 nm is selected for a Gamma. The current data will be lost if the wavelength is changed.

INTEGRATION : Enter integration time in seconds.

- This sets the integration time for which the result is measured.

The current data will be lost if the integration time is changed.

DELAY TIME : Set a time in the range 00.05 to 99.59, using . to separate minutes and seconds.

- This sets the time between pressing RUN and the start of the measurement. The range is from 0 to 99 minutes and 59 seconds.

FACTOR : Enter the factor for the concentration.

Only available in CONC mode, the factor is used to multiply the absorbance result to produce a concentration result. Changing the factor will not cause current results to be lost. They will be recalculated using the new factor.

UNITS : Enter units for concentration using the TEXT ENTRY page.

Only available in CONC mode, used to enter the required description of the concentrations up to 11 alphanumeric characters.

AUTOPRINT : Toggles between ON and OFF.

When **ON,** results are printed automatically after each run.


Before attempting to print at any point during operation of the instrument ensure that the printer is ready to print, ie switched on, on line and supplied with paper. Failure to do so will result in an error condition. Press CLEAR to clear the error message (the system may take a little time to respond). Ensure the printer is ready and retry.

CELL PROG : Switches to the CELL PROG. screen (if fitted).

REF MODE : Toggles the status of the Cell Programmer reference mode (if fitted).

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

FIXED Method Page Function Keys

VIEW RESULTS - Switches to the FIXED RESULTS page.

SAVE METHOD - Brings up the SAVE page from where the method can be saved to the Methods Library.

MENU – Switches to the MENU page.

- Pressing RUN starts a fixed measurement using the current method and switches to the FIXED RESULTS page.
- Pressing ZERO starts a zero using the current method.

Any changes to the WAVELENGTH or INTEGRATION parameters will invalidate the current results.

If AUTOPRINTING is selected, a change to the MODE parameter will invalidate the current results.


FIXED RESULTS Page

- The layout of the page depends on the Mode and λ Select in use.

SINGLE λ	In ABS or %T modes up to 2 columns of results are displayed per page. In CONC mode a single column of results is displayed on each page. Results accumulate on the same page until it is full.
MULTI λ	Two columns of results are displayed per page. Results of each Sample always start on a new page.

- To move up or down pages of results use the Up/Down arrow keys.
- Results are numbered sequentially from 1 to 600. Conc values are from 0.0001 to 9999.9. Any concentration value outside this range is marked "UNDER RANGE" or "OVER RANGE".

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

CLEAR RESULTS - All results are cleared ready for the start of the next batch.

FIXED PAGE - Returns to the FIXED METHOD page.

PRINT LIST - Prints the list on the selected printer.

- Press RUN to take another sample measurement.
- Press ZERO/BASE to zero the instrument at the wavelength(s) specified in the method. (See page 5).

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

QUANT

Instrument and analysis parameters are set up on the QUANT page. Move the cursor to the required parameter using the Up/Down Arrow keys. Change the parameter by pressing the ENTER key.

QUANT METHOD Page

```
          * QUANT *

ID                    :
WAVELENGTH            : 550.0 nm
INTEGRATION           : 1 s
STANDARDS             : 3
REPLICATES            : 1
UNITS                 :
CURVE FIT             : LINEAR

AUTOPRINT             : OFF
CELL PROG             : AUTO
REF MODE              : OFF
```

ID : Enter a description using the TEXT ENTRY screen.

WAVELENGTH : Select a wavelength between 190.0 nm for Gamma or 315.0nm for Delta, and 1100.0nm.

The Deuterium lamp will be switched on if a wavelength less than 325.0 nm is selected for a Gamma. The current data will be lost if the wavelength is changed.

INTEGRATION : Enter integration time in seconds.

This sets the integration time for which the result is measured.
The current data will be lost if the integration time is changed.

STANDARDS : Brings up the Standards Entry Page

Use the up and down arrow keys to move through the list of standards. When the standard to be entered or edited is highlighted, press ENTER to display the Edit pop-up. Use the numeric keys to enter the concentration of the standard and press ENTER when finished. The instrument returns to the Standards Entry page with the highlight on the next standard on the list. Up to 20 standards can be specified.

Changing the standards will cause any current data to be lost.

REPLICATES : Enter number of replicates for each standard.

- Sets the number of times each standard is measured (maximum 3). Each value obtained is used in the calibration.

UNITS : Enter units for concentration using the TEXT ENTRY page.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

CURVE FIT : Select from LINEAR / LINEAR TO 0 / QUADRATIC / QUAD TO 0.

- Selects the curve fit algorithm used in the calibration.

LINEAR	performs a linear calibration. At least two standards are required.
LINEAR TO 0	performs a linear calibration forced through zero.
QUADRATIC	performs a quadratic fit on the data. At least three standards are required.
QUAD TO 0	performs a quadratic fit with the data forced through zero. At least two standards are required.

Changing the curve fit will cause the existing calibration to be recalculated. Any results associated with the previous calibration will be lost.

AUTOPRINT : Toggles between ON and OFF.


When **ON,** results are printed automatically after each run.

Before attempting to print at any point during operation of the instrument ensure that the printer is ready to print, ie switched on, on line and supplied with paper. Failure to do so will result in an error condition. Press CLEAR to clear the error message (the system may take a little time to respond). Ensure the printer is ready and retry.

CELL PROG : Switches to the CELL PROG. screen (if fitted).

REF MODE : Toggles the status of the Cell Programmer Reference Mode (if fitted).

QUANT Page function keys

VIEW CALIB - Switches to the QUANT GRAPH page if a valid calibration exists.

VIEW RESULTS - Switches to the Quant Results page if there are any sample results.

SAVE METHOD - Brings up the SAVE page from where the method can be saved to the Methods Library.


CALIBRATE – Starts the measurements of the standards.

MENU – Switches to the MENU page.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

QUANT STANDARDS Page

- This page lists the standards as defined in the QUANT METHOD. Before the system can be calibrated each standard must have a concentration entered.
- If a calibration has been done then the correlation coefficient and the equation are displayed.
- If a calibration has not been done pressing RUN causes the warning prompt "CANNOT RUN WITHOUT CALIBRATION" to appear, otherwise it takes a sample measurement and switches to the Quant Results screen. Pressing ZERO/BASE starts a zero using the current method.

VIEW CALIB - Switches to the QUANT GRAPH page if a valid calibration exists

VIEW RESULTS - Switches to the QUANT RESULTS page if there are any results in the table.

QUANT PAGE - Switches back to the QUANT METHOD page.

EDIT STD – Available after calibration. Allows each standard to be used, ignored or remeasured.


EDIT CURVE – Available after calibration. Allows curve fit to be changed.

QUANT CALIBRATION

- Press ZERO/BASE to zero the instrument with the current method. (See page 5).
- To calibrate the system press return to the QUANT METHOD page and press CALIBRATE. The QUANT CALIBRATION graph will be displayed and the instrument will prompt for each standard (and replicate) in turn. As the measurements of the standards proceed the datapoints are marked on the graph. When all the standards have been measured the system calculates the equation, rescales the graph then draws and displays the line of best fit on the graph.
- A calibration can be stopped by pressing the STOP function key. The calibration will be aborted and the software will return to the QUANT METHOD page. Any values obtained will be lost.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

- If a calibration has not been done pressing RUN causes the warning prompt "CANNOT RUN WITHOUT CALIBRATION" to appear, otherwise it takes a sample measurement and switches to the Quant Results screen.

VIEW RESULTS - Switches to the Quant Results page if there are results.

QUANT PAGE - Switches back to the Quant Method page.


STANDARDS - Switches to the Quant Standards page.

PRINT GRAPH - Prints the Quant method and calibration graph.

SAVE METHOD - Switches to the SAVE page.

QUANT RESULTS Page

- To view further pages of results use the Up and Down arrow keys.
- Pressing RUN takes another sample measurement and displays the result.
- Results are numbered sequentially and any batch can be of up to 600 samples.

CLEAR RESULTS - Deletes all results from the results table.

QUANT PAGE - Switches back to the Quant Method page.

PRINT LIST - Gives a printout of the results.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

RATIO

- Instrument and analysis parameters are set up on the RATIO METHOD page. Move the cursor to the required parameter using the Up/Down Arrow keys. Change the parameter by pressing the ENTER key.

- Once the method has been set up press ZERO/BASE to zero the instrument for the current method (see page 5) and press RUN. The spectrophotometer will perform a measurement and display the result on the RATIO RESULTS page.

RATIO METHOD Page

```
          * RATIO *
MODE              : ABS

ID                :
λ1                : 550.0 nm
λ2                : 550.0 nm
CORR. RATIO       : YES
CORR. λ           : 550.0 nm
INTEGRATION       : 1 s
DELAY TIME        : 00:00

AUTOPRINT         : OFF
CELL PROG         : AUTO
REF MODE          : OFF
```

MODE : Select from ABS or %T.

ABS Selects Absorbance.
%T Selects % Transmittance.

ID : Enter a description using the TEXT ENTRY screen.

λ1, λ2, CORR. λ : Use the numeric key pad to enter the required wavelength into the pop-up box. Press ENTER when finished.

The Deuterium lamp will be switched on if a wavelength less than 325.0 nm is selected for a Gamma. The current data will be lost if the wavelength is changed.

CORR.RATIO : Toggles between YES and NO with ENTER.

- When set to NO, the formula Abs(λ1) / Abs(λ2) is used to calculate the ratio. When corrected ratios are selected, the formula (Abs(λ1) – Abs(λC)) / (Abs(λ2) – Abs(λC)) is used instead.

INTEGRATION : Enter integration time in seconds.

- This sets the integration time for which the result is measured.

The current data will be lost if the integration time is changed.

DELAY TIME : Set a time in the range 00.05 to 99.59, using . to separate minutes and seconds.

- This sets the time between pressing RUN and the start of the measurement. The range is from 0 to 99 minutes and 59 seconds.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

AUTOPRINT : Toggles between ON and OFF.


When **ON,** results are printed automatically after each run.

Before attempting to print at any point during operation of the instrument ensure that the printer is ready to print, ie switched on, on line and supplied with paper. Failure to do so will result in an error condition. Press CLEAR to clear the error message (the system may take a little time to respond). Ensure the printer is ready and retry.

CELL PROG : Switches to the CELL PROG. screen (if fitted).

REF MODE : Toggles the status of the Cell Programmer reference mode (if fitted).

RATIO Method Page Function Keys

VIEW RESULTS - Switches to the RATIO RESULTS page.

SAVE METHOD - Brings up the SAVE page from where the method can be saved to the Methods Library.

MENU – Switches to the MENU page.

- Pressing RUN starts a ratio measurement using the current method and switches to the RATIO RESULTS page.
- Pressing ZERO starts a zero using the current method.


Any changes to the WAVELENGTHS or INTEGRATION parameters will invalidate the current results.

If AUTOPRINTING is selected, a change to the MODE parameter will invalidate the current results.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

RATIO RESULTS Page

- To move up or down pages of results use the Up/Down arrow keys.
- Results are numbered sequentially from 1 to 600.

CLEAR RESULTS - All results are cleared ready for the start of the next batch.

RATIO PAGE - Returns to the RATIO METHOD page.

PRINT LIST - Prints the list on the selected printer.

- Press RUN to take another sample measurement.
- Press ZERO/BASE to zero the instrument at the wavelength(s) specified in the method. (See page 5).

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

SCAN

- To select Scan highlight the SCAN option on the MENU PAGE and press ENTER. The SCAN Methods page is displayed and from here the instrument and analysis parameters can be set up.


- Move the cursor to the required parameter using the Up/Down Arrow keys. Press ENTER to enable a parameter to be changed.

- Once the method has been set up press ZERO/BASE to perform a baseline scan with the current method (see page 5) and then press RUN. The spectrophotometer will perform the scan and display the result on the SCAN GRAPH page. From here the spectrum can be manipulated and printed.

- The scan speed is set to 600 nm per min., and the data interval is 0.5nm.

SCAN PARAMETERS Page

Note: The current spectrum will be lost if any of the method parameters are changed.

MODE : Select from ABS / %T / I.

ABS	Selects Absorbance.
%T	Selects % transmittance.
I	Selects Intensity mode. This will measure the intensity of the signal in the sample beam.

The current spectrum will be lost if the Scan Mode is changed.

ID : Enter a description using the TEXT ENTRY screen.

START : Select a wavelength between 190.0 nm (315.0 nm for Delta) and 1096.0 nm.

- Selects start wavelength.

The Deuterium lamp will be switched on if a wavelength less than 325.0 nm is selected for a Gamma. The Start wavelength must be at least 4 nm less than Stop wavelength. The current spectrum will be lost if the Start wavelength is changed.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

STOP : Select a wavelength between 194.0 nm (319.0 nm for Delta) and 1100.0 nm.

- Selects stop wavelength.

The Stop wavelength must be at least 4 nm greater than the Start wavelength. The current spectrum will be lost if the Stop wavelength is changed.

PEAK TABLE : Select from OFF / PEAKS / TRACK / RATIO / CORR. RATIO / PEAK HEIGHT.

This selects the type of peak picking done automatically as part of the method. Results are reported on the Peaks Page..

OFF Sets Peak Table to Off. No peaks information is produced as part of the scan.

PEAKS Picks the highest peaks in a spectrum up to a maximum of 10 peaks.

TRACK This function allows the Absorbance (or other mode) values to be reported at up to 10 user selected wavelengths. To enter the Desired wavelengths select PEAK TABLE TRACK then press VIEW GRAPH.

You do not have to have a spectrum present on the graph to enter the selected wavelengths. Press MANIPULATE and select TRACK. For each wavelength move the cursor to the desired position press ENTER. Once all the wavelengths have been entered go back to the SCAN METHOD page and save the method.

RATIO This function allows the ratio λ_1/λ_2 to be automatically calculated at the end of the scan. To enter the desired wavelengths select PEAK TABLE and press ENTER then select RATIO. A pop up box Appears in which to enter the first wavelength. Enter the desired wavelength and press ENTER. Repeat for the second wavelength. Once all the method parameters have been set save the method.

CORR RATIO This function allows the ratio of two wavelengths to be calculated relative to a third wavelength $(\lambda_1-\lambda_3)/(\lambda_2-\lambda_3)$ automatically at the end of a scan. To enter the desired wavelengths select PEAK TABLE and press ENTER then select CORR RATIO. A pop up box appears in which to enter the first wavelength. Enter the desired wavelength and press ENTER. Repeat for the second and correction wavelengths. Once all the method parameters have been set save the method.

PEAK HEIGHT This function allows the height of a peak to be calculated relative To a local baseline rather than y = 0. To enter the desired wavelengths, select PEAK TABLE and press ENTER then select PEAK HEIGHT. A pop-up box appears in which to enter the wavelengths required. λ_1 and λ_3 define the baseline, λ_2 defines the peak. Once all the parameters have been set, save the method.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

GRAPH HIGH : Select from range (GRAPH LOW + 0.01) to 6.00.

Sets the upper graph limits on the SCAN GRAPH page.
GRAPH HIGH must be 0.01 greater than GRAPH LOW.

GRAPH LOW : Select from range -0.3 to (GRAPH HIGH- 0.01).

Sets the lower graph limits on the SCAN GRAPH page.
GRAPH LOW must be 0.01 less than GRAPH HIGH.

AUTOPRINT : Toggles between ON and OFF.


When **ON,** results are printed automatically after each run.

Before attempting to print at any point during operation of the instrument ensure that the printer is ready to print, ie switched on, on line and supplied with paper. Failure to do so will result in an error condition. Press CLEAR to clear the error message (the system may take a little time to respond). Ensure the printer is ready and retry.

CELL PROG : Appears when the Cell Programmer is fitted. Switches to the CELL PROG screen.

REF MODE : Appears when the Cell Programmer is fitted. Toggles the status of the Cell Programmer reference mode.

SCAN PARAMETERS page function keys

VIEW RESULTS - Switches to the Scan Peak Table screen if any of the peak functions have been performed or the Track Table screen if track has been used.

VIEW GRAPH - Switches to the Scan Graph screen.


SAVE METHOD - Brings up the Filename Function screen and then saves the method including ID and track wavelengths if the PEAK TABLE parameter is set to TRACK.

MENU – Switches to the MENU page.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

SCAN GRAPH Page

- This page displays spectra and allows them to be manipulated.

VIEW RESULTS - Switches to the PEAK TABLE page.

SCAN PAGE - Returns to the SCAN page.

CHART OUTPUT – Plots the spectrum on an analogue chart recorder.

PRINT GRAPH - Provides a hardcopy of the results as shown on screen.

MANIPULATE - Displays the Manipulate options.

- Pressing RUN starts a scan using the current method.
- Pressing ZERO/BASE starts a baseline using the current method. (See page 5).

MANIPULATE OPTIONS

MANIPULATE
TRACK RESCALE MODE PEAKS ORIGINAL

TRACK Reports x and y axis values using the tracking cursor.

RESCALE Changes x and y axis scales automatically or manually.

MODE Changes mode. Select from %T or ABS.

PEAKS Finds spectral peaks. Select from PEAKS / RATIO / CORR. RATIO / PK HEIGHT.


ORIGINAL Resets the graph to display the data as originally collected.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

TRACK

- To move the vertical cursor across the screen use the Left and Right Arrow keys. The cursor always moves to a data point regardless of the displayed scales. Pressing ENTER places a marker at the current wavelength. Up to 10 wavelengths can be selected.

TRACK page function keys

- Pressing CLEAR will delete markers in turn, highest number first.
- The x-axis values are listed on the TRACK table page. Further markers can be added to the spectrum at any time; however selecting TRACK will cause any previous PEAK TABLE information to be lost.

VIEW TABLE - Switches to the TRACK TABLE page.

FAST/SLOW - Toggles between two cursor speeds. In FAST mode the cursor jumps 5% of the graph or to the next data point whichever is the greater. In SLOW mode the cursor jumps to the next data point or the next display pixel whichever is the greater. The function key label shows the next speed ie the opposite to the one selected.

CLEAR ALL - Clears all the markers and the TRACK TABLE.

PRINT GRAPH - Provides a hardcopy of the results showing the markers and x and y-axis values.

SCAN GRAPH - Returns to the SCAN GRAPH page.

RESCALE

RESCALE
AUTO
GRAPH HIGH
GRAPH LOW
GRAPH START
GRAPH STOP
PROCEED

- This option gives pop-up menus for changing the graph x and y -axis scales.
- Move the cursor to one of the options and press ENTER to select an operation.

AUTO	Displays the SCAN GRAPH with the x and y-axes rescaled so that the Spectrum fills the screen.
GRAPH HIGH	Pops up a window for entry of the GRAPH HIGH limit.
GRAPH LOW	Pops up a window for entry of the GRAPH LOW limit.
GRAPH START	Pops up a window for entry of the required start wavelength.
GRAPH STOP	Pops up a window for entry of the required stop wavelength.
PROCEED	Used after GRAPH HIGH, GRAPH LOW, GRAPH START or GRAPH STOP to return to the SCAN GRAPH page with the graph rescaled using the new parameters.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

MODE : Select from ABS / %T.

ABS	Selects Absorbance.
%T	Selects % transmittance.

PEAKS

FUNCTION
PEAKS
RATIO
CORR RATIO
PK. HEIGHT

- This option enables the spectrum to be automatically searched for peaks. Move the cursor to one of the options and press ENTER to perform a search. When the search is complete the spectrum is displayed with the peak positions marked. For a peak to be found there must be more than 15 data points between that point and a previous peak.


For RATIO and CORR RATIO, enter the wavelengths as prompted. All results can be viewed by pressing ○○○○● VIEW RESULTS.

PEAKS	Marks the 10 highest peaks.
RATIO	Calculates the ratio λ_1/λ_2.
CORR RATIO	Calculates the ratio $(\lambda_1-\lambda_3)/(\lambda_2-\lambda_3)$.
PK HEIGHT	Calculates the peak maximum relative to a local baseline.

ORIGINAL

- This removes any manipulation and displays the spectrum as originally collected and specified by the scan method.

TRACK TABLE Page

- The list shows the y-axis values of the spectrum for the wavelengths marked during TRACK. The measured values will be ABS, %T, or INTENSITY depending on the current mode.

VIEW GRAPH - Switches to the TRACK page.


PRINT LIST - Prints the list on the selected printer.

CHART OUTPUT – Plots the spectrum on an analogue chart recorder.

SCAN GRAPH - Returns to the SCAN GRAPH page.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

PEAK TABLE Page


- The list shows the positions and values of the peaks as calculated by the function selected in MANIPULATE PEAKS. The measured values will be ABS, %T, or INTENSITY depending on the current mode and are sorted by wavelength.

CHART OUTPUT – Plots the spectrum on an analogue chart recorder.

PRINT LIST - Prints the list on the selected printer.

SCAN GRAPH - Returns to the SCAN GRAPH page.

RATIO TABLE Page

- The page shows the positions and values of the wavelengths and the ratio as selected by the RATIO or CORR RATIO functions.


VIEW GRAPH - Returns to the SCAN GRAPH page.

CHART OUTPUT – Plots the spectrum on an analogue chart recorder.

PRINT LIST - Prints the list on the selected printer.

SCAN GRAPH - Returns to the SCAN GRAPH page.

PEAK HEIGHT Page

- This page shows the position and values of the wavelengths and the peak height as selected by the PEAK HEIGHT function.

VIEW GRAPH - Returns to the SCAN GRAPH page.

CHART OUTPUT – Plots the spectrum on an analogue chart recorder.

PRINT LIST - Prints the list on the selected printer.

SCAN GRAPH - Returns to the SCAN GRAPH page.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

RATE

- To select Rate highlight the RATE option on the MENU PAGE and press ENTER. The RATE Methods page is displayed and from here the instrument and analysis parameters can be set up.

- Move the cursor to the required parameter using the Up/Down Arrow keys. Press ENTER to enable a parameter to be changed.

- Once the method has been set up press RUN. The spectrophotometer will perform the rate and display the result on the RATE GRAPH page. From here the data can be manipulated and printed.

RATE PARAMETERS Page

Note: The current data will be lost if any of the method parameters (except for the ID, Slope, Factor, or Units) are changed.

```
            * RATE
*
ID                    : TEST 1
WAVELENGTH            : 340.0  nm
MEASURE TIME          : 00:30
DELAY TIME            : 00:00
SLOPE                 : POSITIVE
RANGE                 : 0.5
FACTOR                : 1.000
UNITS                 : I/U

AUTOPRINT             : OFF
CELL PROG             : AUTO
REF MODE              : OFF
```

ID : Enter a description using the TEXT ENTRY screen.

WAVELENGTH : Select a wavelength between 190.0nm for Gamma or 315.0nm for Delta, and 1100.0 nm.

The Deuterium lamp will be switched on if a wavelength less than 325.0 nm is selected for a Gamma. Any current data will be lost if the wavelength is changed.

MEASURE TIME : Set a time in the range 00:05 to 99:59.

- This sets the time over which the sample will be measured. The range is from 5 seconds to 99 minutes and 59 seconds in steps of 1 second.

- If the cell programmer is ON then the measurement time sets the time between individual measurements on the first cell (i.e. the time for each cycle).

DELAY TIME : Set a time in the range 00:05 to 99:59.

This sets the time between pressing RUN and the start of the measurement. The range is from 0 to 99 minutes and 59 seconds in steps of 1 second.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

SLOPE : Select from POSITIVE or NEGATIVE.

- Sets the graph to display positive or negative changes in Absorbance. Choose POSITIVE if Absorbance increases with time. Choose NEGATIVE if Absorbance decreases with time.

RANGE : Set a number in the range 0 to 3A.

- This sets the graph y-axis. Enter a number slightly larger than the change in Absorbance expected.

FACTOR : Enter the factor for Activity as a number in the range 0.001 to 9999.999.

UNITS : Enters the units of Activity using the TEXT ENTRY page.

- Enters the required description or units of Activity up to 11 alphanumeric characters.

AUTOPRINT : Toggles between ON and OFF.


When **ON,** results are printed automatically after each run.

Before attempting to print at any point during operation of the instrument ensure that the printer is ready to print, ie switched on, on line and supplied with paper. Failure to do so will result in an error condition. Press CLEAR to clear the error message (the system may take a little time to respond). Ensure the printer is ready and retry.

CELL PROG : Switches to the CELL PROG. screen (if fitted).

REF MODE : Toggles the status of the Cell Programmer reference mode (if fitted).


RATE PARAMETERS Page function keys

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

RATE GRAPH Page

- This page displays the RATE curve and allows it to be manipulated.

VIEW RESULTS - Switches to the RATE RESULTS page.

RATE PAGE - Returns to the RATE PARAMETERS page.

PRINT - Provides a hardcopy of the results i.e. RATE GRAPH and RATE RESULTS.

MANIPULATE - Displays the Manipulate options.

- Pressing RUN starts a measurement using the current method.
- Pressing ZERO/BASE zeros the instrument using the wavelength specified in the current method. (See page 5).

MANIPULATE

MANIPULATE
TRACK RESCALE ABS DISPLAY ORIGINAL ANOTHER CELL

TRACK Sets the start and stop time for the rate calculations.

RESCALE Changes y-axis scale automatically or manually.

ABS DISPLAY Sets the display of absolute absorbance values or absorbances Relative to the first measurement(s) of the run.

ORIGINAL Resets the graph to display the data as originally collected.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com


ANOTHER CELL Only present if the Cell Programmer has been used. Enables the results of another cell from the same run to be displayed.

TRACK

- To move the vertical cursor across the screen use the Left and Right Arrow keys. The cursor always moves to a data point regardless of the displayed scales. Pressing ENTER places a marker at the current time.
- To delete a marker, place the cursor on the marker and press CLEAR.
- The x-axis values are used to recalculate the rate of change of Absorbance between the new start and stop times. Results are listed on the RATE RESULTS page.

Up to four discrete pairs of cursors can be placed on the graph. Arrows are placed on the cursors and results are displayed on the RATE RESULTS page for those parts of the graph indicated by the arrows.

- The minimum interval between TRACK cursors is one second.

VIEW RESULTS - Switches to the RATE RESULTS page.

FAST/SLOW - Toggles between two cursor speeds. In FAST mode the cursor jumps 5% of the graph or to the next data point whichever is the greater. In SLOW mode the cursor jumps to the next data point or the next display pixel whichever is the greater.The function key label shows the next speed ie the opposite to the one selected.

CLEAR ALL - Clears all the markers.

RATE GRAPH - Returns to the RATE GRAPH page.


RESCALE

- This option gives a pop-up menu for changing the graph y -axis scale.
- Move the cursor to one of the options and press ENTER to select an operation.

The rescale options depend on whether Absolute or Relative Absorbance has been selected.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com


Absolute Absorbance:

AUTO - Displays the RATE GRAPH with the y-axis rescaled so that the trace fills the screen.

GRAPH HIGH, GRAPH LOW – Allow the user to set the upper and lower limits for the RATE GRAPH y-axis.

Relative Absorbance:

AUTO - Displays the RATE GRAPH with the y-axis rescaled so that the trace fills the screen.

RANGE - Allows the user to set the upper y-axis value.

ABS DISPLAY

- This option enables the RATE GRAPH to be redisplayed with Absolute Absorbance or Relative Absorbance values.

ORIGINAL

This removes any manipulation and displays the rate graph as originally specified by the rate method.

ANOTHER CELL


- If more than one rate has been run in parallel using the Cell Programmer then this function allows the results from any cell in the run to be selected and displayed. Enter the number of cells you wish to display.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

RATE RESULTS Page

- The Rate Results page displays the Initial and Final Absorbance, Initial and Final Time, the change in absorbance per minute, calculated Activity, Correlation Coefficient of the best fit line and finally the smoothing parameter used.

If the rate curve has been tracked the Initial and Final Absorbance with the Initial and Final Time will reflect those chosen by the two cursors.

VIEW GRAPH - Returns to the RATE GRAPH page.

RATE PAGE - Returns to the RATE PARAMETERS page.

PRINT - Provides a hardcopy of the results i.e. RATE GRAPH and RATE RESULTS.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

PARALLEL RATE MEASUREMENTS

- The 7 Cell Programmer can be used in conjunction with the RATE software to measure between 2 and 7 cells in parallel.
- For the Cell Programmer to be used it must be ON with MODE set to AUTO.
- The RATE METHOD setup and MANIPULATE functions are exactly the same except that the MEASURE TIME now sets the time between each cycle, i.e. the length of time between measurements on the first cell. The number of measurements taken on each cell is set by the CELL CYCLES parameter on the cell programmer page. The total time over which the measurements are made is the result of the time between measurements (MEASURE TIME) and the number of measurements (CELL CYCLES). For example an analysis using 4 cells with MEASURE TIME set to 15 seconds and CELL CYCLES set to 20 would give a total measurement time of 5 minutes.

To set up a parallel rate method do the following:

- Set up the method as normal on the RATE PARAMETERS page.

Note that MEASURE TIME specifies the time between each measurement on the first cell.

- Go to the Cell Programmer page. Press HOME then ACCESSORIES and select CELL PROG.
- Set up the cell programmer method as required.
- Return to the RATE PARAMETERS page, press ZERO/BASE to zero the instrument (if required) then press RUN.
- The analysis will proceed with the rate graph for the first cell being drawn on screen.
- To view the result of any other cell press MANIPULATE on the RATE GRAPH page and select ANOTHER CELL. The cell number currently displayed is shown to the right of the ID: line.
- To print the results press PRINT on either the RATE GRAPH page or the RATE RESULTS page.

Three further print options are available: ALL OVERLAY, ALL SEQUENTIAL, and ONE RESULT.

ALL OVERLAY	Prints the results of all cells in the run together in batches of 4 (3 if using the internal printer). TRACK markers are not included.
ALL SEQUENTIAL	Prints each result in the run separately. TRACK markers are included.
ONE RESULT	Prints the result currently displayed. TRACK markers are included.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

SETUP

Overview of Setup Options

- This section describes how to set up the instrument.
- The main system setup options are available directly from the SETUP function key on the MENU page.

SETUP Page

- From the SETUP page move the cursor to the required option using the Up/Down Arrow keys. Select the option by pressing the ENTER key.

```
     * SETUP *

CLOCK
PRINTER
ENVIRONMENT
WAVELENGTH
INITIALISE
WHITE LIGHT
RECORDER
```

CLOCK : Switches to the CLOCK page.
PRINTER : Switches to the PRINTER page.
ENVIRONMENT : Switches to the ENVIRONMENT page.
WAVELENGTH : Switches to the WAVELENGTH CALIBRATION page(Gamma only).
INITIALISE : Switches to the OPTICAL INITIALISATION page.
WHITE LIGHT : Switches to the WHITE LIGHT Page.
RECORDER : Switches to the RECORDER Page.

CLOCK Page

```
     * CLOCK *

     * TIME *

HOURS  :  16
MINS   :  32



     * DATE *

DAY    :  25
MONTH  :  12
YEAR   :  96
```

- From this page the internal spectrophotometer clock/calendar can be reset.
- To reset the time or date highlight the required parameter and press ENTER. Enter the new value using the number keys and press ENTER.
- Once all the parameters have been changed press ACCEPT.

The date or time will not be changed unless ACCEPT is pressed.

- Pressing CANCEL cancels the edit, leaving the previous values unchanged.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

PRINTERS Page

- This page sets the system to work with the selected printer.
- The printer type is always highlighted. To choose a printer press ENTER to display the list of supported printers, and select using cursor keys. Press ENTER to confirm entry.

PRINTER
EPSON 9 PIN INTERNAL

Printer options	Supported Printers
EPSON 9 PIN	Epson 9 or 24 Pin Dot Matrix, using ESC/P language.
INTERNAL	Supplied thermal printer only

Printers not on the above list that claim Epson 9 pin / 24 pin / ESC/P compatibility should work with the instrument but are not guaranteed to do so and are therefore not supported. If in doubt contact your local Thermo Spectronic approved Customer Support Organisation.

Note: Printers designed to work only in a Windows environment are not compatible with Local Control Software.

Before attempting to print using an external printer at any point during operation of the instrument, ensure that the printer is ready to print. Failure to do so will result in an error condition. Press CLEAR to clear the error message. Then rectify the problem with the printer, and try again.

ENVIRONMENT page

This page is used to select the language used for the software, the use of the beep, date format, etc.

```
    * ENVIRONMENT *

LANGUAGE              : ENGLISH
SOUND                 : OFF
DATE FORMAT           : dd/MM/yy

AUTOPRINT RESULTS     : OFF
PROTECT METHODS       : OFF
CHANGE PASSWORD       :
HIDDEN PASSWORD       : OFF
```

LANGUAGE : Select from the list. The language used immediately changes to the one selected.

SOUND : Turns the warning beeper ON or OFF. If set to OFF then the only indication of any error is the screen message.

DATE FORMAT: Defaults to dd/MM/yy, but will toggle with ENTER to MM/dd/yy.

AUTOPRINT RESULTS : Toggles between ON and OFF.

When **ON,** results are printed automatically after each run.

Before attempting to print at any point during operation of the instrument ensure that the printer is ready to print, ie switched on, on line and supplied with paper.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Failure to do so will result in an error condition. Press Clear to clear the error message (the system may take a little time to respond). Ensure the printer is ready and retry.

PROTECT METHODS : Defaults to OFF. When ON methods are protected against unauthorised alteration. Changing between ON and OFF is password protected.

- To turn on Protect Methods, highlight it in the menu and press ENTER. This brings up the text entry page for the password to be entered.
- The default password is “ADMIN”. Note that the password is case sensitive, and in this password all letters are upper case.
- When the correct password has been entered, press ACCEPT to return to the Environment Page. Protect Methods will now be ON. To return to OFF select and enter password as above.

CHANGE PASSWORD : Enables the password to be changed.

- It is strongly recommended that a new password is entered as soon as possible.
- Note that passwords are case sensitive.
- To change the password, select Change Password to bring up the Text Entry Page, enter the existing password and press ACCEPT. The message ENTER NEW PASSWORD appears at the bottom of the Text Entry Screen. Delete the old password by pressing the CLEAR key, or using the right hand function key to delete one character at a time, and enter the new one. Press ACCEPT when finished. You will then be prompted to enter the new password a second time for verification. When this has been done correctly the system returns to the Environment Page with the password changed.
- An error at any stage in this procedure will result in the system returning to the Environment Page with the old password unchanged.

HIDDEN PASSWORD : Enables the password to be displayed or hidden. Changing between ON and OFF is password protected.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

WAVELENGTH CALIBRATION Page (Gamma only)

```
* WAVELENGTH CALIBRATION *

CALIBRATION    :   USING D2 LAMP
```

- From this page the Wavelength calibration of the instrument can be adjusted.

WARNING: DO NOT ATTEMPT TO RECALIBRATE THE INSTRUMENT UNLESS YOU ARE ABSOLUTELY SURE YOU NEED TO DO SO. IF IN ANY DOUBT CONTACT THERMO SPECTRONIC APPROVED CUSTOMER SUPPORT.

ENSURE THE SAMPLE BEAM IS CLEAR BEFORE ATTEMPTING ANY WAVELENGTH CALIBRATION.

- This option will optimise the wavelength calibration of the instrument by measuring the Deuterium lamp emission line at 656.1nm and fitting the calibration accordingly. It should only be attempted if for any reason the instrument no longer achieves its quoted wavelength accuracy specification.

The calibration will take at least 10 minutes.

- Ensure Deuterium Lamp is ON.
- Ensure instrument is fully warmed up and the beams are clear.
- From WAVELENGTH CALIBRATION page press CALIBRATE.
- Once calibration is complete switch the instrument OFF.
- Switch instrument and Deuterium lamp ON and allow to warm up.
- Run a new Default Baseline.

OPTICAL INITIALISATION Page

- This page is used to reset the instrument and define its initialisation and default baseline. These procedures ensure the optimum performance of the spectrophotometer.

INITIALISATION TYPE: When selected this displays a pop-up menu to choose between initialisation of optics or baseline.

OPTICS - During initialisation the instrument performs some simple hardware checks, calculates various data tables and measures the dark current. The filter wheel is then initialised before the instrument drives to the default wavelength and performs an autozero.

BASELINE - This re-measures the default baseline. Ensure that both lamps are on and that the spectrophotometer is fully warmed up. This process will take about one hour.

The default baseline should be re-measured whenever one of the source lamps is changed or if the instrument is working at temperatures significantly different from 25°C, or if the wavelength calibration is altered.

INITIALISE WITH D2 : (Gamma only) This sets the instrument to initialise with or without the Deuterium lamp on. If set to ON then the instrument will automatically strike the Deuterium lamp during initialisation.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

WHITE LIGHT Page

- The **WHITE LIGHT** feature is used to facilitate alignment of optical accessories in the sample compartment.
- When the **INITIALISE** function key is pressed the instrument will align the grating so that the zero order diffraction passes through the sample compartment. This provides a beam of white light which can be seen when a white card or similar is placed in the light path.
- When alignment is completed, pressing the **STOP** function key returns the grating to its normal position and pressing the **SETUP PAGE** function key returns to the Setup Page.

RECORDER Page

- The Chart High and Chart Low parameters set the 0 to 1V full scale deflection (fsd) on the analogue chart recorder output for each of the available measurement modes.

```
        * RECORDER
*

CHART HIGH (ABS)      : 3.000
CHART LOW (ABS)       : -0.200

CHART HIGH (%T)       : 200.0
CHART LOW (%T)        : 0.1

CHART HIGH (I)        : 99.9999
CHART LOW (I)         : 0.0000
```


- On startup, these limits are set to the maximum measurement ranges (as shown above).
- To reset the limits highlight the required parameter and press ENTER. Enter the new value using the number keys and press ENTER.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

LAMPS Page - Gamma

- The lamp functions are available directly from the LAMPS key on the SETUP and MENU pages.
- This page shows the status of the Tungsten Halogen and Deuterium lamps whether ON, OFF or FAILED, and their appropriate energy levels. It also allows the lamp hours to be reset and the Deuterium lamp to be switched on or off.
- The HOURS parameter states the number of hours that the lamp has been in use.

The Tungsten Halogen lamp should be replaced after 2000 hours. The Deuterium lamp should be replaced after 1000 hours.
Whenever a lamp is changed then the hours parameter should be reset to zero.

RESET W HRS - Resets hours and measures energy at the appropriate wavelength. **ALLOW THE LAMP AT LEAST 10 MINUTES TO WARM UP BEFORE RESETTING ITS HOURS.**

RESET D2 HRS - Resets hours and measures energy at the appropriate wavelength. **ALLOW THE LAMP AT LEAST 10 MINUTES TO WARM UP BEFORE RESETTING ITS HOURS.**

W ENERGY - **ALLOW THE LAMP AT LEAST 10 MINUTES TO WARM UP BEFORE MEASURING ITS ENERGY. CLEAR BOTH SAMPLE AND REFERENCE BEAMS BEFORE MEASURING LAMP ENERGIES.**

D2 ENERGY - **ALLOW THE LAMP AT LEAST 10 MINUTES TO WARM UP BEFORE MEASURING ITS ENERGY. CLEAR BOTH SAMPLE AND REFERENCE BEAMS BEFORE MEASURING LAMP ENERGIES.**

SWITCH D2 -Turns the Deuterium lamp ON or OFF depending on its current status.


Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

LAMPS Page - Delta

- The lamp functions are available directly from the LAMPS key on the SETUP and MENU pages.
- This page shows the status of the Tungsten Halogen lamp, whether ON, OFF or FAILED, and the energy level. It also allows the lamp hours to be reset.
- The HOURS parameter states the number of hours that the lamp has been in use.

The Tungsten Halogen lamp should be replaced after 2000 hours.

Whenever a lamp is changed then the hours parameter should be reset to zero.

RESET W HRS - Resets hours and measures energy at the appropriate wavelength. **ALLOW THE LAMP AT LEAST 10 MINUTES TO WARM UP BEFORE RESETTING ITS HOURS.**

W ENERGY - **ALLOW THE LAMP AT LEAST 10 MINUTES TO WARM UP BEFORE MEASURING ITS ENERGY. CLEAR BOTH SAMPLE AND REFERENCE BEAMS BEFORE MEASURING LAMP ENERGIES.**

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com


CELL PROGRAMMER

- The 7 Cell Programmer is fitted as standard and enables up to seven samples to be presented for measurement sequentially.

- This section describes the operation of the accessory with the software and also describes its removal / refit.

When the Cell Programmer is installed then a status box appears on the top right hand side of the screen, to the left of the Absorbance/Wavelength box, on all the pages (except HOME) indicating the the presence of the accessory and its setting.

If this status box is displayed, then the Cell Programmer can be moved manually using the < > arrow keys. To step foward/backward many positions, repeated pressing of the key will be 'remembered' by the system.

- This page allows the Cell Programmer to be set up for the required analysis. To reach this page highlight the CELL PROGRAMMER option from the FIXED, QUANT, RATIO, SCAN, or RATE method pages, and press ENTER.

- Alternatively, press ACCESSORIES from the MENU page then select CELL PROG and press ENTER.

CELL POS. : Used to change the current cell position, using <> arrow keys. This change is reflected in the Status box.

MODE : Selects from MANUAL / RUN&STEP / AUTO / OFF.

MANUAL Cell position is only changed by use of the <> key on the keyboard.

RUN&STEP Allows measurement on the current cell then automatically moves to the next cell ready for the next measurement.

AUTO Performs a measurement on each of the cells in turn.

OFF Turns off the Cell Programmer. The instrument now behaves as if it has a single cell only.

The number of cells automatically measured in sequence (maximum of 7) will be dependent on the value of the LAST CELL parameter

REF. MODE : Toggles between ON and OFF with Enter.

- When ON, in all the above modes CELL 1 is assigned as the REFERENCE position and a zero or baseline will be performed.

LAST CELL : Sets the number of cells to be used in the range 1 to 7.

CELL CYCLES : Sets the number of cycles (up to 300). For example, if set to 4 then each cell will be measured in turn according to the application method being used four times. MODE must be set to AUTO.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

SPEED : Brings up a pop-up list to enable HIGH, MEDIUM or LOW rotation speed to be selected.

FUNCTION KEYS

INITIALISE Resets the accessory and places cell 1 in the sample beam.

CELL PROGRAMMER - REMOVAL AND REFIT

This procedure is essential for the fitting of any of the alternative cell holders available for long pathlength cells, single cell thermostatting, or the MiniSipper.

If the Cell Programmer basket is removed / replaced at any time, the 'Refit' procedure must always be repeated.


REMOVAL:

- Holding the basket firmly with one hand, unscrew the central screw anti-clockwise until the basket is released.
- **REPLACE THE COVER ON THE OPTOSENSOR**, and secure the cell holder.

REFIT:

- Remove the cover from the optosensor, and remove the cell holder.
- Indentify the position of the 'keyway' on the motor shaft.
- View the underside of the basket to locate the position of the moulded 'key'.
- Re-locate the basket on the shaft ensuring a positive location of the key in the keyway, and tightened the central screw clockwise.
- **From the CELL PROG page use the INITIALISE function key to correctly align the basket with the spectrophotometer.**

This last action is important to ensure the correct operation of the system. Whilst there is significant resistance to manual movement of the motor, it is suggested that 'as routine' initialisation of the Cell Programmer should be performed as a check before any measurements are performed.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

SUPERSIPPER

- The SuperSipper is an optional accessory that enables samples to be drawn into a flowcell of the user's choice for automatic measurement. After the measurement is complete the sample may be sent to waste or returned to its original vessel. A continuous pumping mode allows the system to be washed through when required, e.g. between applications.

- This section describes the operation of the SuperSipper with the Local Control software. Full details of the installation and operation of the accessory are described in the SuperSipper Installation and Maintenance Manual (9499 230 29611) supplied with the SuperSipper.

- To operate the SuperSipper, install as described in the above manual. Present the sample to the SuperSipper and press the switch. The required sample volume will be drawn into the tubing. When the system beeps remove the sample and the required air gap will be drawn in. Once the measurement has taken place press the switch again. The sample will be pumped out of the flowcell either to waste or returned to the sample vessel.

When the sipper is connected a status box is displayed on the right hand side of all the method screens that indicates the presence of the SuperSipper and its status.

SIPPER Page

```
        * SIPPER *

SIPPER              : OFF
MODE                : SIP
AIR GAP             : 50 cm
SAMPLE VOL.         : 1.000 ml
SAMPLE              : WASTE
LOW VOL.            : OFF
```

- This page allows the SuperSipper to be set up for the required analysis. The method set on this page will be saved with any data produced by the software.

- To reach this page press ACCESSORIES from the MENU page then select SIPPER and press ENTER.

To change a parameter highlight the required value using the Up/Down Arrow keys and press ENTER.

SIPPER : Selects from ON, OFF or STANDBY.

In STANDBY mode, the sipper pumps a small volume approximately every 30 minutes. This is to change the point at which pressure is applied to the Sipper tubing, thus preventing the formation of a permanent kink. The first 6 movements are in the Return direction, and the next 6 are in the Waste direction. The total volume pumped is sufficiently small as to ensure that any sample present remains in the tubing. When the sipper is operated normally, its clock is reset and the Standby process re-starts.

MODE : Selects from SIP / SIP&RUN / CONTINUOUS / AUTOSAM.

SIP Sets the system to fill the flowcell. If SAMPLE is set to RETURN

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

then alternate switch presses will fill and empty the flowcell. In this mode instrument operation is completely independent of the SuperSipper

SIP&RUN Sets the system to fill the flowcell and automatically perform a measurement. If SAMPLE is set to RETURN then alternate switch presses will fill and empty the flowcell. The current method used to produce the result (e.g. if FIXED is current then the sample will be measured using the FIXED METHOD as set).

CONTINUOUS Sets the system to pump continuously to waste. Alternate switch presses will start and stop pumping. In this mode instrument operation is completely independent of the SuperSipper.

AUTOSAM Sets the system to work with the Gilson 221XL and 222XL Autosamplers. Refer to the Autosampler Interface manual for further details.

AIR GAP **:** Enter value between 0 and 500 cm.

Sets the gap between the trailing meniscus of the current sample and the leading meniscus of the next sample. The gap is measured to the nearest centimetre.

For best results set the airgap no less than 8cm from the flowcell.

SAMPLE VOL **:** Enter a value between 0.2 and 9.999 ml.

Sets the volume of sample to be pumped.

SAMPLE **:** Selects from WASTE or RETURN.

WASTE After measurement the sample is pumped through the flowcell To waste by the act of pumping the next sample.

RETURN After measurement the pump direction is reversed and the sample is returned to the sample vessel.

LOW VOL **:** Toggles ON or OFF.

Automatically adjusts the pumping time to maintain the correct air gap for narrow uptake tubing.

OFF Use standard internal diameter (1.1mm) uptake tube.
ON Use narrow internal diameter (0.8mm) uptake tube

FUNCTION KEYS

VIEW CALIB Goes to the SIPPER CALIBRATION page.

CALIBRATE Starts the Sipper calibration procedure.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Sipper Calibration

- This calibrates the SuperSipper to take account of variations in pump and uptake tubing and sample viscosities. A volume is set and using the appropriate solvent and tubing several sips are performed. The actual volume sipped is entered and a calculation done to produce a calibration factor. This factor is then used to adjust the pumping time to ensure that the correct volume of sample is always used.

- Details of the calibration used are displayed on the SIPPER CALIBRATION page.

CALIBRATE SIPPER Page

```
*                    SIPPER *
CALIBRATE

NOMINAL VOL.         : 1.000 ml
NO. SIPS DONE        : 5
```

- Using the solvent and tubing which will be used for the sample solutions offer a measuring cylinder filled to the highest gradation to the sipper uptake tube and press the switch plate.

The sipper will pump a sample and the spectrophotometer will issue a beep. Withdraw the measuring cylinder and the sipper will pump the air gap. The values used for sample volume and air gap are those set on the SIPPER page.

- Repeat this process for a number of cycles up to a maximum of 10 then press ENTER.

- Measure the total volume taken from the measuring cylinder and enter this. The calibration will be displayed.

FUNCTION KEY

SIPPER PAGE	Returns to the SIPPER page and abandons the calibration.

SIPPER CALIBRATION Page

```
*                    CALIBRATION *
SIPPER               16:47
        25/10/96

NOMINAL VOL.         : 1.000 ml
NO. SIPS DONE        : 5
TOTAL VOL SIPPED     : 5.100 ml
TUBING CAL           : 1.020
```

- This page displays the current sipper calibration.

- To alter the calibration, press CALIBRATE

FUNCTIONS KEYS

SIPPER PAGE	Returns to the SIPPER page.
CALIBRATE	Starts the calibration procedure.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

MINISIPPER

- The MiniSipper is an optional accessory that enables samples to be drawn into a flowcell of the user's choice for automatic measurement. After the measurement is complete the sample is sent to Waste. A continuous pumping mode allows the system to be washed through when required, e.g. between applications.

- This section describes the operation of the MiniSipper with the Local Control software. Full details of the installation and operation of the accessory are described in the MiniSipper Installation and Maintenance Manual (9499 230 45111) supplied with the MiniSipper.

- To operate the MiniSipper install as described in the above manual. Present the sample to the MiniSipper and press the switch. The required sample volume will be drawn into the tubing. When the system beeps remove the sample and the required air gap will be drawn in. Once the measurement has taken place press the switch again. The sample will be pumped out of the flowcell to waste.

When the sipper is connected a status box is displayed on the right hand side of all the method screens that indicates the presence of the MiniSipper and its status.

SIPPER Page

```
        * SIPPER *

SIPPER              : OFF
MODE                : SIP
AIR GAP             : 50 cm
SAMPLE VOL.         : 1.000 ml
SAMPLE              : WASTE
LOW VOL.            : OFF
```

- This page allows the MiniSipper to be set up for the required analysis. The method set on this page will be saved with any data produced by the software.

- To reach this page press ACCESSORIES from the MENU page then select SIPPER and press ENTER.

To change a parameter highlight the required value using the Up/Down Arrow keys and press ENTER.

SIPPER : Toggle ON or OFF.
Switches the MiniSipper ON or OFF.

MODE : Selects from SIP / SIP&RUN / CONTINUOUS.

SIP Sets the system to fill the flowcell.

SIP&RUN Sets the system to fill the flowcell and automatically perform a measurement. The current method is used to produce the result (e.g. if FIXED is current then the sample will be measured using the FIXED METHOD as set).

CONTINUOUS Sets the system to pump continuously to waste. Alternate switch presses will start and stop pumping. The instrument will not process any key presses while the MiniSipper is pumping in continuous mode.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

AIR GAP : Enter value between 0 and 500 cm.

Sets the gap between the trailing meniscus of the current sample and the leading meniscus of the next sample. The gap is measured to the nearest centimetre.

For best results set the airgap no less than 8cm from the flowcell.

SAMPLE VOL : Enter a value between 0.5 and 9.999 ml.

Sets the volume of sample to be pumped.

SAMPLE : Set to WASTE.

After measurement the sample is pumped through the flowcell to waste by the act of pumping the next sample.

LOW VOL : Set to OFF.

A standard internal diameter (4.0mm) uptake tube is used.

FUNCTION KEYS

VIEW CALIB Goes to the SIPPER CALIBRATION page.

CALIBRATE Starts the Sipper calibration procedure.

Sipper Calibration

- This calibrates the MiniSipper to take account of variations in pump and uptake tubing and sample viscosities. A volume is set and several sips are performed using the appropriate solvent. The actual volume sipped is entered and a calculation done to produce a calibration factor. This factor is then used to adjust the pumping time to ensure that the correct volume of sample is always used.
- Details of the calibration used are displayed on the SIPPER CALIBRATION page.

CALIBRATE SIPPER Page

```
      * CALIBRATE   SIPPER *

NOMINAL VOL.        : 1.000 ml
NO. SIPS DONE       : 5
```

- Using the solvent which will be used for the sample solutions, offer a measuring cylinder filled to the highest gradation to the sipper uptake tube and press the switch.

The sipper will pump a sample and the spectrophotometer will issue a beep. Withdraw the measuring cylinder and the sipper will pump the air gap. The values used for sample volume and air gap are those set on the SIPPER page.


Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

- Repeat this process for a number of cycles up to a maximum of 10 then press ENTER.
- Measure the total volume taken from the measuring cylinder and enter this. The calibration will be displayed.

FUNCTION KEY

SIPPER PAGE	Returns to the SIPPER page and abandons the calibration.

SIPPER CALIBRATION Page


- This page displays the current sipper calibration.
- To alter the calibration press CALIBRATE

FUNCTIONS KEYS

SIPPER PAGE	Returns to the SIPPER page.
CALIBRATE	Starts the calibration procedure.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

INTERNAL PRINTER

- The Internal Printer is a factory-fitted thermal-head printer. It is supplied with a roll of 11.2cm wide single-ply thermal paper. Paper rolls may be ordered separately under the following part number: 4401 161 00391
- For long-term storage, it is recommended that the paper is kept away from light and at room temperature.
- A red warning stripe is printed on the paper within 1 metre of the end of the roll. It is recommended that the roll is replaced as soon as the warning stripe is visible to avoid possible paper jams.
- The printer is fitted with a Line Feed Button. Press this once to switch the printer off-line and feed the paper though the printer. Press the button a second time to put the printer back on-line. Do not press the Line Feed Button during printing as this will cause a internal error condition which can only be rectified by restarting the instrument.

To replace the paper roll:

- Cut the paper at the end of the new roll in a “V” to make a point at the end of the paper.
- Carefully feed the printer paper (from the bottom of the roll, i.e. shiny side down) into the back of the printer using the Line Feed Button.
- Keep the Line Feed Button depressed until the paper emerges “squarely” through the top slot of the printer housing.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

MAINTENANCE

- The information given in this section deals only with those parts of maintenance or service which can be safely carried out by the user. Work other than that detailed should be carried out by a service engineer.

- **ALWAYS ENSURE THAT THE SAMPLE BEAM IS CLEAR BEFORE SWITCHING ON THE INSTRUMENT. Failure to do so will produce abnormal results.**

- Low lamp energy values can be caused by leaving cells in the sample beam during energy measurement. **ALWAYS CHECK THAT THE SAMPLE BEAM IS CLEAR BEFORE MEASURING LAMP ENERGIES.**

- Abnormal results will be produced if a sample is left in the beam when ZERO/BASE is pressed. **ALWAYS ENSURE THAT THE SAMPLE IS REMOVED AND THAT THE SAMPLE BEAM IS CLEAR OR CONTAINS THE APPROPRIATE ZERO REFERENCE BEFORE ZEROING THE INSTRUMENT OR PERFORMING A BASELINE SCAN.**

- **VERY OFTEN POOR INSTRUMENT PERFORMANCE OR FAILURE CAN BE ATTRIBUTED TO SIMPLE FAILURE OF THE TUNGSTEN LAMP - THEREFORE REPLACE (AS BELOW) USING THE SPARE LAMP PROVIDED BEFORE SEEKING FURTHER ASSISTANCE.**

- If any fault occurs (including the above lamp failure), these are reported by the system as an 'Error condition', and an 'Exxxx' number is generated. Descriptive text is also included with this message.

 DETAILED BELOW ARE THE ERROR CODES PRODUCED IF:

 1. The tungsten lamp fails or is poorly aligned
 E3010 E3011 E3104 E3015 E3030

 2. The deuterium lamp fails (Gamma only)
 E3003 E3004 E3005 E3006 E3007 E3008 E3009 E3012 E3013
 E3022 E3029 E3044 E3045

 3. The beam is blocked on initialisation
 E3027 E3056

 4. The cell programmer is stalled in use (or fails to initialise)
 E3001 E3002 E3054 E3055 E3082 E3084

 5. The sample compartment is open
 E3053 E3062 E3068 E3069 E3071

- A comprehensive list of these codes is available in the Service Manual for this product. Generation of a code not related to replacement of either the tungsten or deuterium lamps usually requires you to contact your local Thermo Spectronic approved Customer Support Organisation.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

ROUTINE MAINTENANCE

Very little maintenance is required to keep the spectrophotometer in good working condition. The interior should be kept as dust free as possible and the sample compartment cleaned regularly; wipe off spilt chemicals immediately.

Replacement sample compartment liners are available under the following part number:

Helios Single Beam : 9423 UV9 7000E

Cleaning Instrument Exterior

The exterior of the instrument can be cleaned periodically as follows:

CAUTION: Do not allow moisture to leak into the instrument.


- Switch off the spectrophotometer and disconnect from the mains supply.
- Using a lint free cloth dampened with a weak solution of detergent and water, wipe the exterior surface of the instrument as necessary.
- Wipe over with a cloth dampened with plain water.
- Dry the surface with another cloth.


Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

REMOVAL AND REPLACEMENT OF TUNGSTEN HALOGEN LAMP

WARNING: Switch off and disconnect the spectrophotometer from the mains and allow the lamp to cool for at least 15 minutes before proceeding.

- Remove the back corner cover by turning the fastener one quarter turn anti-clockwise and slide the cover up to remove.
- Now remove vertically upwards the metal lamp cover.

- If fitted, remove the spring clip
- Hold the W lamp and pull upwards to remove.

When fitting the new tungsten halogen lamp, avoid handling the silica envelope. Finger marks become burnt on and cannot be removed after the lamp is switched on. As this can affect the output characteristics, handle only the base of the lamp. If the silica envelope does become contaminated, clean with a powerful degreasing solvent such as absolute alcohol before the lamp is switched on.

- Use the new lamp's protective sleeve, a polythene bag or a piece of tissue paper wrapped around the lamp and insert the pins into the socket.
- Replace the spring clip.

These lamps are manufactured to very high tolerances, but to ensure optimum energy throughput, align the lamp filament exactly as shown in the diagram, with the white line on the lamp base facing towards the front of the instrument.

- Refit both the metal and rear covers.
- Reconnect the spectrophotometer to the mains supply and switch on.
- Lamp hours and energy must be reset from the controlling software.


Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

REMOVAL AND REPLACEMENT OF DEUTERIUM LAMP (Gamma only)

WARNING: **(1) Switch off and disconnect the spectrophotometer from the mains supply and allow the lamp to cool for at least 15 minutes before proceeding.**

(2) UV radiation from a Deuterium lamp can be harmful to the skin and eyes. Always view the lamp through protective glasses that will absorb UV radiation. Avoid looking directly at the Deuterium arc. Do not expose the skin to direct or reflected UV radiation.

- Set the power switch to off and disconnect the spectrophotometer from the mains supply.
- Remove the back corner cover by turning the fastener one quarter turn anti-clockwise and sliding the cover up to remove.

- Make sure the lamp has cooled, disconnect the lamp at the in-line connector, Using the key provided loosen the three locating screws, rotate lamp assembly anti-clockwise and lift lamp out.

When fitting the Deuterium lamp avoid handling the silica envelope. Finger marks become burnt on and cannot be removed after the lamp is switched on. This can affect the light output characteristics. Handle only the base of the lamp or the mounting plate. If the silica envelope becomes contaminated, clean with a powerful degreasing solvent such as absolute alcohol before the lamp is switched on.

- Take the new lamp, identify the notch in the mounting plate. Locate the lamp such that the notch points towards the lamp change mirror. Tighten locating screws down with the key provided.
- Re-connect the new lamp at the in-line connector.
- Refit the rear cover.
- Reconnect the spectrophotometer to the mains supply and switch on. Allow half an hour for warm up time.
- Lamp hours must be reset by the controlling software.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Fault Finding Guide

Problem	Symptom	Possible Cause
Instrument dead	The fans are not running	1. Ensure the mains lead is firmly pressed home. Some leads are a very tight fit in the three pin IEC connector on the rear of the instrument. 2. If a switched outlet is being used, ensure it is on. 3. Check the fuse in the plug. 4. Try another mains lead. 5. Try another mains outlet. 6. Ensure the power switch on the instrument has been fully operated. 7. If all of the above are OK the power supply may have failed. Contact your local agent.
Display Blank	No text on display, a slight glow can be seen around the edges of the display in subdued lighting.	1. The contrast control is incorrectly set. Go to the MENU page, the right and left arrow keys adjust the contrast. If nothing happens it may be that error messages are being displayed. Press the CLEAR key 5 times at 10 second intervals, to clear the error list, and try again. The display contrast is affected by temperature and may need adjusting from day to day, or as the instrument warms up. 2. The instrument has been put into remote control without allowing it to initialize first. The display will return next time that the instrument is switched on. 3. An unsuccessful attempt has been made to update the software. Connect the instrument to a PC and repeat the software upgrade process.
Any performance problem, Failure to initialise or error 1053	No light can be seen from the ventilation slots in the lamphouse cover.	1. The tungsten lamp has failed, fit a new lamp. A small number of tungsten lamps do fail very rapidly. If this happens fit a new lamp and contact your agent for a free replacement. It is likely that the spare lamp shipped with the instrument will be from the same batch as the one in the instrument.
Lamp energy low Error 1053, Drift at visible wavelengths	Instrument noisy in the visible region.	1. Check the tungsten lamp is correctly fitted and in good condition. The envelope should not show any signs of blackening or opacity. If in doubt

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

		replace with a new one. A small number of tungsten lamps do rapidly degrade. If this happens fit a new lamp and contact your agent for a free replacement.
Lamp energy low Error 1053. Noisy signal.	Cell programmer partially blocking the beam.	1. If the carousel is fitted after using another cell holder, go to the Cell Programmer page and key "Initialise". 2. Make sure the carousel is not jamming or being fouled on tubing. 3. Make sure there are no cells in the beam.
Fails to initialise.	Error 3093	1. Check that the sample compartment lid is correctly shut. 2. Check that the accessory panel, found at the left hand side of the sample compartment, is in place. 3. Position the instrument so that it is not directly in strong sunlight. 4. If using a water circulator, ensure that black water tubing is used.
Deuterium lamp energy low.	Lamp energy is indicated as low but performance seems OK.	1. Lamp energy was measured with a cell in the beam. Remove all cells and re-measure the lamp energy.
Deuterium lamp energy low.	Performance is poor in the UV region.	1. Plastic cells that do not transmit in the UV are being used. 2. The deuterium lamp may need replacing.
Fails to initialise.	Error 3027	1. Check that the beams are clear and retry. 2. Record all the error messages. 3027 is a general failure message, it is usually preceded by more specific error codes. 3. Check that the tungsten lamp is working. 4. Check that the cell programmer is not blocking the beam. 5. If the problem persists connect the instrument to a PC running a terminal programme, see below for method. Switch the instrument on whilst holding down the RUN key. Debug data will be sent to the PC. Collect this and Email it to customer support for further help.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

CONNECTING HELIOS TO A PC

Terminal programmes.


You will need software that allows you to write to and read from the computer's serial port. The communications parameters are:

Baud rate 9600
Data bits 8
Stop bits 1
No parity
Flow control off.
Local Echo on

HyperTerminal


You can also use HyperTerminal that comes with Windows 95. To establish a connection, start HyperTerminal and select "New Connection" from the file menu. In the dialogue that appears, type in Helios for the name and choose an icon. Key OK. In the next dialogue box select "Direct to COM1" and key OK.

Fill in the COM 1 properties dialogue box as shown below.

Then from the FILE menu select PROPERTIES and then click on the SETTINGS tab. Fill in the dialogue box as shown below. Click on OK.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Once communication has been established you will need to save the debug data to a file. Do this using the TRANSFERS/CAPTURE TEXT as shown below.

Connecting Helios to a PC

1) Use a null modem cable (part number 4013 172 82111) to connect the RS232C port of Helios to a free COM port on the PC.

 Here is the pin out data for the RS232 Cable

Pin		Pin
9	to	9
2		3
3		2
4		8
8		4
6		1 and 7
1 and 7		6
5		5
Earth		Earth

 Switch on the PC and set up a terminal programme as described above.


Collecting Data

Once you have established control of the instrument and switched on the debug software you will need to collect the data that is returned and save it to file. This can be done by several methods, the next section details one of them.


Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Using Hyperterminal to collect data


1) From the menu bar select TRANSFERS and then CAPTURE TEXT.

2) The CAPTURE TEXT dialogue box appears. Type in a path and file name then key START.

3) Once the instrument has stopped sending data, from the menu bar Select TRANSFER then CAPTURE TEXT then key STOP to terminate the transfer and save the file.

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Heλios Gamma & Delta

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9499 230 45811
Édition N°1a
001201

Thermo Spectronic

A Thermo Electron business

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

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Artisan Technology Group is an independent supplier of quality pre-owned equipment

Gold-standard solutions

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

Learn more!

Visit us at **artisantg.com** for more info on price quotes, drivers, technical specifications, manuals, and documentation.

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

We're here to make your life easier. How can we help you today?

(217) 352-9330 | sales@artisantg.com | artisantg.com

ARTISAN® TECHNOLOGY GROUP