ARTISAN®
TECHNOLOGY GROUP

Your **definitive** source for quality pre-owned equipment.

Artisan Technology Group
(217) 352-9330 | sales@artisantg.com | artisantg.com

Full-service, independent repair center
with experienced engineers and technicians on staff.

We buy your excess, underutilized, and idle equipment
along with credit for buybacks and trade-ins.

Custom engineering
so your equipment works exactly as you specify.


- Critical and expedited services
- In stock / Ready-to-ship
- Leasing / Rentals / Demos
- ITAR-certified secure asset solutions

Expert team | Trust guarantee | 100% satisfaction

All trademarks, brand names, and brands appearing herein are the property of their respective owners.

Find the *Wheaton / Socorex Acura* at our website: ***Click HERE***

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DESCRIPTION

Adjustable volume multichannel micropipette for use with detachable polypropylene tips, for an easy and continuous work with standard 96-well microplates, 9 mm spacing.

The sample is drawn into the polypropylene tip without any contact with the plunger assembly. The integrated tip ejector is located on the handpiece allowing the automatical release of the tips.

Simultaneous delivery of 8 or 12 samples from 5 µl to 200 µl by means to two pipette sizes.

Each instrument is individually calibrated and comes along with its control certificate.

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HOLDING THE MULTICHANNEL MICROPIPETTE

Recommended way of holding the pipette is shown on fig. 1: place fingerrest on the phalanx of your forefinger with the ejector knob turned towards yourself.

The multichannel casing revolves on the handpiece and allows the selection of the best working position.

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VOLUME SETTING (fig. 2)

Set desired volume by rotating graduated knob No 31. Graduation mark must be even with vernier ring, vertical line on zero position. By moving vertical line clockwise along vernier graduation, volume will increase by 1 µl per division.

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Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com


Fig. 1 Operation


Fig. 2 Setting of the volume


Fig. 3 Calibration

TIPS AND FILTERS

For consistent results and proper operation of the ejector, use only our original polypropylene tips, or their equivalent.

These are supplied on 96 pieces racks, with the 9 mm spacing of the barrels of the multichannel micropipette.

Place the instrument on a row of tips, press firmly to ensure a solid fit of the tips on each channel to avoid leakage.

High-sitting tips (collar longer than regular tips) can be used after easy replacement of ejector head No 3 with alternative HST-Ejector head available from your distributor.

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PIPETTING OPERATION (fig. 1)

Phase 1 Depress plunger knob to the intermediate stop ending calibrated stroke before immersing tips into one of the various reservoirs.

Phase 2 Immerse tips approximately 5 mm and *smoothly* release plunger knob, thus drawing up simultaneously liquid into the 8 or 12 tips.
Withdraw tips from sample solution, sliding them along wall of the reservoir.

Phase 3 Place tips on the wells of the microplate and slowly depress plunger to first stop, thus emptying tips; after approximately one second, depress knob to lower discharge position: the excess air volume produced by depressing the plunger beyond the calibrated stroke blows out possible fluid particles which may cling to the tips.
Remove the multichannel pipette from the wells of the microplate, with plunger knob in its lowest position.
Release knob slowly.
To eject used tips, press ejector knob E.
Apply new tips for every liquid; change tips if droplets are retained after blow out. For better accuracy, fill new tip once before sampling.

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Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

MAINTENANCE

To check tightness of the multichannel pipette, draw liquid into the tips. Droplet forming at the extremity of one of the tips means this channel is leaking. O-ring of the corresponding barrel has to be changed.

Fluid getting accidentally into the multichannel pipette may clog one of the barrels or cause the plungers to drag.

In order to clean the plungers or to change the O-rings, remove the two screws No 4.

Cleaning the plungers: depress the plunger knob No 31: the plungers will protrude from the casing. Clean them, grease slightly and wipe with a clean cloth.

Changing the O-rings: to have access to the barrels No 6, pull out the two ejector shafts No 12 off the ejector head No 3 by pressing together the fingers of the ejector shaft: the ejector head will come out.
Remove the two screws No 5 together with spacers No 18. Take out the leaking barrel and change the two O-rings No 9 and 2 on the pipette 851-200 and the O-ring No 2 on the pipette 851-050.

Depress the plunger knob No 31 when tightening the screws No 4 and 5.

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CALIBRATION

Each Acura© has been calibrated permanently at the factory and individually controlled. The Q.C. certificate accompayning the instrument includes control data, serial number and operator identification. Calibration is carried out at a constant temperature (+/– 0,5° C) comprised between 20 and 25° C with Aqua bidest, according to DIN 12650, part 6.

Under normal conditions, a pipette does not need to be recalibrated. After replacing a part, which could influence the measure, the calibration should be checked, preferably by gravimetry according to the above mentioned norm; an analytical balance is needed. Ten measurements should be performed after checking the seal tightness of the instrument.

In-lab calibration (fig. 3): Whenever required, for instance when pipetting viscous or non-aqueous liquids, calibration adjustments can be performed by laboratory personnel. Reach calibration screw located inside plunger button by removing comfort pad a). Unlock side allen screw b), then delicately turn calibration screw c) clockwise to decrease volume, and vice versa. Re-tighten locking screw b). Check resulting volumes after each fine tuning. Cover push-button with comfort pad again.

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WARRANTY

Your multichannel micropipette is guaranteed for a period of two years against any factory defect. Return it to the dealer from whom you purchased it if the instrument presents a defect which you cannot eliminate by applying the ordinary maintenance procedure as described. The identification is facilitated by its serial number.

Decontaminate instrument prior to returning it.

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11/3/00

Parts List

artistically crafted by linda lee steve smith

Wheaton Science Products Acura Multi-Channel Micropipette 851 (A-Series)						
ITEM #	DESCRIPTION	5-50uL	5-50uL	50-200uL	50-200uL	NOTES
		8 POS.	12 POS.	8 POS.	12 POS.	
	COMPLETE UNIT	851240	851241	851242	851243	
1	PLUNGER ASSEMBLY	I053035	I053037	I053036	I053038	
2	O-RING	I051249	I051249	I051646	I051646	
3	EJECTOR HEAD	I054126	I054127	I054126	I054127	
4	SCREW	I053039	I053039	I053039	I053039	
5&18	SCREW & SPACER	I053040	I053040	I053040	I053040	
6	BARREL	851244	851244	851245	851245	
8	BOTTOM PLATE	I053770	I053771	I053770	I053771	
9	O-RING	N/A	N/A	I051248	I051248	
10	BARREL HOLDER	I053941	I053943	I053942	I053944	
11	SPRING	I053043	I053043	I053043	I053043	
12	EJECTOR SHAFT	I054148	I054148	I054148	I054148	
13	SOCKET	I051559	I051559	I051559	I051559	
14	SCREW	I051560	I051560	I051560	I051560	
15	EJECTOR PUSH BUTTON	I051561	I051561	I051561	I051561	
16	EJECTOR SLEEVE	I053047	I053047	I053047	I053047	
17	HANDLE	I053048	I053048	I053049	I053049	
21	CIRCLIP	I052398	I052398	I052398	I052398	
22	WASHER	I053052	I053052	I053052	I053052	
23	LOWER STOP	I053053	I053053	I053053	I053053	
24	OVERSHOOT SPRING	I053054	I053054	I053054	I053054	
25	SPRING	I053055	I053055	I053055	I053056	
26	OVERSHOOT STOP	I053057	I053057	I053057	I053057	
27	SLEEVE	I053058	I053058	I053058	I053058	
28	CASING	I053059	I053060	I053059	I053060	
29	LOCKING RING	I059190	I059190	I059190	I059190	
30	YELLOW DISK - ALUM.	I050471	I050471	I050471	I050471	S/N <0841000
30	YELLOW DISK - RUBBER	I054521	I054521	I054521	I054521	S/N >8441000
31	MICROMETER	I052580	I052580	I052581	I052581	S/N >0141000
31	MICROMETER ASSY.					
	INCL. # 30,33,34,35,40	I058506	I058506	I053063	I053063	S/N <0141000
31	MICROMETER ASSY.					
	INCL. # 30, 35, 40	I054524	I054524	I054525	I054525	S/N >0841001
32	UPPER STOP	I053064	I053064	I053064	I053064	
33	MICROMETER SCREW	I051630	I051630	I051631	I063061	S/N >0141000
34	SCREW	I053065	I053065	I053065	I053065	
35	CALIBRATION SCREW	I053066	I053066	I053066	I053066	
36	RATCHET	I053067	I053067	I053067	I053067	
37	SPRING	I052588	I052588	I052588	I052588	
38	CIRCLIP	I050475	I050475	I050475	I050475	
39	CAP	I054961	I054961	I054961	I054961	
40	LOCKING SCREW	I053071	I053071	I053071	I053071	
56	VERNIER	I053072	I053072	I053073	I053073	

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
11/3/00

Parts List

artistically crafted by linda lee steve smith


851 Acura Multi-Channel Micropipette

11/3/00

Parts List

artistically crafted by linda lee steve smith

Artisan Technology Group is an independent supplier of quality pre-owned equipment

Gold-standard solutions

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

Learn more!

Visit us at **artisantg.com** for more info on price quotes, drivers, technical specifications, manuals, and documentation.

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

We're here to make your life easier. How can we help you today?

(217) 352-9330 | sales@artisantg.com | artisantg.com

ARTISAN® TECHNOLOGY GROUP