Berger SFC Thermal Control Module



Limited Availability
Used and in Excellent Condition

Open Web Page

https://www.artisantg.com/51452-1

All trademarks, brandnames, and brands appearing herein are the property of their respective owners.

- Critical and expedited services
- In stock / Ready-to-ship

- · We buy your excess, underutilized, and idle equipment
- · Full-service, independent repair center

ARTISAN'
TECHNOLOGY GROUP

Your **definitive** source for quality pre-owned equipment.

Artisan Technology Group

(217) 352-9330 | sales@artisantg.com | artisantg.com

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

AutoChem

Optimize your SFC for

accurate temperature control

■ BERGER SFC Analytix

Thermal Control Module



Multiple valving options

The Thermal Control Module's stirred air bath design assures accurate temperature control with your analytical SFC. The controlled temperature ranges from ambient +7° C to 150° C. Cryogenic cooling using liquid CO₂ allows subambient operation extending the temperature range down to -50° C.

The oven provides an accurately controlled stable column environment. For column switching, solvent switching, or back flush applications the oven can be customized with additional automated valves.

The cabinet with its spacious column racks allows you to fit up to 6 analytical columns in sizes up to 30 cm.

Optionally up to two analog input boards can be mounted for interfacing with external detectors.

METTLER TOLEDO

Specifications:	
Operating range	7° C above ambient to 150° C ±1% (to -50° C with cryogenic cooling)
Internal dimension	11 cm deep, 15 cm wide, 36 cm high
Options	Manual injection valve (optional-10 µl external sample loop injection valve) Analog Input Board (for external signal input) Column Switching Valve (up to 6 columns) Solvent Swtiching Valve (up to 6 solvents)
External dimensions	58 cm deep, 20 cm wide, 58 cm high, 21 kg weight
Frequency	50 or 60 Hz ±5%

Mettler-Toledo AutoChem, Inc. 8223 Cloverleaf Drive, Suite 120 Millersville, MD 21108, USA

Toll Free 866 333 MTAC (866 333 6822)

Fax 410 987 2626

Internet http://www.mtautochem.com

METTLER TOLEDO

© Copyright 2004. 0401_TCM-ds/1m

Artisan Technology Group is an independent supplier of quality pre-owned equipment

Gold-standard solutions

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

Learn more!

Visit us at artisantg.com for more info on price quotes, drivers, technical specifications, manuals, and documentation.

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

We're here to make your life easier. How can we help you today? (217) 352-9330 | sales@artisantg.com | artisantg.com

