

\$1295.00

In Stock **Qtv Available: 3 Used and in Excellent Condition** 

**Open Web Page** 

https://www.artisantg.com/53365-1

All trademarks, brandnames, and brands appearing herein are the property of their respective owners.

- Critical and expedited services
- In stock / Ready-to-ship

- · We buy your excess, underutilized, and idle equipment
- · Full-service, independent repair center



Your definitive source for quality pre-owned equipment.

**Artisan Technology Group** 

(217) 352-9330 | sales@artisantg.com | artisantg.com

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

### **ACCESSORIES: Steri-Cult 200 and Forced Draft Incubators**

		Model Numbers					
Stock #	Description	3033	3035	3029	3039	3860	3862
190087	Caster Dolly	•	•				
190236	Caster Dolly					•	•
224136	Stainless Steel Shelf Kit*	•	•				
224138	Stainless Steel Shelf Kit*			•	•		
224144	Copper Shelf					•	•
224145	Electropolished Stainless Steel Shelf Kit*					•	•
224159	Electropolished Stainless Steel Shelf Kit*	•	•				
190035	Duplex Outlet Kit, 120 V (FI** interior rear wall)	•				•	
190036	Duplex Outlet Kit, 230 V (FI** interior rear wall)		•				•
505099	Duplex Outlet Kit,120 V (FI** interior rear wall)			•			
505094	Duplex Outlet Kit, 230 V (FI** interior rear wall)				•		
513028	Stacking Adapter (For Two Units)			•	•		
523028	Floor Stand (to Stack Units)	•	•	•	•	•	•
760045	Disposable HEPA Filter (99.97%) Assembly	•	•			•	•
770001	Disposable Microbiological Gas Filter, 99.97%	•	•	•	•	•	•

<sup>\*</sup> Channels Included

#### **GAS GUARD**

Forma Scientific's Model 3050 Gas Guard takes the guesswork out of gas usage. It monitors any non-hazardous gas supply such as  $\mathrm{CO}_2$  or  $\mathrm{N}_2$ , sounds an alarm, and automatically switches from one cylinder to the other when the supply is exhausted. The alarm reset switch will silence the alarm. Dual voltage provided with selector switch.

Complete with remote alarm contacts.

Forma offers an optional built-in Gas Guard system on water jacketed incubators.

3050

Exterior Dimensions	7.13"W X 5.13"H X 9.25"F-B (18.1cm x 13.1cm x 23.5cm)			
Gas Connections	Inlet - two 1/4" barb fittings; Outlet - one 1/4" barb fitting			
Inlet Pressure	15 PSIG (1.0 bar)			
Construction	Heavy-gauge metal powder coated for a durable easily			
	maintained surface.			
Electrical	90-125V/180-250V, 1 PH, 50/60 HZ, 0.3/0.15 FLA			
Shipping Weight	13 lbs. (5.9 Kg)			

<sup>\*\*</sup>Factory Installed

### **GENERAL ACCESSORIES**

	Stock/Model#	Description
GAS ANALYZER KITS	155021	CO <sub>2</sub> Fyrite Kit 0-20%
	155022	Duplex Fyrite Kit CO <sub>2</sub> , 0-20%; O <sub>2</sub> , 0-21%
	155023	O <sub>2</sub> Fyrite Kit 0-60%
	155024	O <sub>2</sub> Fyrite Kit 0-21%
FYRITE FLUID	220051	Replacement Fluid for CO <sub>2</sub> Fyrite Kit, Catalog #'s 155021 and
		155022; three bottles per carton
	155009	O <sub>2</sub> Fyrite Fluid 21% and 60% for catalog #'s 155022, 155023 &
		155024; three bottles per carton
HYGROMETER &		
THERMOMETERS	155010	Hygrometer 4" Dial, 0-100% RH
	285722	Glass Thermometer, 0°C to +100°C
	853227	LCD Pocket Digital Thermometer -90°C to + 60°C
	292001	Accessory Probe for use with 853227
	940246	Glass Thermometer -40°C to +30°C
REGULATORS	961027	Two-Stage Regulator, $N_2$
	965010	Two-Stage Regulator, CO <sub>2</sub>
DISINFECTANTS	170001	O-SYL, six (6) 1 gallon bottles
	170003	Roccal II, four (4) 1 gallon bottles
MONITORS & ALARMS	1535	Programmable Monitor/Alarm System, 90-125 VAC 50/60 HZ,
		1PH, 0.7 FLA or 180-250 VAC, 50/60 HZ, 1PH, 0.35 FLA
	Modules	Maximum 26 VDC, fused in Model 1535
ACCESSORIES FOR	190456	Cryogenic Converter Module. RTD temperature probe to
1535 MONITOR/		RS-485. Includes an RTD probe with 6 ft. lead and one (1)
ALARM SYSTEM		set of dry contacts (temperature range -200°C to -90°C).
	190457	Converter Module. RTD temperature probe to RS-485.
		Includes an RTD probe with 6 ft. lead and one (1) set of dry
		contacts (temperature range -90°C to +70°C).
	190458	Converter Module. RS-232 to RS-485. Monitors one (1)
	100.170	RS-232 port and one (1) set of dry contacts.
	190459	Converter Module. Analog to RS-485. Monitors four (4)
	1000.17	analog signals (0-1V) and two (2) sets of dry contacts.
INCUBATOR	400047	Sensaphone Telephone Dialing System, 110 VAC, 60HZ.
CONTROLLERS	3050	Gas Guard, 90-125V/180-250V, 50/60 HZ
CONTROLLERS	3057	CH/P Portamatic <sup>™</sup> , Portable CO <sub>2</sub> Controller, 120V, 60HZ CH/P Portamatic <sup>™</sup> , Portable CO <sub>2</sub> Controller, 230V, 50HZ
CELL ROLL SYSTEM	3063 4862	15 Position Cell Roller 120V 60HZ
OLLE NOLL STOTEM	4868	15 Position Cell Roller 230V 50HZ
	190049	5 Position Add-On Tier for Cell Roller
CELLROLL		110mm X 285mm bottle (4 per case)
PRODUCTION APPARATUS	475560 475580	110mm X 285mm bottle (4 per case)  110mm X 480mm bottle (4 per case)
	475580	
OTHER ACCESSORIES	190043	Modular Incubator Chamber  Extra Tray for Modular Incubator Chamber
	190238	Extra Tray for Modular Incubator Chamber

# Artisan Technology Group is an independent supplier of quality pre-owned equipment

### **Gold-standard solutions**

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

# We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

#### Learn more!

Visit us at artisantg.com for more info on price quotes, drivers, technical specifications, manuals, and documentation.

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

We're here to make your life easier. How can we help you today? (217) 352-9330 | sales@artisantg.com | artisantg.com

