



Artisan Technology Group is your source for quality new and certified-used/pre-owned equipment

- FAST SHIPPING AND DELIVERY
- TENS OF THOUSANDS OF IN-STOCK ITEMS
- EQUIPMENT DEMOS
- HUNDREDS OF MANUFACTURERS SUPPORTED
- LEASING/MONTHLY RENTALS
- ITAR CERTIFIED SECURE ASSET SOLUTIONS

SERVICE CENTER REPAIRS

Experienced engineers and technicians on staff at our full-service, in-house repair center

*InstraView*SM REMOTE INSPECTION

Remotely inspect equipment before purchasing with our interactive website at www.instraview.com ↗

WE BUY USED EQUIPMENT

Sell your excess, underutilized, and idle used equipment. We also offer credit for buy-backs and trade-ins. www.artisanng.com/WeBuyEquipment ↗


LOOKING FOR MORE INFORMATION?

Visit us on the web at www.artisanng.com ↗ for more information on price quotations, drivers, technical specifications, manuals, and documentation

Contact us: (888) 88-SOURCE | sales@artisanng.com | www.artisanng.com

Outlook® 100 Safety Infusion System

Building a Foundation for Safety



Easy-to-use device helps to prevent IV medication errors. As the basis for our smart Extendable Platform Technology, it lets you plan for the future today.

B. Ready.

B | BRAUN
SHARING EXPERTISE



"It's important to me that an IV infusion system provides safeguards to reduce medication errors, while also recognizing the importance of my clinical judgment."


Prevent IV Medication Errors

- Promotes patient safety by helping to avoid accidental under- or overdosing using audible and visual alerts
- Provides flexibility for each care area—ability to customize dose limits to your hospital drug formulary
- Promotes clinician choice with dose limit override when limits are exceeded during manual programming and titration
- Visually indicates dose selection with DoseGuard™ software


This value is outside of the allowable dosage limits. Continue?
YES NO
Cursor ← → & Press Enter

The Proper Balance of Safety and Ease of Use

- Easy to learn
 - Manual operation RATE, VOLUME, RUN
 - Automatic Dose Rate Calculator
 - Large, highly visible keys
- Reduce nuisance air-in-line alarms with air-trapping cassette
- Stackable modular design reduces bedside cord clutter and errors due to channel confusion
- Comprehensive line of needle-free administration sets including UV light resistant tubing that protects drug efficacy



B. Sure.



Set-Based, Free-Flow protection meets JCAHO Goals

Smaller-sized decimals help prevent errors

Automatic Dose Rate Calculator for added patient safety

Dose Limit Indicator for visual confirmation and titration safety

Quickclamp™ Assembly for easy addition/removal of pumps

Unique Decimal Keys for easy identification

Laptop Style Battery Compartment for quick replacement


"It's important to me that I can access an alert log, and that safety features, like dose limits, are integral to the operation of the pump."



Outlook® 100 Safety Infusion System

SPECIFICATIONS

Mechanism:	Volumetric, Positive Pressure Displacement Reservoir
Width:	30.5 cm (12 in)
Height:	11.4 cm (4.5 in)
Depth:	22.2 cm (8.75 in)
Weight:	5 kg (11 lbs)
Power:	120 VAC, 3A, 60 Hz USA excluding accessories
Grounding Resistance:	Meets UL Standard 2601-1
Leakage current:	Meets or Exceeds UL Standard 2601-1
Battery:	12 V sealed lead-acid battery
Battery Operating Time:	3 hours at 125 mL/hr
Rate Range:	0.1 to 999.9 mL/hr in 0.1 mL increments Micro Mode: 0.1 to 99.9 mL/hr in 0.1 mL increments
Volume to be Infused Range:	0.1 to 9999.9 mL in 0.1 mL increments
Variable Occlusion Pressure:	75, 100, 200, 300, 400, 500, and 750 mmHg
KVO Rate:	3 mL/hr or selected rate whichever is lower
Volumetric Accuracy:	± 5%
Alarms:	Air-in-Line, Bag Empty, Door Open, Downstream Occlusion, Hold time Exceeded, Low Battery, Plug-in, System Error, Upstream Occlusion, Check Set, DoseGuard™ Alert
Memory:	Permanent data retention after the pump is turned off unless alternative Data Retention Defaults have been selected in the Biomed Options
Tubing:	Horizon® Pump Sets
Modes/Features:	Dose Mode with Drug Menu and DoseGuard™ dose limits, Primary and Piggyback Modes, Set Time and Volume, Schedule Next Run, Profile Mode for TPN, Program Mode, q(x) Mode, Site Trend, Variable Occlusion Pressure, Micro Mode, Epidural Mode



Promote Safe IV Medication Delivery.

Preventing medication errors begins with confirming the right patient with the right drug and the right dose along with the right route, the right time, and the right clinician.



Extend your Outlook on Safety.

Schedule a consultation with a
B. Braun Infusion Specialist by calling
1-800-227-2862.

B. Safe.

B | BRAUN
SHARING EXPERTISE

www.bbraunusa.com



Artisan Technology Group is your source for quality new and certified-used/pre-owned equipment

- FAST SHIPPING AND DELIVERY
- TENS OF THOUSANDS OF IN-STOCK ITEMS
- EQUIPMENT DEMOS
- HUNDREDS OF MANUFACTURERS SUPPORTED
- LEASING/MONTHLY RENTALS
- ITAR CERTIFIED SECURE ASSET SOLUTIONS

SERVICE CENTER REPAIRS

Experienced engineers and technicians on staff at our full-service, in-house repair center

*InstraView*SM REMOTE INSPECTION

Remotely inspect equipment before purchasing with our interactive website at www.instraview.com ↗

WE BUY USED EQUIPMENT

Sell your excess, underutilized, and idle used equipment. We also offer credit for buy-backs and trade-ins. www.artisanng.com/WeBuyEquipment ↗

LOOKING FOR MORE INFORMATION?

Visit us on the web at www.artisanng.com ↗ for more information on price quotations, drivers, technical specifications, manuals, and documentation

Contact us: (888) 88-SOURCE | sales@artisanng.com | www.artisanng.com