

Piramoon F155-8x50c  
**Fixed Angle Rotor**



**In Stock**

**New From Surplus Stock**

**Open Web Page**

<https://www.artisanng.com/66509-3>

All trademarks, brandnames, and brands appearing herein are the property of their respective owners.



Your **definitive** source  
for quality pre-owned  
equipment.

**Artisan Technology Group**

(217) 352-9330 | [sales@artisanng.com](mailto:sales@artisanng.com) | [artisanng.com](http://artisanng.com)

- Critical and expedited services
- In stock / Ready-to-ship

- We buy your excess, underutilized, and idle equipment
- Full-service, independent repair center

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.



**F21S-8x50y**



## Rotor Specifications

<b>Catalog No.</b>	096-084018
<b>Model</b>	F21S-8x50y
<b>Max. Speed</b>	21000
<b>Max. g-force</b>	52356
<b>k factor</b>	669
<b>Rotor Mass</b>	4.23 kg (9.34 lbs.)
<b>Liquid Containment</b>	Yes
<b>Will Replace</b>	SS-34
<b>Will Spin In</b>	RC2, RC5, RC6, RC26, Evolution RC
<b>Rmax</b>	106
<b>Rmin</b>	33
<b>Tube Angle</b>	34
<b>Application</b>	This rotor can be used for subcellular fractions, bacteria pelleting, lysate clarification, cell membranes, sucrose gradient isolation of mitochondria, endoplasmic reticulum, and gradient separation of cells with Percoll, Ficoll, Optiprep, or Nycodenz solutions.








## Replacement Parts

<b>Refresher Kit</b>	021-084018
<b>Extraction Tool</b>	Not Applicable
<b>Lid O-rings</b>	022-084018
<b>Lid Assembly</b>	099-084018

## Bottles

Part No.	Description
 3139-0050	50 mL Nalgene Oak Ridge Polypropylene Bottle with Sealing Cap.
 3138-0050	50 mL Nalgene Oak Ridge Polycarbonate Bottle with Sealing Cap.

## Adapters

	Adapter Part No.	Description	UOM
	010-1068	50 mL to 8.0/10 mL	pk/2
	010-1127	50 mL to 4x1.0 mL	pk/2
	010-0382	50 mL to 16 mL	pk/2
	010-1123	50 mL to 15 mL Conical	pk/2
	010-1128	50 mL to 3.0 mL	pk/2
	010-0167	50 mL to 30 mL	pk/2
	010-1306	50 mL to 10 mL	pk/2

# Artisan Technology Group is an independent supplier of quality pre-owned equipment

## Gold-standard solutions

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

## We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

## Learn more!

Visit us at [artisanTG.com](https://www.artisanTG.com) for more info on price quotes, drivers, technical specifications, manuals, and documentation.

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

**We're here to make your life easier. How can we help you today?**

(217) 352-9330 | [sales@artisanTG.com](mailto:sales@artisanTG.com) | [artisanTG.com](https://www.artisanTG.com)

