

ARTISAN®
TECHNOLOGY GROUP

Your **definitive** source for quality pre-owned equipment.

Artisan Technology Group
(217) 352-9330 | sales@artisantg.com | artisantg.com

Full-service, independent repair center
with experienced engineers and technicians on staff.

We buy your excess, underutilized, and idle equipment
along with credit for buybacks and trade-ins.

Custom engineering
so your equipment works exactly as you specify.

- Critical and expedited services
- In stock / Ready-to-ship
- Leasing / Rentals / Demos
- ITAR-certified secure asset solutions

Expert team | Trust guarantee | 100% satisfaction

All trademarks, brand names, and brands appearing herein are the property of their respective owners.

Visit our website - *Click HERE*



Chiral Columns

Normal Phase HPLC Columns

3μm - coating type | 5μm - coating type | 10μm/20μm - coating type | Feature

3μm - coating type (Normal phase)

Product name	Product type	Internal Diameter(mm)	Length (mm)	Particle size (μm)	P/N
CHIRALPAK® AY-3	Guard cartridge (×3)	4	10	3	47511
	Analytical column	4.6	50	3	47522
	Analytical column	4.6	150	3	47524
	Analytical column	4.6	250	3	47525
	Analytical column	2.1	150	3	47594
	Analytical column	2.1	250	3	47595
CHIRALPAK® AD-3	Guard cartridge (×3)	4	10	3	19511
	Analytical column	4.6	50	3	19522
	Analytical column	4.6	150	3	19524
	Analytical column	4.6	250	3	19525
	Analytical column	2.1	150	3	19594
	Analytical column	2.1	250	3	19595
CHIRALPAK® AS-3	Guard cartridge (×3)	4	10	3	20511
	Analytical column	4.6	50	3	20522

Cookieの使用について

当ウェブサイトでは、より快適にウェブサイトをご利用いただくため、Cookieを利用しております。当ウェブサイトをご覧いただく際は、Cookieの使用に同意ください。なお、このまま当ウェブサイトをご覧いただいた場合、Cookieの使用にご同意いただいたものとさせていただきます。

詳細はこちら

同意する

Product name	Product type	Internal Diameter(mm)	Length (mm)	Particle size (µm)	P/N
	Analytical column	2.1	250	3	20595
CHIRALCEL® OZ-3	Guard cartridge (×3)	4	10	3	42511
	Analytical column	4.6	50	3	42522
	Analytical column	4.6	150	3	42524
	Analytical column	4.6	250	3	42525
	Analytical column	2.1	150	3	42594
	Analytical column	2.1	250	3	42595
CHIRALCEL® OD-3	Guard cartridge (×3)	4	10	3	14511
	Analytical column	4.6	50	3	14522
	Analytical column	4.6	150	3	14524
	Analytical column	4.6	250	3	14525
	Analytical column	2.1	150	3	14594
	Analytical column	2.1	250	3	14595
CHIRALCEL® OJ-3	Guard cartridge (×3)	4	10	3	17511
	Analytical column	4.6	50	3	17522
	Analytical column	4.6	150	3	17524
	Analytical column	4.6	250	3	17525
	Analytical column	2.1	150	3	17594
	Analytical column	2.1	250	3	17595
	Guard cartridge holder	4	10	-	00011

5µm - coating type (Normal phase)

Product name	Product type	Internal Diameter(mm)	Length (mm)	Particle size (µm)	P/N
CHIRALPAK® AY-H	Guard cartridge (×3)	4	10	5	47311
	Analytical column	4.6	150	5	47324
	Analytical column	4.6	250	5	47325
	Guard column for semi-prep	10	20	5	47337
	Semi-prep	10	250	5	47335
	Semi-prep	20	250	5	47345
	Analytical column	2.1	150	5	47394
	Analytical column	2.1	250	5	47395
CHIRALPAK® AD-H	Guard cartridge (×3)	4	10	5	19311
	Analytical column	4.6	150	5	19324

Product name	Product type	Internal Diameter(mm)	Length (mm)	Particle size (μm)	P/N
	Analytical column	4.6	250	5	19325
	Guard column for semi-prep	10	20	5	19337
	Semi-prep	10	250	5	19335
	Semi-prep	20	250	5	19345
	Analytical column	2.1	150	5	19394
	Analytical column	2.1	250	5	19395
CHIRALPAK® AS-H	Guard cartridge (×3)	4	10	5	20311
	Analytical column	4.6	150	5	20324
	Analytical column	4.6	250	5	20325
	Guard column for semi-prep	10	20	5	20337
	Semi-prep	10	250	5	20335
	Semi-prep	20	250	5	20345
	Analytical column	2.1	150	5	20394
	Analytical column	2.1	250	5	20395
CHIRALPAK® AZ-H	Guard cartridge (×3)	4	10	5	61311
	Analytical column	4.6	150	5	61324
	Analytical column	4.6	250	5	61325
	Guard column for semi-prep	10	20	5	61337
	Semi-prep	10	250	5	61335
	Semi-prep	20	250	5	61345
	Analytical column	2.1	150	5	61394
	Analytical column	2.1	250	5	61395
CHIRALCEL® OZ-H	Guard cartridge (×3)	4	10	5	42311
	Analytical column	4.6	150	5	42324
	Analytical column	4.6	250	5	42325
	Guard column for semi-prep	10	20	5	42337
	Semi-prep	10	250	5	42335
	Semi-prep	20	250	5	42345
	Analytical column	2.1	150	5	42394
	Analytical column	2.1	250	5	42395
CHIRALCEL® OD-H	Guard cartridge (×3)	4	10	5	14311
	Analytical column	4.6	150	5	14324
	Analytical column	4.6	250	5	14325
	Guard column for semi-prep	10	20	5	14337

Product name	Product type	Internal Diameter(mm)	Length (mm)	Particle size (μm)	P/N
	Semi-prep	10	250	5	14335
	Semi-prep	20	250	5	14345
	Analytical column	2.1	150	5	14394
	Analytical column	2.1	250	5	14395
CHIRALCEL® OJ-H	Guard cartridge (×3)	4	10	5	17311
	Analytical column	4.6	150	5	17324
	Analytical column	4.6	250	5	17325
	Guard column for semi-prep	10	20	5	17337
	Semi-prep	10	250	5	17335
	Semi-prep	20	250	5	17345
	Analytical column	2.1	150	5	17394
	Analytical column	2.1	250	5	17395
CHIRALCEL® OB-H	Guard cartridge (×3)	4	10	5	12311
	Analytical column	4.6	150	5	12324
	Analytical column	4.6	250	5	12325
	Analytical column	2.1	150	5	12394
	Analytical column	2.1	250	5	12395
CHIRALCEL® OC-H	Guard cartridge (×3)	4	10	5	13311
	Analytical column	4.6	150	5	13324
	Analytical column	4.6	250	5	13325
	Analytical column	2.1	150	5	13394
	Analytical column	2.1	250	5	13395
CHIRALCEL® OX-H	Guard cartridge (×3)	4	10	5	63311
	Analytical column	4.6	150	5	63324
	Analytical column	4.6	250	5	63325
	Guard column for semi-prep	10	20	5	63337
	Semi-prep	10	250	5	63335
	Semi-prep	20	250	5	63345
	Analytical column	2.1	150	5	63394
	Analytical column	2.1	250	5	63395
CHIRALPAK® AZ-H	Guard cartridge (×3)	4	10	5	61311
	Analytical column	4.6	150	5	61324
	Analytical column	4.6	250	5	61325
	Guard column for semi-prep	10	20	5	61337

Product name	Product type	Internal Diameter(mm)	Length (mm)	Particle size (µm)	P/N
	Semi-prep	10	250	5	61335
	Semi-prep	20	250	5	61345
	Analytical column	2.1	150	5	61394
	Analytical column	2.1	250	5	61395
	Guard cartridge holder	4	10	-	00011

10µm/20µm - coating type (Normal phase)

Product name	Product type	Internal Diameter(mm)	Length (mm)	Particle size (µm)	P/N
CHIRALPAK® AD	Guard column for analysis	4.6	50	10	19022
	Analytical column	4.6	250	10	19025
	Guard column for semi-prep	10	50	10	19032
	Semi-prep	10	250	10	19035
	Guard column for semi-prep	20	50	10	19042
	Semi-prep	20	250	10	19045
	Guard column for preparative	50	100	20	19253
	Preparative column	50	500	20	19256
CHIRALPAK® AS	Guard column for analysis	4.6	50	10	20022
	Analytical column	4.6	250	10	20025
	Guard column for semi-prep	10	50	10	20032
	Semi-prep	10	250	10	20035
	Guard column for semi-prep	20	50	10	20042
	Semi-prep	20	250	10	20045
	Guard column for preparative	50	100	20	20253
	Preparative column	50	500	20	20256
CHIRALCEL® OA	Guard column for analysis	4.6	50	10	11022
	Analytical column	4.6	250	10	11025
CHIRALCEL® OB	Guard column for analysis	4.6	50	10	12022
	Analytical column	4.6	250	10	12025
CHIRALCEL® OC	Guard column for analysis	4.6	50	10	13022
	Analytical column	4.6	250	10	13025
CHIRALCEL® OD	Guard column for analysis	4.6	50	10	14022
	Analytical column	4.6	250	10	14025
	Guard column for semi-prep	10	50	10	14032

Product name	Product type	Internal Diameter(mm)	Length (mm)	Particle size (μm)	P/N
	Semi-prep	10	250	10	14035
	Guard column for semi-prep	20	50	10	14042
	Semi-prep	20	250	10	14045
	Guard column for preparative	50	100	20	14253
	Preparative column	50	500	20	14256
CHIRALCEL® OF	Guard column for analysis	4.6	50	10	15022
	Analytical column	4.6	250	10	15025
	Guard column for semi-prep	10	50	10	15032
	Semi-prep	10	250	10	15035
	Guard column for semi-prep	20	50	10	15042
	Semi-prep	20	250	10	15045
CHIRALCEL® OG	Guard column for analysis	4.6	50	10	16022
	Analytical column	4.6	250	10	16025
CHIRALCEL® OJ	Guard column for analysis	4.6	50	10	17022
	Analytical column	4.6	250	10	17025
	Guard column for semi-prep	10	50	10	17032
	Semi-prep	10	250	10	17035
	Guard column for semi-prep	20	50	10	17042
	Semi-prep	20	250	10	17045
	Guard column for preparative	50	100	20	17253
	Preparative column	50	500	20	17256
CHIRALCEL® OK	Guard column for analysis	4.6	50	10	18022
	Analytical column	4.6	250	10	18025

Feature

◯Polyssacharide-based chiral stationary phases (CSPs) are most widely used for chiral separation in the world. More than 80% of enantiomers can be separated with our polysaccharide-based chiral columns.
◯We prepare various size of columns from Semi-micro to Preparative type.
◯We provide 3 micron series of various type of Normal Phase chiral columns.

Product list

Immobilized type →
Normal Phase →
Reversed-Phase →
AGP/CBH/HSA →
CROWNPAK →
Qseries →
Other →
Accessories →
SFC →







Worldwide Network READ MORE →

Chiral Technologies Worldwide serves our customers with full technical support through our Group companies and distributors at your area.
Please visit the WEB page of our Group companies and distributor at your area.

CHIRAL TECHNOLOGIES INC
CHIRAL TECHNOLOGIES EUROPE
DAICEL CHIRAL TECHNOLOGIES (China)
DAICEL CHIRAL TECHNOLOGIES (India)

© 2021 Daicel Corporation. ALL RIGHTS RESERVED.

Artisan Technology Group is an independent supplier of quality pre-owned equipment

Gold-standard solutions

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

Learn more!

Visit us at **artisantg.com** for more info on price quotes, drivers, technical specifications, manuals, and documentation.

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

We're here to make your life easier. How can we help you today?

(217) 352-9330 | sales@artisantg.com | artisantg.com

ARTISAN® TECHNOLOGY GROUP