

ARTISAN®
TECHNOLOGY GROUP

Your **definitive** source for quality pre-owned equipment.

Artisan Technology Group
(217) 352-9330 | sales@artisantg.com | artisantg.com

Full-service, independent repair center
with experienced engineers and technicians on staff.

We buy your excess, underutilized, and idle equipment
along with credit for buybacks and trade-ins.

Custom engineering
so your equipment works exactly as you specify.

- Critical and expedited services
- In stock / Ready-to-ship
- Leasing / Rentals / Demos
- ITAR-certified secure asset solutions

Expert team | Trust guarantee | 100% satisfaction

All trademarks, brand names, and brands appearing herein are the property of their respective owners.

Visit our website - *Click HERE*


Specifications

Carl Zeiss, Inc.
One Zeiss Drive
Thornwood, NY 10594
(914) 747-1800

ZEISS

West Germany

Photomicroscope III



Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Contents of this brochure

Photomicroscope, basic equipment . . . . . Page 3

Accessories for film labeling . . . . . Page 3

Illuminating systems for transmitted or reflected light . . . . . Page 3

Optional bulbs, sockets and power supplies for 100W lamp housing . . . . . Page 3

Optical equipment for transmitted light . . . . . Page 4

Accessories for transmitted light systems . . . . . Page 4

Accessories for fluorescence microscopy . . . . . Page 5

Optical equipment for reflected light . . . . . Page 5, 6

Accessories for vertical illuminators . . . . . Page 6

Body tubes . . . . . Page 6

Accessories for TV . . . . . Page 6

Specimen stages . . . . . Page 7,8

Multiple beam interference equipment . . . . . Page 9

Mirau interference equipment . . . . . Page 9


Projection equipment . . . . . Page 10

Microhardness tester . . . . . Page 10


Attachment camera MC 63 . . . . . Page 11

Computer flash . . . . . Page 11


Separate information is available on the following categories.




Camera Systems
Contents:
Attachment cameras
Adapters for direct attachment of SLR cameras
Electronic Flash




Light Filters
Contents:
Contrast filters
Interference filters
Heat filters
Neutral filters



Microscope Illuminators
Contents:
Illuminators for transmitted light
Illuminators for reflected light
Replacement bulbs and burners



Optical Systems for Reflected Light
Contents:
Objectives
Eyepieces
Micrometer discs and other measuring accessories
Interference contrast equipment




Optical Systems for Transmitted Light
Contents:
Objectives
Eyepieces
Micrometer discs and other measuring accessories
Interference contrast equipment

(Continued on page 3)

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Basic equipment




	47 21 90	Photomicroscope stand with tube head and Optovar, built-in set of gray filters, 1 conversion filter and green filter VG9
	46 70 57-9904	Diaphragm insert for microscope base
	47 20 83	Power supply 100-240V,50/60 Hz, 190VA
	47 90 71	Housing for power supply
	38 00 74-0160	Power cable, 110V
	47 30 12-9902	Inclined binocular body 45°
(2)	46 40 48-9902	Widefield focusing eyepieces KPL 10x/20
	47 20 26-9901	Film cassette with counter
	47 15 40	Clamp-on stage carrier
	47 34 15	Square mechanical stage with graduations and low-position coaxial controls for right hand use
	47 34 48	Specimen holder for stage
	47 93 01	Plastic dust cover
		Photomicroscope III, basic equipment

Accessories for film labeling

47 20 27-9901	Film casette for film labeling
47 20 91	Data recording strips (25 per pkg.)
47 20 92	Sheet with image scale for 47 20 27-9901

Illuminating systems for transmitted or reflected light



	46 72 57	Lamp housing 12V, 60W with collector
	46 80 15	Lamp socket 60/1
	46 70 40-9903	Connecting tube with relay lens
(3)	38 00 18-2520	Bulbs 12V, 60W
		12V, 60W illuminator complete



	46 72 59-9901	Lamp housing 100 without collector
	46 72 73-9901	Two lens collector 100
	46 80 19-9901	Lamp socket 100
	46 70 40-9903	Connecting tube with relay lens
(2)	38 00 59-1660	Halogen bulbs 12V, 100W
		12V, 100W halogen illuminator complete

Optional bulbs, sockets and power supplies for 100W lamp housing

38 16 19	Mercury burner HBO 50
46 80 32-9902	Lamp socket for mercury burner
46 72 74-9902	Three-lens collector 100
91 00 70	Power supply for mercury burner

91 01 43	Mercury burner HBO 50 W/L1 DC
91 01 38	Lamp socket for HBO 50 DC
46 72 74-9902	Three-lens collector 100
91 01 35	Power supply for HBO 50 DC

38 00 18-4060	Mercury burner HBO 100W
46 80 30-9999	Lamp socket for HBO 100W
46 72 74-9902	Three-lens collector 100
91 02 35	DC power supply for HBO 100 with built-in fan and high voltage ignition

91 01 42	Xenon burner XBO 75 W2/DC
46 80 31-9999	Lamp socket for XBO 75 DC
46 72 76	Four-lens collector
91 02 34	Power supply for XBO 75 DC, stabilized


Cooling attachment for lamphousing 100W

46 72 95	Fan attachment for lamphousing 100W
39 25 81-9002	Fan transformer
	Fan attachment complete

(Continued on page 4)

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Optical equipment for transmitted light



47 15 58-9901	Attachable condenser carrier
47 08 61-9901	Centerable carrier with auxiliary lens
47 31 59	Revolving nosepiece 5x with Telan lens
46 52 67	Achromatic-aplanatic condenser 0.32-1.4 Pol
46 52 58	Condenser head 1.4
46 01 10-9906	Planachromat 2.5/0.08
46 03 10	Planachromat 6.3/0.16
46 05 10	Planachromat 16/0.35
46 07 10	Planachromat 40/0.65
46 19 10	Planachromat 100/1.25 oil

Complete optical equipment for brightfield with planachromats

Photomicroscope III with 12V, 60W illuminator and optical equipment for brightfield with planachromats

47 15 58-9901	Attachable condenser carrier
47 08 61-9901	Centerable carrier with auxiliary lens
47 31 59	Revolving nosepiece 5x with Telan lens
46 52 77-9901	Achromatic-aplanatic phase/darkfield condenser VZ (1.4)
46 29 29	Phase contrast mounting medium
46 03 10	Planachromat 6.3/0.16 (for brightfield only)
46 05 21	Neofluar 16/0.40 Ph
46 06 21	Neofluar 25/0.60 Ph
46 07 21	Neofluar 40/0.75 Ph
46 19 21-9903	Neofluar 100/1.30 Ph oil

Complete optical equipment for brightfield and phase contrast with neofluars

Photomicroscope III with 12V, 60W illuminator and optical equipment for brightfield and phase contrast with neofluars

47 15 58-9901	Attachable condenser carrier
47 08 65-9902	Centerable carrier with auxiliary lens and rotatable polarizer
47 31 56	Revolving nosepiece 5x with intermediate rings for DIC
46 52 85	Achromatic-aplanatic phase interference contrast condenser VZ (1.4)
47 36 63-9903	Simple analyzer slider with movable quartz white plate
47 37 04	Quartz plate red 1
46 03 10	Planachromat 6.3/0.16 (not for DIC)
46 05 10	Planachromat 16/0.35
47 45 51	DIC slider for Planachromat 16
46 07 10	Planachromat 40/0.65
47 45 71	DIC slider for Planachromat 40
46 19 10	Planachromat 100/1.25 oil
47 45 91	DIC slider for Planachromat 100

Complete optical equipment for interference contrast

Photomicroscope III with 12V, 100W halogen illuminator and optical equipment for brightfield and interference contrast

Accessories for transmitted light systems

Nosepieces

47 31 69	Revolving nosepiece 7x with Telan lens
47 31 55	Revolving nosepiece 5x for widefield observation with Telan lens 0.63

(Continued on page 5)


Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Accessories for fluorescence microscopy



46 62 49-9904	Vertical illuminator IIIRS for Epi excitation without dichroic reflectors and filters
48 77 02	Exciter-barrier filter and reflector combination, UV exciter (Catecholamine)
48 77 05	Exciter-barrier filter and reflector combination, Violet exciter (Quinacrine Mustard)
48 77 16	Exciter-barrier filter and reflector combination, Blue exciter (FITC)
48 77 14	Exciter-barrier filter and reflector combination, Green exciter (Rhodamine)
46 04 20	Neofluar 10/0.30
46 15 26	Plan-Neofluar 16/0.50 water/glycerin Ph2
46 16 26	Plan-Neofluar 25/0.80 oil/water/glycerin Ph2 with correction mount
46 17 26	Plan-Neofluar 40/0.9 oil/water/glycerin Ph3
46 18 36	Plan-Neofluar 63/1.25 oil
46 70 34	Incident insert for lateral illumination
46 72 59-9901	Lamp housing 100
46 72 74-9902	Three lens collector 100
46 80 32-9902	Lamp socket for HBO 50
(2) 38 16 19	Mercury burners HBO 50
91 00 70	Power supply for mercury burner
	Complete optical equipment for EPI-flourescence microscopy

Optical equipment for reflected light




46 62 44	Vertical illuminator IIIC/45 with quadruple revolving nosepiece
46 62 60-9903	Plane glass reflector HPL
46 62 62-9903	Darkfield reflector D
46 03 69-9904	Epiplan 8/0.20 HD
46 05 69-9906	Epiplan 16/0.35 HD
46 07 69	Epiplan 40/0.85 HD
46 08 69-9901	Epiplan 80/0.95 HD
(2) 46 72 52	Filter holders for 32 mm dia. filters
47 20 75	Insert with aperture diaphragm for vertical illumination
(4) 46 29 89-9999	Intermediate rings 12mm
	Complete optical equipment for brightfield and darkfield
	Photomicroscope III with 12V, 60W illuminator and optical equipment for brightfield and darkfield

46 62 44	Vertical illuminator IIIC/45 with quadruple revolving nosepiece
46 62 65	Plane glass reflector H-PI Pol
46 62 62-9903	Darkfield reflector D
(2) 46 72 52	Filter holders for 32 mm dia. filters
47 20 75	Insert for vertical illumination with aperture diaphragm
(2) 46 29 89-9999	Intermediate rings 12 mm
46 03 69-9904	Epiplan 8/0.20 HD
46 21 23	LD-Epiplan 16/0.30 Pol
46 29 15	Protection cap
47 44 63-9902	Interference attachment for LD-Epiplan 16
46 29 96	Intermediate ring W0.8 x 1/36 on 24 x 0.75
46 21 24-9901	LD-Epiplan 40/0.6 Pol
46 29 16	Protection cap
47 44 64-9902	Interference attachment for LD-Epiplan 40
46 29 96	Intermediate ring W0.8 x 1/36 on 24 x 0.75
46 08 69-9901	Epiplan 80/0.95 HD
47 36 16	Epi polarizer with scale
47 36 63-9903	Simple analyzer slider with movable quartz white plate
	Complete optical equipment for interference contrast after Nomarski (CNRS) brightfield and darkfield
	Photomicroscope III with 12V, 100W halogen illuminator and optical equipment for brightfield, darkfield and Nomarski interference contrast


(Continued on page 6)

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Optional vertical illuminators



	46 62 37	Vertical illuminator IIC
	46 62 61	Prism reflector HPr Pol
	46 62 60-9903	Plane glass reflector HPL
(2)	46 72 52	Filter holders for 32 mm dia. filters
	46 78 50-9901	Conversion filter CB 12
	46 78 26	Neutral filter NG9
	47 36 16	Epi polarizer with scale
	47 36 62-9903	Rotatable analyzer 360°, with movable quartz white plate
	47 20 71	Insert with aperture diaphragm for vertical illumination
	46 20 01	Epiplan 4/0.10 Pol
	46 20 02	Epiplan 8/0.20 Pol
	46 20 03	Epiplan 16/0.35 Pol
	46 20 04	Epiplan 40/0.85 Pol
	46 20 05-9903	Epiplan 100/1.25 Pol
(5)	46 62 56	Centering change rings
		Complete optical equipment for reflected polarized light



	46 62 46	Vertical illuminator IIIC/33 with quadruple revolving nosepiece
	46 62 60-9903	Plane glass reflector HPL
	46 62 62-9903	Darkfield reflector D
(2)	46 72 52	Filter holders for 32 mm dia. filters
	47 20 71	Insert with aperture diaphragm for vertical illumination
		Vertical illuminator IIIC/33 for brightfield and darkfield without HD objectives
	46 62 45	Vertical illuminator IIIC/33 for widefield with quadruple revolving nosepiece
	46 62 60-9903	Plane glass reflector HPL
	46 62 62-9903	Darkfield reflector D
(2)	46 72 52	Filter holders for 32 mm dia. filters
	47 20 71	Insert for vertical illumination with aperture diaphragm
		Vertical illuminator IIIC/33 for widefield

Accessories for vertical illuminators


For vertical illuminator IIC, IIIC/33 and IIIC/45

46 62 60-9903	Plane glass reflector HPL
46 62 62-9903	Darkfield reflector D
46 62 63	Reflector HD including swing-in annular diaphragm
46 62 65	Plane glass reflector HPL-Pol

Body tubes

47 30 00	Inclined monocular tube, 25 mm dia.
47 30 34	Inclined binocular tube Pol

Accessories for TV




47 79 01	TV tube, long
47 79 21	C-mount adapter for TV cameras

(Continued on page 7)

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Specimen stages




47 15 40 Clamp-on stage carrier .....



Square stage

47 34 20-9904 Square stage with 2 stage clips 128 x 143 mm (Simplest form of specimen support) .....

41 33 02 Attachable mechanical stage 25/75 with graduation and low position controls (For 47 34 20 only) .....



Plain mechanical stage

41 34 22-9901 Simple mechanical stage with graduation travel range 25 x 75mm .....




Large mechanical stage

47 34 15 Square mechanical stage with graduations 160 x 140 mm, low coaxial drive for specimen movement, travel range 50 x 75 mm for right hand operation .....

47 34 16 Square mechanical stag same as above for left hand operation .....

47 34 48 Specimen holder .....



Mechanical stage with coaxial controls both sides

47 34 28 Square mechanical stage, graduated 137 x 136mm travel range 50 x 75mm .....



Rotating autolevelling circular glide stage

47 33 19-9904 Autolevelling stage for polished specimens .....

47 33 00 Centering device for stage .....



Rotating and centering circular glide stage

47 34 54 Diameter 150 mm. Displacement range all around 16mm from center (32mm). Swiftly moving specimens can be easily followed. The stage can be centered and rotation can be clamped .....



Rotating and centering mechanical stage

47 34 56 Diameter 140mm, range 50 x 75mm ungraduated .....

47 34 57 Same as above but graduated .....



Simple Pol rotary stage

47 34 64 Diameter 115m. With 2 nonii on the graduation for read-off of rotary angle to 0.1° exactly .....



Ball bearing Pol rotary stage

47 34 66 Diameter 150mm. Read-off of rotary angle to 0.1 exactly. Stage can be clamped in any position. 45 click-stop for marking initial azimuth. Drive mechanism for fine adjustment of rotary angle .....




Attachable Pol mechanical stage

47 33 25-9903 Range 30 x 40mm .....

(Continued on page 8)

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com



Screening stages with motorized drive

137 x 136mm
Range 50 x 75mm
Drive rate from 0 to max. 450µm/sec continuously
Choice of speed and screening direction manually or with foot switch

47 34 89	Screening stage with motorized drive in both X and Y direction.
47 74 73	Control panel for X and Y drive
38 00 74-0160	Power cable 110V

Optional

47 74 70	Foot control pedal for 47 74 73

Heating and cooling devices

Temperature range — 30°C ± 70°C


47 80 13	Heating and cooling chamber

or for heating only:

91 02 24	Power supply 6V/12V, 125W with power cable

Heating and cooling chamber, Cat. #47 80 13, can be attached to the microscope stage in place of the attachable mechanical stage with:

47 80 76	Adapter for mechanical stages 47 34 23/24
47 80 77	Adapter for mechanical stages 47 34 25/26



47 80 15	Heating and cooling stage

or for heating only:

91 02 24	Power supply 6V/12V, 125W with power cable

Heating and cooling stage, Cat. #47 80 15, can be directly mounted in centering device, Cat. #47 33 00 on clamp-on stage carrier, Cat. #47 15 40




Optional

47 80 50	Nitrogen chamber Z for Cat. #47 80 13
47 80 51	Nitrogen chamber T for Cat. #47 80 15

Both nitrogen chambers can be used to avoid condensation at cooling temperatures.

Scanning stages



47 34 81	Scanning stage, 0.5 µm
47 34 41	Specimen holder, right with vernier
* 47 34 42	Specimen holder, left

Scanning stage 0.5 µm, 25 × 75mm travel

47 34 82	Scanning stage, 10µm
47 34 41	Specimen holder, right with vernier
* 47 34 42	Specimen holder, left

Scanning stage 10 µm, 25 × 75mm travel

47 34 86	Scanning stage 100 µm
47 34 41	Specimen holder, right with vernier
* 47 34 42	Specimen holder, left

Scanning stage 100 µm 25 × 75mm travel

*** Optional**

47 34 43	Specimen holder, left with spring finger

For control of stage

Requires computer and interfacing components for operation. Consult your Zeiss instrument specialist for details.

(Continued on page 9)

Multiple beam interference equipment




47 44 20	Multiple beam interference equipment in case
47 33 19-9904	Autoleveling stage for polished specimens
47 33 00	Centering device for stage
47 15 40	Clamp-on stage carrier
47 44 40	Reference plate R = 0.04
47 44 41	Reference plate R = 0.15
47 44 42	Reference plate R = 0.35
47 44 43	Reference plate R = 0.80
46 78 08	Interference line filter PIL 546 nm ±2 nm
	Multiple beam interference equipment (For IIC Vertical Illuminator)
46 70 42	Connecting tube
46 72 59-9901	Lamp housing 100W without collector
46 72 74-9902	Three lens collector 100
46 80 32-9902	Lamp socket for HBO 50
(2) 38 16 19	Mercury burners HBO 50
91 00 70	Power supply for mercury burner
	Mercury burner HBO 50 for reflected light
	Multiple beam interference equipment complete with mercury burner HBO 50


Also required for multiple beam interference attachment:

47 44 63-9902	Interference attachment for LD 16
46 21 23	LD-Epiplan 16/0.30 Pol
47 44 64-9902	Interference attachment for LD 40
46 21 24-9901	LD-Epiplan 40/0.60 Pol
(2) 46 62 58	Fixed change rings W 0.8"

Mirau interference equipment




47 44 25	Double beam interference head
47 44 26	Interference attachment for LD-Epiplan 16
47 44 27	Interference attachment for LD-Epiplan 40
46 78 08	Interference line filter PIL 546nm ± 2nm
46 21 23	LD-Epiplan 16/0.30 Pol
46 21 24-9901	LD-Epiplan 40/0.60 Pol



Intermediate rings


For Vertical Illuminator IIIC/45:

46 29 88	Intermediate ring, extension 23mm, I/E W0.8"
46 29 96	Intermediate ring, extension 0mm, I W0.8"/EM24



For Vertical Illuminator IIIC/33:

46 29 87	Centerable intermediate ring, extension 12mm, I/E W0.8"
46 29 96	Intermediate ring, extension 0mm, IW0.8"/EM24




For Vertical Illuminator IIC:

46 62 58	Fixed change ring W0.8"
46 29 87	Centerable intermediate ring, extension 12mm, I/E W0.8"

(Continued on page 10)

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

For projection on 6″ screen and 4 x 5 photography



47 30 85 47 30 81

47 30 85	Projection attachment for Universal and Photomicroscope
47 30 81	Screen insert
	Projection attachment complete
47 30 85	Projection attachment for Universal and Photomicroscope
47 30 87	Glarex projection screen
39 25 72	Transformer for Glarex
38 00 74-0160	Power cable 110V
	High resolution projection attachment
47 30 20	Straight monocular tube 25 mm dia.
47 78 50	Deflecting prism
46 49 10	Eyepiece clamping ring 25 mm dia.
46 35 78	Projection eyepiece f = 80 mm
	Projection prism for microprojection eyepieces





Microhardness tester

47 40 41	Microhardness Tester IB including galvanometer unit, monocular draw tube, measuring eyepiece, polarizer, analyzer, interference contrast slider, vertical illuminator, Bertrand lens, object micrometer, balance with weight and tweezers, mirror, stage with centering device and stage carrier ultra-fine adjustment and coupling piece, cover plate, case
	Required in addition:
47 33 25-9903	Attachable mechanical stage Pol with graduation
47 40 45	Vickers pyramid with Neofluar 63/0.9
	Microhardness Tester 1B, complete
Optional	
47 40 46	Knoop pyramid with Neofluar 63/0.9

(Continued on page 11)

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

47 60 72-9901



47 74 25-9903

47 60 12-9901

Attachment camera MC 63

47 60 12-9901	Special body with shutter M
47 60 13-9901	Photo diode sensor ME
47 74 25-9903	Exposure control panel
38 00 74-0160	Power cable 110V
47 60 21	Photo reticule MC for visual eyepiece

Attachment camera MC63 without camera back



47 74 20-9901

Optional control panel

47 74 20-9901 Exposure control panel for MC63 A



47 60 72-9901 and 47 60 79-9901



47 60 57

Camera backs


47 60 72-9901	Camera back 35mm for MC 63
47 60 74	Camera back M35F for micro flash
47 60 79-9901	Motor drive with cable for 47 60 72-9901
47 60 57	Camera back 4" x 5"

Photoeyepieces


46 40 40	Photoeyepiece S-Kpl 10x/20 for Plan optics only
46 42 21	Photoeyepiece Kpl 16x/12M

Phototube

47 30 24	Phototube M for Universal and Photomicroscope stands



Photomicroscope III



Computer flash

46 70 45	Connecting tube
46 72 50	Lamphousing with collector (frosted)
46 80 10-9904	Lamp socket 6V, 15W
38 00 18-1740	Bulb 6V, 15W
46 70 20	Double collector for flash
46 80 46	Flash bulb in slider
39 29 03-9901	Flash generator (110-220V)
38 00 78-4060	Synchronizing cable (included in 39 29 03-9901)
39 79 02-8003	Sensor cable
38 00 74-0160	Power cable, 110V

Computer flash for Photomicroscope III complete

For lamphousing 60 or 100:

42 70 09	Dovetail for flash collector
46 70 42	Connecting tube for 12V illuminators

SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE

Artisan Technology Group - Quality Instrumentation ... Guaranteed | (888) 88-SOURCE | www.artisantg.com

Artisan Technology Group is an independent supplier of quality pre-owned equipment

Gold-standard solutions

Extend the life of your critical industrial, commercial, and military systems with our superior service and support.

We buy equipment

Planning to upgrade your current equipment? Have surplus equipment taking up shelf space? We'll give it a new home.

Learn more!

Visit us at **artisantg.com** for more info on price quotes, drivers, technical specifications, manuals, and documentation.

Artisan Scientific Corporation dba Artisan Technology Group is not an affiliate, representative, or authorized distributor for any manufacturer listed herein.

We're here to make your life easier. How can we help you today?

(217) 352-9330 | sales@artisantg.com | artisantg.com

ARTISAN® TECHNOLOGY GROUP